


**Florida Department of Agriculture and Consumer Services  
Division of Agricultural Environmental Services**



**Fertilizer Laboratory Methods Manual**

**Bureau of Agricultural Environmental Laboratories (BAEL)**

**Rev 01/31/25**

**Rule 5E-1.014, F.A.C.**

# Table of Contents

| Fertilizer Method                         | BAEL Method              | Document Title                                                                             | Page Numbers |
|-------------------------------------------|--------------------------|--------------------------------------------------------------------------------------------|--------------|
| N/A                                       | PREP 450                 | Riffling and Splitting of Fertilizer Samples                                               | 1-5          |
| N/A                                       | PREP 451                 | Grinding of Fertilizer Samples                                                             | 6-10         |
| N/A                                       | PREP 452                 | Dilution of Fertilizer Samples                                                             | 11-12        |
| N/A                                       | PREP 453                 | Weighing of Fertilizer Samples                                                             | 13-16        |
| N/A                                       | PREP 453<br>Attachment A | Fertilizer Laboratory Standard and Batch Size                                              | 17           |
| FM604                                     | PREP 454                 | Sample Preparation for Water Soluble Potash and Magnesium Analysis                         | 18-20        |
| FM803                                     | PREP 455                 | Sample Preparation for Water Soluble Boron Analysis                                        | 21-23        |
| FM822<br>FM825<br>FM860<br>FM875<br>FM902 | PREP 456                 | Sample Preparation for Total Micronutrients and Heavy Metal Analysis                       | 24-26        |
| FM505                                     | PREP 458                 | Sample Preparation for Available Phosphorus Analysis                                       | 27-31        |
| FM521                                     | PREP 459                 | Sample Preparation for Total Phosphorus Analysis                                           | 32-36        |
| N/A                                       | PREP 461                 | Sample Preparation for Available Phosphorus and Soluble Potash in Fertilizer               | 37-40        |
| FM 830<br>FM 831<br>FM832<br>FM833        | PREP 462                 | Preparation for Soluble Cu, Fe, Mn, and Zn                                                 | 41-46        |
| FM 840<br>FM 841                          | PREP 464                 | Sample preparation for Chelated Fe, Mn, Cu, Zn                                             | 47-50        |
| FM430                                     | PREP 466                 | Sample Preparation for Nitrate, Urea and Ammoniacal Analysis using Flow Injection Analyzer | 51-54        |
| FM425<br>FM430                            | METHOD 500               | Sample Preparation and Analysis of Direct Water Soluble N                                  | 54-63        |
| FM445                                     | METHOD 501               | Sample Preparation and Analysis of Direct Insoluble N                                      | 64-71        |
| FM411                                     | METHOD 503               | Analysis of Total Nitrogen using LECO Combustion Analyzer                                  | 72-76        |
| FM420                                     | METHOD 504               | Sample Preparation and Analysis of Direct Ammoniacal N with Urea Absent                    | 77-83        |
| FM415                                     | METHOD 505               | Sample Preparation and Analysis of Direct Ammoniacal N with Urea Present                   | 84-90        |
| FM435                                     | METHOD 506               | Sample Preparation and Analysis of Direct Ammoniacal N with Urea and Insoluble N Present   | 91-97        |
| FM822<br>FM825<br>FM860<br>FM875<br>FM902 | METHOD 509               | Analysis of Total Metals and Heavy Metals using ICP-OES                                    | 98-106       |

# Table of Contents

| Fertilizer Method                                                                            | BAEL Method                | Document Title                                                                                                                                    | Page Numbers |
|----------------------------------------------------------------------------------------------|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| <b>FM611</b>                                                                                 | METHOD 510                 | Sample Preparation and Analysis of Chloride in Fertilizer                                                                                         | 107-111      |
| <b>FM604</b><br><b>FM803</b>                                                                 | METHOD 513                 | Analysis of Water Soluble Potash, Water Soluble Magnesium, and Water Soluble Boron using Inductively Coupled Plasma Optical Emission Spectrometry | 112-119      |
| <b>N/A</b>                                                                                   | METHOD 514                 | Fertilizer Preparation and Analysis for Potash by STPB using an Autotitrator                                                                      | 120-127      |
| <b>FM452</b><br><b>FM462</b><br><b>FM472</b>                                                 | METHOD 515                 | Analysis of NAU using Flow Injection Analyzer                                                                                                     | 128-134      |
| <b>N/A</b>                                                                                   | METHOD 515<br>Attachment A | Nitrate Manifold Diagram                                                                                                                          | 135          |
| <b>N/A</b>                                                                                   | METHOD 515<br>Attachment B | Ammonia Manifold Diagram                                                                                                                          | 136          |
| <b>N/A</b>                                                                                   | METHOD 515<br>Attachment C | Urea/Ammonia Manifold Diagram                                                                                                                     | 137-138      |
| <b>N/A</b>                                                                                   | METHOD 516                 | Analysis of Available Phosphorous and Soluble Potash in Fertilizers using ICP-OES                                                                 | 139-146      |
| <b>FM805</b>                                                                                 | METHOD 520                 | Sample Preparation and Analysis of Combined Sulfur                                                                                                | 147-151      |
| <b>FM806</b>                                                                                 | METHOD 521                 | Sample Preparation and Analysis of Free Sulfur                                                                                                    | 152-157      |
|                                                                                              | METHOD 521<br>Attachment A | Soxhlet Apparatus Assembly                                                                                                                        | 158-160      |
| <b>FM830</b><br><b>FM831</b><br><b>FM832</b><br><b>FM833</b><br><b>FM840</b><br><b>FM841</b> | METHOD 523                 | Analysis of Soluble and Chelate Metals by AA                                                                                                      | 161-166      |
| <b>FM415</b><br><b>FM420</b><br><b>FM435</b>                                                 | METHOD 530                 | Automated Direct Distillation of Ammoniacal N in Fertilizer                                                                                       | 167-175      |
| <b>FM405</b><br><b>FM425</b><br><b>FM440</b><br><b>FM445</b>                                 | METHOD 531                 | Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer                                                                                   | 176-190      |

Note: The page numbers listed here are for navigation within this compilation and do not correspond to the internal pagination of each individual SOP document.

# ***STANDARD OPERATING PROCEDURE*** **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

## **PREP 450** ***Riffling and Splitting of Fertilizer Sample***

Version: 3.0

Effective Date: 01/27/2021

Page 1 of 5

### **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for creating a homogenous portion of the dry fertilizer sample. This instruction allows for a portion that is representative of the entire sample but still easy to handle.

### **2. SCOPE**

2.1. The procedure delineated below is applicable to Fertilizer samples prepared in the Bureau of Agricultural Environmental Laboratories (BAEL). This procedure applies to the converting of solid fertilizer into uniformly consistent analytical samples while preserving the integrity of the original material. The laboratory prepares three portions that represent the original sample as a whole.

### **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves, and laboratory coat is required.

### **4. OUTLINE of PROCEDURE**

- 4.1. 9.1 Sample Type Determination
- 4.2. 9.2 Riffling "Normal" Sample
- 4.3. 9.3 Splitting "Normal" Sample
- 4.4. 9.4 Riffling and Splitting Sample with Dolomite or Gypsum
- 4.5. 9.5 Splitting of Custom Blend Sample
- 4.6. 9.6 Comminuting Fertilizer Spike Sample
- 4.7. 9.7 Problems during Splitting
- 4.8. 9.8 Cleaning of Riffling and Splitting Equipment

### **5. INTERFERENCE**

5.1. N/A

### **6. APPARATUS and MATERIAL**

- 6.1. Bottle, plastic
- 6.2. Butcher paper
- 6.3. Camera
- 6.4. Cloth
- 6.5. Divider (small approximately 180-cubic inch, large approximately 450-cubic inch)
- 6.6. Downdraft Bench (with hiflow arm and compressed air hose)
- 6.7. Feeder
- 6.8. Spatula
- 6.9. Splitter

### **7. REAGENTS and MEDIA**

7.1. N/A

### **8. ASSOCIATED DOCUMENTS**

8.1. Form F00023, Storage Custody

### **9. SPECIFIC PROCEDURES**

9.1. Sample Type Determination

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

### **PREP 450**

#### *Riffling and Splitting of Fertilizer Sample*

Version: 3.0

Effective Date: 01/27/2021

Page 2 of 5

9.1.1. The laboratory encounters different types of Fertilizer samples. The composition of the individual sample will control the riffling or splitting process to be used. Determine the sample type before beginning the riffling or splitting process.

9.1.2. A “normal” sample is defined as a dry, granular, or pelletized sample that is not dolomite, gypsum, custom blend, or fertilizer spike.

9.1.3. A pelletized dolomite or gypsum sample is treated as a “normal” sample.

#### **9.2. Riffling “Normal” Sample**

9.2.1. The bag containing the submitted sample and three plastic bottles (two 250mL and one 500mL) should be located next to the Downdraft Bench.

9.2.2. Place two trays beneath the smaller divider to catch the split sample. Ensure that the gate is closed.

9.2.3. Only one sample (bags 1-4) should be on the splitting bench during this process.

9.2.4. Match the sample number on the bag with all three bottles.

9.2.5. Take a picture of the bag and all three bottles. Ensure that the numbers are in focus.

9.2.6. Open the bag and pour the entire sample into the center of the feeder to create an elongated cone-shaped pile.

9.2.6.1. If the entire sample is contained in two bags, then both bags must be combined into the feeder at the same time.

9.2.6.2. If the entire sample does not fit into the top feeder, use the larger divider. Different sized dividers are located in the laboratory.

9.2.6.3. Alternatively, sample can be run through the divider to make room for more sample; but all the sample must be run through in the same riffling.

9.2.6.4. If there are particles too large to pass through the divider, break them up manually. If the particles are concrete or rock and not fertilizer, they can be discarded.

9.2.7. After the sample has passed through the divider, dislodge any trapped particles.

9.2.8. Once the sample is collected into the two trays, return both halves of the sample into the top feeder to be riffled again.

9.2.9. Riffle the sample a minimum of three times depending on sample size.

9.2.10. If the large divider is used, discard the right tray and riffle the sample in the left tray through the small divider three more times.

#### **9.3. Splitting “Normal” Sample**

9.3.1. Flip the power switch (located in the back of the feeder) to on.

9.3.2. Set the speed level to 30.

9.3.3. If the sample has trouble flowing through the tray, an adjustment of the vibration can be done within a 25 to 35 range for the speed level. The dial on the front of the machine adjusts the speed.

9.3.4. Flip the power switch (located at the back of the splitter) to on.

9.3.5. Press the power button on the front of the machine, and then press the start button. When the start button lights up, it indicates that the machine is on.

9.3.6. Screw the three labeled plastic bottles to the openings of the splitter.

9.3.7. Pour the entire contents of one tray into the top loader of the feeder. Keep the other tray (containing the other half of the sample) on the bench. Push start on the feeder.

9.3.7.1. If the sample has a pink sticker, place an alternate container under the splitter to collect the excess; and transfer to the designated hazardous waste container upon collection.

# STANDARD OPERATING PROCEDURE BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

## PREP 450 *Riffling and Splitting of Fertilizer Sample*

Version: 3.0

Effective Date: 01/27/2021

Page 3 of 5

9.3.8. If all three bottles are not receiving sample, this indicates a problem. Suggestions on how to deal with the problem are addressed below.

9.3.9. After the entire tray of sample passes through the splitter, unscrew the bottles and cap them.

9.3.10. Place each bottle in the designated location.

| Bottle                          | Location                         |
|---------------------------------|----------------------------------|
| 250mL bottle with white label   | grinding area                    |
| 500mL bottle with white label   | storage area for retained sample |
| 250mL bottle with the red label | storage area for normal sample   |

9.3.11. Discard the sample in the other tray into the designated disposal bucket.

9.3.12. A yellow sticker indicates that calcium nitrate ( $\text{CaNO}_3$ ) is present in the sample. If the bottles are identified with a yellow sticker, the 250mL bottle with the white label must be placed into the freezer until it is ready for grinding.

### 9.4. Riffling and Splitting Sample with Dolomite or Gypsum

9.4.1. Dolomite or gypsum samples are identified as such on the sticker.

9.4.2. Place butcher paper on the table.

9.4.3. Pour the entire sample in the center of the butcher paper creating a cone shape.

9.4.4. Mix the sample thoroughly, with a spatula, at least four times.

9.4.5. Divide the sample into four approximately equal portions. Dispose of two opposing portions and combine the remaining two portions.

9.4.6. Repeat the splitting process above until just enough remains to fill the three bottles. Of the remaining two portions, place one portion into the 500mL bottle, and divide the other portion between the two 250mL bottles.

9.4.7. Place bottle in the designated location.

| Bottle                                    | Location                                                     |
|-------------------------------------------|--------------------------------------------------------------|
| 250mL bottle labeled for moisture / sieve | weighing area<br><i>sample will be used in moisture test</i> |
| 500mL bottle with white label             | storage area for retained sample                             |
| 250mL bottle labeled for oven             | oven<br><i>sample will be ground after drying</i>            |

### 9.5. Splitting of Custom Blend Sample

9.5.1. Custom Blend samples are indicated as such by a sticker and accompanied by five plastic bottles (four 250mL and one 500mL).

9.5.2. Riffle the sample as per the riffling routine of a normal sample.

9.5.2.1. Process the first tray into the three bottles as per the splitting routine of a normal sample but

9.5.2.2. Do not discard the second tray and

9.5.2.3. Detach only two of the 250mL bottles and cap them.

9.5.3. Attach the remaining two 250mL bottles and process the second tray through the splitter.

9.5.3.1. Transfer the second set of bottles to the Pesticide Formulation Storage location. Use form F00023, Storage Custody, to indicate transfer.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**PREP 450**  
*Riffling and Splitting of Fertilizer Sample*

Version: 3.0

Effective Date: 01/27/2021

Page 4 of 5

**9.6. Comminuting Fertilizer Spike Sample**

9.6.1. Place the Fertilizer Spike sample in a clean bag.

9.6.1.1. Seal the bag by twisting.

9.6.1.2. Crush the sample to make it as homogeneous as possible. A hammer or maul may be used.

9.6.2. Transfer the sample to the grinding area

**9.7. Problems during Splitting**

9.7.1. If the bottle does not receive the sample, check that the rotating arm is spinning inside the splitter.

9.7.2. If a problem occurs during the sample's run through the machine, all bottles should be emptied, and the machine returned to its starting position.

9.7.3. When the second tray of sample is used, a clean designated disposal container should be placed under the splitter to catch all discarded sample. This will ensure sample purity in case the sample needs to be split a third time.

**9.8. Cleaning of Riffling and Splitting Equipment**

9.8.1. The feeder, splitter and divider shall be cleaned between each sample. Use the compressed air hose to blow out any trapped particles in the feeder, splitter or divider.

9.8.2. After all the samples are finished, wipe any accessible metal surfaces with a water-dampened cloth.

9.8.3. The following items are rinsed with water in a sink after all samples are finished: the divider, divider trays, bottom half of the splitter, tray of the vibratory feeder, and the funnel of the vibratory feeder.

**10. WASTE MANAGEMENT**

10.1. Fertilizer containing pesticide (indicated on the bottle by a pink sticker) should be retained separately and discarded as hazardous waste.

10.2. All other solid fertilizer is collected and may be used by employees or groundskeeping.

**11. REFERENCES**

11.1. SOP for Preparing Fertilizer Samples Appendix A, 6/26/2002

**12. VERSION HISTORY**

| Version | Date       | Description                                                                                                            | Author      | Editor          |
|---------|------------|------------------------------------------------------------------------------------------------------------------------|-------------|-----------------|
| 1.0     | 07.29.2015 | Original Version                                                                                                       | Sarah Smith |                 |
| 2.0     | 11.30.2018 | Changed to new format<br>Added 9.1.2, 9.1.3, 9.2.3, 9.3.7.1,<br>9.8.3, Section 10, and 11.1<br>Minor grammatical fixes |             | Taleya Williams |
| 3.0     | 11.25.2020 | Replaced Elizabeth Schneider with<br>Quentin Cooper in Signature<br>History                                            |             | Quentin Cooper  |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 450  
*Riffling and Splitting of Fertilizer Sample*

Version: 3.0

Effective Date: 01/27/2021

Page 5 of 5

**13. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director

Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director

Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer/ Editor

Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Wet Chemistry and Preparation Manager

Date\_\_\_\_\_

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 451  
*Grinding of Fertilizer Sample*

Version: 3.0

Effective Date: 01/15/2021

Page 1 of 5

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the grinding of Fertilizer samples to be used for further analysis.

## 2. SCOPE

2.1. The procedure delineated below is applicable to Fertilizer samples in the Bureau of Agricultural Environmental Laboratories (BAEL). This procedure applies to the pulverizing of a sample and converts dry solid fertilizer into homogeneous analytical samples while preserving compositional integrity of the original material.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required. Use of downdraft for dust control is required.

## 4. OUTLINE of PROCEDURE

- 4.1. 9.1 Setting Up the Grinding Mill
- 4.2. 9.2 Assembling the Grinding Tools
- 4.3. 9.3 Grinding a Sample
- 4.4. 9.4 Cleaning and Reassembly
- 4.5. 9.5 Laboratory Information Management System (LIMS) Entry

## 5. INTERFERENCE

5.1. N/A

## 6. APPARATUS

- 6.1. Bottle, wide mouth
- 6.2. Bottlebrush, large
- 6.3. Brush, wire brass wooden handle
- 6.4. Dehumidifier
- 6.5. Downdraft Unit
- 6.6. Funnel, medium plastic
- 6.7. Glazed butcher's paper, 24 inches by 36 inches
- 6.8. Humidity Monitor
- 6.9. Mortar and pestle
- 6.10. Retsch ZM 200 equipped with ring sieves (1.0mm, 0.75mm, 0.50mm, 0.25mm) (notched) or equivalent
- 6.11. Retsch ZM 200 Ultra Centrifugal Mill with anti-rebound fitting or equivalent
- 6.12. Shop vacuum with attachments
- 6.13. Sieve, U.S. Standard No. 20
- 6.14. Sieve, U.S. Standard No. 40
- 6.15. Spatula, large

## 7. REAGENTS and MEDIA

7.1. N/A

## 8. ASSOCIATED DOCUMENTS

- 8.1. ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer
- 8.2. Form F00041, Sieving Pulverized Samples

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 451  
*Grinding of Fertilizer Sample*

Version: 3.0

Effective Date: 01/15/2021

Page 2 of 5

## 9. SPECIFIC PROCEDURES

### 9.1. Setting Up the Grinding Mill

9.1.1. At the beginning of each work period, activate the mill by turning the power on using the toggle switch at the lower, left rear of the unit.

**NOTE:** The mill cover will not open without the power on. To unlock the cover, use the pictogram that displays - "open unit."

9.1.2. Use the menu and setting button, to adjust the speed. 14,000 – 16,000 revolutions per minute (rpm) are most frequently used, with available range being 6,000 – 18,000.

9.1.3. The humidity is captured by the SensoScientific system, however the technician should visually check the humidity reading before starting grinding. The humidity in the Sample Room and the integrity of the sample will ultimately dictate the most effective rpm to use. Sample grinding should not take place when Sample Room humidity is above 40%.

### 9.2. Assembling the Grinding Tools

9.2.1. The labyrinth is placed with the torsion notch in the "6 o'clock" position.

9.2.2. The rotary blades are placed over the drive shaft.

9.2.3. Insert the collecting vessel (cassette) with the notch over the torsion lock of the labyrinth.

9.2.4. The ring sieve fits over the rotary blades, notch side down.

9.2.5. The ring sieve is placed with its notch over the torsion lock of the labyrinth (0.25mm for a First grind, Dolomite or Gypsum; 0.50mm for a Second grind; 0.75mm for a Calcium Nitrate and other designated samples).

9.2.6. The cassette cover is then installed.

9.2.7. Close the mill cover, using moderate pressure.

9.2.8. If the grinding tools have been properly assembled, a "ratcheting" sound can be heard.

9.2.9. The unit automatically locks itself.

9.2.10. If the unit will not close, use the pictogram to open the cover, inspect and/or reassemble the components until the unit's cover automatically closes.

**NOTE:** Unit will not close/lock if grinding tools are not correctly assembled.

### 9.3. Grinding a Sample

9.3.1. Samples requiring a First grind will be determined by checking the FERT\_PREP section of the Fertilizer Prep backlog list. These samples will be ground using a 0.25mm ring sieve.

9.3.2. Samples requiring a Second grind will be determined by checking the FERT\_GRIND\_2 section of the Fertilizer Prep backlog list. These samples will be ground using a 0.50mm ring sieve.

9.3.3. Calcium Nitrate samples will be indicated on the sample backlog as "CA\_NO3\_FREEZE\_PREP". These samples will be retrieved from the freezer and ground immediately using a 0.75mm ring sieve.

9.3.4. Once the mill is closed/locked, turn the downdraft bench on.

9.3.5. Select the "start" function on the mill display.

9.3.6. The performance screen is displayed and should be near zero percent.

9.3.7. Obtain the 8-ounce sample bottle.

9.3.8. Using the knob, adjust the rpm to the appropriate speed on the feeder and slowly feed the sample (while turning the bottle) into the filling hopper on the lid housing until all sample has been introduced to the mill.

9.3.9. Listen to the grinder and watch the display for any potential issues.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 451  
*Grinding of Fertilizer Sample*

Version: 3.0

Effective Date: 01/15/2021

Page 3 of 5

9.3.10. During grinding, the 100% range on the display should be avoided because the device is overloaded in this performance range and reproducibility of grind is no longer possible.

9.3.11. When the display indicates that the sample isn't grinding properly, evaluate and adjust the speed and screen size accordingly.

9.3.12. Compressed air should be blown into the filling hopper before you select the stop function, select the "stop" function.

9.3.13. The unit will unlock and begin to open. The lid can then be opened the rest of the way.

9.3.14. Remove the cassette lid. Using the bottle brush, carefully brush the pulverized sample from the lid onto the butcher's paper. Place lid aside.

9.3.15. Next, remove the ring sieve from the rotary blades and carefully brush the entire sample onto the butcher's paper. Place blades aside.

9.3.16. Remove the cassette, filled with pulverized sample from the assembly and empty the contents onto the butcher's paper. Brush the cassette to remove the entire sample. Place cassette aside.

9.3.17. Using a large spatula, thoroughly mix the pulverized sample by using a mix and roll technique. Ensure that all sample remains on the butcher's paper.

9.3.18. Fold the butcher's paper in half; set sample aside.

9.3.19. Quality Control (QC) of the grinding process is done each day that samples are to be ground. This is done for the first sample of the day on a First grind and for the first sample of the day on a Second grind.

9.3.19.1. After grinding the first sample of the day for a First grind, obtain the No. 40 sieve.

9.3.19.2. Place the sieve on top of a second piece of butcher's paper.

9.3.19.3. Retrieve the sample and pour the sample into the sieve; all pulverized sample must pass through the screen.

9.3.19.4. Sample trapped by the mesh must be re-pulverized by (a) grinding using a smaller ring sieve or (b) mortar and pestle.

9.3.19.5. After grinding the first sample of the day needing a Second grind, obtain the No. 20 sieve.

9.3.19.6. Place on top of a second piece of butcher's paper.

9.3.19.7. Retrieve the sample and pour the sample into the sieve; all pulverized sample must pass through the screen.

9.3.19.8. Sample trapped by the mesh must be re-pulverized by (a) grinding using a smaller ring sieve or (b) mortar and pestle.

9.3.19.9. Record the QC results on form F00041, Sieving Pulverized Sample.

9.3.20. Using a medium funnel pour pulverized sample back into 8-ounce bottle.

9.3.21. Place completed samples on a lab cart with a sign "Ready To Be Weighed".

9.4. Cleaning and Reassembly

9.4.1. After each sample, clean the cassette cover, ring sieve, rotor blades, cassette and butcher's paper using appropriate brush and compressed air.

9.4.2. Reassemble the grinding tools.

9.5. Laboratory Information Management System (LIMS) Entry

9.5.1. Log on to LIMS. Refer to ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer.

9.5.2. Select template FERT\_LOGIN for First grind samples, or FERT\_GRIND\_2 for Second grind samples.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 451  
*Grinding of Fertilizer Sample*

Version: 3.0

Effective Date: 01/15/2021

Page 4 of 5

9.5.3. Scan or manually enter the sample ID's for the completed samples.

9.5.4. Enter "C" or appropriate sieve size in the GRIND-1 or GRIND-2 box for each completed sample.

9.5.5. In the GRIND-1 template, check to make sure codes such as CaNO3 and initial dilution codes are completed. If the grind code is grayed out, indicating test codes haven't been entered for the sample, the administrative Environmental Specialist 1 (ES1) will be notified.

9.5.6. The administrative ES1 should be notified to add missing tag codes. If the tag codes are missing, set the samples aside until the codes have been added.

## 10. WASTE MANAGEMENT

10.1. N/A

## 11. REFERENCES

11.1. Preparing "First" and "Miscellaneous" Samples Using Retsch ZM 200 Ultra Centrifugal Mill. FSFL-SOP 303, (01/09/2006).

11.2. Preparing Commercial and Deficient Samples ("Seconds") Using Retsch ZM 200 Ultra Centrifugal Mill. FSFL-SOP 304, (01/09/2006).

11.3. Standard Operating Procedures For Preparing Fertilizer Samples, (06/26/2002).

11.4. Operating Instructions Ultra Centrifugal Mill Type ZM 200 Retsch GmbH, (08/02/12).

## 12. VERSION HISTORY

| Version | Date       | Description                                                                                                                         | Author        | Editor          |
|---------|------------|-------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------|
| 1.0     | 2015.12.22 | Original Version                                                                                                                    | George Taylor |                 |
| 2.0     | 11.30.2018 | Updated outline format<br>Grammatical and verbiage fixes                                                                            |               | Taleya Williams |
| 3.0     | 12.09.2020 | Added or equivalent to items in<br>Section 6 Apparatus.<br>Replaced Elizabeth Schneider with<br>Quentin Cooper in Signature History |               | Taleya Williams |

***STANDARD OPERATING PROCEDURE***  
***BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES***

PREP 451  
*Grinding of Fertilizer Sample*

Version: 3.0

Effective Date: 01/15/2021

Page 5 of 5

**13. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer  
Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Wet Chemistry and Preparation Manager/Editor  
Date\_\_\_\_\_

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 452

### *Dilution of Fertilizer Sample*

Version: 4.0

Effective Date:08/02/2024

Page 1 of 2

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for dilution of liquid fertilizer samples.

## **2. SCOPE**

2.1. The procedure delineated below is applicable to fertilizer samples in the Bureau of Agricultural Environmental Laboratories (BAEL). This procedure designates how to properly dilute liquid fertilizer samples.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## **4. OUTLINE OF PROCEDURE**

- 4.1. 9.1 Check Balance
- 4.2. 9.2 Sample and Container Preparation
- 4.3. 9.3 Sample Dilution and LIMS data Entry

## **5. INTERFERENCE**

5.1. N/A

## **6. APPARATUS and EQUIPMENT**

- 6.1. Balance, top loader
- 6.2. Camera
- 6.3. Container, 1000 ml square Nalgene or other appropriate container
- 6.4. Labels
- 6.5. Laboratory Information Management System (LIMS)

## **7. REAGENTS and MEDIA**

7.1. Deionized (DI) Water, Type 2

## **8. ASSOCIATED DOCUMENTS**

- 8.1. ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer
- 8.2. Form F00010, Balance Verification
- 8.3. LABOP 122, Maintenance, Operation, and Program of Balances

## **9. SPECIFIC PROCEDURES**

- 9.1. Check Balance
  - 9.1.1. Balance should be clean and level.
  - 9.1.2. Check balance calibration daily prior to first use, and record on form F00010, Balance Verification. Refer to LABOP 122, Maintenance, Operation, and Program of Balances.
- 9.2. Sample and Container Preparation
  - 9.2.1. Obtain a clean plastic 1000mL square container.
  - 9.2.2. Inspect for discolorations, imperfections and integrity before use.
  - 9.2.3. Examine original sample containment prior to opening for signs of leakage.
  - 9.2.4. Should sample appear to be compromised, identify the problem and notify management via form F00021, Nonconformance and Deviation.

# STANDARD OPERATING PROCEDURE BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

PREP 452

## *Dilution of Fertilizer Sample*

Version: 4.0

Effective Date:08/02/2024

Page 2 of 2

- 9.2.5. Remove sample container from shipping package.
- 9.2.6. Write sample number on original container if omitted in field.
- 9.2.7. Take photo of original containment and all new labeled containers for field identification verification purposes.

### 9.3. Sample Dilution and LIMS data entry

- 9.3.1. Remove cap and weigh original sample.
- 9.3.2. Record as "Total Weight" in LIMS
- 9.3.3. **Note: Remember to weigh and record all samples before proceeding to the next step.**
- 9.3.4. Tare new container and transfer the entire sample to the new container.
- 9.3.5. Rinse original container and lid with DI water enough times to sufficiently dissolve any existing crystals. The rinsate is added to the new container.
- 9.3.6. Add sufficient DI water to almost fill the new container.
- 9.3.7. Weigh new container and record as "New Weight" in LIMS.
- 9.3.8. Allow original sample container to dry, then weigh and record as "Bottle Weight" in LIMS.
- 9.3.9. If crystals are still present after dilution contact a manager for further instructions.
- 9.3.10. Place completed samples on a lab cart with a sign "Ready To Be Weighed".
- 9.3.11. Record the "Dilution Factor" on sample container.
- 9.3.12. Refer to ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer.
- 9.3.13. Select template FERT\_LOGIN.
- 9.3.14. Enter "C" in the GRIND-1 field for each completed sample.
  - 9.3.14.1. If a liquid sample requires further analysis after the initial analysis, all subsequent analyses are performed on the original dilution. "Seconds" cannot be prepared on liquids. GRIND-2 has to be "completed" to trigger additional analyses.

## 10. WASTE MANAGEMENT

- 10.1. Original sample containers can be discarded in the trash after process is complete.

## 11. REFERENCES

- 11.1. N/A

## 12. VERSION HISTORY

| Version | Date       | Description                                                                        | Author          | Editor          |
|---------|------------|------------------------------------------------------------------------------------|-----------------|-----------------|
| 1.0     | 03.09.2016 | Original Version                                                                   | Linda Giardelli |                 |
| 2.0     | 11.30.2018 | Updated outline format<br>Added 9.3.8, 9.3.11, 9.3.13.1<br>Minor grammatical fixes |                 | Taleya Williams |
| 3.0     | 11/25/2020 | Replaced Elizabeth Schneider<br>with Quentin Cooper in<br>Signature History        |                 | Quentin Cooper  |
| 4.0     | 08.02.2024 | Removed Signature History                                                          |                 | Eddie Bello     |

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 453

### *Weighing of Fertilizer Sample*

Version: 4.0

Effective Date: 08/02/2024

Page 1 of 4

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for the weighing of fertilizer for sample preparation.

## **2. SCOPE**

2.1. The procedure delineated below is applicable to Fertilizer samples in the Bureau of Agricultural Environmental Laboratories (BAEL). See the Containers and Codes table in section 9 below for included preparation methods.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## **4. OUTLINE of PROCEDURE**

- 4.1. 9.1 Preparation for Weighing
- 4.2. 9.2 Login Laboratory Standards
- 4.3. 9.3 Weighing
- 4.4. Attachment A, Fertilizer Laboratory Standard and Batch Size

## **5. INTERFERENCE**

5.1. Water and humidity can cause issues with this process; historical data has shown that performing these procedures with a humidity less than 40% is required.

## **6. APPARATUS and MATERIAL**

- 6.1. Balance, analytical
- 6.2. Barcode Reader
- 6.3. Beakers, disposable, 50mL
- 6.4. Centrifuge tubes, 50mL
- 6.5. Computer
- 6.6. Erlenmeyer flask, 500mL
- 6.7. Laboratory Information Management System (LIMS)
- 6.8. Pipette, disposable
- 6.9. Spatula
- 6.10. Sticker, orange
- 6.11. Volumetric flask, 200mL
- 6.12. Volumetric flask, 250mL

## **7. REAGENTS and MEDIA**

- 7.1. SRM 200b
- 7.2. SRM 695
- 7.3. SRM 88b

## **8. ASSOCIATED DOCUMENTS**

- 8.1. ADMIN 023 Attachment C, LIMS Sample Login for Multi Sample
- 8.2. ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer
- 8.3. LABOP 122, Maintenance, Operation and Program of Balances

## **9. SPECIFIC PROCEDURES**

- 9.1. Preparation for Weighing

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 453

#### *Weighing of Fertilizer Sample*

Version: 4.0

Effective Date: 08/02/2024

Page 2 of 4

9.1.1. Check that the balance is clean, level, and calibrated. Refer to LABOP 122, Maintenance, Operation, and Program of Balances.

9.1.2. Generate or print out either Fert\_Prep\_Weight or Fert\_Prep\_Weight2 sections from backlog report.

9.1.3. Obtain samples from grinding room.

9.1.4. **No more than 10 samples should be scanned for weighing.** Samples denoted as "Priority" should be weighed first. Priority samples may be included in smaller batches in order to expedite their preparation.

9.1.5. Program requirements for proficiency samples mandate two weights/preps on different days.

9.1.6. Using backlog, generate labels and line up appropriate weighing vessel behind each sample.

9.1.7. The table below shows the weighing vessel for each preparation method.

| Method       | Core Weight Code | Core Suggested Weight Code | Container               |
|--------------|------------------|----------------------------|-------------------------|
| 430 NAU      | 430_NAU_WT       | 430_XXX_WT_SUG             | 250mL volumetric flask  |
| 430 KJL      | 430_KJL_WT       | 430_XXX_WT_SUG             | 250mL volumetric flask  |
| 461          | <b>461_WT</b>    | <b>461_APA_WT_SUG</b>      | 500mL Erlenmeyer flask  |
| 505 APA      | 505_WT           | 505_WT_SUG                 | 500mL Erlenmeyer flask  |
| 521 TP       | 521_WT           | 521_WT_SUG                 | 500mL Erlenmeyer flask  |
| 541 CIP      | 541_WT           | 541_WT_SUG                 | 500mL Erlenmeyer flask  |
| 604 SMGK     | 604_WT           | 604_WT_SUG                 | 50mL disposable beakers |
| 803 SB       | 803_WT           | 803_WT_SUG                 | 50mL centrifuge tube    |
| 822 TOTAL    | 822_WT           | 822_WT_SUG                 | 50mL centrifuge tube    |
| 825 CO/AL    | 825_WT           | 825_WT_SUG                 | 50mL centrifuge tube    |
| 830 SMN      | 830_WT           | 830_WT_SUG                 | 200mL volumetric flask  |
| 831 SCU      | 831_WT           | 831_WT_SUG                 | 200mL volumetric flask  |
| 832 SFE      | 832_WT           | 832_WT_SUG                 | 200mL volumetric flask  |
| 833 SZN      | 833_WT           | 833_WT_SUG                 | 200mL volumetric flask  |
| 840 CFE      | 840_WT_F         | 840_WT_SUG_F               | 200mL volumetric flask  |
| 841 CFE      | 841_WT_F         | 841_WT_SUG_F               | 200mL volumetric flask  |
| 841 CMN      | 841_WT_M         | 841_WT_SUG_M               | 200mL volumetric flask  |
| 841 CCU      | 841_WT_C         | 841_WT_SUG_C               | 200mL volumetric flask  |
| 841 CZN      | 841_WT_Z         | 841_WT_SUG_Z               | 200mL volumetric flask  |
| 860 TCA      | 860_TCA_WT       | 860_WT_SUG_TCA             | 50 mL centrifuge tube   |
| 875 DOL LIME | 875_WT           | 875_WT_SUG                 | 50mL centrifuge tube    |
| 902 HM       | 902_WT           | 902_WT_SUG                 | 50mL centrifuge tube    |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 453

#### *Weighing of Fertilizer Sample*

Version: 4.0

Effective Date: 08/02/2024

Page 3 of 4

9.1.8. Log on to Laboratory Information Management System.

#### 9.2. Login Laboratory Standards

9.2.1. To log in laboratory (lab) standards, refer to ADMIN 023 Attachment C, LIMS Sample Login for Multi Sample.

9.2.2. Print or record the lab standard sample numbers generated.

9.2.3. Suggested weights for lab standards will be listed at the top of that sample's column in the result entry template.

9.2.4. Refer to Attachment A, Fertilizer Laboratory Standard and Batch Size.

#### 9.3. Weighing

9.3.1. In LIMS, open Result Entry using the appropriate entry template:

9.3.1.1. For zero weights, use template "FERT\_PREP\_WEIGHT\_0"

9.3.1.2. For first, second, and referee weights, select "FERT\_PREP\_WEIGHT2"

9.3.2. Scan sample barcodes into appropriate template.

9.3.3. Place a labeled container on the balance and tare.

9.3.4. Confirm that the label on the sample bottle matches the sample number in the LIMS template.

9.3.5. For solid samples:

9.3.5.1. The sample must be mixed thoroughly by shaking and rotating the bottle back and forth to break up clumps before weighing the sample. A spatula may be used to break up the sample before the mixing.

9.3.5.2. Using a spatula, weigh the amount listed in LIMS under the section for suggested weight code. Wipe spatula between samples.

9.3.5.3. Confirm that the Initial\_Liquid\_Dilution in LIMS for a solid sample, including lab standard, is 1. If it is not, change it to 1.

9.3.5.4. If there is another value in the test code and this is not the first time the sample is being weighed (1st/2nd preps) staff should notify the section manager or administrative ESI because there may already be sample results calculated using the incorrect value and changing it here will not update the results in LIMS.

9.3.6. For liquid samples:

9.3.6.1. Shake the sample bottle vigorously, then inspect sample bottle for residual crystallization. If any crystals are found, notify supervisor and postpone weighing until further instructions are given.

9.3.6.2. Using a disposable pipette, weigh double the amount listed in LIMS under the suggested weight code for diluted samples. If the sample is not diluted, do not double the suggested weight.

9.3.6.3. Use a clean pipette for each sample.

9.3.6.4. Confirm that the liquid dilution factor on the bottle matches the liquid dilution code in LIMS.

9.3.6.5. If there is another value in the test code and this is not the first time the sample is being weighed (1st/2nd preps) staff should notify the section manager or administrative ESI because there may already be sample results calculated using the incorrect value and changing it here will not update the results in LIMS.

9.3.7. Wait for the weight reading on the balance to stabilize.

9.3.8. Press **transfer button** to load sample weight reading from balance directly into LIMS. If this function is unavailable, manually enter the weight as displayed by the balance.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 453

#### *Weighing of Fertilizer Sample*

Version: 4.0

Effective Date: 08/02/2024

Page 4 of 4

9.3.9. Weigh out all requested analyses based on suggested weight for one sample at a time.

9.3.9.1. Core test codes may be followed by \_1, \_2, R\_1, R\_2

9.3.9.2. If %DIV/0 is in the suggested weight code, contact supervisor via email for correction of the problem.

9.3.9.3. If more than one suggested weight is shown, weigh the higher suggested weight.

9.3.10. Place each weighed sample vessel in a **bin**.

9.3.11. The laboratory standard must be weighed during the process of weighing for each batch and included in the bin for the corresponding batch.

9.3.12. Completed batches in bins should be placed on their appropriate shelves outside the balance room.

9.3.13. Be sure, at the end of each day, to return original sample bottles to the grinding room and place them in their correct numerical order on shelves. Certain samples, such as commercial and proficiency samples, may have their own specific location.

## 10. WASTE MANAGEMENT

10.1. N/A

## 11. REFERENCES

11.1. N/A

## 12. VERSION HISTORY

| Version | Date       | Description                                                                                                                                                                                                                                                                                            | Author        | Editor          |
|---------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------|
| 1.0     | 01.13.2016 | Original Version                                                                                                                                                                                                                                                                                       | George Taylor |                 |
| 2.0     | 12.05.2018 | Update outline format<br>Added 5.1, 9.3.5.4, 9.3.6.5<br>Minor grammatical fixes                                                                                                                                                                                                                        |               | Taleya Williams |
| 3.0     | 11.25.2020 | Replaced Elizabeth Schneider with Quentin Cooper in Signature History                                                                                                                                                                                                                                  |               | Quentin Cooper  |
| 4.0     | 08.02.2024 | Removed Signature History.<br>9.1.4 Changed entire first sentence to a new one limiting the number of samples to "no more than 10".<br>9.1.7. Corrected codes for method 461<br>9.3.8. Changed sentence to now refer to "transfer button".<br>9.3.10. Simplified sentence to just end with "in a bin". |               | Eddie Bello     |

# ATTACHMENT

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 453 Attachment A

#### *Fertilizer Laboratory Standard and Batch Size*

Version: 5.0

Effective Date: 08/02/2024

Page 1 of 1

| BAEL SOP                 | FM Method                                                      | Batch Size* | Location Code          | Standard Material        |
|--------------------------|----------------------------------------------------------------|-------------|------------------------|--------------------------|
| METHOD 502               | 405 TKN                                                        | 24          | QA_FERT_LS_TKN_405     | SRM 695                  |
| METHOD 505<br>METHOD 530 | 415 KJL                                                        | 24          | QA_FERT_LS_DIR_AMM_415 | Ammonium Sulfate Std.    |
| METHOD 504<br>METHOD 530 | 420 KJL                                                        | 24          | QA_FERT_LS_DIR_AMM_420 | Ammonium Sulfate Std.    |
| METHOD 500<br>METHOD 531 | 425 KJL                                                        | 24          | QA_FERT_LS_DIR_AMM_425 | Ammonium Sulfate Std.    |
| PREP 466                 | 430 KJL                                                        | 47          | QA_FERT_LS_NAUL        | SRM 695                  |
| PREP 466                 | 430 NAU                                                        | 47          | QA_FERT_LS_NAUS        | SRM 695                  |
| METHOD 506<br>METHOD 530 | 435 KJL                                                        | 24          | QA_FERT_LS_DIR_AMM_435 | SRM 695                  |
| METHOD 500<br>METHOD 531 | 440 KJL                                                        | 24          | QA_FERT_LS_DIR_AMM_440 | SRM 695                  |
| METHOD 501<br>METHOD 531 | 445 KJL                                                        | 24          | QA_FERT_LS_WIN_445     | SRM 695                  |
| METHOD 461<br>P&K        | N/A                                                            | 21          | QA_FERT_LS_PANDK       | SRM 695                  |
| PREP 454                 | 604 SMGK                                                       | 18          | QA_FERT_LS_SMGK        | SRM 695                  |
| METHOD 510               | 611 CLO                                                        | 20          | QA_FERT_LS_CHLORINE    | SRM 695                  |
| PREP 455                 | 803 SB                                                         | 47          | QA_FERT_LS_SB          | SRM 695                  |
| METHOD 520               | 805 CS<br>860 CS                                               | 21          | QA_FERT_LS_SULFUR_COMB | Ammonium Sulfate Std.    |
| METHOD 520               | 805 CS<br>860 CS                                               | 21          | QA_FERT_LS_SULFUR_COMB | Calcium Thiosulfate Std. |
| METHOD 521               | 806 FS                                                         | 8           | QA_FERT_LS_SULFUR_FREE | Free Sulfur Std.         |
| PREP 456                 | 822 TOTAL<br>825 Co/AL<br>860 TCA<br>875 DOL<br>LIME<br>902 HM | 47          | QA_FERT_LS_TOTAL       | SRM 695                  |
| PREP 462                 | 830 SMN                                                        | 24          | QA_FERT_LS_SMN         | SRM 695                  |
| PREP 462                 | 831 SCU                                                        | 24          | QA_FERT_LS_SCU         | SRM 695                  |
| PREP 462                 | 832 SFE                                                        | 24          | QA_FERT_LS_SFE         | SRM 695                  |
| PREP 462                 | 833 SZN                                                        | 24          | QA_FERT_LS_SZN         | SRM 695                  |
| PREP 464                 | 840 CFE                                                        | 25          | N/A                    |                          |
| PREP 464                 | 841                                                            | 25          | N/A                    |                          |

\*Batch Size includes batch spike.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**PREP 454**

*Sample Preparation for Water Soluble Potash and Magnesium Analysis*  
*[FM 604]*

Version: 3.0

Effective Date: 10/16/2023

Page 1 of 3

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for the preparation of water soluble magnesium (Mg) and water soluble potash ( $K_2O$ ) nutrients in fertilizer samples for analysis.

## **2. SCOPE**

2.1. The procedure delineated below is applicable to Fertilizer samples in the Bureau of Agriculture Environmental Laboratories (BAEL). This preparation is for the analysis of water soluble potash and magnesium by inductively coupled plasma optical emission spectrometry.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## **4. OUTLINE of PROCEDURE**

- 4.1. 10.1 Collection of Weighed Samples
- 4.2. 10.2 Digestion
- 4.3. 10.3 Final Volume
- 4.4. 10.4 Data Entry

## **5. INTERFERENCE**

- 5.1. N/A

## **6. APPARATUS**

- 6.1. Centrifuge tube, disposable (15mL)
- 6.2. Glass Bead
- 6.3. Hot plate
- 6.4. Parafilm
- 6.5. Rack
- 6.6. Volumetric flask (500mL)

## **7. REAGENTS and MEDIA**

- 7.1. Water, Deionized (DI) Type II

## **8. REAGENT PREPARATION GUIDELINES**

- 8.1. N/A

## **9. ASSOCIATED DOCUMENTS**

- 9.1. PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size

## **10. SPECIFIC PROCEDURES**

- 10.1. Collection of Weighed Samples

10.1.1. Weigh 0.5g of sample into a 50mL disposable beaker.

10.1.2. If samples are pre-weighed, Collect weighed samples.

10.1.3. Batch size and Laboratory Control Sample (LCS) material are defined in PREP 453 attachment A.

- 10.2. Digestion

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

### **PREP 454**

#### *Sample Preparation for Water Soluble Potash and Magnesium Analysis [FM 604]*

Version: 3.0

Effective Date: 10/16/2023

Page 2 of 3

10.2.1. Pre-heat hot plate by turning it on.

10.2.2. Record on each volumetric flask the appropriate sample number as listed on weigh vessel label.

10.2.3. Transfer sample from weigh vessel to labeled 500mL volumetric flask with thorough rinsing with DI water.

10.2.4. Add approximately 300mL of DI water and approximately 3-5 glass beads to each volumetric flask.

10.2.5. Place the volumetric flasks onto the hot plate.

10.2.6. Bring samples to a boil and continue boiling for at least one hour. The volumetric flasks **should be consistently monitored to ensure all samples are vigorously boiling. Flasks can** be rotated periodically to ensure proper heat distribution. Do not allow samples to dry out. More DI water can be added if necessary.

10.2.7. Remove flasks from the hotplate and allow cooling to room temperature.

#### **10.3. Final Volume**

10.3.1. Bring each sample to 500mL final volume with DI water and stopper the flask.

**10.3.2.** Mix the contents by **thoroughly** inverting the flask **multiple times and allow to settle**

10.3.3. **Pour extracts into** labeled the 15mL centrifuge tubes.

#### **10.4. Data Entry**

10.4.1. When preparation is complete, make sure all applicable preparation information is recorded in LIMS correctly in the appropriate template.

10.4.2. Refer to ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer, for result entry guidance.

10.4.3. Deliver samples and batch folder to the analysis section, for analysis by METHOD 513.

### **11. WASTE MANAGEMENT**

11.1. Material, marked for disposal, is removed via the following mechanism:

| <b>Material</b>        | <b>Disposal Mechanism</b>                                                |
|------------------------|--------------------------------------------------------------------------|
| Digested Sample        | Poured down the sink.                                                    |
| Glass Beads            | Capture and dispose of in trash. Make sure they do not go down the sink. |
| Used Screening Columns | Dispose of in trash                                                      |

### **12. REFERENCES**

12.1. AOAC Official Method 937.02 Magnesium (Water Soluble) in Fertilizers 16<sup>th</sup> Edition Volume 1

12.2. AOAC Official Method 983.02 Potassium in Fertilizers Flame Photometric Method (Manual or Automated) 16<sup>th</sup> Edition Volume 1

***STANDARD OPERATING PROCEDURE  
BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES***

PREP 454

*Sample Preparation for Water Soluble Potash and Magnesium Analysis  
[FM 604]*

Version: 3.0

Effective Date: 10/16/2023

Page 3 of 3

### 13. VERSION HISTORY

| Version | Date       | Description                                                                                                                                                                                                                                                                               | Author           | Editor          |
|---------|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|-----------------|
| 1.0     | 08.12.2015 | Original Version                                                                                                                                                                                                                                                                          | Joshua Hartstein |                 |
| 2.0     | 04.05.2019 | Updated Outline Format<br>Added section 7.1, 9.1, 10.1.2<br>Clarified 10.1.1, 10.2.6<br>Removed dilution section<br>Renamed 10.4 and combined with 10.5<br>Added Waste Management Section, Section 11                                                                                     |                  | Laketa Levy     |
| 3.0     | 10.16.2023 | Section 10.1.1 and section 10.1.2. added weighing instructions.<br>Section 10.2.6. expanded upon instructions to ensure vigorous boiling.<br>Section 10.3.2. added instructions for the mixing of the sample.<br>Section 10.3.3 added instructions for pouring aliquots into 15 mL tubes. |                  | Taleya Williams |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

PREP 455

### *Sample Preparation for Water Soluble Boron Analysis* *[FM 803]*

Version: 3.0

Effective Date: 01/16/2024

Page 1 of 3

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the preparation of water soluble boron (B) micronutrients in Fertilizer samples.

## 2. SCOPE

2.1. The procedure delineated below is applicable to Fertilizer samples in the Bureau of Agricultural Environmental Laboratories (BAEL). This preparation is for the analysis of water soluble boron by inductively coupled plasma optical emission spectrometry.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## 4. OUTLINE of PROCEDURE

- 4.1. Collection of Weighed Samples
- 4.2. Digestion
- 4.3. Final Volume
- 4.4. Data Entry

## 5. INTERFERENCE

5.1. Glass apparatus can interfere with boron analysis.

## 6. APPARATUS and MATERIAL

- 6.1. Centrifuge tube, non-glass disposable (15mL)
- 6.2. Parafilm
- 6.3. DigiTUBE, 50mL designed for Digester
- 6.4. SCP DigiPREP hot block or equivalent brand
- 6.5. Watch glasses, plastic, disposable

## 7. REAGENTS and MEDIA

7.1. Water, Deionized (DI) Type II

## 8. REAGENT PREPARATION GUIDELINES

8.1. N/A

## 9. ASSOCIATED DOCUMENTS

- 9.1. ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer
- 9.2. PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size

## 10. SPECIFIC PROCEDURES

10.1. Collection of Weighed Samples

10.1.1. Weigh samples per the table below:

| Guarantee | Weight |
|-----------|--------|
| 0.074%    | 1g     |
| 0.075%    | 0.8g   |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 455

#### *Sample Preparation for Water Soluble Boron Analysis [FM 803]*

Version: 3.0

Effective Date: 01/16/2024

Page 2 of 3

10.1.2. If samples are pre-weighed, collect weighed samples.

10.1.3. Batch size and Laboratory Control Sample (LCS) material are defined in PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size.

#### 10.2. Digestion

10.2.1. Pre-heat hot block with temperature setting of  $106 \pm 1^\circ\text{C}$ .

10.2.2. Make sure the tube with the temperature probe is filled at least to the 40mL mark with water or approximately 50% nitric acid. Approximately 50% nitric acid must be used if other samples that are prepped in acid will be digesting at the same time.

10.2.3. Add approximately 30-40mL of DI water to each sample tube.

10.2.4. Place the sample tube onto hot block and place a plastic disposable watch glass on each tube.

10.2.5. Visually verify that the heating block is at  $106 \pm 1^\circ\text{C}$ .

10.2.6. Digest samples for a minimum of 60 minutes at  $106 \pm 1^\circ\text{C}$ . The time can be reset if needed to ensure 60 minutes digestion.

#### 10.3. Final Volume

10.3.1. After digestion is complete, remove the samples from the hot block and allow the samples to cool to room temperature.

10.3.2. Remove the watch glasses and bring each sample up to final volume of 50 mL with DI water.

10.3.3. Cap and invert each sample. Allow samples to settle.

#### 10.4. Data Entry

10.4.1. When preparation is complete, make sure all applicable preparation information is recorded in LIMS correctly in the appropriate template.

10.4.2. Refer to ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer, for result entry guidance.

10.4.3. Deliver samples and batch folder to analysis section for analysis by applicable method.

## 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material               | Disposal Mechanism                            |
|------------------------|-----------------------------------------------|
| Sample not used        | Dispose as non-hazardous waste down the drain |
| Watch Glasses          | Dispose of in trash                           |
| Used Screening Columns | Dispose of in trash                           |

## 12. REFERENCES

12.1. E00142, Boron Analysis, FM 803, (02/06/2003)

12.2. E00143, Hot Block Preparation of Boron Samples for ICP-OES Analysis. FSFL-SOP 513 (12/02/2005)

12.3. E00163, Boron (Acid and Water Soluble) in Fertilizers. AOAC 16<sup>th</sup> Edition Method 982.01, (1995)

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 455  
*Sample Preparation for Water Soluble Boron Analysis*  
*[FM 803]*

Version: 3.0

Effective Date: 01/16/2024

Page 3 of 3

### 13. VERSION HISTORY

| Version | Date       | Description                                                                                                                                                                                                                       | Author           | Editor          |
|---------|------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|-----------------|
| 1.0     | 03.09.2016 | Original Version                                                                                                                                                                                                                  | Joshua Hartstein |                 |
| 2.0     | 04.05.2019 | Updated Outline Format<br>Added section 5.1, 9.2, 10.1.1, 10.1.2, 10.2.6,<br>10.3.1, 10.3.2, 10.3.3<br>Clarified 6.1, 10.2.2<br>Renamed 10.3<br>Renamed 10.4 and combined with 10.5<br>Added Waste Management Section, Section 11 |                  | Taleya Williams |
| 3.0     | 01.16.2024 | Removed 10.3.4-5<br>Added weight table<br>Removed screening column from materials<br>Removed signature page<br>Added external document references                                                                                 |                  | Taleya Williams |

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

### **PREP 456**

#### *Sample Preparation for Total Micronutrients and Heavy Metal Analysis [FM 822, FM 825, FM 860, FM 875, FM 902]*

Version: 5.0

Effective Date: 12/05/2023

Page 1 of 3

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for the preparation of acid soluble micronutrients and heavy metals [Aluminum (Al), Arsenic (As), Cadmium (Cd), Calcium (Ca), Chromium (Cr), Cobalt (Co), Copper (Cu), Iron (Fe), Lead (Pb), Magnesium (Mg), Manganese (Mn), Molybdenum (Mo), Selenium (Se), and Nickel (Ni) and Zinc (Zn)] in fertilizer samples.

## **2. SCOPE**

2.1. The procedure delineated below is applicable to Fertilizer samples in the Bureau of Agricultural Environmental Laboratories (BAEL). Acid soluble micronutrient metals and heavy metals are extracted from fertilizer sample by digesting the sample with nitric and hydrochloric acids using a hot block. The sample extract is then analyzed by inductively coupled plasma optical emission spectrometry.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

3.2. The acids may react with samples to produce noxious nitrogen oxides. This is accelerated with heat. Work in the fume hood to avoid breathing these vapors.

## **4. OUTLINE of PROCEDURE**

- 4.1. 9.1 Preparation and Set up
- 4.2. 9.2 Sample Digestion
- 4.3. 9.3 Data Entry

## **5. INTERFERENCE**

- 5.1. N/A

## **6. APPARATUS and MATERIAL**

- 6.1. Autodiluter
- 6.2. Centrifuge tube, disposable glass (15mL)
- 6.3. Heating block, DigiPREP MS or equivalent
- 6.4. DigiTUBE, plastic 50mL designed for Digester or equivalent
- 6.5. Parafilm
- 6.6. Tube rack
- 6.7. Watch glass (50mL), disposable

## **7. REAGENTS and MEDIA**

- 7.1. Hydrochloric acid (HCl), trace metal grade or equivalent
- 7.2. Nitric acid (HNO<sub>3</sub>), trace metal grade or equivalent
- 7.3. Water, Deionized (DI) Type I

## **8. ASSOCIATED DOCUMENTS**

- 8.1. ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer
- 8.2. Form 00026, Dispenser/Diluter Verification
- 8.3. LABOP 120 Attachment C, Maintenance and Operation of Hot Block
- 8.4. PREP 453, Weighing of Fertilizer Sample
- 8.5. PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 456

#### *Sample Preparation for Total Micronutrients and Heavy Metal Analysis [FM 822, FM 825, FM 860, FM 875, FM 902]*

Version: 5.0

Effective Date: 12/05/2023

Page 2 of 3

## 9. SPECIFIC PROCEDURES

### 9.1. Preparation and Set up

9.1.1. Make sure hood is on and the ribbon indicates flow.

9.1.2. Turn on heating block and set the temperature to  $106 \pm 1^{\circ}\text{C}$  and set the heating time for 90 minutes. Please refer to LABOP 120 Attachment C, Maintenance and Operation of Hot Block for instructions of use of hot block.

9.1.3. Check the DigiPROBE storage well container to ensure there are no cracks, holes, or that the container is not melted. If there are any cracks or holes, replace with new 50 mL DigiTUBE. If the tube shows signs of melting, the hot block should not be used.

9.1.4. Make sure the probe storage tube is filled with approximately 30 mL  $\text{HNO}_3$  and 10 mL HCl to at least the 40 mL mark.

9.1.5. Collect weighed samples. Samples are weighed according to PREP 453, Weighing Fertilizer.

9.1.6. Batch size and Laboratory Control Sample (LCS) material are defined in PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size.

### 9.2. Sample Digestion

9.2.1. Weigh 0.5g of fertilizer into a 50mL digestion tube.

9.2.2. **Under a fume hood:** Add 4.5mL of trace metal grade  $\text{HNO}_3$  allowing samples to react before adding 1.5mL of trace metal grade HCl to each sample.

9.2.3. Place the samples in the hot block and cover each sample tube with a disposable watch glass.

9.2.4. Visually verify that the heating block is at  $106 \pm 1^{\circ}\text{C}$ .

9.2.5. Digest samples for a minimum of 60 minutes (without boiling) at  $106 \pm 1^{\circ}\text{C}$ . The time can be reset if needed to ensure 60 minutes digestion.

9.2.6. After digestion is complete, remove the samples from the hot block. Leave the watch glasses in place and put the samples in a tube rack. Leave the tube rack in the hood and allow the samples to cool to room temperature.

9.2.7. Remove and rinse the watch glass with DI water into the sample and bring the sample up to 50mL with DI water.

9.2.8. Cap and mix each sample. Allow samples to sit for particulates to settle.

### 9.3. Data Entry

9.3.1. When preparation is complete, make sure all applicable preparation information is recorded in LIMS correctly in the appropriate template.

9.3.2. Refer to ADMIN 023 Attachment Q, Manual Result Entry for Fertilizer, for result entry guidance.

9.3.3. Deliver samples and batch folder to analysis section for analysis by applicable METHOD 509.

## 10. WASTE MANAGEMENT

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 456

#### *Sample Preparation for Total Micronutrients and Heavy Metal Analysis* [FM 822, FM 825, FM 860, FM 875, FM 902]

Version: 5.0

Effective Date: 12/05/2023

Page 3 of 3

10.1. Material, marked for disposal, is removed via the following mechanism:

| Material         | Disposal Mechanism                 |
|------------------|------------------------------------|
| Screening Column | Dispose of in trash                |
| Acid Waste       | Neutralized and disposed down sink |

## 11. REFERENCES

11.1. Secondary/Micronutrient Analysis. FM 822, (05/20/2002) [External 00144](#)

11.2. Preparation and Analysis of Acid Soluble (Total) Micronutrients by ICP-OES FSFL-SOP 504, (04/15/2009) [External 00145](#)

11.3. *Official Methods of Analysis of AOAC INTERNATIONAL* (2019) 21st Ed., AOAC INTERNATIONAL, Rockville, MD, USA, **Method 937.02**, [www.eoma.aoac.org](http://www.eoma.aoac.org) [External 00135](#)

11.4. *Official Methods of Analysis of AOAC INTERNATIONAL* (2019) 21st Ed., AOAC INTERNATIONAL, Rockville, MD, USA, **Method 2017.02**, [www.eoma.aoac.org](http://www.eoma.aoac.org) [External 00135](#)

## 12. VERSION HISTORY

| Version | Date       | Description                                                                                                                                                                                                                                                                                            | Author        | Editor                               |
|---------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|--------------------------------------|
| 1.0     | 06.01.2016 | Original Version                                                                                                                                                                                                                                                                                       | George Taylor |                                      |
| 2.0     | 04.19.2019 | Updated Outline Format<br>Combined with PREP 457 and added applicable information<br>Reformatted section 9.1, 9.4, to make consistent with other PREP SOPs<br>Added Waste Management, section 10<br>Removed different preparations for 0s,1s, and 2s<br>Grammatical fixes and clarification throughout |               | Taleya Williams /<br>Madeline Funaro |
| 3.0     | 2020.09.18 | Updated Scope<br>Updated Health and Safety<br>Updated references to include AOAC<br>Updated to include procedures for the digestion of fertilizer samples with Nitric and Hydrochloric acid                                                                                                            |               | William Meeks Jr.                    |
| 4.0     | 2020.10.26 | Modified temperature settings for heat blocks through out to reflect standard method.                                                                                                                                                                                                                  |               | William Meeks, Jr.                   |
| 5.0     | 2023.12.05 | Removed signature page<br>Removed steps related to final filtration of samples<br>Added external document numbers                                                                                                                                                                                      |               | Taleya Williams                      |

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 458

*Sample Preparation for Available Phosphorus Analysis [FM 505]*

Version: 2.0

Effective Date: 10/21/2020

Page 1 of 5

### **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for the preparation of available phosphorous (APA) nutrients in fertilizer samples.

### **2. SCOPE**

2.1. The procedure delineated below is applicable to fertilizer samples in the Bureau of Agricultural Environmental Laboratories (BAEL). This is a double extraction for the determination of available phosphoric acid in fertilizer samples that can be hydrolyzed to orthophosphate and have a range of 0 to 20%  $P_2O_5$ . The preparation is achieved by extracting the phosphate from the fertilizer sample with neutral ammonium citrate solution at 65°C. The range may be increased or decreased by changing the sample weight and/or the use of dilutions.

2.2. This method is applicable to fertilizer samples containing any amount of calcium present in the sample. The calcium is extracted from the sample with water at 65°C and prevented from forming insoluble calcium/magnesium phosphate complexes by the addition of nitric acid to the flask prior to the citrate extraction. This procedure must also be used when triple super phosphate or super phosphate are present in the sample.

### **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

### **4. OUTLINE of PROCEDURE**

- 4.1. 10.1 Collection of Weighed Samples
- 4.2. 10.2 Calcium Extraction
- 4.3. 10.3 Citrate Extraction
- 4.4. 10.4 Final Volume and Filtration
- 4.5. 10.5 Data Entry into Laboratory Information Management System
- 4.6. 10.6 Transfer of Extracts

### **5. INTERFERENCE**

5.1. N/A

### **6. APPARATUS and MATERIAL**

- 6.1. Centrifuge tube, disposable (15mL)
- 6.2. Filter disk, 7 cm Whatman No. 1 or equivalent
- 6.3. Flask, 500mL Kohlrausch (class "A")
- 6.4. Flask, Erlenmeyer (250mL)
- 6.5. Funnel, Buchner 7cm
- 6.6. Label maker
- 6.7. Laboratory Information Management System (LIMS)
- 6.8. Printer
- 6.9. Rack, tube
- 6.10. Screening column
- 6.11. Shaker Bath, Constant Temperature Circulating
- 6.12. Stopper, Two Hole for 500mL Kohlrausch Flask
- 6.13. Stopper, One Hole for 250mL Erlenmeyer flask
- 6.14. Thistle Tube, Polyethylene
- 6.15. Motorized Stirrer with ~1½ foot reach
- 6.16. Timer
- 6.17. pH meter

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**PREP 458**

*Sample Preparation for Available Phosphorus Analysis [FM 505]*

Version: 2.0

Effective Date: 10/21/2020

Page 2 of 5

## **7. REAGENTS and MEDIA**

- 7.1. Ammonium citrate dibasic, ACS grade or equivalent
- 7.2. Ammonium hydroxide, ACS grade or equivalent
- 7.3. Citric acid monohydrate, ACS grade or equivalent
- 7.4. Nitric acid, ACS grade or equivalent
- 7.5. Water, Deionized (DI) Type I

## **8. REAGENT PREPARATION GUIDELINES**

- 8.1. Nitric Acid 1:1
  - 8.1.1. Pour 1.0L of DI water into a 2L pre-labeled container. Place on a large magnetic stir plate. Turn on stir plate. Slowly add 1.0L of nitric acid to the container. Allow to stir for 10 minutes.
- 8.2. Neutral Ammonium Citrate Solution
  - 8.2.1. Add 400g of ammonium citrate dibasic to carboy, add approximately 20L of hot DI water. Add 114mL of ammonium hydroxide to citrate solution. Stir solution and allow to cool completely before testing. Adjust the pH to 7.00 by adding small increments of ammonium hydroxide or citric acid monohydrate.
- 8.3. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015.
- 8.4. The specific ingredient amounts can vary to accommodate the batch size.
- 8.5. Expiration Date: 6 months from date of preparation unless otherwise specified.
- 8.6. Store at room temperature.

## **9. ASSOCIATED DOCUMENTS**

- 9.1. ADMIN 023 Attachment P, Manual Result Entry
- 9.2. ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates
- 9.3. PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size
- 9.4. Paradigm 3, Document Control Software, Improvement Module

## **10. SPECIFIC PROCEDURES**

- 10.1. Collection of Weighed Samples
  - 10.1.1. Pick up pre-weighed sample batches from Completed Sample Weighing shelf.
  - 10.1.2. The batch size should not exceed the limit specified in PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size.
- 10.2. Calcium Extraction
  - 10.2.1. Preheat shaker bath to 65°C.
  - 10.2.2. Add 100mL of hot DI water to Erlenmeyer flask containing the sample, stopper with 1 hole stopper containing thistle tube. Place in 65°C shaker bath and set timer for 1 hour.
  - 10.2.3. Add 10mL of Nitric acid 1:1, to Kohlrausch flasks.
  - 10.2.4. Set up Buchner funnel – Kohlrausch flask – vacuum apparatus.
  - 10.2.5. After one hour, remove Erlenmeyer flask from shaker bath, rinse thistle tube and stopper with hot DI water into Erlenmeyer flask.
  - 10.2.6. Place filter disk into the Buchner funnel and wet with DI water.
  - 10.2.7. Transfer contents of Erlenmeyer flask to wetted filter disk in Buchner funnel – Kohlrausch flask – vacuum apparatus.
  - 10.2.8. Thoroughly rinse Erlenmeyer flask with hot DI water through filter disk in bottom of Buchner Funnel.
  - 10.2.9. Vacuum dry filter disk.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**PREP 458**

*Sample Preparation for Available Phosphorus Analysis [FM 505]*

Version: 2.0

Effective Date: 10/21/2020

Page 3 of 5

10.2.10. Perform one additional rinse of the Buchner funnel to ensure that any remaining sample is collected.

10.2.11. Carefully remove filter disk and place in Erlenmeyer flask.

**10.3. Citrate Extraction**

10.3.1. Add 100mL of hot Neutral Ammonium Citrate Solution, to Erlenmeyer flask.

10.3.2. Seal with 1 hole stopper containing a thistle tube. Place in 65°C shaker bath and set timer for 1 hour.

10.3.3. After a minimum of 90 minutes, remove Erlenmeyer flask from shaker bath, rinse thistle tube and stopper with hot DI water into Erlenmeyer flask.

10.3.4. Place filter disk into the Buchner funnel and wet with DI water

10.3.5. Transfer contents of Erlenmeyer flask to wetted filter disk in Buchner funnel-Kohlrusch flask – vacuum apparatus.

10.3.6. Thoroughly rinse Erlenmeyer flask with hot DI water through filter disk in bottom of Buchner Funnel.

10.3.7. Vacuum dry filter disk.

10.3.8. Perform one additional rinse of the Buchner funnel to ensure that any remaining sample is collected.

10.3.9. Remove the filter disk and discard.

10.3.10. Remove top of funnel.

10.3.11. Rinse the bottom funnel into the Kohlrusch flask.

**10.4. Final Volume and Filtration**

10.4.1. Allow extracts to cool for approximately 30 minutes, then dilute to 500mL with DI water.

10.4.2. Shake well.

**NOTE:** 30-45 seconds of vigorous shaking may be required to mix sample completely.

10.4.3. Filter approximately 12-14mL of sample extract into a disposable tube using a screening column.

**10.5. Data Entry into Laboratory Information Management System (LIMS)**

10.5.1. Using the applicable LIMS result entry template, complete the APA flask volume test codes in the LIMS. Refer to ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates, for a list of Result Entry templates and Attachment P, Manual Result Entry, for result entry instructions.

10.5.2. Print a copy of the LIMS result entry screen to include as the paperwork for the prep batch.

**10.6. Transfer of Extracts**

10.6.1. Transfer all extracts and paperwork to designated area in building 5 using a laboratory approved sample carrier.

**11. WASTE MANAGEMENT**

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material                           | Disposal Mechanism                                   |
|------------------------------------|------------------------------------------------------|
| Remaining Sample Extract from 10.4 | Can be discarded down drain with sufficient flushing |
| Screening Columns and Filter Disks | Discarded in regular waste bin                       |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 458

*Sample Preparation for Available Phosphorus Analysis [FM 505]*

Version: 2.0

Effective Date: 10/21/2020

Page 4 of 5

## 12. REFERENCES

- 12.1. Available Phosphorus (Dual Extraction) – Flow Injection. FM 505, (07/15/2009)  
 12.2. APA (P<sub>2</sub>O<sub>5</sub>) Extraction –Citrate Soluble for Flow Injection. FSFL-SOP 500, (02/03/2006).

## 13. VERSION HISTORY

| Version | Date       | Description                                                                                                                                             | Author      | Editor          |
|---------|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|-----------------|
| 1.0     | 10.12.2016 | Original Version                                                                                                                                        | Sylvia Mehr |                 |
| 2.0     | 04.15.2019 | Updated outline format<br>Revised scope to remove FM504<br>Added Section 8<br>Added Section 11<br>Minor grammatical fixes and clarifications throughout |             | Taleya Williams |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 458

*Sample Preparation for Available Phosphorus Analysis [FM 505]*

Version: 2.0

Effective Date: 10/21/2020

Page 5 of 5

**14. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director

Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director

Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer

Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Wet Chemistry and Preparation Manager / Editor

Date\_\_\_\_\_

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 459

*Sample Preparation for Total Phosphorus Analysis*  
*[FM 521]*

Version: 4.0

Effective Date: 02/01/2024

Page 1 of 5

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the preparation of total phosphorous (TP) nutrients in fertilizer samples.

## 2. SCOPE

2.1. The procedure delineated below is applicable to fertilizer samples in the Bureau of Agricultural Environmental Laboratories (BAEL). This preparation is for the analysis of total phosphorous by flow injection analysis. Phosphorus that can be hydrolyzed to orthophosphate with a concentration range of 0 to 20%  $P_2O_5$  may be determined by this method. The range may be increased by decreasing sample weight and/or the use of dilutions. Laboratory Information Management System (LIMS) suggested weights must be followed. This method is applicable to fertilizer samples containing any amount of calcium.

2.2. Available Phosphorous: This method is also performed on a sample that has failed the normal available extraction methods. The following APA source materials will almost always fail the citrate soluble APA extractions: bone meal, basic slag, ammonium polyphosphate, and phosphorous acid (not to be confused with phosphoric acid).

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

3.2. Acid digestions must be performed under a fume hood.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Weigh Sample
- 4.2. 10.2 Handling Pre-Weighed Samples
- 4.3. 10.3 Acid Digestion
- 4.4. 10.4 Citrate Extraction
- 4.5. 10.5 Final Volume
- 4.6. 10.6 Data Entry into Laboratory Information Management System
- 4.7. 10.7 Transfer of Extracts

## 5. INTERFERENCE

5.1. N/A

## 6. APPARATUS and MATERIAL

- 6.1. Centrifuge tube, disposable (15mL)
- 6.2. Erlenmeyer flask, 250mL
- 6.3. Filter disk, 7cm Whatman No. 1 or equivalent
- 6.4. Flask, 500mL Kohlrausch (class "A")
- 6.5. Funnel, Buchner 7cm
- 6.6. Hot plate
- 6.7. Laboratory Information Management System (LIMS)
- 6.8. pH meter

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 459

*Sample Preparation for Total Phosphorus Analysis*  
*[FM 521]*

Version: 4.0

Effective Date: 02/01/2024

Page 2 of 5

- 6.9. Rack, tube
- 6.10. Shaker bath, constant temperature circulating
- 6.11. Stopper, one hole for 250mL Erlenmeyer flask
- 6.12. Stopper, two hole for 500mL Kohlraush Flask
- 6.13. Thistle tube, polyethylene
- 6.14. Timer

## **7. REAGENTS and MEDIA**

- 7.1. Ammonium citrate dibasic, ACS grade or equivalent
- 7.2. Ammonium hydroxide, ACS grade or equivalent
- 7.3. Citric acid monohydrate, ACS grade or equivalent
- 7.4. Water, deionized (DI) Type II
- 7.5. Concentrated Hydrochloric acid, (HCl) ACS grade or equivalent
- 7.6. Concentrated Nitric acid, (HNO<sub>3</sub>) ACS grade or equivalent

## **8. REAGENT PREPARATION GUIDELINES**

### **8.1. 1:1 Nitric Acid**

8.1.1. Directions: Using a clean 1000mL graduated cylinder, measure 1000mL of DI Water and pour into the 2L storage container. Using a clean 1000mL graduated cylinder measure 1000mL of concentrated HNO<sub>3</sub>. Slowly add the HNO<sub>3</sub> to the DI water, cap the container and mix thoroughly using stir bar for approximately 10 minutes. Diluting the HNO<sub>3</sub> can cause heat; be careful not to add the acid too fast.

### **8.2. Neutral Ammonium Citrate Solution**

Directions: Add 600g of ammonium citrate dibasic to a 30L vessel. Add approximately 24L of hot DI water to the vessel and turn on stir bar. Add 100mL of ammonium hydroxide to the citrate solution. Allow solution to stir for approximately 2 hours. Adjust the pH to 7.00 by adding small increments of ammonium hydroxide or citric acid monohydrate.

8.3. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015

8.4. The specific ingredient amounts can vary to accommodate the batch size.

8.5. Expiration Date: 12 months from date of preparation

8.6. Store at room temperature

## **9. ASSOCIATED DOCUMENTS**

- 9.1. ADMIN 023 Attachment P, Manual Result Entry
- 9.2. ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates
- 9.3. PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 459

*Sample Preparation for Total Phosphorus Analysis*  
*[FM 521]*

Version: 4.0

Effective Date: 02/01/2024

Page 3 of 5

## 10. SPECIFIC PROCEDURES

### 10.1. Weigh Sample

10.1.1. Weigh samples as per the table below:

| Guarantee | Weight |
|-----------|--------|
| <5%       | 3g     |
| 5-10%     | 2g     |
| 10-20%    | 1g     |
| 20-30%    | 0.5g   |
| >30%      | 0.25g  |

### 10.2. Handling Pre-Weighed Samples

10.2.1. If samples are pre-weighed, collection of Weighed Samples

10.2.2. Pick up pre-weighed sample batches from Completed Sample Weighing shelf.

10.2.3. The batch size should not exceed the limit specified in PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size.

### 10.3. Acid Digestion

10.3.1. Into 250mL Erlenmeyer flask, with weighed sample, add 20mL of Nitric Acid 1:1 and 5mL of concentrated HCl.

10.3.2. Place the sample with solution on hot plate, boil approximately 15 minutes or until solution is clear. Do not let the sample go dry.

10.3.3. Allow the sample to cool, then proceed to 10.4, Citrate Extraction.

### 10.4. Citrate Extraction

10.4.1. Preheat the shaker bath to 65°C.

10.4.2. Add 100mL of hot (approximately 65°C) Neutral Ammonium Citrate Solution to Erlenmeyer flask.

10.4.3. Seal with 1 hole stopper containing a thistle tube. Place in 65°C shaker bath and set timer for 15 minutes.

10.4.4. After 15 minutes remove Erlenmeyer flask from shaker bath, rinse thistle tube and stopper with hot DI water into Erlenmeyer flask.

10.4.5. Place filter disk into the Buchner funnel and wet with DI water.

10.4.6. Transfer contents of Erlenmeyer flask onto wetted filter disk in Buchner funnel-Kohlrausch flask – vacuum apparatus.

10.4.7. Thoroughly rinse Erlenmeyer flask with hot DI water through filter disk.

10.4.8. Vacuum dry filter disk.

10.4.9. Perform one additional rinse of the Buchner funnel to ensure that any remaining sample is collected.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 459

*Sample Preparation for Total Phosphorus Analysis*  
*[FM 521]*

Version: 4.0

Effective Date: 02/01/2024

Page 4 of 5

10.4.10. Remove the filter disk and discard.

10.4.11. Remove top of funnel.

10.4.12. Rinse the bottom of the funnel into the Kohlrausch flask.

10.5. Final Volume

10.5.1. Allow extracts to cool to room temperature and dilute to 500mL.

10.5.2. Mix the contents by thoroughly inverting the flask multiple times and allow to settle.

10.5.3. Pour approximately 12-14mL of sample extract into a disposable centrifuge tube.

10.6. Data Entry into Laboratory Information Management System

10.6.1. Using the applicable Laboratory Information Management System (LIMS) result entry template, complete the APA flask volume test codes in the LIMS. Refer to ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates, for the list of result entry templates and Attachment P, Manual Result Entry, for result entry instructions.

10.6.2. Print a copy of the LIMS result entry screen to include as the paperwork for the prep batch.

10.7. Transfer of Extracts

10.7.1. Transfer all extracts and paperwork to designated area in building 5 using a laboratory approved sample carrier.

10.7.2. Record initials and date on the paperwork.

## 11. WASTE MANAGEMENT

11.1. Remaining sample extract from 10.5 can be discarded down the drain with sufficient flushing using tap water.

## 12. REFERENCES

12.1. E00164, Total Phosphorus -Flow Injection. FM 521, (07/15/2009)

12.2. E00146, Total P<sub>2</sub>O<sub>5</sub> Sample Preparation, FSFL-SOP 501, (02/03/2006)

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 459  
*Sample Preparation for Total Phosphorus Analysis*  
*[FM 521]*

Version: 4.0

Effective Date: 02/01/2024

Page 5 of 5

**13. VERSION HISTORY**

| <b>Version</b> | <b>Date</b> | <b>Description</b>                                                                                                                                                                                                                                                                                                                                      | <b>Author</b> | <b>Editor</b>   |
|----------------|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------|
| 1.0            | 10.12.2016  | Original Version                                                                                                                                                                                                                                                                                                                                        | Sylvia Mehr   |                 |
| 2.0            | 12.05.2018  | Updated Outline Format<br>Moved Reagents back into SOP (Section 8)<br>Minor grammatical fixes and clarifications                                                                                                                                                                                                                                        |               | Taleya Williams |
| 3.0            | 1/13/2021   | Clarified language in section 2.2<br>Removed section 10.3 Acid Digestion for Deficient Samples Containing Basic Slag, Iron or Aluminum Phosphate, or Liquid Samples and renumbered sections accordingly.<br>Removed Reagent 01 and Reagent 02 through to align with other SOPs<br>Replaced Elizabeth Schneider with Quentin Cooper in Signature History |               | Taleya Williams |
| 4.0            | 02.01.2024  | Removed signature page<br>Removed steps related to final filtration<br>Added weight table<br>Added external document numbers                                                                                                                                                                                                                            |               | Taleya Williams |

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

### **PREP 461**

#### *Sample Preparation for Available Phosphorus and Soluble Potash in Fertilizer*

Version: 2.0

Effective Date: 10/16/2023

Page 1 of 4

## **1. PURPOSE**

1.1 This standard operating procedure (SOP) provides instruction for the preparation of available phosphorous (APA) and soluble potash ( $K_2O$ ) nutrients in fertilizer samples.

## **2. SCOPE**

2.1 The procedure delineated below is applicable to fertilizer samples in the Bureau of Agricultural Environmental Laboratories (BAEL). APA ( $P_2O_5$ ) and  $K_2O$  are extracted from fertilizer with a citrate-EDTA while shaking at 65°C. After the extract is brought to volume with deionized water, the extract is analyzed by ICP-OES.

2.2 This method is adapted from AOAC Method 2015.18. The method described in this SOP deviates from Method 2015.18 by bringing the samples to a final volume of 500 mL instead of 250 mL, as stated in Method 2015.18. Thus, the citrate-EDTA of the sample is diluted prior to analysis, eliminating the need to dilute the sample inline by increasing the tubing size of the buffer/internal standard peristaltic tube and reducing the tubing size of the sample peristaltic tube as described in Method 2015.18. Sample weight was increased from 0.5 g to 0.75 g to achieve required detection limit.

## **3. HEALTH and SAFETY**

3.1 A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## **4. OUTLINE of PROCEDURE**

- 4.1. 10.1 Extraction
- 4.2. 10.2 Data Entry into Laboratory Information Management System
- 4.3. 10.3 Transfer of Extracts

## **5. INTERFERENCE**

5.1. Fertilizers containing phosphite will also be extracted by this method.

5.2. Due to the susceptibility of the citrate-EDTA solution to microbial degradation, the citrate-EDTA solution has a 2-week expiration when stored in the dark and a 16-hour expiration for sample extracts. An inhouse longevity study demonstrated microbial degradation could be mitigated by utilizing refrigeration. The expiration dates found in this SOP have been extended based on this study.

## **6. APPARATUS and MATERIAL**

- 6.1. Centrifuge tube, disposable (15-mL)
- 6.2. Flask, 500 mL Kohlrausch (class "A") or Volumetric Flask
- 6.3. Flask, Erlenmeyer (250 mL)
- 6.4. Label maker
- 6.5. Magnetic stir bar
- 6.6. pH meter
- 6.7. Rack, tube
- 6.8. Shaker Bath, Constant Temperature Circulating
- 6.9. Water Bath

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**PREP 461**

*Sample Preparation for Available Phosphorus and Soluble Potash in  
Fertilizer*

Version: 2.0

Effective Date: 10/16/2023

Page 2 of 4

- 6.10. Beaker
- 6.11. Balance Stir Plate
- 6.12. Stopper, rubber for 500 mL Kohlrausch Flask
- 6.13. Thistle Tube, Polyethylene
- 6.14. Time

## **7. REAGENTS and MEDIA**

- 7.1. onium citrate dibasic, ACS grade or equivalent
- 7.2. DTA, disodium salt, dihydrate, >99%
- 7.3. Ammonium hydroxide, ACS grade or equivalent
- 7.4. Nitric acid, ACS grade or equivalent
- 7.5. Water, Deionized (DI) Type I

## **8. REAGENT PREPARATION GUIDELINES**

### **8.1. Citrate-EDTA Extraction Solution**

8.1.1. Directions: Add 25 g disodium EDTA and 50 g dibasic ammonium citrate to a 2L volumetric flask containing ~1500 mL DI water. Adjust to near neutral by adding 30 mL of a solution of ammonium hydroxide-water (1+1, v/v), under a hood. (Volumes are scalable as needed.)

8.1.2. Using a pH electrode while stirring, adjust the final pH to 7.00 using a pH by adding the ammonium hydroxide-water (1+1, v/v) solution, or 1:1 Nitric Acid dropwise. After a stable pH of :7.00 ±0.02 has been attained, dilute to volume with DI water and mix.

8.1.3. Citrate-EDTA solution has a **2-months expiration** when stored refrigerated.

## **9. ASSOCIATED DOCUMENTS**

- 9.1. ADMIN 023 Attachment P, Manual Result Entry
- 9.2. DMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates
- 9.3. PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size

## **10. SPECIFIC PROCEDURES**

### **10.1. Extraction**

10.1.1. Fill shaker bath until water level covers the flask holders.

10.1.2. Turn on shaker bath and set temperature to 65° C.

Note: It takes roughly 90 min to heat to 65° C from room temperature.

10.1.3. Fill the water bath approximately two-thirds full.

Note:

10.1.4. Turn on water bath and set temperature 65° C.

10.1.5. Heat citrate-EDTA to reach 65 ±2°C

Note: It takes at least 30 min for the citrate-EDTA to reach 65 ±2°C.

10.1.6. Weigh 0.75 g of fertilizer into a separate 250-mL Erlenmeyer flask. See PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size for batch size and Laboratory Standard information.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**PREP 461**

*Sample Preparation for Available Phosphorus and Soluble Potash in Fertilizer*

Version: 2.0

Effective Date: 10/16/2023

Page 3 of 4

10.1.7. Add 100 mL of the 65 ±2°C preheated citrate-EDTA extraction solution to the 250-mL Erlenmeyer flask containing the sample, stopper with a rubber stopper with thistle tube and place in the 65 ±2°C shaker bath. Shake at 200 reciprocations/min. for 1 hour.

10.1.8. After one hour, remove the 250-mL Erlenmeyer flask from the shaker bath, **rinse rubber stopper with thistle tube with DI water into the 250-mL Erlenmeyer flask,** and cool to room temperature.

10.1.9. Transfer the extract into a 500-mL Kohlrausch flask **or volumetric flask,** rinsing the flask with DI water a minimum of 3 times to ensure a complete transfer. Dilute to volume with DI water; stopper and mix by inverting 10 to 12 times.

10.1.10. Transfer approximately 15 mL of the extract into a disposable 15-mL centrifuge tube. Filter or centrifuge any extract containing suspended debris.

10.1.11. Extracts must be stored refrigerated and are to be analyzed within **3 days.** **Note:** **Storage location is in Building 5, Refrigerator 6**

**10.2. Data Entry into Laboratory Information Management System (LIMS)**

10.2.1. Using the applicable LIMS result entry template, complete the P&K flask volume test codes in the LIMS. Refer to ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates, for a list of Result Entry templates and Attachment P, Manual Result Entry, for result entry instructions.

10.2.2. Print a copy of the LIMS result entry screen to include as the paperwork for the prep batch.

**10.3. Transfer of Extracts**

10.3.1. Transfer all samples and paperwork to the designated area in building 5 using a laboratory-approved sample carrier.

**11. REFERENCES**

11.1 *Official Methods of Analysis of AOAC INTERNATIONAL* (2019) 21st Ed., AOAC INTERNATIONAL, Rockville, MD, USA, **Method 2015.18**, [www.eoma.aoac.org](http://www.eoma.aoac.org), **External Document 00135.**

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

|                                                                                                        |                            |             |
|--------------------------------------------------------------------------------------------------------|----------------------------|-------------|
| <b>PREP 461</b><br><i>Sample Preparation for Available Phosphorus and Soluble Potash in Fertilizer</i> |                            |             |
| Version: 2.0                                                                                           | Effective Date: 10/16/2023 | Page 4 of 4 |

**12. VERSION HISTORY**

| Version | Date       | Description                                                                                                                                                                                                                                                                                                                                                                   | Author        | Editor                              |
|---------|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------------------------------------|
| 1.0     | 2021.09.17 | Original Version                                                                                                                                                                                                                                                                                                                                                              | William Meeks |                                     |
| 2.0     | 10.16.2023 | Section 6 updated Apparatus and Material list.<br>Section 8.1.2. added option to use Nitric Acid to adjust pH.<br>Added section 10.1.4. and 10.1.5. and related note.<br>Added instructions for the rubber stopper with thistle to section 10.1.8.<br>Added note to section 10.1.11.<br>Added external document number to reference in section 11.1<br>Removed signature page |               | Taleya Williams<br>Jonnie Carpenter |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 462

#### *Preparation for Soluble Cu, Fe, Mn, and Zn* *[FM 830] [FM 831] [FM 832] [FM 833]*

Version: 2.0

Effective Date: 05/10/2023

Page 1 of 6


## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the determination of soluble iron, manganese, copper, and zinc in fertilizers by Atomic Absorption analysis.

## 2. SCOPE

2.1. The procedure delineated below is applicable to fertilizer that is analyzed in the Bureau of Agricultural Environmental Laboratories (BAEL).

2.2. Soluble metal samples are extracted at acidic pH. Depending on the analyzed metal, the sample is solubilized in a solution of a specific pH, at which metal oxides are slightly soluble and are not a source of significant positive bias (ferric oxides, 1-2% solubility; manganese oxides, 1-1.5% solubility; copper oxides, 4-5% solubility). The respective extraction solutions also ensure that all the metal sulfate sources are fully recoverable. In samples containing high phosphate (%P<sub>2</sub>O<sub>5</sub> guarantee ≥ 10), from monoammonium, diammonium or polyammonium phosphate sources, the soluble metals can be precipitated as metal ammonium phosphate. The most probable reaction between the soluble micronutrient and the phosphate salt is a recombination reaction which forms a stable (reciprocal) salt pair. This reaction proceeds in an irreversible manner to near completion, in which the salt pair is of much lower solubility than the original compounds.



This reaction has been documented in Lehr (1972), where he listed several such reactions that may take place during the manufacture or storage of fertilizers. Many reactions require only water which is attracted to the hygroscopic fertilizer ingredients. Thus, samples must be protected from moisture, before extraction, and must be analyzed as soon as possible following extraction. After appropriate dilution, samples are analyzed by atomic absorption (A.A.) spectroscopy.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Collection of Weighed Samples
- 4.2. 10.2 Extraction of Samples
- 4.3. 10.3 Transfer of Samples
- 4.4. 10.4 Data Entry into Laboratory Information Management System (LIMS)
- 4.5. 10.5 Delivery

## 5. INTERFERENCE

5.1. Polyphosphate that is incompletely solubilized and superphosphates interfere with SCu and SZn extraction.

## 6. APPARATUS and MATERIAL

- 6.1. Graduated Centrifuge Tube, 15mL
- 6.2. Syringe with Luer-Lok tips (without needle), 20mL
- 6.3. Analytical Balance (accuracy to 0.0001g)
- 6.4. Aqueous 0.45µm nonsterile Syringe Filters 25mm disc

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**PREP 462**

*Preparation for Soluble Cu, Fe, Mn, and Zn*  
*[FM 830] [FM 831] [FM 832] [FM 833]*

Version: 2.0

Effective Date: 05/10/2023

Page 2 of 6

- 6.5. Bin to accommodate 25 200mL volumetric flasks
- 6.6. Filter screen column, 200 to 300µm or equivalent
- 6.7. Glass pipets (Class A)
- 6.8. Green rubber stoppers, size 3
- 6.9. Magnetic Stir Bar (approximately 3 inches)
- 6.10. Nalgene Container, 20L
- 6.11. Parafilm
- 6.12. pH meter (accuracy to 0.01 pH units)
- 6.13. 100mL graduated cylinder
- 6.14. Volumetric flasks Class A 200mL
- 6.15. Centrifuge

## **7. REAGENTS and MEDIA**

- 7.1. Sodium Hydroxide Solution, 1N, purchased
- 7.2. 2,2-Bis(hydroxymethyl)-2,2',2''-nitrilotriethanol, C<sub>8</sub>H<sub>19</sub>NO<sub>5</sub>, 99%+ (BIS-TRIS)
- 7.3. Solution for pH meter calibration, pH 4, Certified
- 7.4. Solution for pH meter calibration, pH 7, Certified
- 7.5. Hydrochloric acid (HCl) Certified A.C.S. grade or equivalent
- 7.6. Potassium Hydrogen Phthalate (KHP) Certified A.C.S grade or equivalent
- 7.7. Sodium Hydroxide (NaOH), pellets
- 7.8. SRM695 (NIST)
- 7.9. Sulfuric Acid, Certified A.C.S grade or equivalent
- 7.10. Water, Deionized (DI) Type I

## **8. REAGENT PREPARATION GUIDELINES**

- 8.1. Soluble Iron Buffer Solution
  - 8.1.1. Directions: 200mL of Sulfuric Acid, Certified A.C.S grade or equivalent in a 20L calibrated reservoir containing 10L of DI Water. Fill up to the 20L mark and mix well.
- 8.2. 1N Sodium Hydroxide Solution
  - 8.2.1. Directions: Add 8g of NaOH pellets to approximately 100mL of DI Water in a 200mL volumetric flask. Cool to room temperature. Bring to volume with DI Water and mix well. A commercially made 1N NaOH solution can be used.
- 8.3. Soluble Manganese Buffer Solution (0.05M KHP)
  - 8.3.1. Preparation of 1L: Dissolve 10g of KHP in 800mL of DI Water. Adjust pH (4.95 to 5.0) by adding 20mL of 1.0N Sodium Hydroxide and bring to volume with DI Water; mix well.
  - 8.3.2. Preparation of 20L: Dissolve 100g KHP in approximately 3L of DI Water in a 4L Beaker with handle. Transfer to the 20-L reservoir containing 10L DI Water. Repeat with 100g of KHP again, and transfer to the 20-L container. Fill to the 19L mark. Add 390mL of 1.0N Sodium Hydroxide, mix well and ensure that pH is adjusted (4.95-5.0). Fill with DI Water to the 20L mark; mix well.
- 8.4. Soluble Zn/Cu Buffer Solution (0.24M BIS-TRIS)
  - 8.4.1. Preparation of 1L: Dissolve 50g of 2,2-Bis(hydroxymethyl)-2,2',2''-nitrilotriethanol in 700mL of DI water. Add 4mL of concentrated HCl and mix well. Verify the final pH of the buffer solution (5.95-6). Adjust final pH by adding 1:1 HCl, or 1:1 NaOH, if needed. Bring solution to 1L with DI Water and mix well.
  - 8.4.2. Preparation of 20L: Dissolve 500g of BIS-TRIS in approximately 3L of DI Water. Transfer to the 20L reservoir containing 10L of DI Water. Repeat with 500g of TRI-BIS. Mixing well, add 250 mL of 1:1 HCL. Adjust final pH to 6.0 by adding concentrated HCL or 1:1 NaOH, if needed, and bring solution to 20L with DI water. Mix well.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 462

#### *Preparation for Soluble Cu, Fe, Mn, and Zn* *[FM 830] [FM 831] [FM 832] [FM 833]*

Version: 2.0

Effective Date: 05/10/2023

Page 3 of 6

#### 8.5. 1:1 Hydrochloric Acid/Water (v/v) Solution

8.5.1. Directions: Add 250mL hydrochloric acid to approximately 200mL of DI Water in a 500mL volumetric flask. Bring to volume with DI Water, and mix well

**Caution:** Always add the acid to the water. Adding water to the acid causes the acid to spatter!

#### 8.6. 1:1 Sodium Hydroxide/water (w/v) Solution

8.6.1. Directions: Add 250g of NaOH to approximately 200mL DI Water in a 500mL volumetric flask, bring to volume DI Water, and mix well.

8.7. Label the container appropriately, Label L003 indicates required information. Add relevant information to the Reagent Preparation, Form F00015.

8.8. The specific ingredient amounts can vary to accommodate the batch size.

8.9. Expiration Date: 12 months from date of preparation, unless otherwise specified.

8.10. Store at room temperature

### 9. ASSOCIATED DOCUMENTS

9.1. ADMIN 023, Attachment AA, Fertilizer QC Loc Code and Result Entry Templates.

9.2. F00015, Reagent Preparation

9.3. L003, Reagent

9.4. PREP 453, Weighing of Fertilizer Sample

9.5. PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size

### 10. SPECIFIC PROCEDURES

#### 10.1. Collection of Weighed Samples

10.1.1. Samples are weighed according to PREP 453, Weighing of Fertilizer Sample.

10.1.2. The recommended sample weights are found in the following tables:

| SCu, SFe, SMn Suggested Weight |                    |
|--------------------------------|--------------------|
| % Guarantee                    | Weight (g)         |
| <0.008                         | 5                  |
| 0.008 - 0.08                   | 0.04 / % Guarantee |
| 0.08 - 0.8                     | 0.2 / % Guarantee  |
| >0.8                           | 0.25               |

| SZn Suggested Weight |                    |
|----------------------|--------------------|
| % Guarantee          | Weight (g)         |
| <0.008               | 5                  |
| 0.008 - 0.16         | 0.04 / % Guarantee |
| >0.16                | 0.25               |

Where,

weight factor = (target ppm) (extraction vol) / (% conversion factor)

0.04 = (2 ppm) (200 mL) / (10000)

0.2 = (10 ppm) (200 mL) / (10000)

10.1.3. Samples are weighed into 200mL volumetric flasks.

10.1.4. Weigh 0.5g SRM695 (NIST) to be used as spike sample.

10.1.5. The batch size should not exceed the limit specified in PREP 453, Attachment A, Fertilizer Laboratory Standard and Batch Size.

10.1.6. Collect batch samples from the properly identified bin(s).

#### 10.2. Extraction of Samples

10.2.1. Using a graduated cylinder, add 100mL of the appropriate buffer to the weighed samples:

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 462

#### *Preparation for Soluble Cu, Fe, Mn, and Zn* *[FM 830] [FM 831] [FM 832] [FM 833]*

Version: 2.0

Effective Date: 05/10/2023

Page 4 of 6

- 10.2.1.1. Soluble Iron Buffer, for Fe samples
- 10.2.1.2. Soluble Manganese Buffer Solution (0.05M KHP), for Mn samples
- 10.2.1.3. Soluble Zn/Cu Buffer Solution (0.24M BIS-TRIS), for Soluble Copper/Zinc **Samples**
- 10.2.2. Bring to volume with DI water.
- 10.2.3. Shake the flask for a minimum of 30 seconds, vigorously, and allow the sample to stand for a least 45 minutes.
- 10.2.4. Filter **all cloudy samples or** liquid samples containing “free-floating” particles and undissolved particles **using an aqueous 25mm, 0.45µm non-sterile Syringe Filter .**
- 10.2.5. Verify pH of buffer solution.
- 10.3. Transfer of Samples
  - 10.3.1. Assemble a set of 15mL Centrifuge Tubes, one per sample.
  - 10.3.2. Place the sample tubes on a suitable test tube rack, in a numerical sequence.
  - 10.3.3. Samples can be centrifuged, or syringe filtered.**
  - 10.3.4. Label sample tubes with the sample ID number.
  - 10.3.5. Pour approximately 14mL of each sample into the test tubes. This allows for a gap between the top of the tube and the sample surface that prevents spills.
  - 10.3.6. Cover the samples with parafilm and make a seal to avoid spills or cross contamination.
  - 10.3.7. Samples must be analyzed on the same day as extracted or must be filtered immediately following preparation. Samples extracted for SZn must be analyzed within one hour or filtered immediately following preparation.**
- 10.4. Data Entry into Laboratory Information Management System (LIMS)
  - 10.4.1. Log into Laboratory Information Management System (LIMS).
  - 10.4.2. Double click the “Results” folder and single click “Results Entry”.
  - 10.4.3. Refer to ADMIN 023, Attachment AA, Fertilizer QC Loc Code and Result Entry Templates, for specific result entry template names for each metal:
    - 10.4.3.1. Soluble Iron: FERT\_AA\_SOLUBLE\_IRON
    - 10.4.3.2. Soluble Manganese: FERT\_AA\_SOLUBLE\_MANGANESE
    - 10.4.3.3. Soluble Copper: FERT\_AA\_SOLUBLE\_COPPER
    - 10.4.3.4. Soluble Zinc: FERT\_AA\_SOLUBLE\_ZINC
  - 10.4.4. For each sample, enter the flask volume and save.
  - 10.4.5. Reopen the template to verify that the data was saved.
  - 10.4.6. Print and close the template.
  - 10.4.7. Refer to ADMIN 023, Attachment AA, Fertilizer QC Loc Code and Result Entry Templates, for specific result entry template names. For example: select the:
    - 10.4.7.1. Soluble Iron: FERT\_PRINT\_AA\_SFE
    - 10.4.7.2. Soluble Manganese: FERT\_PRINT\_SMN
    - 10.4.7.3. Soluble Copper: FERT\_PRINT\_SCU
    - 10.4.7.4. Soluble Zinc: FERT\_PRINT\_SZN
  - 10.4.8. Print and close the template.
  - 10.4.9. Indicate the QC sample number on the printouts, initial and date.
- 10.5. Delivery
  - 10.5.1. Deliver all samples and paperwork to the designated area, using a laboratory-approved sample carrier.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

PREP 462

*Preparation for Soluble Cu, Fe, Mn, and Zn*  
*[FM 830] [FM 831] [FM 832] [FM 833]*

Version: 2.0

Effective Date: 05/10/2023

Page 5 of 6

### 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material        | Disposal Mechanism                |
|-----------------|-----------------------------------|
| Excess Extracts | Dispose down drain while flushing |
| Excess Buffers  | Dispose down drain while flushing |

### 12. REFERENCES

- 12.1. Soluble Iron FM-832 - QC Corporation – Baltimore, M.D. (Method # AWWA B402-90) FSFL-SOP-508
- 12.2. Soluble Manganese FM-830 – AOAC 16<sup>th</sup> Edition, Method 972.03 FSFL-SOP-508
- 12.3. Soluble Copper FM-831 – FSFL SOP-508
- 12.4. Soluble Zinc FM-833 – AOAC 16<sup>th</sup> Edition, Method 972.03 FSFL SOP-508
- 12.5. Preparation and Analysis of Soluble and Chelated Micro Nutrient Samples By AA FSFL-SOP 509.

### 13. VERSION HISTORY

| Version | Date                     | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                              | Author                          | Editor                          |
|---------|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|---------------------------------|
| 1.0     | 11.01.2019               | Original Version                                                                                                                                                                                                                                                                                                                                                                                                                                                         | Jack Andreu<br>Christian Amason |                                 |
| 2.0     | 02.11.2021<br>02.22.2023 | <p>Added section 5.1.<br/> Added section 6.15.<br/> Expanded upon instructions in section 8.4.2.<br/> Added recommended sample weight table to section 10.1.2.<br/> Section 10.1.2.3. corrected buffer to samples<br/> Section 10.2.4. expanded upon criteria for filtering samples.<br/> Section 10.3.3 updated step to align with currently laboratory activities.<br/> Section 10.3.7. added instructions for samples extracted for SZn.<br/> Added section 12.5.</p> |                                 | Zachary Tower<br>Patrizia Lemma |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 462

*Preparation for Soluble Cu, Fe, Mn, and Zn*  
*[FM 830] [FM 831] [FM 832] [FM 833]*

Version: 2.0

Effective Date: 05/10/2023

Page 6 of 6

**14. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer  
Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Sample Preparation and Wet Chemistry Manager  
Date\_\_\_\_\_

\_\_\_\_\_  
Zachary Tower  
Hemp Analysis Manager/ Editor  
Date\_\_\_\_\_

\_\_\_\_\_  
Patricia Lemma  
Chemist Specialist/ Editor  
Date\_\_\_\_\_

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 464

#### *Sample preparation for Chelated Fe, Mn, Cu, Zn [FM 840] [FM 841]*

Version: 2.0

Effective Date: 10/16/2023

Page 1 of 4

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the preparation of samples for the analysis of Chelated Iron, Manganese, Copper and Zinc. These samples subsequently undergo Atomic Absorption analysis.

## 2. SCOPE

2.1. The procedure delineated below is applicable to fertilizer that is analyzed in the Bureau of Agricultural Environmental Laboratories (BAEL).

2.2. In this method, a sodium acetate buffer solution is added to samples to adjust their pH between 3 and 4. The acidic extraction solution compensates for the basic condition brought about when large amounts of dolomite or limestone are used as fillers. It also prevents irreversible ligand exchange reactions of the iron chelates with copper, manganese, zinc and other metal ions, and formation of insoluble iron salts. In the absence of a complexing agent, Ferric ( $\text{Fe}^{+3}$ ) iron forms insoluble salts at  $\text{pH} > 3$ . The chelated metal is quantitated by atomic absorption spectroscopy analysis.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Collection of Weighed Samples
- 4.2. 10.2 Extraction of Samples
- 4.3. 10.3 Transfer of Samples
- 4.4. 10.4 Data Entry into Laboratory Information Management System (LIMS)
- 4.5. 10.5 Delivery

## 5. INTERFERENCE

- 5.1. N/A

## 6. APPARATUS and MATERIAL

- 6.1. Graduated Polystyrene Centrifuge Tube, 15mL
- 6.2. Syringe with Luer-Lok tips (without needle), 20mL
- 6.3. Analytical Balance (accuracy to 0.0001g)
- 6.4. Aqueous 0.45 $\mu\text{m}$  nonsterile Syringe Filters, 25mm disc
- 6.5. Bin to accommodate up to 25 200mL volumetric flasks
- 6.6. Glass pipets (Class A)
- 6.7. Green rubber stoppers, size 3
- 6.8. Magnetic Stir Bar, approximately 3 inches
- 6.9. Nalgene Container, 20L
- 6.10. Parafilm
- 6.11. pH meter (accuracy to 0.01 pH units)
- 6.12. Volumetric flasks Class A, 100, 200, and 1000mL

## 7. REAGENTS and MEDIA

- 7.1. Diammonium phosphate (DAP),  $(\text{NH}_4)_2\text{HPO}_4$ , solid
- 7.2. Glacial Acetic Acid ( $\text{CH}_3\text{COOH}$ ), **Caution:** Avoid breathing vapors and skin contact. Use in a fume hood and wear protective equipment.
- 7.3. Sodium Acetate ( $\text{C}_2\text{H}_3\text{NaO}_2$ )
- 7.4. 5% Sodium hypochlorite solution, ( $\text{NaClO}$ )

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**PREP 464**

*Sample preparation for Chelated Fe, Mn, Cu, Zn [FM 840] [FM 841]*

Version: 2.0

Effective Date: 10/16/2023

Page 2 of 4

7.5. Water, Deionized (DI) Type I

## 8. REAGENT PREPARATION GUIDELINES

### 8.1. DAP Solution

8.1.1. Directions: Dissolve 10g of DAP in approximately 500mL of DI Water in a 1000mL volumetric flask, and bring up to volume with DI water; mix well.

### 8.2. 0.1% Sodium Hypochlorite Solution

8.2.1. Directions: Add 20.0mL of 5% sodium hypochlorite to a 1000mL volumetric flask containing approximately 500mL of DI water, and bring up to volume; mix well.

### 8.3. Chelate Buffer

8.3.1. Preparation of 1L: Dissolve 1g of sodium acetate and 25mL of glacial acetic acid in 700mL of DI water; bring to a final volume of 1L with DI water.

8.3.2. Preparation of 20L: Add 20g of sodium acetate and 500mL of acetic acid to a calibrated container; bring up to a final volume of 20L with DI water; mix well.

8.4. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015.

8.5. The specific ingredient amounts can vary to accommodate the batch size.

8.6. Expiration Date: 12 months from date of preparation, unless otherwise specified.

8.7. Store at room temperature.

## 9. ASSOCIATED DOCUMENTS

9.1. ADMIN 023, Attachment AA, Fertilizer QC Loc Code and Result Entry Templates

9.2. F00015, Reagent Preparation

9.3. L003, Reagent

9.4. PREP 453, Weighing of Fertilizer Sample

9.5. PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size

## 10. SPECIFIC PROCEDURES

### 10.1. Collection of Weighed Samples

10.1.1. Samples are weighed according to PREP 453, Weighing of Fertilizer Sample.

10.1.2. Follow Laboratory Information Management System (LIMS) weight recommendations.

10.1.3. Weigh samples into 200mL volumetric flasks **per the table below:**

| Guarantee (Guar) | Weight   |
|------------------|----------|
| <0.008%          | 5g       |
| 0.008-0.08%      | 0.4/Guar |
| 0.08-0.8%        | 0.2/Guar |
| >8%              | 0.25g    |

10.1.4. The batch size should not exceed the limit specified in Prep 453, Attachment A, Fertilizer Laboratory Standard and Batch Size.

10.1.5. Collect batch samples from the properly identified bin(s).

10.1.6. Follow Laboratory Information Management System (LIMS) weight recommendations.

### 10.2. Extraction of Samples

10.2.1. Add 100mL of Chelate Buffer to the samples and bring to volume with DI water.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 464

#### *Sample preparation for Chelated Fe, Mn, Cu, Zn [FM 840] [FM 841]*

Version: 2.0

Effective Date: 10/16/2023

Page 3 of 4

10.2.2. Shake sample flasks for a minimum of 30 seconds, vigorously, and allow samples to stand for at least 45 minutes.

10.2.3. Filter liquid samples containing “free floating” particles and undissolved particles using a screen column.

10.2.4. Filter cloudy samples using an aqueous 25mm, 0.45µm non-sterile syringe filter.

10.2.5. If Chelated Iron results are excessively high (% found > 3x % guarantee), samples must be re-prepared and analyzed.

10.2.5.1. Weigh sample in a 200mL volumetric flask. Add 100mL of Chelate Buffer, 15mL Diammonium phosphate (DAP) Solution and 15mL Sodium Hypochlorite Solution. Bring up to a 200mL final volume and shake for 30 seconds. Allow sample to sit 5-10 minutes to ensure precipitation reaction goes to completion.

#### 10.3. Transfer of Samples

10.3.1. Assemble a set of 15mL Graduated Polystyrene Centrifuge Tubes, one per sample.

10.3.2. Place sample tubes on a suitable test tube rack, in a numerical sequence.

10.3.3. Label sample tubes with the sample ID number.

10.3.4. Pour approximately 14mL of each sample into the test tube. This allows a gap between the top of the tube and sample surface that prevents spills.

10.3.5. Cover samples with parafilm and make a seal to avoid spills or cross contamination.

#### 10.4. Data Entry into Laboratory Information Management System (LIMS)

10.4.1. Log onto Laboratory Information Management System (LIMS).

10.4.2. Double click the “Results” folder and single click “Results Entry”.

10.4.3. Refer to ADMIN 023, Attachment AA, Fertilizer QC Loc Code and Result Entry Templates, for specific result entry template names for each metal:

10.4.3.1. Chelate Iron: FERT\_AA\_CHELATE\_IRON

10.4.3.2. Chelate Manganese: FERT\_AA\_CHELATE\_MANGANESE

10.4.3.3. Chelate Copper: FERT\_AA\_CHELATE\_COPPER

10.4.3.4. Chelate Zinc: FERT\_AA\_CHELATE\_ZINC

10.4.4. For each sample, enter the flask volume and save.

10.4.5. Reopen template to verify that the data was saved.

10.4.6. Print and close the template.

10.4.7. Refer to ADMIN 023, Attachment AA, Fertilizer QC Loc Code and Result Entry Templates, for specific result entry template names. For example: select the:

10.4.7.1. Chelate Iron: FERT\_PRINT\_AA\_CHELATE\_FE

10.4.7.2. Chelate Manganese: FERT\_PRINT\_AA\_MN

10.4.7.3. Chelate Copper: FERT\_PRINT\_AA\_CU

10.4.7.4. Chelate Zinc: FERT\_PRINT\_AA\_ZN

10.4.8. Print and close the template.

10.4.9. Indicate the QC sample number on the printouts, initial and date.

#### 10.5. Delivery

10.5.1. Deliver all samples and paperwork to the designated area, using a laboratory-approved sample carrier.

## 11. WASTE MANAGEMENT

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### PREP 464

*Sample preparation for Chelated Fe, Mn, Cu, Zn [FM 840] [FM 841]*

Version: 2.0

Effective Date: 10/16/2023

Page 4 of 4

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material                   | Disposal Mechanism                |
|----------------------------|-----------------------------------|
| Excess extracts            | Dispose down drain while flushing |
| Excess buffer and reagents | Dispose down drain while flushing |

## 12. REFERENCES

12.1. AOAC 16<sup>th</sup> Edition, Method 938.03

## 13. VERSION HISTORY

| Version | Date       | Description                                                                                                            | Author                          | Editor          |
|---------|------------|------------------------------------------------------------------------------------------------------------------------|---------------------------------|-----------------|
| 1.0     | 11.01.2019 | Original Version                                                                                                       | Jack Andreu<br>Christian Amason |                 |
| 2.0     | 10.16.2023 | Removed reference to filter screen column throughout.<br>Section 10.1.3. added the guarantee/weight determining table. |                                 | Taleya Williams |

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

### **PREP 466**

#### *Sample Preparation for Nitrate, Urea and Ammoniacal Analysis using Flow Injection Analyzer [FM 430]*

Version: 2.0

Effective Date: 09/20/2019

Page 1 of 4

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for the preparation of samples for the individual nitrogen break-down components. These samples subsequently undergo flow injection analysis.

## **2. SCOPE**

2.1. The procedure delineated below is applicable to fertilizer samples in the Bureau of Agriculture Environmental Laboratories (BAEL). This preparation is for the analysis of Nitrate nitrogen, Urea or other water soluble nitrogen, and Ammoniacal nitrogen through the usage of a Flow Injection Analyzer. Fertilizer samples that only contain Nitrate nitrogen, Urea or other water soluble nitrogen, and Ammoniacal nitrogen will be labeled as 430 NAU. Fertilizer samples that contain some or all the analytes in NAU as well as Water Insoluble Nitrogen will be labeled 430 KJL.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## **4. OUTLINE of PROCEDURE**

- 4.1. 10.1 Collection of Weighed Samples
- 4.2. 10.2 Extraction of Samples
- 4.3. 10.3 Transfer of Samples
- 4.4. 10.4 Data Entry into Laboratory Information Management
- 4.5. 10.5 Delivery

## **5. INTERFERENCE**

5.1. N/A

## **6. APPARATUS and MATERIAL**

- 6.1. Bin
- 6.2. Centrifuge tube (15-50mL)
- 6.3. Flask (250mL)
- 6.4. Parafilm
- 6.5. Screening column
- 6.6. Stopper/cork
- 6.7. Test tube, Glass (13-15mL)

## **7. REAGENTS and MEDIA**

- 7.1. Methanol, Optima or equivalent
- 7.2. Water, Deionized (DI) Type I

## **8. REAGENT PREPARATION GUIDELINES**

8.1. N/A

## **9. ASSOCIATED DOCUMENTS**

- 9.1. ADMIN 023, Attachment AA, Fertilizer QC Loc Code and Result Entry Templates.
- 9.2. PREP 453, Weighing of Fertilizer Sample
- 9.3. PREP 453 Attachment A, Fertilizer Laboratory Standard and Batch Size

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

### **PREP 466**

#### *Sample Preparation for Nitrate, Urea and Ammoniacal Analysis using Flow Injection Analyzer [FM 430]*

Version: 2.0

Effective Date: 09/20/2019

Page 2 of 4

## **10. SPECIFIC PROCEDURES**

### **10.1. Collection of Weighed Samples**

- 10.1.1. Samples are weighed according to PREP 453, Weighing of Fertilizer Sample.
- 10.1.2. The batch size should not exceed the limit specified in Prep 453, Attachment A, Fertilizer Laboratory Standard and Batch Size.
- 10.1.3. Collect batch samples from the properly identified bin(s).

### **10.2. Extraction of Samples**

- 10.2.1. For samples 430 KJL, which contain Water Insoluble Nitrogen (WIN), it may be necessary to wet sample with 1 to 2mL of methanol, to prevent clumping before adding water.
- 10.2.2. Bring the samples to volume (i.e. 250mL) using deionized water.
- 10.2.3. Stopper and mix contents by inverting and vigorously shaking the flask and allow to settle.

### **10.3. Transfer of Samples**

- 10.3.1. Assemble two sets of tubes. Every sample should have a glass test tube and a centrifuge tube.
- 10.3.2. Place glass test tubes on numbered rack, starting with the laboratory control sample (LS) then in sample number ascending order.
- 10.3.3. Put screening columns into the glass test tubes, ensuring that there is a filter in the columns.
- 10.3.4. Remove the labels from the flasks and place them on the respective glass test tubes.
- 10.3.5. Make a duplicate label for the centrifuge tubes.
- 10.3.6. Pour approximately 9mL of each sample into the respective test tube in order to leave headspace.
- 10.3.7. Pour sample into centrifuge tube to approximately 15mL.
- 10.3.8. Place parafilm over the test tubes and cap the centrifuge tubes.

### **10.4. Data Entry into Laboratory Information Management System**

- 10.4.1. Log onto Laboratory Information Management System (LIMS).
- 10.4.2. Double click the "Results" folder and single click "Results Entry".
- 10.4.3. Refer to ADMIN 023, Attachment AA, Fertilizer QC Loc Code and Result Entry Templates, for specific result entry template names. For example: select the FERT\_PREP\_NAU\_SHORT for 430 NAU samples or FERT\_PREP\_NAU\_KJL for 430 KJL samples.
- 10.4.4. For each sample, enter the flask volume and save.
- 10.4.5. Reopen template to verify that the data was saved.
- 10.4.6. Print and close the template.
- 10.4.7. Refer to ADMIN 023, Attachment AA, Fertilizer QC Loc Code and Result Entry Templates, for specific result entry template names. For example: select the template FERT\_PRINT\_NAUS for 430 NAU samples or FERT\_PRINT\_NAUL for 430 KJL samples.
- 10.4.8. Print and close the template.
- 10.4.9. Indicate the QC sample number on the printouts, initial and date the printouts.

### **10.5. Delivery**

- 10.5.1. Deliver all samples and paperwork to the designated area, using a laboratory approved sample carrier.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 466

*Sample Preparation for Nitrate, Urea and Ammoniacal Analysis using Flow Injection Analyzer [FM 430]*

Version: 2.0

Effective Date: 09/20/2019

Page 3 of 4

## 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material         | Disposal Mechanism                |
|------------------|-----------------------------------|
| Screening Column | Discard in trash                  |
| Excess Extract   | Discard down drain while flushing |

## 12. REFERENCES

12.1. Preparation for Breakdown Nitrogen Analysis, FSFL-SOP 491 Rev.:1.0 (2006/01/17)

## 13. VERSION HISTORY

| Version | Date       | Description                                                                                     | Author    | Editor          |
|---------|------------|-------------------------------------------------------------------------------------------------|-----------|-----------------|
| 1.0     | 03.29.2017 | Original Version                                                                                | Ian McCoy |                 |
| 2.0     | 09.20.2019 | Update Outline Format<br>Changed Ian Linton to Ian McCoy, Same individual last name was changed |           | Taleya Williams |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

PREP 466

*Sample Preparation for Nitrate, Urea and Ammoniacal Analysis using Flow  
Injection Analyzer [FM 430]*

Version: 2.0

Effective Date: 09/20/2019

Page 4 of 4

**14. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer  
Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Wet Chemistry and Preparation Manager / Editor  
Date\_\_\_\_\_

\_\_\_\_\_  
Ian McCoy  
Author  
Date\_\_\_\_\_

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 500**

*Sample Preparation and Analysis of Direct Water Soluble Nitrogen*  
*[FM 425, FM 440]*

Version: 2.0

Effective Date: 01/16/2024

Page 1 of 9

**1. PURPOSE**

This standard operating procedure (SOP) provides instruction for determination of water-soluble nitrogen.

**2. SCOPE**

The procedure delineated below is applicable to use for the Kjeldahl method in determining water soluble nitrogen by converting nitrate into nitrate gas using ferrous sulfate and distilling the by-product into an acid trap and titrating the excess.

Samples that contain water insoluble nitrogen are associated with the term “440” and sample lacking ammoniacal and water insoluble nitrogen are associated with the term “425”.

**3. HEALTH and SAFETY**

A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required. Additional PPE required for this procedure includes face shield and acid resistant apron.

**4. OUTLINE of PROCEDURE**

- 4.1 10.1 Weighing
- 4.2 10.2 Filtration
- 4.3 10.3 Digestion
- 4.4 10.4 Acid Trap Setup
- 4.5 10.5 Distillation
- 4.6 10.6 Titration
- 4.7 10.7 Result Entry

**5. INTERFERENCE**

N/A

**6. APPARATUS and MATERIAL**

- 6.1 Auto titration system
- 6.2 Beaker, 300-ml
- 6.3 Filter paper
- 6.4 Flask, 250-ml
- 6.5 Flask, round bottom 800-ml
- 6.6 Fume hood
- 6.7 Kjeldahl digestion and distillation unit

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 500**

*Sample Preparation and Analysis of Direct Water Soluble Nitrogen*  
*[FM 425, FM 440]*

Version: 2.0

Effective Date: 01/16/2024

Page 2 of 9

6.8 Stopper

6.9 Weigh container

**7. REAGENTS and MEDIA**

7.1 Ammonium sulfate

7.2 Boiling stones, granular

7.3 Bromocresol green

7.4 Ferrous sulfate

7.5 Methanol

7.6 Methyl red

7.7 Sodium hydroxide, pellets

7.8 Sodium hydroxide, 50%

7.9 SRM-695

7.10 Sulfuric acid, concentrated reagent grade

7.11 Tris-hydroxymethyl aminomethane

7.12 Water, deionized

7.13 Water, tap

**8. REAGENT PREPARATION GUIDELINES**

**8.1 0.5N Sulfuric Acid (Standard Acid)**

8.1.1 Directions: Add approximately 20L of DI Water into a 50L plastic carboy marked at 41L. Slowly add 580mL of sulfuric acid, bring to 41L with DI water. Bubble air through this solution until thoroughly mixed.

8.1.2 Titration of Sulfuric Acid Solution (Option 1): Dry tris-hydroxymethyl aminomethane (THAM) in oven at 105°C for at least 2 hours. Cool it down in the desiccator. Weigh 0.4500g of THAM into each of the six 125mL Erlenmeyer flasks. Add 5 drops of Bromocresol Green to each flask. Add approximately 50mL of DI water to each flask. The color of the solution should be blue. If not blue, add more bromocresol green a drop at a time to the flask. Add 0.5N sulfuric acid solution to 25mL dispensing burette using a pipet. Slowly titrate THAM with sulfuric acid to bromocresol green end point (pH 4.7). Discard the highest and lowest results and average the remaining four results.

8.1.3 Determine the concentration (normality) of each of the six sulfuric acid solutions using the formula:

$$\text{Acid Normality} = \frac{[\text{THAM Weight (g)}] \times 8.2549}{\text{Acid Volume (mL)}}$$

$$\text{NOTE: } 8.2549 = \frac{1000\text{mL/L} \times 2(\text{Conversion from Molarity to Normality})}{121.14\text{g/mol} \times 2(\text{Stoichiometric Ratio})}$$

8.1.4 Titration of Sulfuric Acid Solution (Option 2): Add 4mL of 0.5N Sulfuric Acid (Standard Acid) to seven (7) 300mL beakers then fill them to approximately 250mL with water. Fill a separate 300mL beaker with 200mL of certified 0.5N NaOH; place the line of the titrator into the beaker and purge the syringe of the

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 500**

*Sample Preparation and Analysis of Direct Water Soluble Nitrogen*  
*[FM 425, FM 440]*

Version: 2.0

Effective Date: 01/16/2024

Page 3 of 9

titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.1.4.1 Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Acid 1 – Acid 6.

8.1.5 Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$A_N = \frac{B_N \times B_V}{A_V}$$

Where:  $A_N$ =Acid Normality  
 $A_V$ =Acid Volume  
 $B_N$ =Base Normality  
 $B_V$ =Base Volume

8.1.6 Discard the highest and lowest results and average the remaining results

## 8.2 18.75N Stock Sodium Hydroxide Solution

8.2.1 Directions: Add 1 liter of deionized (DI) water into a 2L volumetric flask with stir bar in it on a stirring plate. Add 1500g of sodium hydroxide pellets in small (~100g) portions with continuous mixing. Bring solution volume to 2L with DI water.

## 8.3 0.2N Sodium Hydroxide (Standard Base)

8.3.1 Directions (Option A): Fill a beaker with 1000ml of DI water. Add stir bar and start mixing. Add 1000ml of 50% Sodium Hydroxide and let it mix for 30 minutes.

Directions (Option B): Fill a 20L glass carboy, marked at 16L, half full of DI water. Add stir bar and start mixing. Add 300mL of stock sodium hydroxide solution to carboy while continuously mixing. Bring solution to a volume of 16L with DI water. Stir solution for approximately 30 minutes.

8.3.2 Titration of Sodium Hydroxide Solution (Option 1): Dry potassium acid phthalate (PAP) in oven at 105°C for at least 2 hours and then allow it to cool in desiccators. Take six 300mL beakers and weigh approximately 0.8000g of PAP into each. Add 150mL of DI water into each of the six beakers. Allow PAP to dissolve. Titrate PAP solution on the QC\_Titrate System.

8.3.3 Titration of PAP Solution on the QC-Titrator: Click on PAP icon. Type set name and click OK. Type the exact weights of PAP. Highlight the first line on time-table and click start. When titration is done, a hard copy is printed automatically.

8.3.4 Calculation of Sodium Hydroxide: Determine the concentration (normality) of the six sodium hydroxide solutions using the formula:

$$\text{Base Normality} = \frac{[\text{PAP Weight (g)}] \times 4.8964}{\text{Base Volume (mL)}}$$

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 500**

*Sample Preparation and Analysis of Direct Water Soluble Nitrogen*  
*[FM 425, FM 440]*

Version: 2.0

Effective Date: 01/16/2024

Page 4 of 9

$$\text{NOTE: } 4.8964 = \frac{1000\text{mL/L} \times 1(\text{Conversion from Molarity to Normality})}{204.23\text{g/mol} \times 1(\text{Stoichiometric Ratio})}$$

8.3.5 Titration of Sodium Hydroxide Solution (Option 2): Add 4mL of 0.35N Sodium Hydroxide to seven (7) 300mL beakers. Fill a separate 300mL beaker with 200mL of the normalized solution from 8.1; place the line of the titrator into the beaker and purge the syringe of the titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.3.5.1 Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Base 1 – Base 6.

8.3.6 Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$B_N = \frac{A_N \times A_V}{B_V}$$

Where:  $A_N$ =Acid Normality  
 $A_V$ =Acid Volume  
 $B_N$ =Base Normality  
 $B_V$ =Base Volume

8.3.7 Discard the highest and lowest results and average the remaining four results. Change base normality on Kjeldahl and Direct Excel spreadsheets for LIMS.

8.4 Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015

8.5 The specific ingredient amounts can vary to accommodate the batch size.

8.6 Expiration Date: 12 months from date of preparation unless otherwise specified.

8.7 Store at room temperature

## **9 ASSOCIATED DOCUMENTS**

9.1 ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Template

9.2 Paradigm 3, Document Control Software, Improvement Module

## **10 SPECIFIC PROCEDURES**

### **10.1 Weighing**

10.1.1 Using ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Template, create a laboratory identification number for the laboratory control sample (LS).

10.1.2 Weigh approximately 0.5 g of the applicable laboratory control sample.

10.1.2.1 SRM-695 for "440"

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 500**

*Sample Preparation and Analysis of Direct Water Soluble Nitrogen*  
*[FM 425, FM 440]*

Version: 2.0

Effective Date: 01/16/2024

Page 5 of 9

10.1.2.2 Ammonium sulfate for “425”.

10.1.3 Weigh samples to approximately 0.9 grams except liquid samples that have been diluted. Diluted liquid samples are approximately doubled in weight (i.e. approximately 1.8 g).

10.1.4 Enter weight into Laboratory Information Management System (LIMS) using the template FERT\_Nitrogen\_KJL. *Be aware of the given acid trap volume because it will be used later in the process.*

10.2 Filtration

10.2.1 Only samples that contain water insoluble nitrogen (“440” sample) require filtration.

10.2.2 Fold the filter paper and place it into funnel that drains into 250-ml flask.

10.2.3 Wet the sides of the filter paper with deionized water to create a seal for the filter.

10.2.4 Pour the sample into filter paper.

10.2.5 Using deionized water, rinse the weigh container clean of sample into filter paper.

10.2.6 If the sample is clumping, use 2 to 5 mL of methanol to break up clumping.

10.2.7 Fill the filter paper with deionized water and allow to drain completely. Repeat the washing and draining process twice, for a total of three times.

10.2.8 Bring the filtrate, in the flask, up to 250-mL mark with deionized water.

10.2.9 After the filtrate is brought to volume, stopper in flask to prevent evaporation.

10.3 Digestion

10.3.1 Add 4 to 8 granular boiling stones to each 800-ml round bottom flask.

10.3.2 Add ¼ teaspoon of ferrous sulfate to each 800-ml round bottom flask.

10.3.3 This part of the process differs depending on whether the sample lacks ammoniacal and water insoluble nitrogen or the sample contains water insoluble nitrogen.

10.3.3.1 For sample lacking ammoniacal and water insoluble nitrogen (“425” sample):

10.3.3.1.1 Wash sample into round bottom flask using deionized water.

10.3.3.1.2 Wash weigh container, at least 3 times, until sample is completely transferred.

10.3.3.2 For sample containing water insoluble nitrogen (“440” sample):

10.3.3.2.1 Mix the contents of the 250-mL flask. To ensure that the filtrated sample is thoroughly mix, shake and invert at least 7 times.

10.3.3.2.2 Add 100 mL of filtrate to each flask using either graduated cylinder or pipette.

10.3.4 Add 20 mL of concentrated sulfuric acid to each flask. Swirl flask to mix solution.

10.3.4.1 Turn on fume hood and rack ventilation systems.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 500**

*Sample Preparation and Analysis of Direct Water Soluble Nitrogen*  
*[FM 425, FM 440]*

Version: 2.0

Effective Date: 01/16/2024

Page 6 of 9

10.3.4.2 Digest sample on the digestion rack.

10.3.4.2.1 Approximately 60 minutes for “425” samples.

10.3.4.2.2 Approximately 75 minutes for “440” samples.

10.3.4.3 If any viscous black\brown material forms on top of round bottom flask, rotate flask on digestion rack to facilitate the digestion of material.

10.3.4.4 After digestion is complete, allow sample to cool to room temperature.

10.3.4.5 Once samples have reached room temperature, fill round bottom flask to approximately half way with deionized water.

#### 10.4 Acid Trap Setup

10.4.1 The preparation of the acid trap, using the acid dispenser, differ depending on the specific sample.

10.4.1.1 For Blank:

10.4.1.1.1 Dispense 4 mL of standard acid into 300-ml beaker.

10.4.1.1.2 Bring to the 250-mL mark with deionized water.

10.4.1.1.3 Placed beaker, with sample, in slot 1 on the auto titrator.

10.4.1.2 For Laboratory Control Sample (used in “440” preparation):

10.4.1.2.1 Dispense 4 mL of standard acid into 300-mL beaker.

10.4.1.2.2 Bring to 50-mL mark with deionized water.

10.4.1.2.3 Add 1-4 drops of methyl red indicator.

10.4.1.3 For Laboratory Control Sample (used in “425” preparation):

10.4.1.3.1 Dispense 18 mL of standard acid into 300-ml beaker.

10.4.1.3.2 Bring to 50-mL mark with deionized water.

10.4.1.3.3 Add 1-4 drops of methyl red indicator.

10.4.1.4 For Sample:

10.4.1.4.1 Dispense the applicable volume of standard acid into the 300-ml beaker. *This volume should be based upon the acid trap volume that was calculated in LIMS earlier in the process.*

10.4.1.4.2 Bring to the 50-mL mark with deionized water.

10.4.1.4.3 Add 1-4 drops of methyl red indicator.

10.4.2 Make extra 4-ml aliquots of standard acid. *This standard acid can be added to distillation acid trap if indicator starts to turn orange/yellow.*

#### 10.5 Distillation

10.5.1 Place acid trap into its corresponding spot on the distillation rack with the tip of the distillation tube visibly submerged in the liquid.

10.5.2 Turn on cooling jackets and fume hood.

10.5.3 Set heaters on high and allow to heat up.

10.5.4 Swirl round bottom flask vigorously ensuring sample is mixed thoroughly.

10.5.5 Slowly add 80 mL of 50% sodium hydroxide to each flask.

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 500**

*Sample Preparation and Analysis of Direct Water Soluble Nitrogen*  
*[FM 425, FM 440]*

Version: 2.0

Effective Date: 01/16/2024

Page 7 of 9

10.5.6 Immediately after sodium hydroxide is added, attach flask to the distillation column and place on heating rack.

10.5.7 Allow distillation to proceed until acid trap volume reaches 250-mL mark.

10.5.8 Remove the acid trap from the distillation column.

10.5.9 Turn off heating racks and fume hood.

10.5.10 Once samples have stopped boiling, turn off cooling jackets.

**10.6 Titration**

10.6.1 Start Pc titrate program.

10.6.2 Once program opens, load “Ammonium” template. NOTE: The icon is two test tubes.

10.6.3 Delete any extra spaces from the template. The amount of spaces needed in the template is one more than the number of samples. If there are 11 samples there should be 12 spaces in template.

10.6.4 Enter sample identification numbers into template.

10.6.5 The blank should be in position 1. Label position 1 with the date and user initials.

10.6.6 The laboratory control sample should be in position 2.

10.6.7 Place acid trap into corresponding locations on auto titration rack.

10.6.8 Select position 1 in the template. NOTE: This is where the auto titrator will start the titration process.

10.6.9 Make sure all acid traps correspond to their location on the template.

10.6.10 Click start and ensure that the titration process begins.

10.6.11 Sign, date and file the instrument’s print out of titration results.

**10.7 Result Entry**

10.7.1 Record the results (from the print out) into the FERT\_NITROGEN\_KJL template in the LIMS. A laboratory notebook may be used as an intermediate location for the results.

10.7.2 Any additional acid used during the distillation must be added to the acid trap value in the LIMS. This updated value will be used in calculations.

10.7.3 The raw data, from the auto titrator, is uploaded to the folder on the “\\tlhaeschrome” drive.

**11 Waste Management**

Material, marked for disposal, is removed via the following mechanism:

| Material                          | Disposal Mechanism                          |
|-----------------------------------|---------------------------------------------|
| Boiling stones                    | Trash                                       |
| Residue (from filtration process) | Trash                                       |
| Sample extract, excess            | Stored in waste container until neutralized |

|                                                                                                                                                 |
|-------------------------------------------------------------------------------------------------------------------------------------------------|
| <p align="center"><b><i>STANDARD OPERATING PROCEDURE</i></b></p> <p align="center"><b>BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES</b></p> |
|-------------------------------------------------------------------------------------------------------------------------------------------------|

|                                                                                                                                        |
|----------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>METHOD 500</b></p> <p><i>Sample Preparation and Analysis of Direct Water Soluble Nitrogen</i></p> <p><i>[FM 425, FM 440]</i></p> |
|----------------------------------------------------------------------------------------------------------------------------------------|

|              |                            |             |
|--------------|----------------------------|-------------|
| Version: 2.0 | Effective Date: 01/16/2024 | Page 8 of 9 |
|--------------|----------------------------|-------------|

|                                    |                                  |
|------------------------------------|----------------------------------|
| Neutralized sample extract, excess | Flushed down the sink with water |
|------------------------------------|----------------------------------|

**12 REFERENCES**

1. E00151, Calibration and use of QC-Titrate™ Titration system, FSFL-SOP 494
2. E00149, Kjeldahl Digestion of Fertilizer Samples, FSFL-SOP 492
3. E00150, Kjeldahl Distillation of Fertilizer Samples, FSFL-SOP 493

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 500  
*Sample Preparation and Analysis of Direct Water Soluble Nitrogen*  
*[FM 425, FM 440]*

|              |                            |             |
|--------------|----------------------------|-------------|
| Version: 2.0 | Effective Date: 01/16/2024 | Page 9 of 9 |
|--------------|----------------------------|-------------|

**13 VERSION HISTORY**

| Version | Date       | Description                                                                                                                         | Author        | Editor                      |
|---------|------------|-------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------------------|
| 1.0     | 2017.05.04 | Original Version                                                                                                                    | Zachary Tower |                             |
| 2.0     | 01.16.2024 | Updated SOP to new format<br>Added reagents.<br>Remove signature page<br>Minor grammatical fixes<br>Added external document numbers |               | Amy Bryant<br>Zachary Tower |

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 501**

*Sample Preparation and Analysis of Water Insoluble Nitrogen*  
*[FM 445]*

Version: 2.0

Effective Date: 01/16/2024

Page 1 of 8

## 1. PURPOSE

This standard operating procedure (SOP) provides instruction for the determination of water insoluble nitrogen in fertilizer.

## 2. SCOPE

The procedure delineated below is applicable to use for the Kjeldahl method for the determination of water insoluble nitrogen in fertilizer samples first by digestion using sulfuric acid and Kjeldahl tablets, then distilling the by-product into an acid trap using a standard acid, and finally titrating the excess using a standard base.

This method does not give the degree to which the water insoluble nitrogen is available.

## 3. HEALTH and SAFETY

A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required. Additional PPE required for this procedure includes face shield and acid resistant apron.

## 4. OUTLINE of PROCEDURE

- 4.1 10.1 Weighing
- 4.2 10.2 Filtration
- 4.3 10.3 Digestion
- 4.4 10.4 Acid Trap Setup
- 4.5 10.5 Distillation
- 4.6 10.6 Titration
- 4.7 10.7 Data Evaluation
- 4.8 10.8 Result Entry
- 0.5 N Sulfuric Acid (Standard Acid)
- 18.75N Stock Sodium Hydroxide
- 0.25 N Sodium Hydroxide (Standard Base)

## 5. INTERFERENCE

N/A

## 6. APPARATUS and MATERIAL

- 1. Auto titration system
- 2. Beaker, 300-ml
- 3. Filter paper
- 4. Flask, 250-ml
- 5. Flask, round bottom flasks 800-ml

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 501**

*Sample Preparation and Analysis of Water Insoluble Nitrogen*  
*[FM 445]*

Version: 2.0

Effective Date: 01/16/2024

Page 2 of 8

6. Fume hood
7. Funnels
8. Kjeldahl digestion and distillation unit
9. Stoppers
10. Weigh container
11. Acid dispenser
12. Class A pipet

## 7. REAGENTS and MEDIA

1. Boiling stones, granular
2. Hollow Glass Beads
3. Bromocresol green
4. FisherTab, Kjeldahl tablets (Kel Tabs)
5. Methanol
6. 0.5% Methyl Red indicator
7. Sodium hydroxide, 50%
8. Sodium hydroxide, 0.5N, certified
9. SRM-695
10. Sulfuric Acid, concentrated reagent grade
11. Water, deionized (type II)
12. Water (Tap)
13. Sodium Hydroxide Pellets NF/FCC

## 8. REAGENT PREPARATION GUIDELINES

### 8.1 0.5N Sulfuric Acid (Standard Acid)

8.1.1 Directions: Add approximately 20L of DI Water into a 50L plastic carboy marked at 41L. Slowly add 580mL of sulfuric acid, bring to 41L with DI water. Bubble air through this solution until thoroughly mixed.

8.1.2 Titration of sulfuric acid (option1): Dry tris-hydroxymethyl aminomethane (THAM) in oven for at least 2 hours. Cool it down in the desiccator. Weigh 0.4500g of THAM into each of the six 125ml Erlenmeyer flask. Add 5 drops of bromocresol green to each flask. Add approximately 50ml of DI water to each flask. The color of the solution should be blue. If not blue, add more bromocresol green, one drop at a time to the flask. Add 0.5N sulfuric acid to 25ml dispensing burette using a pipet. Slowly titrate THAM with sulfuric acid to bromocresol green end point (ph7). Discard the highest and lowest results and average the remaining results. Determine the concentration (normality) of each of the six sulfuric acid solutions using the formula:

$$\text{Acid Normality} = \frac{[\text{THAM weight(g)}] \times 8.2549}{\text{Acid volume}}$$

Acid volume

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 501**

*Sample Preparation and Analysis of Water Insoluble Nitrogen*  
*[FM 445]*

Version: 2.0

Effective Date: 01/16/2024

Page 3 of 8

**NOTE:**  $8.2549 = \frac{1000\text{mL/L} \times 2(\text{Conversion from Molarity to Normality})}{121.14\text{g/mol} \times 2(\text{Stoichiometric Ratio})}$

8.1.3 Titration of Sulfuric Acid Solution (option 2): Add 4mL of Reagent 01 to seven (7) 300mL beakers then fill them to approximately 250mL with water. Fill a separate 300mL beaker with 200mL of certified 0.5N NaOH; place the line of the titrator into the beaker and purge the syringe of the titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.1.3 Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Acid 1 – Acid 6.

8.1.4 Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$A_N = \frac{B_N \times B_V}{A_V}$$

Where:  $A_N$ =Acid Normality  
 $A_V$ =Acid Volume  
 $B_N$ =Base Normality  
 $B_V$ =Base Volume

8.1.5 Discard the highest and lowest results and average the remaining results

8.1.6 Change acid normality on the Kjeldahl and Direct Excel spreadsheets. The new value for acid must be changed in LIMS.

## 8.2 18.75N Stock Sodium Hydroxide

8.2.1 Directions: Add 1L of deionized (DI) water into a 2L volumetric flask with a stir bar in it, on a stir plate. Add 1500g of sodium hydroxide pellets in small (~100g) portions with continuous mixing. Bring solution volume to 2L with DI water.

### 8.2.2 0.25N Sodium Hydroxide (Standard Base)

Directions (Option A): Fill a beaker with 1000ml of DI water. Add stir bar and start mixing. Add 1000ml of 50% Sodium Hydroxide and let it mix for 30 minutes.

Directions (Option B): Fill a 20L glass carboy, marked at 16L half full of DI water. Add stir bar and start mixing. Add 300ml of stock sodium hydroxide solution to carboy while continuously mixing. Bring solution to a volume of 16L with DI water. Stir solution for approximately 30 minutes.

8.2.3 Titration of Sodium Hydroxide Solution (Option 1): Dry potassium acid phthalate (PAP) in oven at 105°C for at least 2 hours and then allow it to cool in desiccators. Take six 300mL beakers and weigh approximately 0.8000g of PAP into each. Add 150mL of DI water into each of the six beakers. Allow PAP to dissolve. Titrate PAP solution on the QC Titrate System.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 501**

*Sample Preparation and Analysis of Water Insoluble Nitrogen*  
*[FM 445]*

Version: 2.0

Effective Date: 01/16/2024

Page 4 of 8

8.2.4 Titration of PAP Solution on the QC-Titrator: Click on PAP icon. Type set name and click OK. Type the exact weights of PAP. Highlight the first line on time-table and click start. When titration is done, a hard copy is printed automatically.

8.2.5 Calculation of Sodium Hydroxide: Determine the concentration (normality) of the six sodium hydroxide solutions using the formula:

$$\text{Base Normality} = \frac{[\text{PAP Weight (g)}] \times 4.8964}{\text{Base Volume (mL)}}$$

$$\text{NOTE: } 4.8964 = \frac{1000\text{mL/L} \times 1(\text{Conversion from Molarity to Normality})}{204.23\text{g/mol} \times 1(\text{Stoichiometric Ratio})}$$

8.2.6 Titration of Sodium Hydroxide Solution (Option 2): Add 4mL of Reagent 03 to seven (7) 300mL beakers. Fill a separate 300mL beaker with 200mL of the normalized solution from 8.1; place the line of the titrator into the beaker and purge the syringe of the titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.2.6.1 Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Base 1 – Base 6.

8.2.7 Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$B_N = \frac{A_N \times A_V}{B_V}$$

Where:  $A_N$ =Acid Normality

$A_V$ =Acid Volume

$B_N$ =Base Normality

$B_V$ =Base Volume

8.2.8 Discard the highest and lowest results and average the remaining results. Change base normality on Kjeldahl and Direct Excel spreadsheets for LIMS.

8.3 Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015.

8.4 The specific ingredient amounts can vary to accommodate the batch size.

8.5 Expiration Date: 12 months from date of preparation unless otherwise specified.

8.6 Store at room temperature.

## 9. ASSOCIATED DOCUMENTS

9.1 ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates

9.2 DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs

9.3 Paradigm 3, Document Control Software, Improvement Module

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 501**

*Sample Preparation and Analysis of Water Insoluble Nitrogen*  
*[FM 445]*

Version: 2.0

Effective Date: 01/16/2024

Page 5 of 8

**9.4 Labop 140 at J, Maintenance, Operation, and Procedure of QC Titrator System**

**10. SPECIFIC PROCEDURES**

**10.1 Weighing**

10.1.1 Using ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Template, create a laboratory identification number for the laboratory control sample (LS).

13.1.2 Weigh approximately 0.5 g of LS, SRM-695

10.1.3 Weigh samples to approximately 0.9 grams except liquid samples that have been diluted. Diluted liquid samples are approximately doubled in weight (i.e. approximately 1.8 g).

10.1.4 Enter weight into Laboratory Information Management System (LIMS) using the template FERT\_Nitrogen\_KJL. *Be aware of the given acid trap volume because it will be used later in the process.*

**13.2 Filtration**

10.2.1 Fold filter paper and place into funnel that drains into 250-ml flask.

10.2.2 Wet sides of the filter paper with deionized water to create a seal for the filter.

10.2.3 Pour sample into filter paper.

10.2.4 Using deionized water, rinse the weigh container clean of sample into filter paper.

10.2.5 If sample is clumping, use 2 to 5 mL of methanol to break up clumping

10.2.6 Fill the filter paper up with deionized water and allow to drain completely.

Repeat the washing and draining process twice, for a total of three times.

10.2.7 Bring the filtrate, in the flask, up to 250-mL mark with deionized water.

10.2.8 After the filtrate is brought to volume, put stopper in flask to prevent evaporation.

**10.3 Digestion**

10.3.1 Add 4 to 8 granular boiling stones **or glass beads**, to each 800 ml round bottoms flask.

10.3.2 Add 2 Kjeldahl tablets to each round bottom flask.

10.3.3 Add filter paper to each round bottom flask making sure that the filter paper reaches the bottom of each flask.

10.3.4 Add 40 ml of concentrated sulfuric acid to each flask. Swirl flask until all filter paper is thoroughly dissolved.

10.3.5 Turn on fume hood and rack ventilation systems.

10.3.6 Digest samples on the digestion rack for approximately 75 minutes.

10.3.7 If any viscous black\brown material forms on top of round bottom flask, rotate flask on digestion rack to facilitate the digestion of material.

10.3.8 After digestion is complete allow samples to cool to room temperature.

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 501**

*Sample Preparation and Analysis of Water Insoluble Nitrogen*  
*[FM 445]*

Version: 2.0

Effective Date: 01/16/2024

Page 6 of 8

10.3.9 Once samples have reached room temperature, fill round bottom flasks half to approximately half way with water from the hose adjacent to digestion rack.

**10.4 Acid Trap setup**

10.4.1 The preparation of the acid trap, using the acid dispenser or class A pipet, differ depending on the specific sample.

10.4.2 For Instrument Blank:

10.4.2.1 Dispense 4 ml of standard acid into 300 ml beaker labeled blank.

10.4.2.2 Bring to the 250-ml mark with deionized water.

10.4.2.3 Place the beaker, with sample, in slot 1 on the auto titrator.

10.4.3 For Laboratory Control Sample

10.4.3.1 Dispense 4 mL of standard acid in 300-mL beaker.

10.4.3.2 Bring to the 50ml mark with deionized water.

10.4.3.3 Add 1-2 drops of Methyl Red indicator. This is the LS.

10.4.4 For Sample:

10.4.4.1 Dispense the applicable volume of standard acid into the 300-ml beaker. *This volume should be based upon the acid trap volume that was calculated in LIMS earlier in the process.*

12.4.4.2 Bring to the 50-mL mark with deionized water.

12.4.4.3 Add 1-2 drops of methyl red indicator.

12.4.4.4 Make extra 4-ml aliquots of standard acid. *This standard acid can be added to distillation acid trap if indicator starts to turn orange/yellow.*

**10.5 Distillation**

10.5.1 Place acid traps into its corresponding spots on the distillation rack with the tip of the distillation tube visibly submerged in the liquid.

10.5.2 Turn on cooling jackets and fume hood.

10.5.3 Set heaters on high and allow to heat up.

10.5.4 Swirl round bottom flasks vigorously to ensuring sample is mixed thoroughly.

10.5.5 Slowly add 80mL of 50% sodium hydroxide to each flask

10.5.6 Immediately after sodium hydroxide is added, attach flask to the distillation column and place on heating rack.

10.5.7 Allow distillation to proceed until acid trap volume reaches 250-mL mark.

**If acid trap starts to turn yellow/orange, add extra 4ml aliquot of standard acid. Record how much acid was added in notebook.**

10.5.8 Remove the acid trap from the distillation column.

10.5.9 Turn off heating racks

10.5.10 Once samples have stopped boiling, turn off cooling jackets **and fume hood.**

**10.6 Titration**

10.6.1 Start Pc titrate program.

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 501**

*Sample Preparation and Analysis of Water Insoluble Nitrogen*  
*[FM 445]*

Version: 2.0

Effective Date: 01/16/2024

Page 7 of 8

10.6.2 Once program opens, load “Ammonium” template. NOTE: The icon is two test tubes.

10.6.3 Delete any extra spaces from the template. The amount of spaces needed in the template is one more than the number of samples. If there are 11 samples there should be 12 spaces in template.

10.6.4 Enter sample identification numbers into template.

10.7 The blank should be in position 1. Label position 1 with the date and user initials.

10.8 The laboratory control sample should be in position 2.

10.9 Place acid trap into corresponding locations on auto titration rack.

10.10 Select position 1 in the template. NOTE: This is where the auto titrator will start the titration process.

10.11 Make sure all acid traps correspond to their location on the template.

10.12 Click start and ensure that the titration process begins.

10.13 Sign, date and file the instrument’s print out of titration results.

**10.14 Data Evaluation**

10.14.1 Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs, to evaluate the data against the established criteria. Evaluations, that supersede those of DATA 051, are addressed in this document.

10.14.2 Evaluate the Lab Control sample results.

10.14.3 To determine if it is within acceptable limits refer to Data 051 Attachment B, Fertilizer Laboratory Control Sample Acceptable Limits.

**10.15 Result Entry**

10.15.1 Record the results (from the print out) into the FERT\_NITROGEN\_KJL template in the LIMS. A laboratory notebook may be used as an intermediate location for the results.

10.15.2 Any additional acid used during the distillation must be added to the acid trap value in the LIMS. This updated value will be used in calculations.

10.15.3 The raw data, from the auto titrator, is uploaded to the folder on the “\\tlhaeschrome” drive.

**11 Waste Management**

11.1 Material, marked for disposal, is removed via the following mechanism:

| Material                              | Disposal Mechanism                                                     |
|---------------------------------------|------------------------------------------------------------------------|
| Boiling stones                        | Trash                                                                  |
| Residue (from the filtration process) | Trash                                                                  |
| Sample extract, excess                | Add to KJL waste, neutralize with acid, and flushed down the sink with |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 501**

*Sample Preparation and Analysis of Water Insoluble Nitrogen*  
*[FM 445]*

Version: 2.0

Effective Date: 01/16/2024

Page 8 of 8

running water.

## 12 REFERENCES

- 12.1 E00151, Calibration and use of QC-Titrate Titration system, FSFL-SOP 494
- 12.2 E00150, Kjeldahl Distillation of Fertilizer Samples, FSFL-SOP 493
- 12.3 E00152, Water Insoluble Nitrogen - Kjeldahl, FM-445
- 12.4 E00153, Breakdown Nitrogen Preparation. FM-430
- 12.5 E00149, Kjeldahl Digestion of Fertilizer Samples, FSFL-SOP 492

## 13 VERSION HISTORY

| Version | Date       | Description                                                                                                                          | Author        | Editor                        |
|---------|------------|--------------------------------------------------------------------------------------------------------------------------------------|---------------|-------------------------------|
| 1.0     | 2017.11.16 | Original Version                                                                                                                     | Zachary Tower |                               |
| 2.0     | 01.16.2024 | Added reagents<br>Updated formatting<br>Minor clarifications throughout<br>Added external document numbers<br>Removed signature page |               | Amy Bryant<br>Katherine Malik |

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 503**  
*Analysis of Total Nitrogen using LECO Combustion Analyzer*  
*[FM 411]*

Version: 2.0

Effective Date: 1/15/2021

Page 1 of 5

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for determination of total nitrogen content in fertilizer samples by the LECO Total Nitrogen Analyzer.

## **2. SCOPE**

2.1. The procedure delineated below is applicable to fertilizer samples in the Bureau of Agricultural Environmental Laboratories (BAEL), requiring the analysis of total nitrogen.

2.2. Nitrogen is released from a fertilizer through combustion at high temperature utilizing ultra-high purity (99.99%) oxygen. The nitrogen is quantitatively measured by a thermal conductivity detector and converted to w/w percent nitrogen in the sample using the appropriate mechanism.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## **4. OUTLINE of PROCEDURE**

- 4.1. 9.1 Sample Handling
- 4.2. 9.2 Calibration Standards
- 4.3. 9.3 Instrument Startup
- 4.4. 9.4 Sample Analysis
- 4.5. 9.5 Troubleshooting the Instrument
- 4.6. 9.6 Exporting Data
- 4.7. 9.7 Instrument Shutdown
- 4.8. 9.8 Data Reduction and Reporting
- 4.9. 9.9 Quality Control

## **5. INTERFERENCE**

5.1. N/A

## **6. APPARATUS and MATERIAL**

- 6.1. Analytical Instrument
  - 6.1.1. LECO TruSpec N Nitrogen Determinator
  - 6.1.2. Sample Carousel
  - 6.1.3. Balance
- 6.2. Tin Foil Cups
- 6.3. Tin Capsules

## **7. REAGENTS and MEDIA**

- 7.1. EDTA (Nitrogen standard)
- 7.2. Sucrose
- 7.3. SRM695
- 7.4. Helium, Ultra Pure (Carrier)
- 7.5. Oxygen (Combustion)
- 7.6. Air, Compressed

## **8. ASSOCIATED DOCUMENTS**

- 8.1. E00059, TruSpecN: Nitrogen Determinator Instruction Manual (March 2006)

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 503  
*Analysis of Total Nitrogen using LECO Combustion Analyzer*  
*[FM 411]*

Version: 2.0

Effective Date: 1/15/2021

Page 2 of 5

8.2. LABOP 140 Attachment G, Maintenance, Operation and Procedure of The LECO Total Combustion Analyzer

8.3. Data Evaluation for Feed, Fertilizer, and Pesticide Residue Programs, DATA 051.

## **9. SPECIFIC PROCEDURES**

### **9.1. Sample Handling**

Weigh ~0.1000g of sample to tin foil cup and add ~0.2000g of sucrose to each. Double the weight for 2X liquid dilutions. Different weights may be necessary for high and low guarantees.

### **9.2. Calibration Standards**

#### **9.2.1. Recommended Standard Levels:**

9.2.1.1. CAL 5, 0.3950-0.4050g of EDTA

9.2.1.2. CAL 4, 0.2950-0.3050g of EDTA

9.2.1.3. CAL 3, 0.1950-0.2050g of EDTA

9.2.1.4. CAL 2, 0.0950-0.1050g of EDTA

9.2.1.5. CAL 1, 0.0450-0.0550g of EDTA

9.2.2. Make sure that the purity percent is updated in the standard if there is a change of purity in the new standard.

### **9.3. Instrument Startup**

9.3.1. Instrument is normally left on. However, when needed, turn power on and make sure computer is on.

9.3.2. Once powered up, open software and allow it to connect to instrument.

9.3.3. Check temperature setpoints and allow to come to temperature. Could take up to 12 hours to stabilize.

9.3.4. Clean the sample carousel and replace the furnace filter

9.3.5. In software, run System Check for leaks.

#### **9.3.6. Analysis Parameters:**

9.3.6.1. Combustion Furnace Temperature 950° C

9.3.6.2. Afterburner Temperature 850° C

#### **9.3.7. Element Parameters: Nitrogen**

9.3.7.1. Minimal Analysis Time 30 seconds

9.3.7.2. Comparator Level 1.00

9.3.7.3. Endline Time 1 second

9.3.7.4. Conversion Factor 1.00

9.3.7.5. Significant Digits 5

9.3.7.6. TC Baseline Delay Time 5 seconds

9.3.7.7. TC Baseline Time 2 seconds

#### **9.3.8. Burn Profile:**

9.3.8.1. Burn Steps 1

9.3.8.2. Time 60 seconds

9.3.8.3. Furnace Flow High

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 503  
*Analysis of Total Nitrogen using LECO Combustion Analyzer*  
*[FM 411]*

Version: 2.0

Effective Date: 1/15/2021

Page 3 of 5

9.3.9. Macro Ballast Parameters: Ballast

9.3.9.1. Equilibrate Time 30 seconds

9.3.9.2. Not Filled Timeout 300 seconds

9.3.10. Aliquot Loop

9.3.10.1. Fill Time 20 seconds

9.3.10.2. Equilibrate Pressure Time 4 seconds

9.4. Sample Analysis

9.4.1. Place the standards and samples (tin foil cup or tin capsule) into sample carousel following the sequence

9.4.1.1. Five blanks must be injected at the beginning of a sample sequence to ensure that the system is free of contaminants.

9.4.1.2. Analyze the 5 point calibration curve.

9.4.1.3. In lieu of full calibration, a drift correction may be used to adjust the original calibration response to match the current instrument condition and response.

9.4.1.4. See DATA 051 for further information on sample sequence, initial calibration, and CCV requirements.

9.4.2. Samples require the proper laboratory information management system (LIMS) codes to allow required uploading.

9.4.2.1. Example; AA12345 411-0

Where: 0 corresponds to the sample weight (0, 1, 2)

9.4.3. Continuing Calibration Check Standard:

9.4.3.1. A CCV must be analyzed after every 10 samples and at the end of the sequence.

9.4.3.2. The CCV may consist of any level used in the initial calibration. Levels near the middle of the curve are recommended.

9.4.3.3. If a CCV response is not within the +/-2% limits, all samples analyzed before and after it must be reanalyzed.

9.4.3.4. If the area reading of any sample exceeds the calibration range, then weigh less sample and reanalyze.

9.5. Troubleshooting the Instrument

9.5.1. If instrument maintenance or troubleshooting is required, the following documents can be used for guidance.

9.5.1.1. E00059, TruSpecN: Nitrogen Determinator Instruction Manual (March 2006)

9.5.1.2. LABOP 140 Attachment G, Maintenance, Operation and Procedure of The LECO Total Combustion Analyzer

9.6. Exporting Data

9.6.1. Highlight the data and export to the flash drive.

9.6.2. Transfer the data from flash drive to the network drive TLHAESCHROME.

9.6.3. Upload the data from the network drive to the LIMS.

9.7. Instrument Shutdown

9.7.1 Refer to: LABOP 140 Attachment G, Maintenance, Operation and Procedure of The LECO Total Combustion Analyzer.

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 503**  
***Analysis of Total Nitrogen using LECO Combustion Analyzer***  
***[FM 411]***

Version: 2.0

Effective Date: 1/15/2021

Page 4 of 5

**9.8. Data Reduction and Reporting**

9.8.1. The Laboratory Control Sample must be evaluated first and if it fails, the associated samples in the batch are not uploaded.

9.8.2. Data for all passing batches is uploaded into LIMS and evaluated. Refer to Feed, Fertilizer, and Pesticide Residue Programs to evaluate the data against the established criteria.

**9.9. Quality Control**

9.9.1. Refer to QUALITY 112, Quality Control, to determine the outcome of Quality Control (QC) samples. Guidelines that supersede those of QUALITY 112 are addressed in this document.

**10. WASTE MANAGEMENT**

10.1. Samples can be discarded in regular waste receptacles.

**11. REFERENCES**

11.1. Nitrogen Analysis. FM-410 (01/29/2003)

11.2. Hazardous Waste Disposal. PLS\_SOP 009, Rev. 1.0 (02/09/2006)

11.3. Laboratory Safety, SAFETY 001

11.4. Using the Laboratory Information Management System. ADMIN 023

11.5. E00060, Nitrogen Determination in Fertilizers, LECO Corporation, 06/09

11.6. AOAC 933.13, Nitrogen (total) in Fertilizer – Combustion Method

**12. VERSION HISTORY**

| Version | Date       | Description                                                                                                                                                  | Author     | Editor           |
|---------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------------------|
| 1.0     | 12/28/2018 | Original Version                                                                                                                                             | Carl Huang |                  |
| 2.0     | 12/02/2020 | Changed CCV frequency and changed Shutdown procedures to refer to LABOP 140 attch G<br>Replaced Elizabeth Schneider with Quentin Cooper in Signature History |            | Christian Amason |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 503  
*Analysis of Total Nitrogen using LECO Combustion Analyzer*  
*[FM 411]*

Version: 2.0

Effective Date: 1/15/2021

Page 5 of 5

**13. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director

Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director

Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer

Date\_\_\_\_\_

\_\_\_\_\_  
Christian Amason  
Instrumental Manager/Editor

Date\_\_\_\_\_

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 504

#### *Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for Sample with Urea Absent [FM 420]*

Version: 2.0

Effective Date: 10/02/2020

Page 1 of 7

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the determination of direct ammoniacal nitrogen in fertilizer that does not contain urea, water soluble organic nitrogen, or water insoluble nitrogen using the Kjeldahl method. This is a direct method and does not undergo digestion.

## 2. SCOPE

2.1. The analytical procedure delineated below is applicable to the direct determination of ammoniacal nitrogen in fertilizer samples that do not contain urea or water insoluble nitrogen by making the sample basic, distilling into an acid trap and titrating excess acid with a standard base.

2.2. Refer to BAEL METHOD 530 for an alternative sample preparation and analysis of direct ammoniacal nitrogen with urea absent using the Buchi.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required. Additional PPE that may be used for this procedure includes face shield and acid resistant apron.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Weighing
- 4.2. 10.2 Acid Trap Setup
- 4.3. 10.3 Distillation
- 4.4. 10.4 Titration
- 4.5. 10.5 Data Evaluation
- 4.6. 10.6 Result Entry

## 5. INTERFERENCE

- 5.1. N/A

## 6. APPARATUS and MATERIAL

- 6.1. Balance
- 6.2. Beaker, 300mL
- 6.3. Acid Dispenser
- 6.4. Class A Pipettes
- 6.5. Auto Titration System
- 6.6. Computer
- 6.7. Kjeldahl Distillation Unit
- 6.8. Flask, round bottom 800mL
- 6.9. Weigh containers

## 7. REAGENTS and MEDIA

- 7.1. Ammonium sulfate
- 7.2. Bromocresol Green, indicator
- 7.3. Boiling stones
- 7.4. Hallow Glass Beads
- 7.5. 0.5% Methyl Red indicator
- 7.6. Sodium hydroxide, 50%
- 7.7. Sodium hydroxide, 0.5N, Certified
- 7.8. Water, Deionized (DI), Type II

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 504**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea Absent [FM 420]*

Version: 2.0

Effective Date: 10/02/2020

Page 2 of 7

**8. REAGENT PREPARATION GUIDELINES**

**8.1. Reagent 01: 0.5N Sulfuric Acid (Standard Acid)**

8.1.1. Directions: Add approximately 20L of DI Water into a 50L plastic carboy marked at 41L. Slowly add 580mL of sulfuric acid, bring to 41L with DI water. Bubble air through this solution until thoroughly mixed.

8.1.2. Titration of Sulfuric Acid Solution (Option 1): Dry tris-hydroxymethyl aminomethane (THAM) in oven at 105°C for at least 2 hours. Cool it down in the desiccator. Weigh 0.4500g of THAM into each of the six 125mL Erlenmeyer flasks. Add 5 drops of Bromocresol Green to each flask. Add approximately 50mL of DI water to each flask. The color of the solution should be blue. If not blue, add more bromocresol green a drop at a time to the flask. Add 0.5N sulfuric acid solution to 25mL dispensing burette using a pipet. Slowly titrate THAM with sulfuric acid to bromocresol green end point (pH 4.7). Discard the highest and lowest results and average the remaining four results.

8.1.3. Determine the concentration (normality) of each of the six sulfuric acid solutions using the formula:

$$\text{Acid Normality} = \frac{[\text{THAM Weight (g)}] \times 8.2549}{\text{Acid Volume (mL)}}$$

$$\text{NOTE: } 8.2549 = \frac{1000\text{mL/L} \times 2(\text{Conversion from Molarity to Normality})}{121.14\text{g/mol} \times 2(\text{Stoichiometric Ratio})}$$

8.1.4. Titration of Sulfuric Acid Solution (Option 2): Add 4mL of Reagent 01 to seven (7) 300mL beakers then fill them to approximately 250mL with water. Fill a separate 300mL beaker with 200mL of certified 0.5N NaOH; place the line of the titrator into the beaker and purge the syringe of the titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.1.4.1. Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Acid 1 – Acid 6.

8.1.5. Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$A_N = \frac{B_N \times B_V}{A_V}$$

Where:  $A_N$ =Acid Normality  
 $A_V$ =Acid Volume  
 $B_N$ =Base Normality  
 $B_V$ =Base Volume

8.1.6. Discard the highest and lowest results and average the remaining results.

**8.2. Reagent 02: 18.75N Stock Sodium Hydroxide Solution**

8.2.1. Directions: Add 1 liter of deionized (DI) water into a 2L volumetric flask with stir bar in it on a stirring plate. Add 1500g of sodium hydroxide pellets in small (~100g) portions with continuous mixing. Bring solution volume to 2L with DI water.

**8.3. Reagent 03: 0.35N Sodium Hydroxide (Standard Base)**

8.3.1. Directions: Fill a 20L glass carboy, marked at 16L, half full of DI water. Add stir bar and start mixing. Add 300mL of stock sodium hydroxide solution to carboy while continuously mixing. Bring solution to a volume of 16L with DI water. Stir solution for approximately 30 minutes.

8.3.2. Titration of Sodium Hydroxide Solution (Option 1): Dry potassium acid phthalate (PAP) in oven at 105°C for at least 2 hours and then allow it to cool in desiccators. Take six 300mL beakers

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 504

#### *Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for Sample with Urea Absent [FM 420]*

Version: 2.0

Effective Date: 10/02/2020

Page 3 of 7

and weigh approximately 0.8000g of PAP into each. Add 150mL of DI water into each of the six beakers. Allow PAP to dissolve. Titrate PAP solution on the QC\_Titrate System.

8.3.3. Titration of PAP Solution on the QC-Titrator: Click on PAP icon. Type set name and click OK. Type the exact weights of PAP. Highlight the first line on time-table and click start. When titration is done, a hard copy is printed automatically.

8.3.4. Calculation of Sodium Hydroxide: Determine the concentration (normality) of the six sodium hydroxide solutions using the formula:

$$\text{Base Normality} = \frac{[\text{PAP Weight (g)}] \times 4.8964}{\text{Base Volume (mL)}}$$

$$\text{NOTE: } 4.8964 = \frac{1000\text{mL/L} \times 1(\text{Conversion from Molarity to Normality})}{204.23\text{g/mol} \times 1(\text{Stoichiometric Ratio})}$$

8.3.5. Titration of Sodium Hydroxide Solution (Option 2): Add 4mL of Reagent 03 to seven (7) 300mL beakers. Fill a separate 300mL beaker with 200mL of the normalized solution from 8.1; place the line of the titrator into the beaker and purge the syringe of the titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.3.5.1. Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Base 1 – Base 6.

8.3.6. Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$B_N = \frac{A_N \times A_V}{B_V}$$

Where:  $A_N$ =Acid Normality  
 $A_V$ =Acid Volume  
 $B_N$ =Base Normality  
 $B_V$ =Base Volume

8.3.7. Discard the highest and lowest results and average the remaining four results. Change base normality on Kjeldahl and Direct Excel spreadsheets for LIMS.

8.4. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015

8.5. The specific ingredient amounts can vary to accommodate the batch size.

8.6. Expiration Date: 12 months from date of preparation unless otherwise specified.

8.7. Store at room temperature

## 9. ASSOCIATED DOCUMENTS

- 9.1. ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates
- 9.2. DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs

## 10. SPECIFIC PROCEDURES

### 10.1. Weighing

10.1.1. Using ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Template, create a laboratory identification number for the laboratory control sample (LS).

10.1.2. Weigh approximately 0.5g of the laboratory control standard (LS), ammonium sulfate.

10.1.3. Weigh samples to approximately 0.5 grams as well except liquid samples that have been diluted. Diluted liquid samples are approximately doubled in weight (i.e. approximately 1.00g).

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 504**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea Absent [FM 420]*

Version: 2.0

Effective Date: 10/02/2020

Page 4 of 7

10.1.4. Enter weight into Laboratory Information Management System (LIMS) using the template FERT\_Nitrogen\_Direct. *Be aware of the given acid trap volume because it will be used later in the process.*

10.1.5. Add 4 to 8 boiling stones or glass beads to each round bottom flask.

10.1.6. Rinse sample with deionized water into round bottom flask from weight container a minimum of 3 times. If sample is still present in weight container, continue to rinse.

10.1.7. Fill round bottom flasks half to approximately half way with tap water.

10.1.8. Turn on the fume hood and both sides of the cooling racks.

10.1.9. Preheat each distillation unit on high.

**10.2. Acid Trap Setup**

10.2.1. Measure approximately 10mL of 50% sodium hydroxide and set aside for each sample.

10.2.2. The preparation of the acid trap, using the acid dispenser or a class A pipette, differ depending on the specific sample.

10.2.2.1. For Instrument Blank:

10.2.2.1.1. Dispense 4mL of standard acid into 300mL beaker labeled blank.

10.2.2.1.2. Bring to the 250mL mark with deionized water.

10.2.2.1.3. Place the beaker, with sample, in slot 1 on the auto titrator.

10.2.2.2. For Laboratory Control Sample

10.2.2.2.1. Dispense 18mL of standard acid in 300mL beaker.

10.2.2.2.2. Bring to the 50mL mark with deionized water.

10.2.2.2.3. Add 1-2 drops of Methyl Red indicator. This is the LS.

10.2.2.3. For Sample:

10.2.2.3.1. Dispense the applicable volume of standard acid into the 300mL beaker. *This volume should be based upon the acid trap volume that was calculated in LIMS earlier in the process.*

10.2.2.3.2. Bring to the 50mL mark with deionized water.

10.2.2.3.3. Add 1-2 drops of methyl red indicator.

10.2.3. Make extra 4mL aliquots of standard acid. This standard acid can be added to distillation acid trap if indicator starts to turn orange/yellow.

**10.3. Distillation**

10.3.1. Place acid traps into its corresponding spots on the distillation rack with the tip of the distillation tube visibly submerged in the liquid.

10.3.2. Turn on cooling jackets and fume hood.

10.3.3. Set heaters on high and allow to heat up.

10.3.4. Swirl round bottom flasks vigorously to ensuring sample is mixed thoroughly.

10.3.5. Quickly add sodium hydroxide to each flask and immediately attach flask to the distillation column and place on heating rack.

10.3.6. Allow distillation to proceed until acid trap volume reaches 250mL.

**NOTE:** If the color starts to change to yellow/orange, add 4mL of the standard acid, write down how much extra acid was added to sample in notebook.

10.3.7. Remove the acid trap from the distillation column.

10.3.8. Turn off heating racks and fume hood once sample stop boiling.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 504

#### *Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for Sample with Urea Absent [FM 420]*

Version: 2.0

Effective Date: 10/02/2020

Page 5 of 7

10.3.9. Once samples have stopped boiling, turn off cooling jackets.

#### 10.4. Titration

10.4.1. Start Pc titrate program.

10.4.2. Once program opens, load "Ammonium" template.

**NOTE:** The icon is two test tubes.

10.4.3. Delete any extra spaces from the template. The amount of spaces needed in the template is one more than the number of samples. If there are 11 samples there should be 12 spaces in template.

10.4.4. Enter sample identification numbers into template.

10.4.4.1. The blank should be in position 1. Label position 1 with the date and user initials.

10.4.4.2. The laboratory control sample should be in position 2.

10.4.5. Place acid trap into corresponding locations on auto titration rack.

10.4.6. Select position 1 in the template.

**NOTE:** This is where the auto titrator will start the titration process.

10.4.7. Make sure all acid traps correspond to their location on the template.

10.4.8. Click start and ensure that the titration process begins.

10.4.9. **Initial,** date and file the instrument's print out of titration results.

#### 10.5. Data Evaluation

10.5.1. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs, to evaluate the data against the established criteria. Evaluations, that supersede those of DATA 051, are addressed in this document.

10.5.1.1. Evaluate the Lab Control sample results.

10.5.2. To determine if it is within acceptable limits refer to Data 051 Attachment B, Fertilizer Laboratory Control Sample Acceptable Limits.

#### 10.6. Result Entry

10.6.1. Record the results (from the print out) into the FERT\_NITROGEN **Direct** template in the LIMS. A laboratory notebook may be used as an intermediate location for the results.

10.6.2. Any additional acid used during the distillation must be added to the acid trap value in the LIMS. This updated value will be used in calculations.

10.6.3. The raw data, from the auto titrator, is uploaded to the folder on the "\\tlhaeschrome" drive.

## 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material                              | Disposal Mechanism                                                                                                     |
|---------------------------------------|------------------------------------------------------------------------------------------------------------------------|
| Boiling stones                        | Trash                                                                                                                  |
| Residue (from the filtration process) | Trash                                                                                                                  |
| Sample extract, excess                | Pour waste into KJL waste container, then neutralize with acid. Flush neutralized waste down drain with running water. |

## 12. REFERENCES

12.1. Calibration and use of QC-Titrate™ Titration system. FSFL-SOP 494, (2006/01/17)

12.2. Nitrogen (urea not present)-Direct FM-420, (2003/02/05)

12.3. Distillation. SOP-493, (2006/01/17)

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 504**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea Absent [FM 420]*

Version: 2.0

Effective Date: 10/02/2020

Page 6 of 7

**13. VERSION HISTORY**

| <b>Version</b> | <b>Date</b> | <b>Description</b>                                                                                                                                                                                                | <b>Author</b> | <b>Editor</b>   |
|----------------|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------|
| 1.0            | 11.16.2017  | Original Version                                                                                                                                                                                                  | Zachary Tower |                 |
| 2.0            | 09.30/2020  | Added section 8.0 Reagent Preparation Guide<br>Updated Scope<br>Updated Apparatus and Materials<br>Updated Reagents and Media<br>Minor clarifications throughout<br>Updated disposal mechanism for sample extract |               | Katherine Malik |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 504

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea Absent [FM 420]*

Version: 2.0

Effective Date: 10/02/2020

Page 7 of 7

**14. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Bureau Chief, Technical Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer  
Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Wet Chemistry and Preparation Manager  
Date\_\_\_\_\_

\_\_\_\_\_  
Katherine Malik  
Editor  
Date\_\_\_\_\_

# **STANDARD OPERATING PROCEDURE**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

### **METHOD 505**

#### *Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for Sample with Urea Present [FM 415]*

Version: 2.0

Effective Date: 10/02/2020

Page 1 of 7

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for the determination of direct ammoniacal nitrogen for samples that may contain urea or water soluble organic nitrogen using the Kjeldahl method. This is a direct method and does not undergo digestion.

## **2. SCOPE**

2.1. The procedure delineated below is applicable to the direct determination of ammoniacal nitrogen in fertilizer samples that contain urea nitrogen by making the sample basic, distilling into an acid trap and titrating excess acid with a standard base.

2.2. Refer to BAEL METHOD 530 for an alternative sample preparation and analysis of direct ammoniacal nitrogen with urea present using the Buchi.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required. Additional PPE that may be used for this procedure includes face shield and acid resistant apron.

## **4. OUTLINE of PROCEDURE**

- 4.1. 10.1 Weighing
- 4.2. 10.2 Acid Trap Setup
- 4.3. 10.3 Distillation
- 4.4. 10.4 Titration
- 4.5. 10.5 Data Evaluation
- 4.6. 10.6 Result Entry

## **5. INTERFERENCE**

- 5.1. N/A

## **6. APPARATUS and MATERIAL**

- 6.1. Acid Dispenser
- 6.2. Auto Titration system
- 6.3. Balance
- 6.4. Beaker, 300mL
- 6.5. Class A Pipettes
- 6.6. Computer
- 6.7. Kjeldahl Distillation unit
- 6.8. Round bottom flask, 800mL
- 6.9. Weigh containers

## **7. REAGENTS and MEDIA**

- 7.1. Ammonium sulfate
- 7.2. Bromocresol Green, indicator
- 7.3. Boiling stones
- 7.4. Hallow Glass Beads
- 7.5. Magnesium Oxide
- 7.6. 0.5% Methyl Red, indicator
- 7.7. Sodium hydroxide, 50%
- 7.8. Sodium hydroxide, 0.5N, Certified

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 505**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea Present [FM 415]*

Version: 2.0

Effective Date: 10/02/2020

Page 2 of 7

7.9. Sulfuric Acid Reagent Grade

7.10. Water, Deionized (DI) Type II

## **8. REAGENT PREPARATION GUIDELINES**

### **8.1. Reagent 01: 0.5N Sulfuric Acid (Standard Acid)**

8.1.1. Directions: Add approximately 20L of DI Water into a 50L plastic carboy marked at 41L. Slowly add 580mL of sulfuric acid, bring to 41L with DI water. Bubble air through this solution until thoroughly mixed.

8.1.2. Titration of Sulfuric Acid Solution (Option 1): Dry tris-hydroxymethyl aminomethane (THAM) in oven at 105°C for at least 2 hours. Cool it down in the desiccator. Weigh 0.4500g of THAM into each of the six 125mL Erlenmeyer flasks. Add 5 drops of Bromocresol Green to each flask. Add approximately 50mL of DI water to each flask. The color of the solution should be blue. If not blue, add more bromocresol green a drop at a time to the flask. Add 0.5N sulfuric acid solution to 25mL dispensing burette using a pipet. Slowly titrate THAM with sulfuric acid to bromocresol green end point (pH 4.7). Discard the highest and lowest results and average the remaining four results.

8.1.3. Determine the concentration (normality) of each of the six sulfuric acid solutions using the formula:

$$\text{Acid Normality} = \frac{[\text{THAM Weight (g)}] \times 8.2549}{\text{Acid Volume (mL)}}$$

$$\text{NOTE: } 8.2549 = \frac{1000\text{mL/L} \times 2(\text{Conversion from Molarity to Normality})}{121.14\text{g/mol} \times 2(\text{Stoichiometric Ratio})}$$

8.1.4. Titration of Sulfuric Acid Solution (Option 2): Add 4mL of Reagent 01 to seven (7) 300mL beakers then fill them to approximately 250mL with water. Fill a separate 300mL beaker with 200mL of certified 0.5N NaOH; place the line of the titrator into the beaker and purge the syringe of the titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.1.4.1. Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Acid 1 – Acid 6.

8.1.5. Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$A_N = \frac{B_N \times B_V}{A_V}$$

Where:  $A_N$ =Acid Normality

$A_V$ =Acid Volume

$B_N$ =Base Normality

$B_V$ =Base Volume

8.1.6. Discard the highest and lowest results and average the remaining results.

### **8.2. Reagent 02: 18.75N Stock Sodium Hydroxide Solution**

8.2.1. Directions: Add 1 liter of deionized (DI) water into a 2L volumetric flask with stir bar in it on a stirring plate. Add 1500g of sodium hydroxide pellets in small (~100g) portions with continuous mixing. Bring solution volume to 2L with DI water.

### **8.3. Reagent 03: 0.35N Sodium Hydroxide (Standard Base)**

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 505

#### *Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for Sample with Urea Present [FM 415]*

Version: 2.0

Effective Date: 10/02/2020

Page 3 of 7

8.3.1. Directions: Fill a 20L glass carboy, marked at 16L, half full of DI water. Add stir bar and start mixing. Add 300mL of stock sodium hydroxide solution to carboy while continuously mixing. Bring solution to a volume of 16L with DI water. Stir solution for approximately 30 minutes.

8.3.2. Titration of Sodium Hydroxide Solution (Option 1): Dry potassium acid phthalate (PAP) in oven at 105°C for at least 2 hours and then allow it to cool in desiccators. Take six 300mL beakers and weigh approximately 0.8000g of PAP into each. Add 150mL of DI water into each of the six beakers. Allow PAP to dissolve. Titrate PAP solution on the QC\_Titrate System.

8.3.3. Titration of PAP Solution on the QC-Titrator: Click on PAP icon. Type set name and click OK. Type the exact weights of PAP. Highlight the first line on time-table and click start. When titration is done, a hard copy is printed automatically.

8.3.4. Calculation of Sodium Hydroxide: Determine the concentration (normality) of the six sodium hydroxide solutions using the formula:

$$\text{Base Normality} = \frac{[\text{PAP Weight (g)}] \times 4.8964}{\text{Base Volume (mL)}}$$

$$\text{NOTE: } 4.8964 = \frac{1000\text{mL/L} \times 1(\text{Conversion from Molarity to Normality})}{204.23\text{g/mol} \times 1(\text{Stoichiometric Ratio})}$$

8.3.5. Titration of Sodium Hydroxide Solution (Option 2): Add 4mL of Reagent 03 to seven (7) 300mL beakers. Fill a separate 300mL beaker with 200mL of the normalized solution from 8.1; place the line of the titrator into the beaker and purge the syringe of the titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.3.5.1. Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Base 1 – Base 6.

8.3.6. Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$B_N = \frac{A_N \times A_V}{B_V}$$

Where:  $A_N$ =Acid Normality  
 $A_V$ =Acid Volume  
 $B_N$ =Base Normality  
 $B_V$ =Base Volume

8.3.7. Discard the highest and lowest results and average the remaining four results. Change base normality on Kjeldahl and Direct Excel spreadsheets for LIMS.

8.4. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015

8.5. The specific ingredient amounts can vary to accommodate the batch size.

8.6. Expiration Date: 12 months from date of preparation unless otherwise specified.

8.7. Store at room temperature

## 9. ASSOCIATED DOCUMENTS

- 9.1. ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates
- 9.2. DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 505**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea Present [FM 415]*

Version: 2.0

Effective Date: 10/02/2020

Page 4 of 7

## 10. SPECIFIC PROCEDURES

### 10.1. Weighing

10.1.1. Using ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Template, create a laboratory identification number for the laboratory control sample (LS).

10.1.2. Weigh approximately 0.5g of laboratory control standard (LS), ammonium sulfate.

10.1.3. Weigh samples to approximately 0.5 grams except liquid samples that have been diluted. Diluted liquid samples are approximately doubled in weight (i.e. approximately 1.00g).

10.1.4. Enter weight into Laboratory Information Management System (LIMS) using the template FERT\_Nitrogen\_Direct. Be aware of the given acid trap volume because it will be used later in the process.

10.1.5. Add 4 to 8 boiling stones or glass beads to each round bottom flask

10.1.6. Rinse sample with DI water into round bottom flask from weigh container a minimum of 3 times. If sample is still present in weigh container, continue to rinse.

10.1.7. With the water hose adjacent to the distillation unit fill each round bottom to approximately half way.

10.1.8. Turn on the fume hood and both sides of the cooling racks.

10.1.9. Preheat each distillation unit on high.

### 10.2. Acid Trap Setup

10.2.1. Weigh approximately 6.5 grams of magnesium Oxide and set aside for each sample.

10.2.2. The preparation of the acid trap, using the acid dispenser or class A pipette, differ depending on the specific sample.

10.2.2.1. For Instrument Blank:

10.2.2.1.1. Dispense 4mL of standard acid into 300mL beaker labeled blank.

10.2.2.1.2. Bring to the 250mL mark with deionized water.

10.2.2.1.3. Place the beaker, with sample, in slot 1 on the auto titrator.

10.2.2.2. For Laboratory Control Sample

10.2.2.2.1. Dispense 18mL of standard acid in 300mL beaker.

10.2.2.2.2. Bring to the 50mL mark with deionized water.

10.2.2.2.3. Add 1-2 drops of Methyl Red indicator. This is the LS.

10.2.2.3. For Sample:

10.2.2.3.1. Dispense the applicable volume of standard acid into the 300-ml beaker. This volume should be based upon the acid trap volume that was calculated in LIMS earlier in the process.

10.2.2.3.2. Bring to the 50mL mark with deionized water.

10.2.2.3.3. Add 1-2 drops of methyl red indicator.

10.2.3. Make extra 4mL aliquots of standard acid. This standard acid can be added to distillation acid trap if indicator starts to turn orange/yellow.

### 10.3. Distillation

10.3.1. Place acid traps into their corresponding spots on the distillation rack with the tip of the distillation tube visibly submerged in the liquid.

10.3.2. Swirl round bottom flasks vigorously to ensure sample is mixed thoroughly.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 505**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea Present [FM 415]*

Version: 2.0

Effective Date: 10/02/2020

Page 5 of 7

10.3.3. Quickly add magnesium oxide to each flask and immediately attach flask to the distillation column and place on heating rack.

**NOTE:** Step must be performed quickly to avoid loss of ammonium gas.

10.3.4. Allow distillation to proceed until acid trap volume reaches 250mL.

**NOTE:** If the color starts to change to yellow/orange, add 4mL of the standard acid. Record how much extra acid is added in notebook.

10.3.5. Remove the acid trap from the distillation column and turn off heating racks, and once samples stop boiling, turn off fume hood and cooling jackets.

10.4. Titration

10.4.1. Start Pc titrate program.

10.4.2. Once program opens, load "Ammonium" template.

**NOTE:** The icon is two test tubes.

10.4.3. Delete any extra spaces from the template. The amount of spaces needed in the template is one more than the number of samples. If there are 11 samples there should be 12 spaces in template.

10.4.4. Enter sample identification numbers into template.

10.4.4.1. The blank should be in position 1. Label position 1 with the date and user initials.

10.4.4.2. The laboratory control sample should be in position 2.

10.4.5. Place acid trap into corresponding locations on auto titration rack.

10.4.6. Select position 1 in the template.

**NOTE:** This is where the auto titrator will start the titration process.

10.4.7. Make sure all acid traps correspond to their location on the template.

10.4.8. Click start and ensure that the titration process begins.

10.4.9. Initial, date and file the instrument's print out of titration results.

10.5. Data Evaluation

10.5.1. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs, to evaluate the data against the established criteria. Evaluations, that supersede those of DATA 051, are addressed in this document.

10.5.1.1. Evaluate the Lab Control sample results.

10.5.2. To determine if it is within acceptable limits refer to Data 051 Attachment B, Fertilizer Laboratory Control Sample Acceptable Limits.

10.6. Result Entry

10.6.1. Record the results (from the print out) into the FERT\_NITROGEN\_Direct template in the LIMS. A laboratory notebook may be used as an intermediate location for the results.

10.6.2. Any additional acid used during the distillation must be added to the acid trap value in the LIMS. This updated value will be used in calculations.

10.6.3. The raw data, from the auto titrator, is uploaded to the folder on the "\\tlhaeschrome" drive.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 505**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea Present [FM 415]*

Version: 2.0

Effective Date: 10/02/2020

Page 6 of 7

## 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material                              | Disposal Mechanism                                                              |
|---------------------------------------|---------------------------------------------------------------------------------|
| Boiling stones                        | Trash                                                                           |
| Residue (from the filtration process) | Trash                                                                           |
| Sample extract, excess                | Pour waste into KJL waste container. Flush waste down drain with running water. |

## 12. REFERENCES

- 12.1. Calibration and use of QC-Titrate™ Titration system. FSFL-SOP 494, (2006/01/17)
- 12.2. Breakdown Nitrogen Preparation. FM-430, (2001/11/12)
- 12.3. Ammoniacal Nitrogen (urea present)-Direct. FM-415

## 13. VERSION HISTORY

| Version | Date       | Description                                                                                                                                                                                                                                                                                | Author        | Editor          |
|---------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------|
| 1.0     | 11.16.2017 | Original Version                                                                                                                                                                                                                                                                           | Zachary Tower |                 |
| 2.0     | 09.30.2020 | Added section 8.0 Reagent Preparation Guide<br>Updated Scope<br>Updated Apparatus and Materials<br>Updated Reagents and Media<br>Minor clarifications throughout<br>Updated disposal mechanism for sample extract<br>Replaced Elizabeth Schneider with Quentin Cooper in Signature History |               | Katherine Malik |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 505**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea Present [FM 415]*

Version: 2.0

Effective Date: 10/02/2020

Page 7 of 7

**14. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Bureau Chief, Technical Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer  
Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Wet Chemistry and Preparation Manager  
Date\_\_\_\_\_

\_\_\_\_\_  
Katherine Malik  
Editor  
Date\_\_\_\_\_

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 506**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea and Water Insoluble Nitrogen Present [FM 435]*

Version: 2.0

Effective Date: 10/02/2020

Page 1 of 7

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for the preparation and analysis of direct ammoniacal nitrogen using the Kjeldahl method. This is a direct method and does not undergo digestion.

## **2. SCOPE**

2.1. The analytical procedure delineated below is applicable to the determination of direct ammoniacal nitrogen in fertilizer samples that contain urea and water insoluble nitrogen by making the sample basic, distilling into an acid trap and titrating excess acid with a standard base.

2.2. Refer to BAEL METHOD 530 for an alternative sample preparation and analysis of direct ammoniacal nitrogen with urea and water insoluble present using the Buchi.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required. Additional PPE which can be used for this procedure includes face shield and acid resistant apron.

## **4. OUTLINE of PROCEDURE**

- 4.1. 10.1 Weighing and Filtering
- 4.2. 10.2 Acid Trap Setup
- 4.3. 10.3 Distillation
- 4.4. 10.4 Titration
- 4.5. 10.5 Data Evaluation
- 4.6. 10.6 Result Entry

## **5. INTERFERENCE**

- 5.1. N/A

## **6. APPARATUS and MATERIAL**

- 6.1. Acid Dispenser
- 6.2. Auto Titration system
- 6.3. Balance
- 6.4. Beaker, 300mL
- 6.5. Class A pipettes
- 6.6. Computer
- 6.7. Flask, round bottom 800mL
- 6.8. Flask, 250mL
- 6.9. Filter paper, 125mm
- 6.10. Funnel
- 6.11. Kjeldahl Distillation unit
- 6.12. Pipette or Graduated cylinder, 100mL
- 6.13. Stoppers
- 6.14. Weigh containers

## **7. REAGENTS and MEDIA**

- 7.1. Granular Boiling stones
- 7.2. Hallow glass beads

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 506**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea and Water Insoluble Nitrogen Present [FM 435]*

Version: 2.0

Effective Date: 10/02/2020

Page 2 of 7

- 7.3. Magnesium Oxide
- 7.4. Methanol
- 7.5. 0.5% Methyl Red indicator
- 7.6. Sodium hydroxide, 50%
- 7.7. Sodium hydroxide, 0.5N Certified
- 7.8. SRM-695
- 7.9. Sulfuric Acid Reagent Grade
- 7.10. Water, Deionized (DI) Type II

## **8. REAGENT PREPARATION GUIDELINES**

### **8.1. Reagent 01: 0.5N Sulfuric Acid (Standard Acid)**

8.1.1. Directions: Add approximately 20L of DI Water into a 50L plastic carboy marked at 41L. Slowly add 580mL of sulfuric acid, bring to 41L with DI water. Bubble air through this solution until thoroughly mixed.

8.1.2. Titration of Sulfuric Acid Solution (Option 1): Dry tris-hydroxymethyl aminomethane (THAM) in oven at 105°C for at least 2 hours. Cool it down in the desiccator. Weigh 0.4500g of THAM into each of the six 125mL Erlenmeyer flasks. Add 5 drops of Bromocresol Green to each flask. Add approximately 50mL of DI water to each flask. The color of the solution should be blue. If not blue, add more bromocresol green a drop at a time to the flask. Add 0.5N sulfuric acid solution to 25mL dispensing burette using a pipet. Slowly titrate THAM with sulfuric acid to bromocresol green end point (pH 4.7). Discard the highest and lowest results and average the remaining four results.

8.1.3. Determine the concentration (normality) of each of the six sulfuric acid solutions using the formula:

$$\text{Acid Normality} = \frac{[\text{THAM Weight (g)}] \times 8.2549}{\text{Acid Volume (mL)}}$$

$$\text{NOTE: } 8.2549 = \frac{1000\text{mL/L} \times 2(\text{Conversion from Molarity to Normality})}{121.14\text{g/mol} \times 2(\text{Stoichiometric Ratio})}$$

8.1.4. Titration of Sulfuric Acid Solution (Option 2): Add 4mL of Reagent 01 to seven (7) 300mL beakers then fill them to approximately 250mL with water. Fill a separate 300mL beaker with 200mL of certified 0.5N NaOH; place the line of the titrator into the beaker and purge the syringe of the titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.1.4.1. Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Acid 1 – Acid 6.

8.1.5. Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$A_N = \frac{B_N \times B_V}{A_V}$$

Where:  $A_N$ =Acid Normality  
 $A_V$ =Acid Volume  
 $B_N$ =Base Normality  
 $B_V$ =Base Volume

8.1.6. Discard the highest and lowest results and average the remaining results.

### **8.2. Reagent 02: 18.75N Stock Sodium Hydroxide Solution**

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 506**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea and Water Insoluble Nitrogen Present [FM 435]*

Version: 2.0

Effective Date: 10/02/2020

Page 3 of 7

8.2.1. Directions: Add 1 liter of deionized (DI) water into a 2L volumetric flask with stir bar in it on a stirring plate. Add 1500g of sodium hydroxide pellets in small (~100g) portions with continuous mixing. Bring solution volume to 2L with DI water.

**8.3. Reagent 03: 0.25N Sodium Hydroxide (Standard Base)**

8.3.1. Directions: Fill a 20L glass carboy, marked at 16L, half full of DI water. Add stir bar and start mixing. Add 300mL of stock sodium hydroxide solution to carboy while continuously mixing. Bring solution to a volume of 16L with DI water. Stir solution for approximately 30 minutes.

8.3.2. Titration of Sodium Hydroxide Solution (Option 1): Dry potassium acid phthalate (PAP) in oven at 105°C for at least 2 hours and then allow it to cool in desiccators. Take six 300mL beakers and weigh approximately 0.8000g of PAP into each. Add 150mL of DI water into each of the six beakers. Allow PAP to dissolve. Titrate PAP solution on the QC\_Titrate System.

8.3.3. Titration of PAP Solution on the QC-Titrator: Click on PAP icon. Type set name and click OK. Type the exact weights of PAP. Highlight the first line on time-table and click start. When titration is done, a hard copy is printed automatically.

8.3.4. Calculation of Sodium Hydroxide: Determine the concentration (normality) of the six sodium hydroxide solutions using the formula:

$$\text{Base Normality} = \frac{[\text{PAP Weight (g)}] \times 4.8964}{\text{Base Volume (mL)}}$$

$$\text{NOTE: } 4.8964 = \frac{1000\text{mL/L} \times 1(\text{Conversion from Molarity to Normality})}{204.23\text{g/mol} \times 1(\text{Stoichiometric Ratio})}$$

8.3.5. Titration of Sodium Hydroxide Solution (Option 2): Add 4mL of Reagent 03 to seven (7) 300mL beakers. Fill a separate 300mL beaker with 200mL of the normalized solution from 8.1; place the line of the titrator into the beaker and purge the syringe of the titrator a minimum of two times into a waste container. Load and run the "Ammonium" template. When titration is complete initial, date, and file results.

8.3.5.1. Enter the sample identification numbers into the template. The Blank will always be in Position 1. The rest of the beakers should be labeled Base 1 – Base 6.

8.3.6. Determine the concentration (normality) of each of the seven sulfuric acid solutions using the formula:

$$B_N = \frac{A_N \times A_V}{B_V}$$

Where:  $A_N$ =Acid Normality  
 $A_V$ =Acid Volume  
 $B_N$ =Base Normality  
 $B_V$ =Base Volume

8.3.7. Discard the highest and lowest results and average the remaining four results. Change base normality on Kjeldahl and Direct Excel spreadsheets for LIMS.

8.4. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015

8.5. The specific ingredient amounts can vary to accommodate the batch size.

8.6. Expiration Date: 12 months from date of preparation unless otherwise specified.

8.7. Store at room temperature

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 506**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea and Water Insoluble Nitrogen Present [FM 435]*

Version: 2.0

Effective Date: 10/02/2020

Page 4 of 7

## 9. ASSOCIATED DOCUMENTS

- 9.1. ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Templates
- 9.2. DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs

## 10. SPECIFIC PROCEDURES

### 10.1. Weighing and Filtering

- 10.1.1. Label weigh containers making the first container the laboratory control sample. Ensure that the beakers, flask, distillation and auto titration sample slots are labeled corresponding to the weigh containers.
- 10.1.2. Using ADMIN 023 Attachment AA, Fertilizer QC Loc Codes and Result Entry Template, create a laboratory identification number for the laboratory control sample (LS).
- 10.1.3. Weigh approximately 0.5g of laboratory control standard (LS), SRM-695. Measure all solid samples to approximately 0.9 grams. Liquid samples should be doubled in weight to approximately 1.8 grams if they have been diluted.
- 10.1.4. Enter weight into Laboratory Information management system (LIMS) using the template FERT\_Nitrogen\_KJL and record the given acid trap volume in the laboratory notebook to use later.
- 10.1.5. Fold filters and place into funnels which drain into 250mL volumetric flasks. Wet sides of the filter paper with deionized water to create a seal for the filter.
- 10.1.6. Pour weighed samples into filter paper using deionized water to rinse the weight container clean of sample into filter paper.
- 10.1.7. If sample is clumping, use 2 to 5mL of methanol to break up clumping
- 10.1.8. Fill the filter paper up with DI water and allow to drain completely between washings. Repeat washing step 2 more times for a total of 3 washes.
- 10.1.9. Bring the filtrate receiving flask up to 250mL volume mark with DI water.
- 10.1.10. Put stopper in after flask is brought to volume to prevent evaporation.
- 10.1.11. Add 4 to 8 boiling stones or glass beads, to each round bottom flask
- 10.1.12. With filtrate from the filtering process measure out approximately 100mL of filtrate into a clean 100mL graduated cylinder or pipette and place into 800mL round bottom. (Make sure to invert filtrate bottle at least 7 times and shake thoroughly to insure sample is uniform before measuring).
- 10.1.13. With the water hose adjacent to the distillation unit, fill each round bottom to approximately half way.
- 10.1.14. Turn on the fume hood and both side of the cooling racks.
- 10.1.15. Preheat each distillation unit on high.

### 10.2. Acid Trap Setup

- 10.2.1. Measure approximately 6.5 grams of magnesium oxide and set aside for each sample.
- 10.2.2. The preparation of the acid trap, using the acid dispenser or class A pipette, differ depending on the specific sample.
  - 10.2.2.1. For Instrument Blank:
    - 10.2.2.1.1. Dispense 4mL of standard acid into 300mL beaker labeled blank.
    - 10.2.2.1.2. Bring to the 250mL mark with deionized water.
    - 10.2.2.1.3. Place the beaker, with sample, in slot 1 on the auto titrator.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 506**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea and Water Insoluble Nitrogen Present [FM 435]*

Version: 2.0

Effective Date: 10/02/2020

Page 5 of 7

10.2.2.2. For Laboratory Control Sample

10.2.2.2.1. Dispense 4mL of standard acid in 300mL beaker.

10.2.2.2.2. Bring to the 50mL mark with deionized water.

10.2.2.2.3. Add 1-2 drops of Methyl Red indicator. This is the LS.

10.2.2.3. For Sample:

10.2.2.3.1. Dispense the applicable volume of standard acid into the 300mL beaker.

**NOTE:** This volume should be based upon the acid trap volume that was calculated in LIMS earlier in the process.

10.2.2.3.2. Bring to the 50mL mark with deionized water.

10.2.2.3.3. Add 1-2 drops of methyl red indicator.

10.2.3. Make extra 4mL aliquots of standard acid. This standard acid can be added to distillation acid trap if indicator starts to turn orange/yellow.

10.3. Distillation

10.3.1. Place acid traps into their corresponding spots on the distillation rack with the tip of the distillation tube visibly submerged in the liquid.

10.3.2. Make sure both cooling jackets are turned on and working and exhaust fan is turned on.

10.3.3. Set heaters on high and allow time to heat up.

10.3.4. Swirl round bottom flasks vigorously insuring sample is mixed thoroughly and then add magnesium oxide to each flask.

10.3.5. Immediately after magnesium oxide is added attach flask to the distillation column and place on heating rack.

10.3.6. Allow distillation to proceed until acid trap volume reaches 250mL. Once volume is reached remove the acid trap for the distillation column and turn off heating racks. When samples stop boiling turn off exhaust vent and cooling jackets.

10.3.7. If acid trap starts to turn orange/yellow add extra 4mL of acid. Record how much extra acid is added in notebook.

10.4. Titration

10.4.1. Start "Pc titrate" program.

10.4.2. Once program opens, load "Ammonium" template.

**NOTE:** Icon is two test tubes.

10.4.3. In template, ensure there are only the amount of spaces plus one loaded in the template (e.g. if you have 11 samples there should be 12 spaces in template). Delete any extra spaces.

10.4.4. Position number 1 should always be the instrument blank. In the "Ammonium" template, label position 1 with the date and initials.

10.4.5. Enter LS in position 2 of the template.

10.4.6. Enter sample ID numbers into template.

10.4.7. Place acid trap beakers into corresponding locations on auto titration rack.

10.4.8. Select position 1 in the template (The selected position should be highlighted blue. This is where the auto titrator will start titration).

10.4.9. Click start and ensure that the titration proceed. Initial, date, and file readout.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 506**

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea and Water Insoluble Nitrogen Present [FM 435]*

Version: 2.0

Effective Date: 10/02/2020

Page 6 of 7

**10.5. Data Evaluation**

10.5.1. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs, to evaluate the data against the established criteria. Evaluations, that supersede those of DATA 051, are addressed in this document.

10.5.1.1. Evaluate the Lab Control sample results.

10.5.2. To determine if it is within acceptable limits refer to DATA 051 Attachment B, Fertilizer Laboratory Control Sample Acceptable Limits.

**10.6. Result Entry**

10.6.1. Record the results (from the print out) into the FERT\_NITROGEN\_KJL template in the LIMS. A laboratory notebook may be used as an intermediate location for the results.

10.6.2. Any additional acid used during the distillation must be added to the acid trap value in the LIMS. This updated value will be used in calculations.

10.6.3. The raw data, from the auto titrator, is uploaded to the folder on the “\\tlhaeschrome” drive.

**11. WASTE MANAGEMENT**

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material                              | Disposal Mechanism                                                             |
|---------------------------------------|--------------------------------------------------------------------------------|
| Boiling stones                        | Trash                                                                          |
| Residue (from the filtration process) | Trash                                                                          |
| Sample extract, excess                | Pour waste into KJL waste container. Flush waste down drain with running water |

**REFERENCES**

11.2. FSFL-SOP 494 Calibration and use of QC-Titrate™ Titration system (2006/01/17)

11.3. FM-430 Breakdown Nitrogen Preparation (2001/11/12)

11.4. FM-435 Ammoniacal nitrogen (urea Present) Kjeldahl 02/04/2003

**12. VERSION HISTORY**

| Version | Date       | Description                                                                                                                                                                                                    | Author        | Editor          |
|---------|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------|
| 1.0     | 11.16.2017 | Original Version                                                                                                                                                                                               | Zachary Tower |                 |
| 2.0     | 09.30.2020 | Added section 8.0 Reagent Preparation Guide<br>Updated Scope<br>Updated Apparatus and Materials<br>Updated Reagents and Media<br>Minor clarifications throughout Updated disposal mechanism for sample extract |               | Katherine Malik |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 506

*Sample Preparation and Analysis of Direct Ammoniacal Nitrogen for  
Sample with Urea and Water Insoluble Nitrogen Present [FM 435]*

Version: 2.0

Effective Date: 10/02/2020

Page 7 of 7

**13. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Bureau Chief, Technical Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer  
Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Wet Chemistry and Preparation Manager  
Date\_\_\_\_\_

\_\_\_\_\_  
Katherine Malik  
Editor  
Date\_\_\_\_\_

# STANDARD OPERATING PROCEDURE BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

## METHOD 509

### *Analysis of Total Metals and Heavy Metals using ICP-OES* *[FM822, FM825, FM860, FM875, FM902]*

Version: 4.0

Effective Date: 10/16/2023

Page 1 of 9

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for determination of acid soluble metals including total micronutrients and heavy metals (Aluminum, Arsenic, Cadmium, Calcium, Chromium, Cobalt, Copper, Iron, Lead, Magnesium, Manganese, Molybdenum, Nickel, Selenium, and Zinc) in mixed or pure material fertilizer samples in both liquid and dry matrices by Inductively Coupled Plasma – Optical Emission Spectrometry (ICP-OES).

## 2. SCOPE

2.1. The procedure delineated below is applicable to ICP-OES determination of acid soluble total micronutrients and heavy metals (Aluminum, Arsenic, Cadmium, Calcium, Chromium, Cobalt, Copper, Iron, Lead, Magnesium, Manganese, Molybdenum, Nickel, Selenium, and Zinc) in fertilizer. It is achieved by measuring the intensity of the electromagnetic radiation produced by the analyte in the inductively coupled plasma. A quantitative determination of the amount of analyte present can be made at the specific wavelength by measuring the excited atoms and ions emitted. The intensity of this emission is indicative of the concentration of the analyte within the sample. Acid match standards solutions are used to calibrate the ICP-OES. Cesium is used as an ionization buffer to minimize easily ionized element effects.

2.2. This method is restricted to use by or under the supervision of a trained analyst. Each analyst must demonstrate the ability to generate acceptable results with this method.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

3.2. Torch and interface regions of the ICP can be very hot and caution is needed when working around these areas if the ICP has recently been running. The ICP generates very high voltages. Specific safeties are in place on the instrument to protect the operator from accidental exposure to electrical shock; these safeties should not be subverted. Extreme caution should be used when performing any maintenance that would expose the operator to high voltages. Protective covers shall remain in place at all times except when special maintenance and/or repairs require their removal.

3.3. Exhaust from combustion fumes and vapors from the ICP torch housing can produce toxic vapors and may pose a health hazard if not adequately vented to the outside.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Standards
- 4.2. 10.2 Instrument and Computer Startup
- 4.3. 10.3 Recommended Instrument Parameters
- 4.4. 10.4 Sample Handling
- 4.5. 10.5 Sample Analysis
- 4.6. 10.6 Data Evaluation
- 4.7. 10.7 Quality Control
- 4.8. 10.8 Export Data
- 4.9. 10.9 Instrument Shutdown

## 5. INTERFERENCE

5.1. Spectral interferences are caused by background emissions from continuous or recombination phenomena, stray light from the line emission of high concentration elements, overlap of a spectral line from another element, ionization interferences, or unresolved overlap of molecular band spectra. Interferences may be overcome by dilution of the digestate, application of interelement correction

# STANDARD OPERATING PROCEDURE BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

## METHOD 509

### *Analysis of Total Metals and Heavy Metals using ICP-OES [FM822, FM825, FM860, FM875, FM902]*

Version: 4.0

Effective Date: 10/16/2023

Page 2 of 9

factors, use of alternate wavelengths, or addition of cesium chloride as an ionization buffer to minimize easily ionized element effects. A target of at least 2200mg/kg cesium chloride solution must be presented at the nebulizer.

5.2. Yttrium is found in low levels in some phosphate ore sources and should not be used as an internal standard.

5.3. Samples containing Fe will interfere with the following wavelengths: As 193.696, Cd 214.439, Cd 226.502, Pb 220.353, Pb 217.000 and Se 196.026. Cr can interfere with As 188.980. Spectral deconvolution must be utilized if Fe and Cr are present in samples at levels of interference for the forementioned elements. Multicomponent Spectral Fitting (MSF) is utilized in this method according to the procedure below:

**5.3.1. MSF Modeling:** Spectra must be collected from a matrix blank (B), a standard of the analyte of interest (A), and a standard of the interfering analyte(s) (I) using the same analysis method that will be used to normally analyze samples.

**Note:** The Peak Algorithm must be set Peak Area or Peak Height within the analysis method. The MSF algorithm is not used during this time.

5.3.2. Spectral profile is independent of concentration, however, to optimize the MSF model, the concentration of the analyte should be approximately 100 times its detection limit at the affected wavelength and the concentration of the interfering analyte should approximate maximum concentration that could be present.

E.g., ~100 times the detection limit of Pb would be approximately 10 ppm. In this example, 10 ppm Pb standard was selected. A 1000 ppm Fe standard was selected to represent the maximum concentration of the interfering Fe (10% Fe in the sample).

**Note:** All Analytes must be in the correct matrix (same acid ratio as the extracts).

5.3.3. After the analysis of the blank, analyte, and interferent, build the MSF model using the Examine MSF window. Label the blank (B), analyte (A) and interferent (I) on the wavelength of interest (e.g., 213.617 nm). Once the MSF file is saved, Select MSF as the Peak Algorithm within the analysis method and specify which MSF file to be used on the Spectral Corrections tab.

5.3.4. Confirm the new MSF model is function as intended by reanalyzing the interferent and a known sample or standard such as a LS. The resulting interferent should not be detected at the wavelength of interest.

5.3.5. If a wavelength calibration is performed, the MSF must be recreated.

**Note:** Wavelength calibration is performed during the annual PM.

## 6. APPARATUS and MATERIAL

6.1. Analytical Instruments:

6.1.1. Dual view ICP-OES, Perkin Elmer Optima 5300 DV and 8300 DV or equivalent.

6.1.2. Autosamplers; SC4 DX and SC2 DX or equivalent

6.1.3. Peristaltic Pump capable of handling three introduction lines

6.1.4. Chiller: PolyScience or equivalent

6.2. Computer system capable of operating integration software which interfaces with the ICP and acquires and stores all the required data.

6.3. Pump and tubings: black/black, grey/grey, orange/green, red/red

# STANDARD OPERATING PROCEDURE BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

## METHOD 509

### *Analysis of Total Metals and Heavy Metals using ICP-OES [FM822, FM825, FM860, FM875, FM902]*

Version: 4.0

Effective Date: 10/16/2023

Page 3 of 9

- 6.4. Volumetric flasks
- 6.5. Volumetric pipettes

## 7. REAGENTS and MEDIA

- 7.1. Cesium Chloride (CsCl), Trace metal grade
- 7.2. Hydrochloric Acid (HCl) concentrated, Trace metal grade
- 7.3. Nitric Acid (HNO<sub>3</sub>) concentrated, Trace metal grade
- 7.4. Triton X-100
- 7.5. Water, Deionized (DI) Type I

## 8. REAGENT PREPARATION GUIDELINES

- 8.1. Carrier Solution
  - 8.1.1. Directions: Add 200mL deionized water in 1L volumetric flask, add 90mL of HNO<sub>3</sub> and 30mL HCl and swirl to mix. Fill the flask with deionized water to the mark, invert to mix.  
Optional: 1-2 drops of Triton X-100 may be added. (Triton X -100 may cause beading in some spray chambers and should be omitted in these cases.)
- 8.2. Wash Solution
  - 8.2.1. Directions: Add approximately 5L DI water to a 15L carboy, add 307mL HCl, swirl to mix. Add DI water to the 15L mark, cap the carboy, and shake to mix.  
Optional: 7.5mL of Triton X-100 may be added. (Triton X -100 may cause beading in some spray chambers and should be omitted in these cases.)
- 8.3. Internal Standard / Ionization Buffer
  - 8.3.1. Directions: Add approximately 1500mL DI water to a 2000mL volumetric flask. Add 24g CsCl and swirl to mix. Add 10 mL of 1000 ppm Scandium standard solution and 40mL of 1000 ppm Beryllium standard solution (the final concentration of the Internal Standards are 5.0 ppm Sc and 20 ppm Be). Add 80mL of HNO<sub>3</sub> and swirl. Cap and invert several times to mix.  
Optional: 1-2 drops of Triton X-100 may be added. (Triton X -100 may cause beading in some spray chambers and should be omitted in these cases.)
  - 8.3.2. Expiration is 6 months from preparation date.
- 8.4. 9% Nitric Acid: 3% Hydrochloric Acid
  - 8.4.1. Directions: Add 200mL deionized water in 1L volumetric flask, add 90mL of HNO<sub>3</sub> and 30mL HCl and swirl to mix. Fill the flask with deionized water to the mark, invert to mix.
- 8.5. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015.
- 8.6. The specific ingredient amounts can vary to accommodate the batch size.
- 8.7. Expiration Date: 12 months from date of preparation unless otherwise specified.
- 8.8. Store at room temperature.

## 9. ASSOCIATED DOCUMENTS

- 9.1. ADMIN 023, Using the Laboratory Information Management System
- 9.2. ADMIN 023 Attachment M, Instrument Sequence Generation for Fertilizer
- 9.3. ADMIN 023 Attachment O, Automatic Result Entry for Fertilizer
- 9.4. DATA 050, Data Control
- 9.5. DATA 051, Data Evaluation for Chemical Analysis
- 9.6. DATA 051 Attachment C, Fertilizer Continuing Calibration Verification Acceptance Limits
- 9.7. DATA 051 Attachment A, Fertilizer Analysis Requirements
- 9.8. LABOP 140 Maintenance, Operation, and Program of ICPOES, Attachment I
- 9.9. QUALITY 112, Quality Control

# STANDARD OPERATING PROCEDURE BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

## METHOD 509

### *Analysis of Total Metals and Heavy Metals using ICP-OES [FM822, FM825, FM860, FM875, FM902]*

Version: 4.0

Effective Date: 10/16/2023

Page 4 of 9

## 10. SPECIFIC PROCEDURES

### 10.1. Standards

10.1.1. Primary standards and internal standards (1000ppm, 5000ppm or 10,000ppm) shall be prepared from neat material or purchased through certified vendors that are NIST traceable.

10.1.2. Working standards are to be used for ICP calibration and unknown sample identification and quantification. They are made by diluting with 9% HNO<sub>3</sub>: 3% HCl to obtain the suggested concentrations.

#### 10.1.3. Recommended Standard Levels:

10.1.3.1. Calibration Standards: The tables below contain the recommended calibration curve levels. Other levels may be used providing they are within the instrument's linear range.

**NOTE:** Continuing Calibration Verification Standards (CCV's) are evaluated for each curve separately.

##### 10.1.3.1.1. Recommended Standards Levels:

| Standard Level (ppm) Low Curve |      |     |     |    |    |
|--------------------------------|------|-----|-----|----|----|
| Analyte                        | 1    | 2   | 3   | 4  | 5  |
| As                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Ca                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Cd                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Co                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Cr                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Cu                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Fe                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Mg                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Mn                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Mo                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Ni                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Pb                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Se                             | 0.10 | 0.5 | 5.0 | 10 | 20 |
| Zn                             | 0.10 | 0.5 | 5.0 | 10 | 20 |

| Standard Level (ppm) High Curve |    |     |     |     |     |
|---------------------------------|----|-----|-----|-----|-----|
| Analyte                         | 6  | 7   | 8   | 9   | 10  |
| Ca                              | 50 | 100 | 300 | 500 | 600 |
| Fe                              | 50 | 100 | 300 | 500 | 600 |
| Mg                              | 50 | 100 | 300 | 500 | 600 |
| Mn                              | 50 | 100 | 300 | 500 | 600 |
| Zn                              | 50 | 100 | 300 | 500 | 600 |

### 10.2. Instrument and Computer Startup

10.2.1. See LABOP 140 Maintenance, Operation, and Program of ICPOES, Attachment I for details.

10.2.2. Turn on the plasma

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 509

#### *Analysis of Total Metals and Heavy Metals using ICP-OES* [FM822, FM825, FM860, FM875, FM902]

Version: 4.0

Effective Date: 10/16/2023

Page 5 of 9

10.2.2.1. Allow the system to stabilize for a minimum of 30 minutes.

#### 10.3. Recommended Instrument Parameters

10.3.1. The following parameters may vary for different ICP makes and models and therefore are listed as recommended and may require minor modifications in order to optimize performance.

##### 10.3.2. For Perkin Elmer Optima 5300 DV

|                       |              |
|-----------------------|--------------|
| Plasma Flow           | 15 L/min     |
| Auxiliary Flow        | 0.2 L/min    |
| Nebulizer Flow        | 0.65 L/min   |
| RF Power              | 1450 watts   |
| Autosampler Flow Rate | 0.8 mL/min   |
| Carrier Tubing        | Black/Black  |
| Internal Std. Tubing  | Orange/Green |
| Waste Tubing          | Red/Red      |
| Wash                  | 0.8 mL/min   |
| Spray Chamber         | Cyclonic     |
| Loop Size             | 1.0 mL       |
| Replicate Reads       | 3            |

##### 10.3.3. For Perkin Elmer Optima 8300 DV

|                       |              |
|-----------------------|--------------|
| Plasma Flow           | 15 L/min     |
| Auxiliary Flow        | 0.2 L/min    |
| Nebulizer Flow        | 0.65 L/min   |
| RF Power              | 1450 watts   |
| Autosampler Flow Rate | 0.44 mL/min  |
| Carrier Tubing        | Black/Black  |
| Internal Std. Tubing  | Orange/Green |
| Waste Tubing          | Grey/Grey    |
| Wash                  | 0.44 mL/min  |
| Spray Chamber         | Cyclonic     |
| Loop Size             | 1.5 mL       |
| Replicate Reads       | 3            |

##### 10.3.4 Recommended Wavelengths and Plasma View:

| Element | Wavelength nm | Plasma View  |
|---------|---------------|--------------|
| As      | 188.979*      | Axial        |
| As      | 193.696       | Axial        |
| Be      | 313.107       | Axial/Radial |
| Be      | 234.861       | Axial/Radial |
| Ca      | 317.933*      | Radial       |
| Ca      | 315.887       | Radial       |
| Cd      | 228.802*      | Axial        |
| Cd      | 214.44        | Axial        |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 509

#### *Analysis of Total Metals and Heavy Metals using ICP-OES* [FM822, FM825, FM860, FM875, FM902]

Version: 4.0

Effective Date: 10/16/2023

Page 6 of 9

|    |          |              |
|----|----------|--------------|
| Co | 228.616* | Axial        |
| Co | 238.892  | Axial        |
| Cr | 267.716* | Axial        |
| Cr | 205.56   | Axial        |
| Cu | 327.393* | Radial       |
| Cu | 324.752  | Radial       |
| Fe | 238.204* | Radial       |
| Fe | 239.562  | Radial       |
| Mg | 285.213* | Radial       |
| Mg | 279.077  | Radial       |
| Mn | 257.610* | Radial       |
| Mn | 259.372  | Radial       |
| Mo | 202.031* | Axial        |
| Mo | 203.845  | Axial        |
| Ni | 231.604* | Axial        |
| Ni | 221.648  | Axial        |
| Pb | 220.353* | Axial        |
| Pb | 217      | Axial        |
| Sc | 361.383  | Axial/Radial |
| Se | 196.026* | Axial        |
| Se | 203.985  | Axial        |
| Zn | 206.200* | Radial       |
| Zn | 213.857  | Radial       |
| Y  | 371.029* | Axial        |

#### 10.4. Sample Handling

10.4.1. Sample extracts are filtered during the preparation procedure. Extracts that contain excessive particulates or cloudy extracts, should be re-filtered or centrifuged before analysis.

#### 10.5. Sample Analysis

10.5.1. Use the Sample Information File (SIF) to create the sequence.

10.5.2. An instrument blank must be injected at the beginning of a sample sequence to ensure that the system is free of contaminants.

10.5.3. After the initial calibration, there must be a continuing calibration check standard analyzed after every 10 samples and at the end of the run.

10.5.4. The Laboratory Control Samples (LCS) should be analyzed with their associated samples, not grouped together at the beginning of the sequence.

10.5.5. The sequence requirements are found in, DATA 051 Attachment A, Fertilizer Analysis Requirements.

10.5.6. The sequence must end with a CCV.

10.5.7. Samples require the proper laboratory information management system (LIMS) codes to allow required uploading. Refer to ADMIN 023 Attachment M, Instrument Sequence Generation for Fertilizer, for additional information.

10.5.7.1. Example: AA99999-0AT

# STANDARD OPERATING PROCEDURE BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

## METHOD 509

### *Analysis of Total Metals and Heavy Metals using ICP-OES [FM822, FM825, FM860, FM875, FM902]*

Version: 4.0

Effective Date: 10/16/2023

Page 7 of 9

10.5.7.1.1. 0 corresponds to the sample weight (0, 1, 2.)

10.5.7.1.2. A corresponds to the different injections due to dilutions (A=1st injection, B=2nd injection, C=3rd injection.)

10.5.7.1.3. T or H corresponds to analysis type (T=Totals and H=Heavy Metals)

10.5.8. Dilution factors must be entered into the instrument sequence.

#### 10.6. Data Evaluation

10.6.1. Refer to DATA 051, Data Evaluation for Chemical Analysis, to evaluate the data against the established criteria.

#### 10.7. Quality Control

##### 10.7.1. Initial Calibration:

10.7.1.1. The calibration is based on plotting the ratio of the analyte response to the internal standard response as a function of the analyte concentration of the standards.

10.7.1.2. The correlation coefficient of the calibration curve R value must be  $\geq 0.999$  or the  $R^2$  value must be  $\geq 0.998$ .

10.7.1.3. Refer to DATA 051, Data Evaluation for Chemical Analysis; for additional calibration procedures.

##### 10.7.2. Continuing Calibration Check Standard:

10.7.2.1. Once the initial calibration has passed, a continuing calibration check standard must be analyzed after every 10 samples and at the end of the sequence. CCVs must be run for each separate curve.

10.7.2.2. CCV acceptance criteria are listed below. Addition information may be found in DATA 051 Attachment C, Fertilizer Continuing Calibration Verification Acceptance Limits.

| Analyte | CCV Limit (%) |
|---------|---------------|
| As      | $\pm 10$      |
| Ca      | $\pm 4$       |
| Cd      | $\pm 10$      |
| Co      | $\pm 10$      |
| Cr      | $\pm 10$      |
| Cu      | $\pm 8$       |
| Fe      | $\pm 8$       |
| Mg      | $\pm 4$       |
| Mn      | $\pm 8$       |
| Mo      | $\pm 4$       |
| Ni      | $\pm 10$      |
| Pb      | $\pm 10$      |
| Se      | $\pm 10$      |
| Zn      | $\pm 8$       |

# STANDARD OPERATING PROCEDURE BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

## METHOD 509

### *Analysis of Total Metals and Heavy Metals using ICP-OES [FM822, FM825, FM860, FM875, FM902]*

Version: 4.0

Effective Date: 10/16/2023

Page 8 of 9

10.7.2.3. The CCV may consist of any concentration level used in the initial calibration. Levels near the middle of the curve are recommended.

10.7.2.4. If a CCV response is not within the limits, all samples analyzed before and after it must be reanalyzed.

10.7.2.5. If upon reanalysis of the CCV, the response is still not within the above acceptance criteria, then a new initial calibration must be performed. It is recommended that routine instrument maintenance be performed prior to the new curve analysis.

10.7.3. If the ppm reading exceeds the calibration range, then the sample must be diluted and reanalyzed prior to quantification. The sample should be diluted with 9% HNO<sub>3</sub> and 3% HCl to maintain the appropriate acid proportion.

10.7.4. Internal Standard criteria for samples is 80% - 120% recovery.

10.7.5. %RSD for the sample replicates shall not exceed 10%.

10.7.6. Refer to QUALITY 112, Quality Control, to determine the outcome of Quality Control (QC) samples and for a description of quality control (QC) samples required in the batch.

#### 10.8. Export Data

10.8.1. Data Manager exports data to the P:Drive. Refer to DATA 050, Data Control.

10.8.2. Upload data from the P:Drive to the LIMS. Refer to ADMIN 023 Attachment O, Automatic Result Entry for Fertilizer.

#### 10.9. Instrument Shutdown

10.9.1. See LABOP 140 Attachment I: Maintenance, Operation, and Program of ICPOES, for details.

## 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material                                                                                                       | Disposal Mechanism                                                                                                                                                                                                                                                                                                 |
|----------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sample digests                                                                                                 | Sample digests are good to analyze for up to 6 months and may be disposed of after analysis unless there is limited sample volume; then they may need to be kept for up to 90 days after laboratory report has been issued.<br>Non-hazardous sample digests are disposed of down a sink while flushing with water. |
| Instrument Waste                                                                                               | Disposed of down a sink while flushing with water.                                                                                                                                                                                                                                                                 |
| Non-hazardous, calibration standards, tune standards, stock standards and internal standards.                  | Disposed of down a sink while flushing with water.                                                                                                                                                                                                                                                                 |
| Samples or standards containing antimony, arsenic, cadmium, chromium, molybdenum, lead, thallium above 100 ppm | Labeled as hazardous waste and are disposed of into hazardous waste container labeled for metals.                                                                                                                                                                                                                  |

# STANDARD OPERATING PROCEDURE BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

## METHOD 509

### *Analysis of Total Metals and Heavy Metals using ICP-OES [FM822, FM825, FM860, FM875, FM902]*

Version: 4.0

Effective Date: 10/16/2023

Page 9 of 9

## 12. REFERENCES

12.1. Preparation and Analysis of Acid Soluble (Total) Micro Nutrients by ICP-OES, FSFL-SOP 504 (04/15/2009).

12.2. U.S. EPA. 1994. "Method 200.7: Determination of Metals and Trace Elements in Water and Wastes by Inductively Coupled Plasma-Atomic Emission Spectrometry," Revision 4.4. Cincinnati, OH

**External 00131**

12.3. Analysis of Water-Soluble Mg and Potash, BAEL\_SOP 201 FSFL-SOP 518 (04/15/2009).

12.4. Operation and Maintenance of Perkin Elmer Inductively Coupled Plasma – Optical Emission Spectrometry 3300 Dual View System, FSFL-SOP 525 (16/07/2007).

12.5. *Official Methods of Analysis of AOAC INTERNATIONAL* (2019) 21st Ed., AOAC INTERNATIONAL, Rockville, MD, USA, Method **2017.02**, **External 00135**.

## 13. VERSION HISTORY

| Version | Date       | Description                                                                                                                                                                                                                                                                                                                                                                                                        | Author        | Editor            |
|---------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------------------|
| 1.0     | 04.04.2017 | Original Version                                                                                                                                                                                                                                                                                                                                                                                                   | Naser Tassobi |                   |
| 2.0     | 05.17.2019 | Updated Outline Format, Updated Title, Updated Outline, Added Section 8, Added required documents to Associated Documents, Added section 10.6 Updated Curve Tables, Verified and Updated wavelength table, Moved maintenance, operation, and procedure of the instrument to LABOP 140 attachment I Added Section 11, Waste Managements                                                                             |               | Christian Amason  |
| 3.0     | 09.30.2020 | Updated Purpose, Scope and Health and Safety<br>Updated Interference Section<br>Updated Curve Tables<br>Updated Reagents and Media<br>Updated Reagent Preparation Guidelines<br>Updated title for DATA 051 through out<br>Updated Instrument and Computer Startup<br>Updated Recommended Instrument Parameters<br>Updated Recommended Wavelengths and Plasma View<br>Updated Quality Control<br>Updated References |               | William Meeks Jr. |
| 4.0     | 10.16.2023 | Added external document number to external references<br>Added section 5.3 to document interferences.                                                                                                                                                                                                                                                                                                              |               | William Meeks Jr. |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 510

#### *Sample Preparation and Analysis of Chloride in Fertilizer* *[FM 611]*

Version: 3.0

Effective Date:02/08/2024


Page 1 of 5

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the determination of chloride in fertilizer using an Automatic Titrator.

## 2. SCOPE

2.1. The procedure delineated below is applicable to determining chlorine in fertilizer samples. Chloride is extracted from the fertilizer by boiling in water. The sample is then titrated to an end point with silver nitrate. The end point of the titration is reached when all of the chloride ions are precipitated.



## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Sample Preparation
- 4.2. 10.2 Standardization of Silver Nitrate
- 4.3. 10.3 Analysis
- 4.4. 10.4 Calculations

## 5. INTERFERENCE

5.1. AgCl precipitate will deposit on the electrode. The electrode's response will slow down and may become erratic. Care must be taken to clean the electrode by wiping the electrode with a wet paper tissue and/or by using a 0.1% NH<sub>3</sub> wash. Silver, in the presence of an ammonia solution will become Ag[(NH<sub>3</sub>)<sub>2</sub>]<sup>+</sup>, silver(I)ammine ion, assisting with the AgCl cleaning.

## 6. APPARATUS and MATERIAL

- 6.1. Mettler Toledo Excellence Titrator T5 or equivalent
- 6.2. Mettler Toledo DMi141-SC silver ring electrode or equivalent
- 6.3. 100mL plastic beakers
- 6.4. 1L glass (amber) container
- 6.5. 50mL digestion tube
- 6.6. Eppendorf pipette or equivalent

## 7. REAGENTS and MEDIA

- 7.1. Type I-DI water
- 7.2. 0.1% Ammonia wash solution
- 7.3. Electrolyte Solution KNO<sub>3</sub> 1 mol/L (potassium nitrate)
- 7.4. Silver Nitrate Solution N/10 (protect from light)
- 7.5. Standard 0.1M NaCl, or KCl

## 8. REAGENT PREPARATION GUIDELINES

- 8.1. 0.1% Ammonia wash solution
  - 8.1.1. ~1mL Ammonia Hydroxide to 1L DI water.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 510**  
*Sample Preparation and Analysis of Chloride in Fertilizer*  
*[FM 611]*

Version: 3.0

Effective Date: 02/08/2024

Page 2 of 5

8.2. Purchased Electrolyte Solution  $\text{KNO}_3$  1 mol/L (potassium nitrate)

8.3. Purchased Silver Nitrate Solution N/10 (protect from light)

8.4. Standard 0.1M NaCl, or 0.1M KCl Standard Solution (Purchased or made from salt)

**8.4.1. 0.1M KCl Standard Solution**

8.4.1.1. Dry KCl at  $106^\circ\text{C}$  for approximately 4 hours, weigh 7.455g of dried KCl into 1L volumetric flask, bring to volume with DI water.

**8.4.2. 0.1M NaCl Standard Solution**

8.4.2.1. Dry NaCl at  $106^\circ\text{C}$  for approximately 4 hours, weigh 5.844g of dried NaCl into 1L volumetric flask, bring to volume with DI water.

8.5. Label the containers appropriately, label L003 indicates required information. Add relevant information to Reagent Preparation, form F00015

8.6. The specific ingredient amounts can vary to accommodate the batch size.

8.7. Expiration Date: 12 months from date of preparation, unless specifically stated in the reagent directions or from manufacturer.

8.8. Store at room temperature, unless specifically stated in the reagent directions

## 9. ASSOCIATED DOCUMENTS

9.1. LABOP 121, Maintenance, Operation and Program of Balance

9.2. Label L003, Reagent Label

9.3. Form F00015, Reagent Preparation Form

## 10. SPECIFIC PROCEDURES

### 10.1. Sample Preparation

#### 10.1.1. Weighing sample

10.1.1.1. Label 50mL tubes with the sample number

10.1.1.2. Weigh 0.20g of sample in labeled 50mL digestion tube.

10.1.1.3. Add approximately 30mL of DI water (Type I DI water is required)

10.1.1.4. Digest 60 minutes on hot block at  $\sim 106^\circ\text{C}$

10.1.1.5. Let cool to room temperature

10.1.1.6. Transfer the sample from the 50mL digestion tube to a 100mL beaker. Ensure that the sample has been completely transferred, by rinsing the digestion tube 3 times with deionized water to the beaker.

**NOTE:** Use of small rinse volumes is recommended, so that the total volume does not exceed 60mL.

### 10.2. Standardization of Silver Nitrate

10.2.1. Pipette three 2mL aliquots of the 0.1M NaCl (or 0.1M KCl) standard solution into separate 100mL beakers.

10.2.2. Bring each to approximately 40mL with DI water.

10.2.3. Prime the  $\text{AgNO}_3$  burette.

10.2.4. Open the Mettler LabX Software.

10.2.5. Titrate with 0.1M  $\text{AgNO}_3$ , using the instrument method titled "Standardization 0.1M  $\text{AgNO}_3$  with Cl Standard"

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 510

#### Sample Preparation and Analysis of Chloride in Fertilizer

#### [FM 611]

Version: 3.0

Effective Date:02/08/2024

Page 3 of 5

10.2.6. Enter the volume of the NaCl (or KCl) standard solution into the corresponding sample position.

10.2.7. At the end of the run, the 0.1M AgNO<sub>3</sub> titrant is automatically standardized using the titer factor found above.

#### 10.3. Analysis

10.3.1. Prime the AgNO<sub>3</sub> burette.

10.3.2. Open the Mettler LabX Software.

10.3.3. Titrate samples with 0.1M AgNO<sub>3</sub> using the instrument method titled "BAEL Method 510".

10.3.4. Enter the sample ID and weight into the corresponding sample position.

10.3.5. Enter % found into the LIMS.

#### 10.4. Calculations

##### 10.4.1. Standardization of AgNO<sub>3</sub>

##### 10.4.1.1. Determination of Titer with Liquid Cl

$$\text{Titer} = \frac{C_{(\text{Cl})} \times V_{(\text{Cl})}}{C_{(\text{AgNO}_3)} \times V_{\text{EQ}}}$$

| Symbol                          | Definition                                         | Mettler Symbol |
|---------------------------------|----------------------------------------------------|----------------|
| C <sub>(Cl)</sub>               | Concentration of Cl standard (mol/L)               | cst            |
| V <sub>(Cl)</sub>               | Volume of Cl standard (mL)                         | m              |
| C <sub>(AgNO<sub>3</sub>)</sub> | Concentration of AgNO <sub>3</sub> titrant (mol/L) | c              |
| V <sub>EQ</sub>                 | Volume of titrant consumed (mL)                    | VEQ            |

##### 10.4.1.2. Determination of Titer with Solid KCl\*

$$\text{Titer} = \frac{m_{(\text{KCl})}}{C_{(\text{AgNO}_3)} \times V_{\text{EQ}} \times \text{MW}_{(\text{KCl})}}$$

| Symbol                          | Definition                                         | Mettler Symbol |
|---------------------------------|----------------------------------------------------|----------------|
| m <sub>(KCl)</sub>              | mass of KCl standard (mg)                          | m              |
| C <sub>(AgNO<sub>3</sub>)</sub> | Concentration of AgNO <sub>3</sub> titrant (mol/L) | c              |
| V <sub>EQ</sub>                 | Volume of titrant consumed (mL)                    | VEQ            |
| MW <sub>(KCl)</sub>             | molecular weight of KCl (g/mol)                    | M              |

\*If using other solid standards, substitute the appropriate standard in place of KCl. e.g. using a solid NaCl standard, replace KCl concentration and molecular weight for that of NaCl.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 510

#### *Sample Preparation and Analysis of Chloride in Fertilizer*

#### *[FM 611]*

Version: 3.0

Effective Date:02/08/2024

Page 4 of 5

#### 10.4.2. % Cl Sample Calculations

$$\% \text{ Cl} = \frac{C_{(\text{AgNO}_3)} \times V_{\text{EQ}} \times \text{Titer} \times \text{MW}_{(\text{Cl})} \times 100}{m_{(\text{sample})} \times 1000}$$

| Symbol                    | Definition                                                                                                                                                                      | Mettler Symbol |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------|
| $C_{(\text{AgNO}_3)}$     | Concentration of $\text{AgNO}_3$ titrant (mol/L)                                                                                                                                | c              |
| $V_{\text{EQ}}$           | Volume of titrant consumed (mL)                                                                                                                                                 | VEQ            |
| $m_{(\text{sample})}$     | mass of sample (g)*                                                                                                                                                             | m              |
|                           | *When aliquot of a sample extract is used<br>$m_{(\text{sample})} = \frac{\text{weight}_{(\text{g})} \times \text{aliquot}_{(\text{mL})}}{\text{extraction vol}_{(\text{mL})}}$ |                |
| Titer                     | Titer correction factor from standardization                                                                                                                                    | TITER          |
| $\text{MW}_{(\text{Cl})}$ | molecular weight of Cl (35.4527g/mol)                                                                                                                                           | M              |

#### 10.4.3. Mettler Titrator Formula

$$\% \text{ Cl} = \frac{Q \times C}{m}$$

| Mettler Symbol | Definition                                                             |
|----------------|------------------------------------------------------------------------|
| Q              | mmol of titrant consumed ( $\text{VEQ} \times c \times \text{TITER}$ ) |
| C              | constant (M/10)                                                        |

## 11. WASTE MANAGEMENT

Sample waste from the titration will contain silver chloride precipitate along with some free silver ions.  $\text{NaCl}$  can be added to this waste at a concentration of 0.1% (1g  $\text{NaCl}$  / 1 L waste, then invert bottle to ensure reaction), to precipitate any remaining silver in solution. The  $\text{AgCl}$  precipitate can then be filtered and submitted to hazardous waste. The aqueous portion may be disposed of as non-hazardous waste.

## 12. REFERENCES

- 12.1. E00154, *Chlorine (water soluble) in Fertilizer*. Method 928.02A, AOAC 15<sup>th</sup> Edition.
- 12.2. E00165, Potash Analysis, Chlorine. FM-611, December 9, 2009
- 12.3. E00166, Operation and Maintenance of Mettler Titrator. FSFL-SOP 528, Rev 1.0

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 510  
*Sample Preparation and Analysis of Chloride in Fertilizer*  
*[FM 611]*

Version: 3.0

Effective Date:02/08/2024

Page 5 of 5

**13. VERSION HISTORY**

| Version | Date       | Description                                                                                                                                 | Author        | Editor                                            |
|---------|------------|---------------------------------------------------------------------------------------------------------------------------------------------|---------------|---------------------------------------------------|
| 1.0     | 10.17.2018 | Original Version                                                                                                                            | William Meeks |                                                   |
| 2.0     | 12.09.2020 | Corrected typo in section 1<br>Corrected numbering in section 4<br>Replaced Elizabeth Schneider with<br>Quentin Cooper in Signature History |               | Taleya<br>Williams                                |
| 3.0     | 02/08/2024 | Added option to purchase KCl<br>standard solution<br>Added instructions for creating NaCl<br>standard solution                              |               | William Meeks<br>Jr.<br>Patrizia<br>Lemma, Ph. D. |

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 513**

*Analysis of Water Soluble Potash, Water Soluble Magnesium, and Water Soluble Boron using Inductively Coupled Plasma Optical Emission Spectrometry*  
*[FM604, FM803]*

Version: 4.0

Effective Date: 10/01/2021

Page 1 of 8

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for the determination of water soluble potash ( $K_2O$ ), water soluble Magnesium (Mg), and water soluble Boron (B) in mixed or pure material fertilizer samples in both liquid and dry matrices by Inductively Coupled Plasma – Optical Emission Spectrometry (ICP-OES).

## **2. SCOPE**

2.1. The procedure delineated below is applicable to an ICP-OES determination of water soluble potash, water soluble magnesium, and water soluble boron in fertilizer and is achieved by measuring the intensity of the electromagnetic radiation produced by the analyte in the inductively coupled plasma. A quantitative determination of the amount of analyte present can be made at the specific wavelength by measuring the excited atoms and ions emitted. The intensity of this emission is indicative of the concentration of the analyte within the sample. The determined amount of potassium (K) is converted to potash ( $K_2O$  equivalent) and reported as such.

2.2. Fertilizer samples analyzed for Mg and  $K_2O$  are prepared by digesting the sample according to PREP 454, Sample Preparation for Water Soluble Potash and Magnesium Analysis.

2.3. Fertilizer samples analyzed for B are prepared by digesting the sample according to PREP 455, Sample Preparation for Water Soluble Boron Analysis.

2.4. This method is restricted to use by or under the supervision of a trained analyst. Each analyst must demonstrate the ability to generate acceptable results with this method. Refer to ADMIN 016, Training and Competency, for demonstration of competence procedures.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## **4. OUTLINE of PROCEDURE**

- 4.1. 10.1 Sample Handling
- 4.2. 10.2 Standards
- 4.3. 10.3 Instrument Startup
- 4.4. 10.4 Sample Analysis
- 4.5. 10.5 Export Data
- 4.6. 10.6 Instrument Shutdown
- 4.7. 10.7 Data Reduction and Reporting
- 4.8. 10.8 Quality Control
- 4.9. 10.9 Instrument Blank

## **5. INTERFERENCE**

5.1. Spectral interferences are caused by background emission from continuous or recombination phenomena, stray light from the line emission of high concentration elements, overlap of a spectral line from another element, ionization interferences (for K), or unresolved overlap of molecular band spectra. Interferences may be overcome by dilution of the digestate, application of interelement correction factors, use of alternate wavelengths, or addition of cesium chloride ionization buffer (for K). A target of at least 2200mg/kg cesium chloride solution must be presented at the nebulizer.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 513**

*Analysis of Water Soluble Potash, Water Soluble Magnesium, and Water Soluble Boron using Inductively Coupled Plasma Optical Emission Spectrometry*  
*[FM604, FM803]*

Version: 4.0

Effective Date: 10/01/2021

Page 2 of 8

## **6. APPARATUS and MATERIAL**

### **6.1. Analytical Instruments**

- 6.1.1. Dual view ICP-OES
- 6.1.2. Autosampler
- 6.1.3. Peristaltic Pump capable of handling three introduction lines
- 6.1.4. Chiller
- 6.1.5. Integration software

### **6.2. ICP autosampler tubes**

- 6.3. Pump tubing: black/black, orange/green, red/red
- 6.4. Volumetric flasks
- 6.5. Volumetric pipettes
- 6.6. 0.45µm Syringe Filters

## **7. REAGENTS and MEDIA**

- 7.1. Cesium Chloride (CsCl), Trace metal grade, solid
- 7.2. Hydrochloric Acid (HCl) concentrated, Trace meal grade
- 7.3. Nitric Acid (HNO<sub>3</sub>) concentrated, Trace metal grade
- 7.4. Triton X-100
- 7.5. Scandium Solution, 1000 ppm
- 7.6. Water, Deionized (DI), Type I

## **8. REAGENT PREPARATION GUIDELINES**

### **8.1. Carrier Solution**

8.1.1. Directions: Fill a 1L volumetric flask almost to volume with DI water. Add 1-2 drops Triton X-100, then fill to the mark with DI water. Some spray chambers do not work with the Triton X. It may be left out in those cases. Invert flask several times to mix.

### **8.2. Wash Solution**

8.2.1. Directions: Fill the 15L carboy with approximately 14 liters of DI water. Add approximately 7.5mL Triton X-100 (optional) and 307mL HCl and fill to the mark with DI water. Some spray chambers do not work with the Triton X. It may be left out in those cases. Swirl to mix.

### **8.3. Internal Standard / Ionization Buffer**

8.3.1. Directions: Add approximately 1500mL DI water to a 2000mL volumetric flask. Add 24g CsCl and swirl to mix. Add 10mL of 1000ppm Scandium Internal Standard Solution (the final concentration of scandium Internal Standard is 5.0ppm). Add 80mL of HNO<sub>3</sub> and swirl to mix. Add 1-2 drops Triton X-100 (optional). Some spray chambers do not work with the Triton X. It may be left out in those cases. Fill to volume with DI water. Cap and invert to mix.

8.3.2. Expiration is 6 months from preparation date.

8.4. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015.

8.5. The specific ingredient amounts can vary to accommodate the batch size.

8.6. Expiration Date: 12 months from date of preparation unless otherwise indicated.

8.7. Store at room temperature.

## **9. ASSOCIATED DOCUMENTS**

- 9.1. ADMIN 023 Attachment M Instrument Sequence Generation for Fertilizer
- 9.2. ADMIN 023 Attachment O, Automatic Result Entry for Fertilizer

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 513**

*Analysis of Water Soluble Potash, Water Soluble Magnesium, and Water Soluble Boron using Inductively Coupled Plasma Optical Emission Spectrometry*  
*[FM604, FM803]*

Version: 4.0

Effective Date: 10/01/2021

Page 3 of 8

- 9.3. DATA 051, Data Evaluation for Feed, Fertilizer, and Pesticide Programs
- 9.4. DATA 051 Attachment A, Fertilizer Analysis Requirements
- 9.5. DATA 051 Attachment B, Fertilizer Laboratory Control Sample Acceptance Limits
- 9.6. DATA 051 Attachment C, Continuing Calibration Verification Acceptance Limits for AA and ICP
- 9.7. LABOP 140 Attachment I: Maintenance, Operation and Program of ICPOES
- 9.8. PREP 454, Sample Preparation for Water Soluble Potash and Magnesium Analysis
- 9.9. PREP 455, Sample Preparation for Water Soluble Boron Analysis
- 9.10. QUALITY 112, Quality Control
- 9.11. SAFETY 002, Waste Management

## **10. SPECIFIC PROCEDURES**

### **10.1. Sample Handling**

- 10.1.1. Sample extracts containing particulates should be filtered before analysis.
- 10.1.2. Aqueous sample digestates may be held for up to two weeks after digestion.

### **10.2. Standards**

- 10.2.1. Primary standards and internal standards (1000ppm or 10,000ppm) should be prepared from neat material or purchased through certified vendors that are NIST traceable.
- 10.2.2. Working standards are to be used for ICP calibration and unknown sample identification and quantification. They are made by diluting into DI water to obtain the following concentrations.
- 10.2.3. Recommended Standard Levels:

| <b>Standard Level</b> | <b>Concentration (ppm) K</b> | <b>Concentration (ppm) Mg</b> | <b>Concentration (ppm) B</b> |
|-----------------------|------------------------------|-------------------------------|------------------------------|
| <b>Low</b>            | <b>Curve</b>                 |                               |                              |
| 1                     | -                            | 0.5                           | 0.5                          |
| 2                     | 5                            | 2.5                           | 2.5                          |
| 3                     | 10                           | 5                             | 5                            |
| 4                     | 20                           | 10                            | 10                           |
| 5                     | 50                           | -                             | -                            |
| <b>High</b>           | <b>Curve</b>                 |                               |                              |
| 4                     | -                            | 10                            | 10                           |
| 5                     | 50                           | 25                            | 25                           |
| 6                     | 100                          | 50                            | 50                           |
| 7                     | 200                          | 100                           | 100                          |
| 8                     | 400                          | 200                           | 200                          |
| 9                     | 600                          | -                             | -                            |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 513

*Analysis of Water Soluble Potash, Water Soluble Magnesium, and Water Soluble Boron using Inductively Coupled Plasma Optical Emission Spectrometry  
[FM604, FM803]*

Version: 4.0

Effective Date: 10/01/2021

Page 4 of 8

#### 10.3. Instrument Startup

10.3.1. See LABOP 140 Attachment I: Maintenance, Operation and Program of ICPOES

10.3.2.

| Element             | Wavelength nm | Plasma View   |
|---------------------|---------------|---------------|
| K                   | 766.490*      | Radial        |
| K                   | 404.721       | Radial        |
| B                   | 249.677*      | Axial         |
| B                   | 249.772       | Axial         |
| Mg                  | 285.213*      | Axial/Radial  |
| Mg                  | 279.077       | Axial/ Radial |
| Sc                  | 361.383*      | Axial/Radial  |
| *Primary Wavelength |               |               |

10.3.3.

|                            | PE Optima 5300 | PE Optima 8300 |
|----------------------------|----------------|----------------|
| Plasma Flow                | 15 L/min       | 12 L/min       |
| Auxiliary Flow             | 0.2 L/min      | 0.2 L/min      |
| Nebulizer Flow             | 0.65 L/min     | 0.68 L/min     |
| RF Power                   | 1450 watts     | 1450 watts     |
| Peristaltic Pump Flow Rate | 0.8 mL/min     | 0.44 mL/min    |
| Carrier Tubing             | Black/Black    | Black/Black    |
| Internal Std. Tubing       | Orange/Green   | Orange/Green   |
| Waste Tubing               | Red/Red        | Red/Red        |

#### 10.4. Sample Analysis

10.4.1. Use the Sample Information File (SIF) to create the sequence.

10.4.2. An instrument blank must be injected at the beginning of a sample sequence to ensure that the system is free of contaminants.

10.4.3. After the initial calibration, there must be a continuing calibration check standard analyzed after every **ten** samples and at the end of the run.

10.4.4. The Laboratory Control Samples (LCS) must be analyzed with their associated samples, not grouped together at the beginning of the sequence.

10.4.5. The sequence requirements are found in, DATA 051 Attachment A, Fertilizer Analysis Requirements.

10.4.6. The sequence must end with a CCV. The ending CCV should be at a different concentration than the mid-sequence CCVs.

10.4.7. Samples require the proper laboratory information management system (LIMS) codes to allow required uploading. Refer to ADMIN 023 Attachment M Instrument Sequence Generation for Fertilizer for additional information.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 513**

*Analysis of Water Soluble Potash, Water Soluble Magnesium, and Water Soluble Boron using Inductively Coupled Plasma Optical Emission Spectrometry*  
*[FM604, FM803]*

Version: 4.0

Effective Date: 10/01/2021

Page 5 of 8

10.4.7.1. Example: AA99999-0AS

10.4.7.1.1. 0 corresponds to the sample weight (0, 1, 2.)

10.4.7.1.2. A corresponds to the different injections due to dilutions (A=1st injection, B=2nd injection, C=3rd injection.)

10.4.7.1.3. S corresponds to analysis type (S=soluble)

10.4.8. Dilution factors must be entered into the instrument sequence.

10.4.9. Initial Calibration:

10.4.9.1. The calibration curve is based on response area versus concentration in parts per million (ppm) of analyte.

10.4.9.2. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs; for calibration procedures and acceptance criteria.

10.4.9.3. %RSD for the sample replicates shall not exceed 10%.

10.4.10. Continuing Calibration Check Standard:

10.4.10.1. Refer to Data 051 Attachment C, Fertilizer Continuing Calibration Verification Acceptance Limits, for acceptance criteria.

10.4.10.2. Once the initial calibration has passed, a continuing calibration check standard must be analyzed after every 10 samples and at the end of the sequence. CCVs must be run for each separate curve.

10.4.10.3. The CCV may consist of any concentration level used in the initial calibration. Levels near the middle of the curve are recommended.

10.4.10.4. If a CCV response is not within the limits listed in DATA 051 Attachment C, Continuing Calibration Verification Acceptance Limits for AA and ICP, all samples analyzed before and after it must be reanalyzed.

10.4.10.5. If upon reanalysis of the CCV, the response is still not within the above acceptance criteria, then a new initial calibration must be performed. It is recommended that routine instrument maintenance be performed prior to the new curve analysis.

10.4.10.6. The Internal Standard for CCVs shall not exceed 80%-120% recovery.

10.4.10.7. %RSD for the sample replicates shall not exceed 10%.

10.4.11. If the ppm reading exceeds the calibration range, then the sample will be diluted and reanalyzed prior to quantification. Depending upon whether the ppm reading is below or above the calibration curve, a higher dilution of the sample or a larger amount weighed for the sample will be required.

10.4.12. Internal Standard criteria for samples is 80%-120% recovery.

10.4.13. %RSD for the sample replicates shall not exceed 10%.

10.5. Export Data

10.5.1. Data Manager exports data to the network drive TLHAESCHROME.

10.5.2. Upload data from the network drive to the LIMS.

10.6. Instrument Shutdown

10.6.1. See LABOP 140 Attachment I: Maintenance, Operation and Program of ICPOES

10.7. Data Reduction and Reporting

10.7.1. The Laboratory Control Sample must be evaluated first and if it fails, the associated samples in the batch are not uploaded.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 513**

*Analysis of Water Soluble Potash, Water Soluble Magnesium, and Water Soluble Boron using Inductively Coupled Plasma Optical Emission Spectrometry  
[FM604, FM803]*

Version: 4.0

Effective Date: 10/01/2021

Page 6 of 8

10.7.2. Data for all passing batches is uploaded into LIMS and evaluated. Refer to DATA 051, Data Evaluation for Feed, Fertilizer, and Pesticide Programs and ADMIN 023 Attachment O, Automatic Result Entry for Fertilizer for additional information.

**10.8. Quality Control**

**10.8.1. Laboratory Control Sample**

10.8.1.1. This standard is weighed with a set of samples and is processed in the same manner as the samples. The percent found must be within the limits listed in DATA 051 Attachment B, Fertilizer Laboratory Control Sample Acceptance Limits.

**10.9. Instrument Blank**

10.9.1. This blank consists of DI water and is used to ensure the instrument is free from contamination. The response of the instrument blank must be below the Method Detection Limit (MDL). Refer to QUALITY 112, Quality Control for additional information.

**11. WASTE MANAGEMENT**

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material         | Disposal Mechanism                                            |
|------------------|---------------------------------------------------------------|
| Aqueous extracts | Disposed of by dumping down the sink and flushing with water  |
| Instrument waste | Disposed of by dumping down the sink and flushing with water. |

**12. REFERENCES**

- 12.1. Analysis of Water Soluble Mg and Potash. FSFL-SOP 518, Rev. 3.0 (04/15/2009)
- 12.2. Hazardous Waste Disposal. PLS\_SOP 009, Rev. 1.0 (02/09/2006)
- 12.3. Laboratory Safety. BAEL\_SOP 200, Rev. 1.0 (09/12/2014)
- 12.4. Method for the Preparation of Water Soluble Magnesium and Potash. F\_SOP 204, Draft (09/01/2013)
- 12.5. Operation and Maintenance of Perkin Elmer Inductively Coupled Plasma–Optical Emission Spectrometry 5300 Dual View System. FSFL-SOP 526, Draft (01/14/2009)
- 12.6. Potash Analysis AA. FM-601, Rev. 1.0 (04/03/2002)
- 12.7. Using the Laboratory Information Management System. ADMIN 023, Rev. 1.0 (09/11/2015)
- 12.8. Determination of Metals and Trace Metals in Water and Wastes by Inductively Coupled Plasma-Atomic Emission Spectrometry. EPA Method 200.7, Rev. 4.4
- 12.9. Determination of Trace Elements in Drinking Water by Axially Viewed Inductively Coupled Plasma-Atomic Emission Spectrometry. EPA Method 200.5, Rev. 4.2
- 12.10. U.S. Code of Federal Regulations. 40 CFR, Parts 405-471.

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 513**

*Analysis of Water Soluble Potash, Water Soluble Magnesium, and Water Soluble Boron using Inductively Coupled Plasma Optical Emission Spectrometry*  
*[FM604, FM803]*

Version: 4.0

Effective Date: 10/01/2021

Page 7 of 8

**13. VERSION HISTORY**

| <b>Version</b> | <b>Date</b> | <b>Description</b>                                                                                                                                                                                                                                                                                                                                   | <b>Author</b>                  | <b>Editor</b>    |
|----------------|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------|------------------|
| 1.0            | 01.13.2016  | Original Version                                                                                                                                                                                                                                                                                                                                     | William Meeks<br>Naser Tassobi |                  |
| 2.0            | 03.29.2019  | Updated Outline Format<br>Moved troubleshooting and maintenance information to LABOP 140 Attachment I<br>Added tables for method specific instrument parameters<br>Added Section 8, reagent information,<br>Added 10.1.2<br>Added Section 11 Waste Management<br>Added required associated documents, in section 9 and in section 10 of the document |                                | Christian Amason |
| 3.0            | 09.10.2019  | Added soluble boron references to SOP and title.<br>Clarified 8.2.1 and 8.3.1<br>Added 9.8, 9.9, 9.11, 10.4.9.3, 10.4.12, 10.4.10.6, 10.4.10.7, 10.4.13<br>Minor grammatical fixes                                                                                                                                                                   |                                | Christian Amason |
| 4.0            | 03.11.2021  | Edited to reflect QC requirements spelled out in data 051                                                                                                                                                                                                                                                                                            |                                | Christian Amason |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 513**

*Analysis of Water Soluble Potash, Water Soluble Magnesium, and Water Soluble Boron using Inductively Coupled Plasma Optical Emission Spectrometry*  
*[FM604, FM803]*

Version: 4.0

Effective Date: 10/01/2021

Page 8 of 8

**14. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer  
Date\_\_\_\_\_

\_\_\_\_\_  
Christian Amason  
Instrumental Manager / Editor  
Date\_\_\_\_\_

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 514**

*Fertilizer Preparation and Analysis of Potash by Volumetric Sodium  
Tetraphenylboron using an Autotitrator*

Version: 2.0

Effective Date: 12/05/2023


Page 1 of 8

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the determination of potash in fertilizer by the volumetric sodium tetraphenylboron (STPB) method, using an Automatic Titrator.

## 2. SCOPE

2.1. The procedure delineated below is applicable to determining potassium in fertilizer samples. Potassium is extracted from the fertilizer with water and ammonium oxalate. STPB is added to the extract in excess. As potassium reacts with the STPB, it forms a precipitate. The excess STPB is then back-titrated with benzalkonium chloride (BAC). Clayton yellow is used as the end-point indicator.



2.2. The end-point is detected using an autotitrator with a DP5 phototrode set at 555nm. Due to the particulates generated during this titration, the mirror must be removed from the DP5 to allow for the reflected light to be detected from the particulates. As the amount of particulates increases, the amount of reflected light leads to an increase in signal. Once the end-point is reached, the sample turns pink, no longer reflects light, and produces a drop in signal.

2.3. This method is not suitable for samples that contain high levels of organic matter.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

3.2. A respirator or dust mask should be used when handling Aluminum Hydroxide. Samples must be boiled under the fume hood. Formaldehyde should be used under the fume hood.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Standardization of BAC solution
- 4.2. 10.2 Standardization of STPB solution
- 4.3. 10.3 Sample Preparation
- 4.4. 10.4 QC Preparation
- 4.5. 10.5 STPB Reaction
- 4.6. 10.6 Autotitration
- 4.7. 10.7 Post Titration Cleanup
- 4.8. 10.8 Calculations
- 4.9. 10.9 Sample Determination of %K<sub>2</sub>O

## 5. INTERFERENCE

5.1. Ammonia interference during the potassium precipitation is avoided by complexing the ammonia with formaldehyde under alkaline conditions before the potassium is precipitated.

5.2. Build-up of precipitate on the phototrode must be removed periodically, after every analyzed batch or every 9 analyzed samples.

5.3. High organic content causes interference.

## 6. APPARATUS and MATERIAL

- 6.1. Mettler Toledo Excellence Titrator T5 or equivalent
- 6.2. Mettler Toledo DP5 phototrode or equivalent
- 6.3. 100mL plastic beaker for titrator autosampler, per sample
- 6.4. 120mL amber container

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 514**

*Fertilizer Preparation and Analysis of Potash by Volumetric Sodium  
Tetraphenylboron using an Autotitrator*

Version: 2.0

Effective Date: 12/05/2023

Page 2 of 8

- 6.5. 50mL volumetric flask, per sample
- 6.6. 250mL volumetric flask, per sample
- 6.7. 250mL Erlenmeyer flask, per sample
- 6.8. Glass funnel, per sample
- 6.9. Filter paper, folded, Whatman no. 2V, 150 mm, per sample
- 6.10. Filter paper, Whatman no. 4 and no. 42
- 6.11. Funnel rack
- 6.12. Mechanical pipette or Class A Volumetric pipette
- 6.13. Hot plate
- 6.14. Smooth boiling stones, glass rods or equivalent (per sample)
- 6.15. Transfer pipette, per sample
- 6.16. Beaker
- 6.17. 3 reagent reservoirs
- 6.18. 2L beaker
- 6.19. 2L filtration flask
- 6.20. Stir bars

## 7. REAGENTS and MEDIA

- 7.1. Water, Deionized (DI) Type I
- 7.2. Triton X-100
- 7.3. Formaldehyde 37% - HCHO
- 7.4. Sodium hydroxide 20%
- 7.5. Aluminum Hydroxide -  $\text{Al}(\text{OH})_3$
- 7.6. Ammonium Oxalate -  $(\text{NH}_4)_2\text{C}_2\text{O}_4$
- 7.7. Potassium phosphate monobasic –  $\text{KH}_2\text{PO}_4$ 
  - 7.7.1. Dry  $\text{KH}_2\text{PO}_4$  at 106°C for approximately 4 hours
- 7.8. Sodium tetraphenylboron (STPB)
- 7.9. 17% Benzalkonium chloride (BAC)

## 8. REAGENT PREPARATION GUIDELINES

- 8.1. 4% Ammonium Oxalate
    - 8.1.1. Weigh 360g of Ammonium Oxalate into a beaker. Transfer content of beaker into the Ammonium Oxalate reservoir. Add DI water to the 9L mark and stir solution overnight.
  - 8.2. ~0.1% Triton X-100 wash solution
    - 8.2.1. Dilute ~ 2mL of Triton X-100 to a volume of ~2L with DI water, and stir.
  - 8.3. 0.04% Clayton yellow indicator solution
    - 8.3.1. Dissolve 0.04g of Clayton yellow in 100mL DI water.
    - 8.3.2. Transfer solution to an amber jar.
  - 8.4. ~0.625% BAC Solution
    - 8.4.1. Add 266mL 17% BAC to the BAC reservoir, fill the reservoir with DI water to the 7L mark, and stir.
    - 8.4.2. Standardize with STPB, section 10.1.
  - 8.5. ~1.2% STPB Solution
    - 8.5.1. Weigh 60g of STPB and transfer to a 2L beaker, add about 1500mL of DI water and stir to dissolve. Add 30g of  $\text{Al}(\text{OH})_3$  and stir.
- Note:** Aluminum hydroxide is added as a preservative. Since it is added in excess, it will not completely dissolve.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 514**

*Fertilizer Preparation and Analysis of Potash by Volumetric Sodium  
Tetraphenylboron using an Autotitrator*

Version: 2.0

Effective Date: 12/05/2023

Page 3 of 8

8.5.2. Filter the solution with rinses using a no. 4 Whatman filter and transfer the filtrate into 2L container. Filter again using a no. 42 Whatman filter. If the solution is still cloudy, re-filter using another no. 42 filter. Add 10mL of 20%NaOH, transfer the filtrate into the STPB reservoir, and fill with DI water to the 5L mark. Stir solution with a stir bar, overnight. Allow solution to stand for 48 hours.

8.5.3. Standardize with  $\text{KH}_2\text{PO}_4$ , section 10.2.

8.6. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015

8.7. The specific ingredient amounts can vary to accommodate the batch size.

8.8. Expiration Date: 12 months from date of preparation

8.9. Store at room temperature

## 9. ASSOCIATED DOCUMENTS

9.1. LABOP 121, Maintenance, Operation and Program of Balance

9.2. LABEL L003, Reagent

9.3. Form F00015, Reagent Preparation

## 10. SPECIFIC PROCEDURES

### 10.1. Standardization of BAC solution

10.1.1. Add exactly 4mL STPB solution to a minimum of 5 100mL plastic beakers.

10.1.2. Add ~20mL of DI water.

10.1.3. Add 2.5mL of 4% ammonium oxalate.

10.1.4. Add 2mL of 20% NaOH.

10.1.5. Add 5mL of 37% formaldehyde.

10.1.6. Add 6 drops or approximately 0.9mL of 0.04% Clayton yellow indicator solution into each 100mL beaker.

10.1.7. Bring to approximately 50mL final volume with DI water.

10.1.8. Prime the BAC burette.

10.1.9. Set DP5 phototrode to 555nm.

10.1.10. Open the Mettler LabX Software.

10.1.11. Titrate STPB with BAC using the instrument method titled "BAC Standardization".

10.1.12. Select "Start," then select "T5 Titrator 2."

10.1.13. Enter the STPB reagent ID and STPB volume into the corresponding sample positions.

10.1.14. BAC Standardization factor will be automatically stored as an "auxiliary value" for BAC in the Mettler LabX software.

10.1.15. The BAC solution is standardized upon preparation. It is re-standardized whenever a new STPB solution is prepared.

### 10.2. Standardization of STPB solution

10.2.1. Weigh 1.0g Potassium Phosphate Monobasic ( $\text{KH}_2\text{PO}_4$ ) standard into a minimum of 5 250mL volumetric flasks.

**NOTE:**  $\text{KH}_2\text{PO}_4$  standard must be dried at 100°C, for at least 2 hours, before being weighed. Leftover can be stored in a desiccator.

10.2.2. Add 50mL of 4% Ammonium Oxalate Solution to the 250-mL flasks.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 514**

*Fertilizer Preparation and Analysis of Potash by Volumetric Sodium  
Tetraphenylboron using an Autotitrator*

Version: 2.0

Effective Date: 12/05/2023

Page 4 of 8

10.2.3. Bring the Potassium standards to volume with DI water. Cap and invert/shake standards to mix.

10.2.4. Aliquot exactly 10.0mL of each standard into a labeled 50mL volumetric flask.

10.2.5. Add 2mL of 20% NaOH to each 50mL volumetric flask.

10.2.6. Under a fume hood, add 5mL of Formaldehyde to each 50mL volumetric flask.

10.2.7. Add exactly 15.0mL STPB solution into each 50mL volumetric flask.

10.2.8. Bring standards to volume with DI water. Cap and invert/shake standards to mix.

10.2.9. Let standards settle for approximately 15 minutes.

10.2.10. Place 2V folded filter paper into the funnels on the funnel apparatus and place labeled flasks on the funnel apparatus to collect filtrate.

10.2.11. Pour each standard into the funnels, allowing the standard to filter into the collection flasks.

10.2.12. Aliquot exactly 25.0mL of each standard into 100mL plastic beakers.

10.2.13. Add 6 drops or approximately 0.9mL of 0.04% Clayton yellow indicator solution into each 100mL beaker.

10.2.14. Bring to approximately 50mL volume.

10.2.15. Prime the BAC burette.

10.2.16. Set DP5 phototrode to 555nm.

10.2.17. Open the Mettler LabX Software.

10.2.18. Titrate standards with BAC using the instrument method titled "STPB Standardization".

10.2.19. Select "Start," then select "T5 Titrator 2."

10.2.20. Enter the standard ID, weight and STPB volume ("Correction Factor") into the corresponding sample position.

10.2.21. STPB factor will be automatically stored as an "auxiliary value" for STPB in the Mettler LabX software.

10.2.22. The STPB solution is standardized upon preparation. It is re-standardized whenever a new BAC solution is prepared, immediately after the standardization of BAC.

**10.3. Sample Preparation**

10.3.1. Samples with guarantees of  $\geq 6\%$   $K_2O$

10.3.1.1. Weigh sample into labeled 250mL volumetric flask according to the suggested weight code, in the LIMS.

10.3.1.1.1. A weight tolerance of  $\pm 0.05g$  is recommended.

10.3.1.2. Add 50mL of 4% Ammonium Oxalate Solution to the 250mL flask containing the sample.

10.3.2. Samples with guarantees of  $< 6\%$   $K_2O$

10.3.2.1. Weigh sample into labeled 50mL volumetric flask according to the suggested weight code, in the LIMS.

10.3.2.2. Add 4mL of 4% Ammonium Oxalate Solution to the 50mL flask containing the sample.

10.3.3. Add 2-3 drops of tributyl citrate.

10.3.4. Add boiling rod (or boiling stones) and boil for 20 minutes using a hot plate.

10.3.5. Note: DI water may be added to foaming samples.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 514**

*Fertilizer Preparation and Analysis of Potash by Volumetric Sodium  
Tetraphenylboron using an Autotitrator*

Version: 2.0

Effective Date: 12/05/2023

Page 5 of 8

10.3.6. Let samples cool to room temperature

10.3.7. Bring samples to volume with DI water. Cap and invert/shake samples to mix.

10.3.8. Let samples stand for 2 hours to settle. (sample solutions may be left overnight)

10.4. QC Preparation

10.4.1. Weigh 1.0 g  $\text{KH}_2\text{PO}_4$  into a 250mL volumetric flask.

10.4.2. If a liquid stock standard is available, an aliquot of the standard equivalent to 0.04g of  $\text{KH}_2\text{PO}_4$  (or 1.0 g x 10/250) can be added into the 50mL volumetric flask, together with 2.5mL of 4% Ammonium Oxalate (step 10.5.1.1.).

10.4.3. QC samples do not require boiling and are carried through the rest of the procedure.

10.5. STPB Reaction

10.5.1. Samples with guarantees of  $\geq 6\%$   $\text{K}_2\text{O}$

10.5.1.1. Aliquot exactly 10.0mL of sample extract into a labeled 50mL volumetric flask.

10.5.2. Samples with guarantees of  $< 6\%$   $\text{K}_2\text{O}$

10.5.2.1. Aliquot exactly 25.0mL of sample extract into a labeled 50mL volumetric flask.

10.5.3. Add 2mL of 20% NaOH to the 50-mL volumetric flask.

10.5.4. Under a fume hood, add 5mL of Formaldehyde to the 50mL volumetric flask.

10.5.5. Add exactly 15.0mL STPB solution into the 50mL volumetric flask.

10.5.6. Bring samples to volume with DI water. Cap and invert/shake samples to mix.

10.5.7. Let samples settle for approximately 15 minutes.

10.5.8. Place 2V folded filter paper into the funnels on the funnel apparatus and place labeled flasks on the funnel apparatus to collect filtrate.

10.5.9. Pour the sample into the funnel, allowing the sample to filter into the collection flask.

10.5.10. Aliquot exactly 25.0mL of sample into a 100mL plastic beaker.

10.5.11. Add 6 drops or approximately 0.9mL of 0.04% Clayton yellow indicator solution into the 100mL beaker.

10.5.12. Bring to approximately 50mL volume.

10.6. Autotitration

10.6.1. Prime the BAC burette.

10.6.2. Set DP5 phototrode to 555nm.

10.6.3. Open the Mettler LabX Software.

10.6.4. If analyzing samples with guarantees  $\geq 6\%$ , titrate samples using the instrument method titled "Potassium in Fertilizer by Volumetric STPB". If analyzing samples with guarantees  $< 6\%$ , use instrument method titled "Potassium in Fertilizer by Volumetric STPB\_low percent."

10.6.5. Select "Start," then select "T5 Titrator 2."

10.6.6. In "Task Name," enter the batch name. The batch name includes the date (yyymmdd), analyst's initials, and number of batch run for the day. Example: 181205jc01.

10.6.7. Enter the sample ID, weight and STPB volume ("Correction Factor") into the corresponding sample position. Start titration.

10.6.8. ~0.1% Triton X-100 may be used as a rinse.

10.6.9. Enter % found into the LIMS.

10.7. Post Titration Cleanup

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 514

#### *Fertilizer Preparation and Analysis of Potash by Volumetric Sodium Tetraphenylboron using an Autotitrator*

Version: 2.0

Effective Date: 12/05/2023

Page 6 of 8

10.7.1. Carefully wipe the phototrode, stirrer and dispenser tube with a Kimwipe to remove particulates accumulated during the titration.

10.7.2. Stirring at high speeds in ~0.1% Triton X-100 will aid in the removal of remaining particulates.

#### 10.8. Calculations

##### 10.8.1. Sample weight determination

10.8.1.1. Samples with K<sub>2</sub>O guarantees ≥6%:

$$\text{Sample Wt (g)} = \frac{0.3 \text{ g K}_2\text{O} * 100\%}{\% \text{ K}_2\text{O guaranteed in sample}}$$

10.8.1.2. Samples with K<sub>2</sub>O guarantees <6%:

$$\text{Sample Wt (g)} = \frac{0.024 \text{ g K}_2\text{O} * 100\%}{\% \text{ K}_2\text{O guaranteed in sample}}$$

##### 10.8.2. Standardization of BAC

$$\text{BAC Factor} = \frac{V_{(\text{STPB})}}{V_{(\text{BAC})}}$$

| Formula Term        | Definition                          | Mettler Symbol |
|---------------------|-------------------------------------|----------------|
| V <sub>(STPB)</sub> | Volume of STPB (mL)                 | m              |
| V <sub>(BAC)</sub>  | Volume of BAC titrant consumed (mL) | VEQ            |
| BAC Factor          | Ratio of STPB to BAC                | H[BAC]         |

##### 10.8.3. Standardization of STPB

$$\text{STPB}_{(\text{EX})} = \frac{V_{(\text{BAC})} \times H_{(\text{BAC})} \times 50\text{mL}}{25\text{mL}}$$

$$\text{STPB Factor} = \frac{34.61\% \times m_{(\text{KH}_2\text{PO}_4)} \times 10\text{mL} \times 25\text{mL} \times 1000}{(V_{(\text{STPB})} - \text{STPB}_{(\text{EX})}) 250\text{mL} \times 50\text{mL} \times 100}$$

| Formula Term                                | Definition                                            | Mettler Symbol |
|---------------------------------------------|-------------------------------------------------------|----------------|
| 34.61%                                      | % K <sub>2</sub> O of KH <sub>2</sub> PO <sub>4</sub> | -              |
| m <sub>(KH<sub>2</sub>PO<sub>4</sub>)</sub> | Weight (g) of KH <sub>2</sub> PO <sub>4</sub>         | m              |
| 250mL                                       | Final volume from extraction step                     | -              |
| 10mL                                        | Aliquot from extraction step                          | -              |
| 50mL                                        | Final Volume after STPB addition                      | -              |
| 25mL                                        | Aliquot from STPB addition step                       | -              |
| 100                                         | Conversion from % to decimal                          | -              |
| 1000                                        | Conversion from g to mg                               | -              |
| V <sub>(STPB)</sub>                         | Volume of STPB (mL)                                   | f              |
| V <sub>(BAC)</sub>                          | Volume of BAC titrant consumed (mL)                   | VEQ            |
| STPB <sub>(EX)</sub>                        | mL of STPB in excess                                  | C              |
| H <sub>(BAC)</sub>                          | Ratio of STPB to BAC                                  | H[BAC]         |
| STPB Factor                                 | mg K <sub>2</sub> O / mL STPB                         | H[STPB]        |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 514

#### *Fertilizer Preparation and Analysis of Potash by Volumetric Sodium Tetraphenylboron using an Autotitrator*

Version: 2.0

Effective Date: 12/05/2023

Page 7 of 8

#### 10.9. Sample Determination of %K<sub>2</sub>O

##### 10.9.1. Samples with K<sub>2</sub>O guarantees ≥6%

$$STPB_{(EX)} = \frac{V_{(BAC)} \times H_{(BAC)} \times 50mL}{25mL}$$

$$\%K_2O = \frac{(V_{(STPB)} - STPB_{(EX)}) 250mL \times 50mL \times F \times 100}{m_{(Sample)} \times 10mL \times 25mL \times 1000}$$

##### 10.9.2. Samples with K<sub>2</sub>O guarantees <6%

$$STPB_{(EX)} = \frac{V_{(BAC)} \times H_{(BAC)} \times 50mL}{25mL}$$

$$\%K_2O = \frac{(V_{(STPB)} - STPB_{(EX)}) 50mL \times 50mL \times F \times 100}{m_{(Sample)} \times 25mL \times 25mL \times 1000}$$

| Formula Term          | Definition                          | Mettler Symbol |
|-----------------------|-------------------------------------|----------------|
| m <sub>(sample)</sub> | Weight (g) of sample                | m              |
| 250mL / 50mL          | Final volume from extraction step   | -              |
| 10mL / 25mL           | Aliquot from extraction step        | -              |
| 50mL                  | Final Volume after STPB addition    | -              |
| 25mL                  | Aliquot from STPB addition step     | -              |
| 100                   | Conversion from decimal to %        | -              |
| 1000                  | Conversion from mg to g             | -              |
| V <sub>(STPB)</sub>   | Volume of STPB (mL)                 | f              |
| V <sub>(BAC)</sub>    | Volume of BAC titrant consumed (mL) | VEQ            |
| STPB <sub>(EX)</sub>  | mL of STPB in excess                | C              |
| H <sub>(BAC)</sub>    | Ratio of STPB to BAC                | H[BAC]         |
| F                     | mg K <sub>2</sub> O / mL STPB       | H[STPB]        |

## 11. WASTE MANAGEMENT

### 11.1. Material, marked for disposal, is removed via the following mechanism

| Material                                                      | Disposal Mechanism                                      |
|---------------------------------------------------------------|---------------------------------------------------------|
| Contents of flask after initial ammonium oxalate boiling step | Hazardous waste                                         |
| Contents of titration cups and leftover reagents              | Hazardous waste                                         |
| Used Filters                                                  | Dry in hood, place in box and submit as Hazardous Waste |

## 12. REFERENCES

- 12.1. *Potassium in Fertilizer, Volumetric Sodium Tetraphenylboron Method I*. Method 958.02, AOAC 15<sup>th</sup> Edition, [External 00138](#).
- 12.2. *Volumetric Sodium Tetraphenylboron Method I*. University of Kentucky, College of Agriculture Division of Regulatory Services, June 26, 2008, [External 00137](#).
- 12.3. Ingram, C.W., Bryant, M., Terry, D. (2008, June). *Comparison of Manual and Automatic Titrametric Analysis of Potash in Fertilizer using AOACI Method 958.02*. Presentation at ASPFCO, New Orleans, LA, [External 00136](#).

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 514**

*Fertilizer Preparation and Analysis of Potash by Volumetric Sodium  
Tetraphenylboron using an Autotitrator*

Version: 2.0

Effective Date: 12/05/2023

Page 8 of 8

**13. VERSION HISTORY**

| Version | Date       | Description                                                                                                                                           | Author                                | Editor                |
|---------|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------|-----------------------|
| 1.0     | 03.01.2019 | Original Version                                                                                                                                      | William Meeks /<br>Dr. Patrizia Lemma |                       |
| 2.0     | 12.05.2023 | Updater reference section to add<br>external document numbers.<br>Removed signature page<br>Reformatted section 8 to align with<br>current practices. |                                       | William Meeks,<br>Jr. |

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 515**

*Analysis of NAU using Flow Injection Analyzer*  
*[FM452, FM472, FM462]*

Version: 4.0

Effective Date: 06/30/2023

Page 1 of 7

## **1. PURPOSE**

1.1. This standard operating procedure (SOP) provides instruction for the analysis of Nitrate ( $\text{NO}_3$ ), Ammonia ( $\text{NH}_4^+$ ) and Urea in fertilizer samples by utilizing a Flow Injection Analyzer (FIA).

## **2. SCOPE**

2.1. The procedure delineated below is applicable to fertilizer samples that are analyzed in the Bureau of Agricultural Environmental Laboratories (BAEL). The method is for the determination of Nitrate-Nitrogen [FM 452], Ammoniacal-Nitrogen [FM 472] and Urea-Nitrogen [FM 462] in DI Water extracts of fertilizer samples.

2.2.  $\text{NO}_3$  -Nitrate is quantitatively reduced to nitrite by passage of the sample through a copperized cadmium column. The nitrite (reduced nitrate plus original nitrite) is then determined by diazotizing with sulfanilamide followed by coupling with N-(1-naphthyl)ethylenediamine dihydrochloride. The resulting water soluble dye has a magenta color which is read at 520 nm.

2.3.  $\text{NH}_4$  -This method is based on the Berthelot reaction. Ammonia reacts with alkaline phenol, then with sodium hypochlorite to form indophenol blue. Sodium nitroprusside (nitroferricyanide) is added to enhance sensitivity. The absorbance of the reaction product is measured at 630 nm, and is directly proportional to the original ammonia concentration/

2.4. Urea -This is a two channel, low flow method, which measures ammonia on one channel (FM-471), and the sum of ammonia and urea on a second channel. The sample is injected into a urease solution. The result on the ammonia channel (FM-471) consists of only the ammonia present in the sample. Urease cleaves urea ( $\text{CH}_4\text{N}_2\text{O}$ ) to  $\text{NH}_3 + \text{CO}_2$  on the second channel. The result consists of the ammonia initially present in the sample, plus the ammonia produced by the action of urease on urea. When ammonia is heated with salicylate and hypochlorite in an alkaline phosphate buffer, a colored product results. Absorbance of the colored product at 660 nm is directly proportional to the ammonia concentration. The color is intensified by the addition of sodium nitroprusside. The method detection limit is 1 mg N/L as  $\text{NH}_4$ . At the end of the analysis, the ammonia channel (FM-471) result is subtracted from the urea + ammonia result to give the percent urea N present in the sample. An ammonia standard is used to calibrate both channels. Urea stock solutions are analyzed to determine the percent conversion of urea to ammonia. Percent conversions of urea to ammonia are generally in the range of 93-98%.

## **3. HEALTH and SAFETY**

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## **4. OUTLINE of PROCEDURE**

- 4.1. Preparation of Reagent
- 4.2. Preparation of Standard
- 4.3. Lachat Operation
- 4.4. Initial Calibration
- 4.5. Laboratory Control Sample
- 4.6. Continuing Calibration Verification
- 4.7. Run Sequence
- 4.8. Printing and Reporting Data
- 4.9. Instrument Shutdown
- 4.10. Quality Control
- 4.11. Attachment A, Nitrate Parameters and Manifold Diagram

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 515**

*Analysis of NAU using Flow Injection Analyzer*  
*[FM452, FM472, FM462]*

Version: 4.0

Effective Date: 06/30/2023

Page 2 of 7

- 4.12. Attachment B, Ammonia Parameters and Manifold Diagram
- 4.13. Attachment C, Urea/Ammonia Parameters and Manifold Diagram

## **5. INTERFERENCE**

- 5.1. Residual chlorine can interfere by oxidizing the reductor column.
- 5.2. Low results for Nitrate would be obtained for samples that contain high concentrations of iron, copper or other metals. In this method, EDTA can be added to the buffer to reduce this interference.
- 5.3. Calcium and magnesium ions may precipitate if present in sufficient concentration. EDTA may be added to the sample in-line in order to prevent this problem
- 5.4. Eliminate any marked variation in acidity or alkalinity among samples because intensity of measured color is pH-dependent. Likewise, ensure that pH of standard ammonia solutions approximates that of samples.
- 5.5. Non-volatile amines such as cysteine, ethanolamine and ethylenediamine may cause a decrease in ammonia sensitivity.
- 5.6. Nitrogen oxides must be eliminated before color is formed to prevent interference.
- 5.7. Iron and dimeric silicates may cause high results, if present.
- 5.8. Sample turbidity, color and particulates must be removed by filtration, prior to analysis.

## **6. APPARATUS and MATERIAL**

- 6.1. Autosampler
- 6.2. Autosampler rack, 15 x 6
- 6.3. Balance, analytical minimum capacity 0.001g
- 6.4. Balance, top loading minimum capacity 0.1g
- 6.5. Beaker
- 6.6. Filter, 13mm IC Acrodisc with 0.45µm syringe filter, or equivalent
- 6.7. Flask, class A volumetric
- 6.8. Flow Injection Analyzer
- 6.9. pH paper
- 6.10. Pipette, Eppendorf or equivalent
- 6.11. Pipette, class A volumetric
- 6.12. Pipette tip
- 6.13. Stir Bar
- 6.14. Syringe, plastic disposable
- 6.15. Test tube, borosilicate culture 13 x100mm or equivalent
- 6.16. Weighing boat, plastic disposable
- 6.17. Dialysis membrane
- 6.18. Cadmium column

## **7. REAGENTS and MEDIA**

- 7.1. Sodium phosphate heptahydrate
- 7.2. Sodium hydroxide
- 7.3. Potassium sodium tartrate tetrahydrate
- 7.4. Sodium salicylate
- 7.5. Sodium nitroprusside
- 7.6. Brij-35 (>95%); Other concentrations may be used.
- 7.7. Urease
- 7.8. Sodium hypochlorite
- 7.9. Ammonium chloride
- 7.10. Hydrochloric acid, concentrated (32%)
- 7.11. Sulfanilamide
- 7.12. N-(1-naphthyl)ethylenediamine dihydrochloride (NED)

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 515**

*Analysis of NAU using Flow Injection Analyzer*  
*[FM452, FM472, FM462]*

Version: 4.0

Effective Date: 06/30/2023

Page 3 of 7

7.13. Phosphoric acid, concentrated (85%)

7.14. Water, Deionized (DI) Type I

## **8. REAGENT PREPARATION GUIDELINES**

### **8.1. 10% w/v Sodium Hydroxide**

8.1.1. In a 100mL flask, add approximately 50mL of deionized water and add 10g of sodium hydroxide (NaOH). Mix thoroughly to dissolve solids and dilute to volume.

**Caution:** the solution will get very hot!

### **8.2. 20% v/v Hydrochloric Acid**

8.2.1. In a 100mL flask, add approximately 60mL of deionized water, add 20mL of concentrated hydrochloric acid (HCl). Dilute to volume and mix by inversion.

**Caution:** the solution will get very hot!

### **8.3. Ammonium Chloride Buffer-Acceptor**

8.3.1. In a 2000mL flask, add approximately 1600mL of deionized water, and add 42.5g of ammonium chloride (NH<sub>4</sub>Cl). Adjust pH to 5.5-6.0 with 10% sodium hydroxide. Dilute to volume and mix with a magnetic stirrer.

### **8.4. Ammonium Chloride Buffer-Acceptor**

8.4.1. In a 2000mL flask, add approximately 1600mL of deionized water. Add 127.5 grams of ammonium chloride (NH<sub>4</sub>Cl). Adjust pH to 5.5-6.0 with 10% sodium hydroxide. Dilute to volume mix with a magnetic stirrer.

### **8.5. Sulfanilamide color reagent**

8.5.1. In a 2000mL flask, add approximately 1200mL of deionized water. Add 80.0g of sulfanilamide [C<sub>6</sub>H<sub>8</sub>N<sub>2</sub>O<sub>2</sub>S] and 2.0g N-(1-naphthyl)ethylenediamine dihydrochloride (NED), stir with a magnetic stirrer. Add 200mL concentrated phosphoric acid (H<sub>3</sub>PO<sub>4</sub>) and dilute to volume.

### **8.6. Buffer Solution**

8.6.1. In a 2000mL flask, add approximately 1200mL of deionized water. Add 40.0g of sodium hydroxide (NaOH), 84.0 g sodium phosphate heptahydrate (Na<sub>2</sub>HPO<sub>4</sub>·7H<sub>2</sub>O) and 100.0 g sodium potassium tartrate tetrahydrate (NaKC<sub>4</sub>H<sub>4</sub>O<sub>6</sub>·4H<sub>2</sub>O), stir with a magnetic stirrer and dilute to volume.

### **8.7. Nitroprusside/Salicylate Color Reagent**

8.7.1. In a 2000mL flask, add approximately 1200mL of deionized water. Add 300.0g of sodium salicylate, 2.0g sodium nitroprusside and 0.8g Brij-35 (>95%), stir with a magnetic stirrer, and dilute to volume. If the Brij-35 is a lower concentration, a larger weight will be needed to meet the necessary final concentration. For example: Brij-35(30%w/v) would need 2.67g weighed.

### **8.8. Hypochlorite Solution**

8.8.1. In a 1000mL flask, add approximately 800mL of deionized water. Add 100mL of approximately 5% sodium hypochlorite, stir with a magnetic stirrer and dilute to volume.

Note: If full peaks are not observed at the higher curve concentrations, add additional hypochlorite.

8.8.2. Expires 5 days after preparation; store in refrigerator

### **8.9. Urease Solution**

8.9.1. In a 500mL flask, add approximately 400mL of deionized water. Add 0.260g urease, stir with a magnetic stirrer and dilute to volume.

8.9.2. Expires 5 days after preparation; store in refrigerator

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 515**

*Analysis of NAU using Flow Injection Analyzer*  
*[FM452, FM472, FM462]*

Version: 4.0

Effective Date: 06/30/2023

Page 4 of 7

8.10. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015

8.11. The specific ingredient amounts can vary to accommodate the batch size.

8.12. Expiration Date: 1 month from date of preparation unless otherwise specified.

8.13. Store at room temperature unless otherwise specified.

## **9. ASSOCIATED DOCUMENTS**

- 9.1. ADMIN 023 Attachment M, Instrument Sequence Generation for Fertilizer
- 9.2. ADMIN 023 Attachment O, Automatic Result Entry for Fertilizer
- 9.3. DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs
- 9.4. DATA 051 Attachment A, Fertilizer Analysis Requirements
- 9.5. DATA 051 Attachment B, Fertilizer Laboratory Control Sample Acceptance Limits
- 9.6. DATA 051 Attachment C, Fertilizer Continuing Calibration Verification Acceptance Limits
- 9.7. Form F00015, Reagent Preparation
- 9.8. Label L003, Reagent
- 9.9. LABOP 140 Attachment F, Maintenance and Operation of Flow Injection Analyzer
- 9.10. QUALITY 112, Quality Control

## **10. SPECIFIC PROCEDURES**

### **10.1. Preparation of Reagent**

- 10.1.1. Reagents are prepared according to the applicable instructions in section 8.

### **10.2. Preparation of Standard**

10.2.1. Stock Standard 2000ppm  $\text{NH}_4\text{-N}$  and 1000ppm  $\text{NO}_3\text{-N}$  in DI Water. In a 500mL volumetric flask, dissolve 4.7170g ammonia sulfate  $(\text{NH}_4)_2\text{SO}_4$  and 3.6110g Potassium Nitrate  $(\text{KNO}_3)$  in 100 ml of DI Water. Dilute to mark with DI Water and invert to mix.

**NOTE:** Stock standard solutions are stable for 1 month. Purchased stock standards should adhere to manufacturer's expiration date.

#### **10.2.2. Intermediate Standards**

- 10.2.2.1. Standard A (600mg  $\text{NH}_4\text{-N/L}$ , 300mg  $\text{NO}_3\text{-N/L}$ )
- 10.2.2.2. Fill a 100mL flask with 50mL of DI Water
- 10.2.2.3. Using pipette, add 30mL of 2000/1000ppm Stock Solution
- 10.2.2.4. Dilute to 100mL with DI water
- 10.2.2.5. Invert to mix

#### **10.2.3. Standard B (300mg $\text{NH}_4\text{-N/L}$ , 150mg $\text{NO}_3\text{-N/L}$ )**

- 10.2.3.1. Fill a 200mL flask with 100mL of DI Water
- 10.2.3.2. Using pipette, add 30mL of 2000/1000ppm Stock Solution
- 10.2.3.3. Dilute to 200mL with DI Water
- 10.2.3.4. Invert to mix

#### **10.2.4. Standard C (150mg $\text{NH}_4\text{-N/L}$ , 75mg $\text{NO}_3\text{-N/L}$ )**

- 10.2.4.1. Fill a 200mL flask with 100mL of DI Water
- 10.2.4.2. Using pipette, add 15mL of 2000/1000ppm Stock Solution
- 10.2.4.3. Dilute to 200mL with DI Water
- 10.2.4.4. Invert to mix

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 515**

*Analysis of NAU using Flow Injection Analyzer*  
*[FM452, FM472, FM462]*

Version: 4.0

Effective Date: 06/30/2023

Page 5 of 7

10.2.5. Standard D (60mg NH<sub>4</sub>-N/L, 30mg NO<sub>3</sub>-N/L)

10.2.5.1. Fill a 100mL flask with 50mL of DI Water

10.2.5.2. Using pipette, add 3mL of 2000/1000ppm Stock Solution

10.2.5.3. Dilute to 100mL with DI Water

10.2.5.4. Invert to mix

10.2.6. Standard E (30mg NH<sub>4</sub>-N/L, 15mg NO<sub>3</sub>-N/L)

10.2.6.1. Fill a 200mL flask with 100mL of DI Water

10.2.6.2. Using pipette, add 3mL of 2000/1000ppm Stock Solution

10.2.6.3. Dilute to 200mL with DI Water

10.2.6.4. Invert to mix

10.3. Lachat Operation

10.3.1. Set up manifolds as shown in Attachment A, Nitrate Manifold Diagram.

10.3.2. Set up manifolds as shown in Attachment B, Ammonia Manifold Diagram.

10.3.3. Set up manifolds as shown in Attachment C, Urea/Ammonia Manifold Diagram.

10.3.4. Place all the lines in DI and switch the pump on to check for any leaks.

10.3.5. Stop the pump.

10.3.6. Take the lines out from DI water and place in the respective reagent containers.

10.3.7. Restart the pump.

10.3.8. Allow 20 to 30 minutes for instrument warm-up time. Monitor the base line until a stable baseline is obtained.

10.4. Initial Calibration

10.4.1. A minimum of five calibration standards must be injected. For analytes that require a second order fit, a minimum of six calibration curve standards is required. One of the standards should be near the Method Detection Limit (MDL) concentration. There should be at least one order of magnitude difference in the curve range.

10.4.2. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs; for calibration procedures.

10.5. Laboratory Control Sample

10.5.1. Refer to Data 051 Attachment A, Fertilizer Analysis Requirements, for acceptance criteria.

10.5.2. Reference material used for the Lab Standard is defined in PREP 453 attachment A.

10.6. Continuing Calibration Verification

10.6.1. Refer to Data 051 Attachment B, Fertilizer Continuing Calibration Verification Acceptance Limits, for acceptance criteria.

10.7. Run Sequence

10.7.1. Refer to Data 051 Attachment C, Fertilizer Analysis Requirements, for sequence generation and ADMIN 023 Attachment M, Instrument Sequence Generation for Fertilizer, for sample nomenclature.

10.8. Printing and Reporting Data

10.8.1. Click on 'Tool' button and select 'Custom Report'

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 515**

*Analysis of NAU using Flow Injection Analyzer*  
*[FM452, FM472, FM462]*

Version: 4.0

Effective Date: 06/30/2023

Page 6 of 7

10.8.2. Print the Report

10.8.3. Be sure the report includes the following

- 10.8.3.1. Cup number
- 10.8.3.2. Sample identification (ID)
- 10.8.3.3. Final concentration
- 10.8.3.4. Dilution factor
- 10.8.3.5. Calibration charts
- 10.8.3.6. Detection date and time

10.8.4. Upload data according to ADMIN 023 Attachment O, Automatic Result Entry for Fertilizer

10.9. Instrument Shutdown

10.9.1. Clean the equipment, auto sampler, sample racks, etc. with a damp cloth.

10.9.2. Refer to LABOP 140 Attachment F, Maintenance and Operation of Flow Injection Analyzer, for shutdown procedures and additional information on maintenance and troubleshooting.

10.10. Quality Control

10.10.1. Refer to DATA 051 Attachments B, Fertilizer Laboratory Control Sample Acceptance Limits, for acceptance criteria.

10.10.2. Refer to QUALITY 112, Quality Control, to determine the outcome of Quality Control (QC) samples.

## **11. WASTE MANAGEMENT**

11.1. Material, marked for disposal, is removed via the following mechanism:

| <b>Material</b>         | <b>Disposal Mechanism</b>                   |
|-------------------------|---------------------------------------------|
| Instrument Effluent     | Neutralize with NaOH and dispose down drain |
| Excess filtered samples | Dispose down drain                          |
| Unused Reagents         | Neutralize with NaOH and dispose down drain |
| Cd Column               | Submit to Hazardous Waste                   |

## **12. REFERENCES**

- 12.1. Determination of Nitrate in Fertilizers by Flow Injection Analysis Colorimetry. QuickChem Method 14-107-04-1-B.
- 12.2. Determination of Ammonia in Fertilizers by Flow Injection Analysis Colorimetry. QuickChem Method 14-107-06-2-C
- 12.3. Determination of Urea in Fertilizers by Flow Injection Analysis Colorimetry. QuickChem Method 14-206-00-2-A
- 12.4. Lachat QC FIA + 8000 Flow Injection Analyzer Manual

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 515**  
*Analysis of NAU using Flow Injection Analyzer*  
*[FM452, FM472, FM462]*

Version: 4.0

Effective Date: 06/30/2023

Page 7 of 7

**13. VERSION HISTORY**

| <b>Version</b> | <b>Date</b> | <b>Description</b>                                                                                                                                                                                                                          | <b>Author</b> | <b>Editor</b>    |
|----------------|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|------------------|
| 1.0            | 03.29.2019  | Original Version                                                                                                                                                                                                                            | Carl Huang    |                  |
| 2.0            | 05.21.2021  | Updated formatting throughout<br>Updated reagent expiration information in sections 8.8.2. and 8.9.2.<br>Updated QAO signature block to replace Elizabeth Schneider with Quentin Cooper.<br>Updated storage specifications in section 8.13. |               | Steven Daley     |
| 3.0            | 07.21.2022  | Section 8.0 adjusted volume of sodium hypochlorite and added note.<br>Section 12.0 Updated References added 12.2 and 12.3.                                                                                                                  |               | Christian Amason |
| 4.0            | 03.31.2023  | Updated sections 7.6. and 8.7.1. to add allowance for different Brij concentrations.<br><br>Removed signature page                                                                                                                          |               | Christian Amason |

# ATTACHMENT

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 515 Attachment A


#### *Nitrate Manifold Diagram*

Version: 1.0

Effective Date: 03/29/2019

Page 1 of 1

### 11.3. NITRATE MANIFOLD DIAGRAM



**Carrier:** DI water

**Manifold Tubing:** 0.8 mm (0.032 in) i.d. This is 5.2  $\mu\text{L}/\text{cm}$ .

**AE Sample Loop:** 8 cm

**QC8000 Sample Loop:** 13.5 cm

**Interference Filter:** 520 nm


**Apparatus:** An injection valve, a 10 mm path length flow cell, and a colorimetric detector module is required.

**4.5:** 70 cm of tubing on a 4.5 cm coil support

**12:** 255 cm of tubing on a 12 cm alternating coil support

**Note 1:** 6 inch pathlength dialysis block and a Type C membrane.

**Note 2:** This is a 2 state switching valve used to place the cadmium column in-line with the manifold.



# ATTACHMENT

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 515 Attachment B


#### *Ammonia Manifold Diagram*

Version: 1.0

Effective Date: 03/29/2019

Page 1 of 1

#### AMMONIA MANIFOLD DIAGRAM



**Carrier:**

**Manifold Tubing:**

**QC8000 Sample Loop:**

**Interference Filter:**

**Apparatus:**

**4.5:**

**12:**

**Note 1:**

**Note 2:**

**Note 3:**

**Note 4:**

**Note 5:**

**D.I. water**

**0.5 mm (0.022 in) i.d. This is 5.2 µL/cm.**

**13 cm**

**660 nm**

An injection valve, a 1 mm path length flow cell, and a colorimetric detector module is required

**70 cm of tubing on a 4.5 cm coil support**

**255 cm of tubing on a 12 cm alternating coil support**

The flow cell is **1 mm** path length.

50 cm back pressure loop, 0.52 mm (0.022 in) i.d.

175 cm of tubing on the heater block.

650 cm of tubing on the heater block.

The same heating module is used for both (see notes 3 and 4) heaters.

# ATTACHMENT

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 515 Attachment C

#### Urea/Ammonia Manifold Diagram

Version: 1.0

Effective Date: 03/29/2019

Page 1 of 2

## FM-462


Urea Nitrogen – Flow Injecti

Probe in Sample Period: 35 s

### Valve Timing:

Load Time: 0 s  
Load Period: 25 s  
Inject Period: 65 s

### UREA/AMMONIA MANIFOLD DIAGRAM



**ATTACHMENT**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 515 Attachment C  
*Urea/Ammonia Manifold Diagram*

Version: 1.0

Effective Date: 03/29/2019

Page 2 of 2

**Carrier:** Urease Solution (Reagent 3)

**Manifold Tubing:** 0.5 mm (0.022 in) i.d. This is 2.5 µL/cm.

**AE Sample Loop:** 7.5 cm x 0.5 mm i.d.

**QC8000 Sample Loop:** 13 cm x 0.5 mm i.d.

**Interference Filter:** 660 nm

---

**Apparatus:** An injection valve, a 1 mm path length flow cell, and a colorimetric detector module is required. The  shows xxx cm (see notes below) of tubing wrapped around the heater block at the specified temperature.

**4.5:** 70 cm of tubing on a 4.5 cm coil support

**12:** 255 cm of tubing on a 12 cm alternating coil support

---

**Note 1:** The flow cell is 1 mm path length

**Note 2:** 50 cm backpressure loop x 0.5 mm i.d.

**Note 3:** 175 cm of tubing on the heater block

**Note 4:** 650 cm of tubing on the heater block

**Note 5:** The same heating module is used for both heaters (see notes 3 and 4)

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 516**

*Analysis of Available Phosphorous and Soluble Potash in Fertilizers  
using ICP-OES*

Version: 2.0

Effective Date: 10/16/2023

Page 1 of 8

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the determination of available phosphorous ( $P_2O_5$ ) and soluble potash ( $K_2O$ ) in mixed or pure material fertilizer samples, in both liquid and dry matrices by Inductively Coupled Plasma – Optical Emission Spectrometry (ICP-OES).

## 2. SCOPE

2.1. The procedure delineated below is applicable to an ICP-OES determination of available phosphorous (APA) and soluble potash in fertilizer. Fertilizer samples are extracted with a citrate-EDTA solution (as per PREP 461, *Sample Preparation for Available Phosphorus and Soluble Potash in Fertilizer*) and analyzed by ICP-OES. The analysis is achieved by measuring the intensity of the electromagnetic radiation produced by the analyte in the inductively coupled plasma. A quantitative determination of the amount of analyte present can be made at the specific wavelength by measuring the excited atoms and ions emitted. The intensity of this emission is indicative of the concentration of the analyte within the sample. The determined amount of potassium (K) is converted to potash ( $K_2O$  equivalent) and reported as such. The determined amount of phosphorous (P) is converted to APA ( $P_2O_5$  equivalent) and reported as such.

2.2. Standards solutions used to calibrate the ICP-OES are made to match the matrix of the final sample extracts. Cesium is used as an ionization buffer to minimize easily ionized element effects.

2.3. This method is adapted from AOAC Method 2015.18. The method described in this SOP deviates from Method 2015.18 by bringing the samples to a final volume of 500 mL instead of 250 mL. Thus, the citrate-EDTA of the sample is diluted prior to analysis, eliminating the need to dilute the sample inline by increasing the tubing size of the buffer/internal standard peristaltic tube and reducing the tubing size of the sample peristaltic tube. Sample weight was increased from 0.5 g to 0.75 g to achieve required detection limit.

2.4. This method is restricted to use by or under the supervision of a trained analyst. Each analyst must demonstrate the ability to generate acceptable results with this method. Refer to ADMIN 016, Training and Competency, for demonstration of competence procedures.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

3.2. Torch and interface regions of the ICP can be very hot and caution is needed when working around these areas if the ICP has recently been running. The ICP generates very high voltages. Specific safeties are in place on the instrument to protect the operator from accidental exposure to electrical shock; these safeties should not be subverted. Extreme caution should be used when performing any maintenance that would expose the operator to high voltages. Protective covers shall remain in place at all times, except when special maintenance and/or repairs require their removal.

3.3. Exhaust from combustion fumes and vapors from the ICP torch housing can produce toxic vapors and may pose a health hazard if not adequately vented to the outside.

## 4. OUTLINE of PROCEDURE

- 10.1 ample Handling
- 10.2 tandards
- 10.3 Instrument Startup
- 10.4 Instrument Parameters

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 516**

*Analysis of Available Phosphorous and Soluble Potash in Fertilizers  
using ICP-OES*

Version: 2.0

Effective Date: 10/16/2023

Page 2 of 8

10.5 ample Analysis  
10.6 port Data  
10.7 Instrument Shutdown  
10.8 Quality Control  
10.9 Data Reduction and Reporting

## 5. INTERFERENCE

5.1. Spectral interferences are caused by background emission from continuous or recombination phenomena, stray light from the line emission of high concentration elements, overlap of a spectral line from another element, ionization interferences (for K), or unresolved overlap of molecular band spectra. Interferences may be overcome by dilution of the digestate, application of interelement correction factors, use of alternate wavelengths, or addition of cesium chloride ionization buffer (for K). A target of at least 2200 mg/kg cesium chloride solution must be presented at the nebulizer.

5.2. Plasma efficiency can be reduced by the carbon in the citrate and EDTA. Therefore, samples must be diluted prior to analysis. This can be achieved by increasing the peristaltic tubing for the buffer/internal standard and decreasing the peristaltic tubing for the sample tubing or diluting the sample to a larger final volume as per PREP 461.

5.3. Samples containing copper (Cu) will interfere with the P wavelengths of 213.617 nm and 214.914 nm due to Cu emission intensity at the wavelengths of 213.598 nm and 214.898 nm. Spectral deconvolution must be utilized if reporting P from the affected wavelengths. **Multicomponent Spectral Fitting (MSF)** is utilized in this method according to the procedure below:

**5.3.1. MSF Modeling:** Spectra must be collected from a matrix blank (B), a standard of the analyte of interest (A), and a standard of the interfering analyte(s) (I) using the same analysis method that will be used to normally analyze samples.

**Note:** The Peak Algorithm must be set Peak Area or Peak Height within the analysis method. The MSF algorithm is not used during this time.

5.3.2. Spectral profile is independent of concentration, however, to optimize the MSF model, the concentration of the analyte should be approximately 100 times its detection limit at the affected wavelength and the concentration of the interfering analyte should approximate maximum concentration that could be present.

E.g., ~100 times the detection limit of P would be approximately 500 ppm. In this example, 400 ppm P standard was selected. A 40 ppm Cu standard was selected to represent the maximum concentration of the interfering Cu (2 or 3% Cu in the sample).

**Note:** All Analytes must be in the same matrix as the samples (final buffer).

5.3.3. After the analysis of the blank, analyte, and interferent, build the MSF model using the Examine MSF window. Label the blank (B), analyte (A) and interferent (I) on the wavelength of interest (e.g., 213.617 nm). Once the MSF file is saved, Select MSF as the Peak Algorithm within the analysis method and specify which MSF file to be used on the Spectral Corrections tab.

5.3.4. Confirm the new MSF model is function as intended by reanalyzing the interferent and a known sample or standard such as a LS. The resulting interferent should not be detected at the wavelength of interest.

5.3.5. If a wavelength calibration is performed, the MSF must be recreated.

**Note:** Wavelength calibration is performed during the annual PM.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 516**

*Analysis of Available Phosphorous and Soluble Potash in Fertilizers  
using ICP-OES*

Version: 2.0

Effective Date: 10/16/2023

Page 3 of 8

5.4. Fertilizers containing phosphite will also be detected by this method.

5.5. Due to the susceptibility of the citrate-EDTA solution to microbial degradation, the citrate-EDTA solution has a 2-week expiration when stored in the dark and a 16-hour expiration for sample extracts. An inhouse longevity study demonstrated microbial degradation could be mitigated by utilizing refrigeration. The expiration dates found in this SOP have been extended based on this study.

5.6. For most fertilizer samples, Scandium (Sc) and Beryllium (Be) can be used as the internal standard for P. Polyphosphates present additional challenges to the plasma that are not detected by Sc; in this case, Be may be used as internal standard for P.

## **6. APPARATUS and MATERIAL**

### **6.1. Analytical Instruments**

6.1.1. Dual view ICP-OES

6.1.2. Autosampler

6.1.3. Peristaltic Pump capable of handling three introduction lines

6.1.4. Chiller

6.1.5. Integration software

### **6.2. ICP autosampler tubes**

6.3. Pump tubing: black/black, orange/green, red/red

6.4. Volumetric flasks

6.5. Volumetric pipettes

6.6. 0.45µm Syringe Filters

## **7. REAGENTS and MEDIA**

7.1. Ammonium citrate dibasic, ACS grade or equivalent

7.2. Ammonium hydroxide, ACS grade or equivalent

7.3. Beryllium Solution, 1000 ppm

7.4. Cesium Chloride (CsCl), Trace metal grade, solid

7.5. EDTA, disodium salt, dihydrate, >99%

7.6. Hydrochloric Acid (HCl) concentrated, Trace metal grade

7.7. Nitric Acid (HNO<sub>3</sub>) concentrated, Trace metal grade

7.8. Scandium Solution, 1000 ppm

7.9. Triton X-100

7.10. Water, Deionized (DI), Type I

## **8. REAGENT PREPARATION GUIDELINES**

### **8.1. Citrate-EDTA Solution**

8.1.1. Directions: Add 25 g disodium EDTA and 50 g dibasic ammonium citrate to a 2-L volumetric flask containing ~1500 mL DI water. Adjust to near neutral by adding 30 mL of a solution of ammonium hydroxide-water (1+1, v/v), under a hood. (Volumes are scalable as needed.)

8.1.2. While stirring and using a pH electrode, adjust the final pH to 7.00 by adding the ammonium hydroxide-water (1+1, v/v) solution dropwise. After a stable pH of 7.00 ±0.02 has been attained, dilute to volume with DI and mix.

8.1.3. Citrate-EDTA solution has a **2-months expiration** when stored refrigerated.

### **8.2. Carrier Solution**

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 516**

*Analysis of Available Phosphorous and Soluble Potash in Fertilizers  
using ICP-OES*

Version: 2.0

Effective Date: 10/16/2023

Page 4 of 8

8.2.1. Directions: Fill a 1L volumetric flask almost to volume with DI water. Add 1-2 drops Triton X-100, then fill to the mark with DI water. Triton X – 100 causes beading in some spray chambers, therefore, it may be left out in those cases. Invert the flask several times to mix.

**8.3. Wash Solution**

8.3.1. Directions: Add approximately 5L DI water to a 15-L carboy, add 307 mL HCl, swirl to mix. Add DI water to the 15-L mark, cap the carboy, and shake to mix.

Optional: 7.5 mL of Triton X-100 may be added. (Triton X – 100 causes beading in some spray chambers, therefore, it may be left out in those cases.).

**8.4. Internal Standard/Ionization Buffer**

Directions: Add approximately 1500 mL DI water to a 2000-mL volumetric flask. Add 24 g CsCl and swirl to mix. Add 10 mL of 1000 ppm Scandium standard solution (40 mL of 1000 ppm Beryllium standard solution may also be added as needed). The final concentration of the Internal Standards are 5.0 ppm Sc (20 ppm Be). Add 80 mL of HNO<sub>3</sub>, cap and invert multiple times to mix.

Optional: 1-2 drops of Triton X-100 may be added. (Triton X – 100 causes beading in some spray chambers, therefore, it may be left out in those cases.).

8.4.1. Expiration is **6 months** from preparation date.

8.5. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015.

8.6. The specific ingredient amounts can vary to accommodate the batch size.

8.7. Expiration Date: 12 months from date of preparation unless otherwise indicated.

8.8. Store at room temperature.

**9. ASSOCIATED DOCUMENTS**

- 9.1. ADMIN 023 Attachment M Instrument Sequence Generation for Fertilizer
- 9.2. ADMIN 023 Attachment O, Automatic Result Entry for Fertilizer
- 9.3. DATA 051, Data Evaluation for Feed, Fertilizer, and Pesticide Programs
- 9.4. DATA 051 Attachment A, Fertilizer Analysis Requirements
- 9.5. DATA 051 Attachment B, Fertilizer Laboratory Control Sample Acceptance Limits
- 9.6. DATA 051 Attachment C, Continuing Calibration Verification Acceptance Limits for AA and ICP
- 9.7. LABOP 140 Attachment I: Maintenance, Operation and Program of ICPOES
- 9.8. PREP 461, Sample Preparation for Available Phosphorus and Water-Soluble Potash in Fertilizer
- 9.9. QUALITY 112, Quality Control
- 9.10. SAFETY 002, Waste Management

**10. SPECIFIC PROCEDURES**

**10.1. Sample Handling**

10.1.1. Sample extracts containing particulates should be filtered or centrifuged before analysis.

10.1.2. Extracts must be stored refrigerated and are to be analyzed within **3 days**.

**10.2. Standards**

10.2.1. Primary standards and internal standards (1000ppm or 10,000ppm) should be prepared from neat material or purchased through certified vendors that are NIST traceable.

**Note:** Purchased standards preserved with acid should not be utilized since the acid will change the pH and matrix of the calibration standards.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 516

#### *Analysis of Available Phosphorous and Soluble Potash in Fertilizers using ICP-OES*

Version: 2.0

Effective Date: 10/16/2023

Page 5 of 8

10.2.2. Working standards are to be used for ICP calibration and unknown sample identification and quantification. They are made by diluting with a solution of citrate-EDTA (that matches the sample extracts) to obtain the desired concentrations.

10.2.3. Working Standards have expiration **2 months** when stored refrigerated.

10.2.4. Recommended Standard Levels:

| Standard Level        | Stock (ppm) | Stock (mL) | Citrate-EDTA Solution (mL) | Final Volume (mL) | K (ppm) | P (ppm) |
|-----------------------|-------------|------------|----------------------------|-------------------|---------|---------|
| Intermediate Standard | 10,000      | 100        | -                          | 1000              | 1,000   | 1,000   |
| 0                     | -           | -          | 200                        | 1000              | -       | -       |
| 1                     | 1,000       | 5          | 200                        | 1000              | 5       | 5       |
| 2                     | 1,000       | 10         | 200                        | 1000              | 10      | 10      |
| 3                     | 1,000       | 20         | 200                        | 1000              | 20      | 20      |
| 4                     | 10,000      | 5          | 200                        | 1000              | 50      | 50      |
| 5                     | 10,000      | 10         | 200                        | 1000              | 100     | 100     |
| 6                     | 10,000      | 20         | 200                        | 1000              | 200     | 200     |
| 7                     | 10,000      | 40         | 200                        | 1000              | 400     | 400     |
| 8                     | 10,000      | 60         | 200                        | 1000              | 600     | 600     |

#### 10.3. Instrument Startup

10.3.1. See LABOP 140 Attachment I: Maintenance, Operation and Program of ICPOES

10.3.2. Turn on the plasma and peristaltic pump.

10.3.2.1. Allow the system to stabilize for a minimum of 30 minutes.

#### 10.4. Recommended Instrument Parameters

10.4.1. The following parameters may vary for different ICP makes and models and therefore require minor modifications in order to optimize performance.

|                            | PE Optima 5300 | PE Optima 8300 |
|----------------------------|----------------|----------------|
| Plasma Flow                | 15 L/min       | 15 L/min       |
| Auxiliary Flow             | 0.2 L/min      | 1.5 L/min      |
| Nebulizer Flow             | 0.65 L/min     | 0.65 L/min     |
| RF Power                   | 1450 watts     | 1450 watts     |
| Peristaltic Pump Flow Rate | 0.8 mL/min     | 0.44 mL/min    |
| Carrier Tubing             | Black/Black    | Black/Black    |
| Internal Std. Tubing       | Orange/Green   | Orange/Green   |
| Waste Tubing               | Red/Red        | Red/Red        |
| Wash                       | 0.8 mL/min     | 0.44 mL/min    |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 516

#### *Analysis of Available Phosphorous and Soluble Potash in Fertilizers using ICP-OES*

Version: 2.0

Effective Date: 10/16/2023

Page 6 of 8

|                 |          |          |
|-----------------|----------|----------|
| Spray Chamber   | Cyclonic | Cyclonic |
| Loop Size       | 1.0 mL   | 1.5 mL   |
| Replicate Reads | 3        | 3        |

#### 10.4.2. Recommended Wavelengths and Plasma View:

| Element             | Wavelength nm | Plasma View  |
|---------------------|---------------|--------------|
| K                   | 766.490*      | Radial       |
| K                   | 404.721       | Axial        |
| P                   | 177.434       | Axial        |
| P                   | 178.222       | Axial        |
| P                   | 213.617*      | Axial        |
| P                   | 214.914       | Axial        |
| Sc                  | 361.383*      | Axial/Radial |
| Be                  | 313.107*      | Axial/Radial |
| *Primary Wavelength |               |              |

#### 10.5. Sample Analysis

10.5.1. Use the Sample Information File (SIF) to create the sequence.

10.5.2. A calibration blank is injected with the calibration standards to correct for background.

10.5.3. After the initial calibration, there must be a continuing calibration check standard analyzed after every 10 samples and at the end of the run.

10.5.4. The Laboratory Control Samples (LCS) must be analyzed with their associated samples, not grouped together at the beginning of the sequence.

10.5.5. The sequence requirements are found in, DATA 051 Attachment A, Fertilizer Analysis Requirements.

10.5.6. The sequence must end with a CCV. The ending CCV should be at a different concentration than the mid-sequence CCVs.

10.5.7. Samples require the proper laboratory information management system (LIMS) codes to allow required uploading. Refer to ADMIN 023 Attachment M Instrument Sequence Generation for Fertilizer for additional information.

10.5.7.1. Example: AA99999-0AP

10.5.7.1.1. 0 corresponds to the sample weight (0, 1, 2.)

10.5.7.1.2. A corresponds to the different injections due to dilutions (A=1st injection, B=2nd injection, C=3rd injection.)

10.5.7.1.3. P corresponds to analysis type (P= P&K Method)

10.5.8. Dilution factors must be entered into the instrument sequence.

#### 10.6. Export Data

10.6.1. Data Manager exports data to the network drive TLHAESCHROME.

10.6.2. Upload data from the network drive to the LIMS.

#### 10.7. Instrument Shutdown

10.7.1. See LABOP 140 Attachment I: Maintenance, Operation and Program of ICPOES

#### 10.8. Quality Control

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 516**

*Analysis of Available Phosphorous and Soluble Potash in Fertilizers  
using ICP-OES*

Version: 2.0

Effective Date: 10/16/2023

Page 7 of 8

10.8.1. Initial Calibration:

10.8.1.1. The calibration is based on plotting the ratio of the analyte response to the internal standard response as a function of the analyte concentration of the standards.

10.8.1.2. The correlation coefficient of the calibration curve R value must be  $\geq 0.999$  or the  $R^2$  value must be  $\geq 0.998$ .

10.8.1.3. The internal standard RSD value in the calibration blank must be less than 2%. If internal standard RSD is greater than 2%, reanalyze the calibration blank until 2% is achieved.

10.8.1.4. Refer to DATA 051, Data Evaluation for Chemical Analysis; for additional calibration procedures.

10.8.2. Continuing Calibration Check Standard:

10.8.2.1. Once the initial calibration has passed, a continuing calibration check standard must be analyzed after every 10 samples and at the end of the sequence. CCVs must be run for each separate curve.

10.8.2.2. CCV acceptance criteria for P and K are  $\pm 6\%$ . Addition information may be found in DATA 051 Attachment C, Fertilizer Continuing Calibration Verification Acceptance Limits.

10.8.2.3. The CCV may consist of any concentration level used in the initial calibration. Levels near the middle of the curve are recommended.

10.8.2.4. If a CCV response is not within the limits, all samples analyzed before and after it must be reanalyzed.

10.8.2.5. If upon reanalysis of the CCV, the response is still not within the above acceptance criteria, then a new initial calibration must be performed. It is recommended that routine instrument maintenance be performed prior to the new curve analysis.

10.8.3. When using the suggested sample weight and volume in PREP 461, any sample reading exceeding the calibration range would exceed the sample's theoretical maximum threshold and must be reanalyzed or re-prepped.

10.8.4. Internal Standard criteria for samples is 80% - 120% recovery.

10.8.5. %RSD for the sample replicates shall not exceed 5%.

10.8.6. Laboratory Control Sample

10.8.6.1. This standard is weighed with a set of samples and is processed in the same manner as the samples. The percent found must be within the limits listed in DATA 051 Attachment B, Fertilizer Laboratory Control Sample Acceptance Limits.

10.8.7. Instrument Blank

10.8.7.1. This blank matches the citrate-EDTA extract matrix (standard 0) and is used to ensure the instrument is free from contamination. The response of the instrument blank must be below the Method Detection Limit (MDL). Refer to QUALITY 112, Quality Control for additional information.

10.9. Data Reduction and Reporting

10.9.1. Calculations for  $P_2O_5$ :

$$P_2O_5\% = \frac{P(\mu\text{g/mL}) \times 500(\text{mL}) \times 141.95}{\text{Weight(g)} \times 30.97 \times 2 \times 10,000}$$

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 516

### *Analysis of Available Phosphorous and Soluble Potash in Fertilizers using ICP-OES*

Version: 2.0

Effective Date: 10/16/2023

Page 8 of 8

Where:

$$P_2O_5\% = \frac{P(\mu\text{g/mL}) \times \text{Final Vol. (mL)} \times P_2O_5(\text{MW})}{\text{Weight (g)} \times P(\text{MW}) \times (\text{mole ratio of } P_2O_5 \text{ to } P) \times (\text{ppm to } \%)}$$

10.9.2. Calculations for K<sub>2</sub>O:

$$K_2O\% = \frac{K(\mu\text{g/mL}) \times 500(\text{mL}) \times 94.2}{\text{Weight (g)} \times 39.09 \times 2 \times 10,000}$$

Where:

$$K_2O\% = \frac{K(\mu\text{g/mL}) \times \text{Final Vol. (mL)} \times K_2O(\text{MW})}{\text{Weight (g)} \times P(\text{MW}) \times (\text{mole ratio of } K_2O \text{ to } K) \times (\text{ppm to } \%)}$$

10.9.3. The Laboratory Control Sample must be evaluated first and if it fails, the associated samples in the batch are not uploaded.

10.9.4. Data for all passing batches is uploaded into LIMS and evaluated. Refer to DATA 051, Data Evaluation for Feed, Fertilizer, and Pesticide Programs and ADMIN 023 Attachment O, Automatic Result Entry for Fertilizer for additional information.

## 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material         | Disposal Mechanism                                            |
|------------------|---------------------------------------------------------------|
| Aqueous extracts | Disposed of by dumping down the sink and flushing with water  |
| Instrument waste | Disposed of by dumping down the sink and flushing with water. |

## 12. REFERENCES

12.1. *Official Methods of Analysis of AOAC INTERNATIONAL* (2019) 21st Ed., AOAC INTERNATIONAL, Rockville, MD, USA, **Method 2015.18**, [www.eoma.aoac.org](http://www.eoma.aoac.org), **External Document 00135**.

## 13. VERSION HISTORY

| Version | Date       | Description                                                                                                                     | Author        | Editor        |
|---------|------------|---------------------------------------------------------------------------------------------------------------------------------|---------------|---------------|
| 1.0     | 09.17.2021 | Original Version                                                                                                                | William Meeks |               |
| 2.0     | 10.16.2023 | Added external document number external reference.<br>Section 5.3 added Multicomponent Spectral Fitting (MSF) procedure to SOP. |               | William Meeks |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 520

#### *Sample Preparation and Analysis of Combined Sulfur* *[FM805]*

Version: 2.0

Effective Date: 01/30/2024

Page 1 of 5

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for preparation and analysis of combined sulfur by gravimetric method in fertilizer samples.

## 2. SCOPE

2.1. The procedure delineated below is applicable to the proper protocol to be followed by all staff in analyzing combined sulfur in fertilizer samples.

2.2. The combined sulfur result is found through a gravimetric method in which the sulfate is precipitated as barium sulfate by the addition of 10% barium chloride in excess. The final precipitate is expressed as percent of combined sulfur (%S).

2.3. The manufacturer's tag should state or specify "combined sulfur."

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required. Additional PPE is included such as oven gloves.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Weighing and LIMS Entry
- 4.2. 10.2 Digestion for Solid Samples
- 4.3. 10.3 Digestion for Liquid Samples
- 4.4. 10.4 Barium Sulfate Precipitation for both Solid and Liquid samples
- 4.5. 10.5 Calculation
- 4.6. 10.6 Data Evaluation

## 5. INTERFERENCE

5.1. N/A

## 6. APPARATUS and MATERIAL

- 6.1. Balance, accuracy to 0.001g
- 6.2. Beaker, 250mL
- 6.3. Bottle, 20L Nalgene plastic or equivalent
- 6.4. Filter paper, Whatman (125mm) 2V folded filter paper or equivalent
- 6.5. Fume hood
- 6.6. Glass fiber pads 2.4cm Whatman 934-AH or equivalent
- 6.7. Glass funnel
- 6.8. Gooch crucibles (40mL) or equivalent
- 6.9. Hotplate
- 6.10. Muffle furnace
- 6.11. Vacuum trap apparatus
- 6.12. Volumetric flasks, 500mL and 1L (class "A")

## 7. REAGENTS and MEDIA

- 7.1. Ammonium Sulfate, certified American Chemical Society (A.C.S.) grade or equivalent
- 7.2. Barium chloride,  $2\text{H}_2\text{O}$ , certified American Chemical Society (A.C.S.) grade or equivalent
- 7.3. Calcium thiosulfate, pure, 30-60% solution in water, or equivalent
- 7.4. Hydrochloric acid, certified A.C.S. grade or equivalent
- 7.5. Hydrogen peroxide, 30% certified A.C.S. grade or equivalent
- 7.6. Sodium hydroxide, NF/FCC Pellets, or equivalent

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 520**

*Sample Preparation and Analysis of Combined Sulfur*  
*[FM805]*

Version: 2.0

Effective Date: 01/30/2024

Page 2 of 5

7.7. Standard Reference Material (SRM) 695 or equivalent

7.8. Water, Deionized (DI) water

## **8. REAGENT PREPARATION GUIDELINES**

8.1. 10% Barium Chloride

8.1.1. Add 1173g of barium chloride ( $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$ ) to approximately 8L of deionized water in a 10L Nalgene plastic bottle. Bring to a final volume of 10L with deionized water and mix.

8.2. 50% Sodium Hydroxide

8.2.1. Slowly add 500g of sodium hydroxide while mixing to approximately 450mL of deionized water. Bring to a final volume of 1L with deionized water and mix.

8.3. 1:1 Hydrochloric Acid

8.3.1 Add 50mL of concentrated hydrochloric acid to approximately 40 mL of deionized water in a 100mL flask and mix. Bring to a final volume of 100mL with deionized water and mix.

8.4. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015

8.5. The specific ingredient amounts can vary to accommodate the batch size.

8.6. Expiration Date: 12 months from date of preparation

8.7. Store at room temperature

## **9. ASSOCIATED DOCUMENTS**

9.1. ADMIN 023, Attachment P, Manual Result Entry

9.2. Data 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs

9.3. DATA 051 Attachment B, Fertilizer Laboratory Control Sample Acceptable Limits

9.4. F00015, Reagent Preparation

9.5. L003, Reagent

9.6. METHOD 196, Determination of Percent Moisture for Fertilizer

9.7. PREP 453, Weighing of Fertilizer Sample

## **10. PROCEDURE**

10.1. Weighing and LIMS Entry

10.1.1. Create a sample list per the Laboratory Information Management System (LIMS).

10.1.2. Locate samples and gather the amount of gooch crucibles and glass fiber pads that is needed to correspond with samples. Add a glass fiber pad to the crucible and dry the gooch crucibles with glass fiber pads at 250°C. Allow them to cool and take them to the weighing room.

10.1.3. Set the crucibles to correspond with the number you assigned to the sample order.

10.1.4. Open the Laboratory Information Management System (LIMS) and enter or scan sample ID numbers in to the FERT\_SULFUR\_COMBINED template. Further LIMS instruction can be found in ADMIN 023 Attachment P, Manual Result Entry.

10.1.5. Gypsum samples being analyzed for combined sulfur must be dried according to METHOD 196.

10.1.6. Using the suggested weight below, weigh directly into a 250mL beaker and enter the weight into the 805\_WT code. Refer to PREP 453, Weighing of Fertilizer Sample, for further instruction.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 520**

*Sample Preparation and Analysis of Combined Sulfur*  
*[FM805]*

Version: 2.0

Effective Date: 01/30/2024

Page 3 of 5

| Guarantee | Weight |
|-----------|--------|
| <2%       | 1.5g   |
| 2-10%     | 1g     |
| 10-12%    | 0.5g   |
| >12%      | 0.4g   |

10.1.6.1. When weighing liquid samples, weigh double the suggested weight.

10.1.7. For solid samples weigh 0.40g of the ammonium sulfate laboratory control sample directly into a 250mL beaker and enter the weight into the 805\_WT code. For liquid samples included in the set also weigh 1.00g of calcium thiosulfate for the laboratory control sample.

10.1.8. Tare the scale and weigh the crucible with the glass pad and record the weight in the LIMS code 805\_TARE\_WT.

**10.2. Digestion for Solid Samples**

10.2.1. Digestion for samples must be completed in a hood.

10.2.2. Add 150mL of deionized water then 10mL of concentrated hydrochloric acid (HCl) to beaker.

10.2.3. Gently boil the sample on a hotplate for approximately 20 minutes.

10.2.4. Filter through Whatman 2V filter paper into another 250mL beaker and wash with hot DI water to obtain a volume of 150mL. Discard filter paper.

**10.3. Digestion for Liquid Samples**

10.3.1. Digestion for samples must be performed in a hood.

10.3.2. Add 50mL of deionized (DI) water, followed by 2mL of 50% sodium hydroxide (NaOH) and then 2mL of 30% hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>).

10.3.3. Cover with watch glass and reflux the sample on a hotplate for approximately 60 minutes.

10.3.4. Slowly add 5mL of 30% H<sub>2</sub>O<sub>2</sub> in 1mL increments into sample.

**Note:** Caution, add 30% H<sub>2</sub>O<sub>2</sub> very slowly to beaker. Reaction is extremely exothermic!

10.3.5. Rinse watch glass into sample and dilute to approximately 175mL with hot DI water.

10.3.6. Add 10mL 1:1 HCl:DI water, boil 5 minutes and filter.

**10.4. Barium Sulfate Precipitation for both Solid and Liquid samples**

10.4.1. Bring filtrate back to a slow boil.

10.4.2. Carefully add 2mL of 10% barium chloride solution, in small increments to beaker to avoid a violent reaction.

10.4.3. Add an additional 15mL of 10% barium chloride solution to beaker, adjust hotplate so that solution is just below boiling for approximately 60 minutes.

10.4.4. A white precipitate (barium sulfate) will form in the beaker after a few minutes

10.4.5. Allow the samples to cool to room temperature and let stand overnight.

10.4.6. Set up the gooch crucible with pad on vacuum apparatus.

10.4.7. Seal pad with water and vacuum filter sample through crucible with hot DI water.

**CAUTION:** Pad seal must not be allowed to leak!

10.4.8. Rinse the beaker several times into crucible with hot DI water.

10.4.9. If necessary, use a glass rod to break up precipitant while rinsing.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 520**

*Sample Preparation and Analysis of Combined Sulfur*  
*[FM805]*

Version: 2.0

Effective Date: 01/30/2024

Page 4 of 5

10.4.10. Remove the crucible from vacuum apparatus.

10.4.11. Dry crucible at 250°C for one hour.

10.4.12. After drying, allow crucible to cool to room temperature.

10.4.13. Weigh the crucible and record the weight in the LIMS code 805\_FINAL\_WT.

10.5. Calculation:

$$\frac{[Final\ Weight_{BaSO_4}(g) - Tare\ Weight_{Crucible\ \&\ Filter}(g)](0.1374_{S/BaSO_4})}{Sample\ Weight\ (g)} \times 100 = \% \text{ Combined Sulfur}$$

10.6. Data Evaluation

10.6.1. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs, to evaluate the data against the established criteria. Evaluations, that supersede those of DATA 051, are addressed in this document.

10.6.2. Evaluate the Lab Control sample results.

10.6.3. To determine if it is within acceptable limits refer to DATA 051 Attachment B, Fertilizer Laboratory Control Sample Acceptable Limits.

10.6.4. If there are multiple sample results, use the result closest to the guarantee if all the results are legal. Use the highest result if all the results are deficient. Refer to DATA 051.

10.6.5. Make sure that the found result is always completed.

## 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material                   | Disposal Mechanism     |
|----------------------------|------------------------|
| Filters and Barium Sulfate | Disposed of in Trash   |
| Filtrate                   | Disposed of down drain |

## 12. REFERENCES

- 12.1. E00154, *Sulfur in Fertilizers, Gravimetric Method*. Method 980.02, AOAC 15<sup>th</sup> Edition
- 12.2. E00155, Combined Sulfur Analysis, FM 805
- 12.3. E00156, Analysis of Combined Sulfur, FSFL-SOP 520

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 520**

*Sample Preparation and Analysis of Combined Sulfur*  
*[FM805]*

Version: 2.0

Effective Date: 01/30/2024

Page 5 of 5

**13. VERSION HISTORY**

| Version | Date       | Description                                                                                                                                                                                                                                                                                                                                                                                                                                           | Author            | Editor                                                 |
|---------|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|--------------------------------------------------------|
| 1.0     | 04.14.2021 | Original Version                                                                                                                                                                                                                                                                                                                                                                                                                                      | William Meeks, Jr |                                                        |
| 2.0     | 01/30/2024 | Sections 6.3 and 6.8 added 'or equivalent'.<br>Section 7.2 added ·2H <sub>2</sub> O to Barium Chloride.<br>Sections 8.2 and 8.3 adjusted volume of water.<br>Section 10.1.7 clarified which laboratory control sample to use for liquid and solid. samples.<br>Section 10.2.2 and 10.3.2 clarified language in step.<br>Section 10.4.3 clarified time for boiling.<br>Added weight table<br>Added external document numbers<br>Removed signature page |                   | Taleya Williams<br>Chris Gismondi<br>William Meeks, Jr |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 521

#### *Sample Preparation and Analysis of Free Sulfur [FM 806]*

Version: 1.0

Effective Date: 05/17/2022

Page 1 of 6

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the determination of free (elemental) sulfur in fertilizer samples by gravimetric method.

## 2. SCOPE

2.1. The procedure delineated below is applicable to fertilizer samples that contain free sulfur component. Free sulfur is extracted by Soxhlet using carbon disulfide (CS<sub>2</sub>.) After solvent removal, the percent sulfur is determined after drying by weight difference.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

3.2. **Inhalation Hazard:** Carbon disulfide and carbon tetrachloride pose dangerous inhalation hazards and must be used with proper ventilation such as within a fume hood.

3.3. **Flammable:** Due to carbon disulfide's high flammability, care must be taken while working with or around heat. Do NOT store carbon disulfide next to heat sources.

3.4. **Carcinogenicity:** Carbon tetrachloride is classified as a suspected carcinogen. The use of PPE such as lab coat, gloves, and safety glasses are mandatory. Carbon tetrachloride shall only be used within a fume hood. Any spills must be cleaned immediately.

3.5. Ensure that a closed system is not created while setting up the distillation apparatus. Pressure build-up from heating organic solvents, such as carbon disulfide, can potentially result in an explosion.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Weighing and LIMS Entry
- 4.2. 10.2 Sample Washing to Remove Urea
- 4.3. 10.3 Drying samples
- 4.4. 10.4 Soxhlet Extraction
- 4.5. 10.5 Solvent Removal
- 4.6. 10.6 Final Weighing
- 4.7. 10.7 Calculation
- 4.8. 10.8 Data Evaluation
- 4.9. 10.9 Quality Control Sample

## 5. INTERFERENCE

5.1. Sulfur-coated urea tends to give elevated results in free sulfur analysis. It leaves a brownish, oily residue in the boiling flask after distillation. A hot water rinse of the sample is used to remove urea.

5.2. Moisture will interfere with the extraction process.

5.3. The use of stopcock grease may be dissolved by the carbon disulfide solution and may lead to positive bias interferences.

## 6. APPARATUS and MATERIAL

- 6.1. Analytical balance, minimum capacity 0.0001g
- 6.2. 125mL boiling flask with 24/40 glass joint
- 6.3. 30 mm Soxhlet extraction tube with 24/40 glass joint
- 6.4. 30 mm extraction condenser with 24/40 glass joint
- 6.5. Cellulose extraction thimbles (25 mm X 80 mm), Whatman Cat. # 2800258 or equivalent

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 521

#### *Sample Preparation and Analysis of Free Sulfur [FM 806]*

Version: 1.0

Effective Date: 05/17/2022

Page 2 of 6

- 6.6. Drying oven
- 6.7. Elbow joint with 24/40 glass joint
- 6.8. Fume hood
- 6.9. Water bath with Thermolyne hot plate or equivalent

## 7. REAGENTS and MEDIA

- 7.1. Water, Deionized
- 7.2. Carbon Disulfide CS<sub>2</sub>, Reagent Grade A.C.S.
- 7.3. Carbon Tetrachloride CCl<sub>4</sub>, 99% Anhydrous

## 8. REAGENT PREPARATION GUIDELINES

- 8.1. Carbon Disulfide Extraction Solution (CS<sub>2</sub>:CCl<sub>4</sub>, 60:40)
  - 8.1.1. **Under a hood** and while wearing a respirator, add 400mL of carbon tetrachloride to 600mL of carbon disulfide. Mix well.
- 8.2. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015
- 8.3. The specific ingredient amounts can vary to accommodate the batch size.
- 8.4. Expiration Date: 12 months from date of preparation
- 8.5. Store at room temperature

## 9. ASSOCIATED DOCUMENTS

- 9.1. ADMIN 023 Attachment C, LIMS Sample Login for Multi Sample
- 9.2. DATA 051
- 9.3. F00015, Reagent Preparation
- 9.4. L003, Reagent
- 9.5. QUALITY 112
- 9.6. PREP 453, Weighing of Fertilizer Sample

## 10. SPECIFIC PROCEDURES

### 10.1. Weighing and LIMS Entry

#### 10.1.1. Preparation

- 10.1.1.1. Create a sample list per backlog
- 10.1.1.2. Locate samples and gather all glassware needed. Check 125mL-boiling flasks and Soxhlet extraction tube for small damages and cracks.
- 10.1.1.3. Use LIMS, according to ADMIN 023 Attachment C, LIMS Sample Login for Multi Sample, to create a sample number for a lab standard.
- 10.1.1.4. Label all glassware to correspond to sample numbers. Set the thimbles to correspond with the sample numbers.
- 10.1.1.5. Reference PREP 453, Weighing of Fertilizer Sample for further instruction on preparation for weighing.

#### 10.1.2. Weighing

- 10.1.2.1. Open the LIMS result entry template Fert\_Sulfur\_Free.
- 10.1.2.2. Place a cellulose extraction thimble in a holder (E.g., beaker or gooch crucible) and tare on the balance.
- 10.1.2.3. In LIMS each sample has a suggested weight.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 521**

*Sample Preparation and Analysis of Free Sulfur [FM 806]*

Version: 1.0

Effective Date: 05/17/2022

Page 3 of 6

10.1.2.4. Weigh the sample directly into the cellulose extraction thimble and record as the sample weight (SW) in LIMS. A holder (E.g., beaker, gooch crucible, etc.) may be used to assist in weighing.

10.1.2.5. Weigh the labeled 125mL-boiling flasks. Record this weight as the tare weight (TW) in LIMS.

10.1.2.6. Save and exit LIMS.

10.1.2.7. Return samples to sample storage area.

**10.2. Sample Washing to Remove Urea**

10.2.1. Place a gooch crucible in the vacuum apparatus.

10.2.2. Place the thimble with the sample in the gooch crucible.

10.2.3. Rinse the sample with approximately 20mL of hot DI water using vacuum. This should be done 2 more times for a total of 3 times.

Note: DI water should be brought to a boil and used immediately.

10.2.4. Remove the cellulous thimbles from the vacuum and place them in their respective holders.

**10.3. Drying samples**

10.3.1. Place all samples in the drying oven for 1- 2 hours at 105°C, or until cellulose thimble is dry.

**Note:** Moisture left in the cellulose thimble will interfere with extraction process.

10.3.2. Remove samples from oven and cool

**10.4. Soxhlet Extraction**

10.4.1. **The following steps of the Soxhlet extraction must take place in a fume hood.**

10.4.2. See METHOD 521 attachment A for the Soxhlet assembly setup. Ensure that the hood is on and in working order.

10.4.3. Add approximately 65mL of carbon disulfide extraction solution to the 125mL-boiling flask.

10.4.4. Attach the Soxhlet extraction tube to the 125mL-boiling flask and place the thimble with sample into the Soxhlet extraction tube.

10.4.5. Add approximately 15mL of carbon disulfide extraction solution to the Soxhlet extraction tube.

10.4.6. Attach the condenser to the Soxhlet extraction tube. Ensure that water is flowing through the condenser.

10.4.7. Turn the extraction heater or hot plate on high and reflux for 2 hours.

10.4.8. After refluxing is complete turn off heater/hot block.

10.4.9. **Allow flask and heating block to cool to room temperature.**

10.4.10. Remove the condenser assembly.

10.4.11. Use forceps to lift thimble and drain carbon disulfide-extraction solution. Allow thimble to drip dry in the Soxhlet extraction tube.

10.4.12. After draining is complete, remove cellulose thimble from Soxhlet extraction tube. Set it aside under the hood to allow CS<sub>2</sub> to completely evaporate before disposal.

10.4.13. Separate 125mL-boiling flask from the Soxhlet extraction tube.

10.4.14. Pour remaining carbon disulfide-extraction solution from Soxhlet extraction tube into the 125mL-boiling flask.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 521

#### *Sample Preparation and Analysis of Free Sulfur [FM 806]*

Version: 1.0

Effective Date: 05/17/2022

Page 4 of 6

#### 10.5. Solvent Removal

10.5.1. Attach a 125mL-boiling flask to a glass elbow.

10.5.2. Place 125mL-boiling flask into a hot water bath and attach the other end of the elbow to the cooling condenser.

10.5.3. Distill and collect carbon disulfide extraction solution for disposal.

10.5.4. Remove the 125mL-boiling flask from the hot water bath and disconnect the glass elbow.

10.5.5. Place flask in drying oven for 1 hour at approximately 105°C or until dry and cool to room temperature.

10.5.6. The yellowish residue in the flask is the free (elemental) sulfur.

#### 10.6. Final Weighing

10.6.1. Open LIMS result entry template "Fert\_Free\_Sulfur."

10.6.2. Weigh the 125mL-boiling flask containing the residue (sulfur). This weight will be entered as the final weight (FW).

#### 10.7. Calculation

10.7.1. Calculate the Free Sulfur result using the formula below:

$$\%S \text{ (Free Sulfur)} = (FW - TW) / (SW) * 100$$

S = Free Sulfur

FW = Final weight of 125mL-boiling flask with sulfur residue

TW = Tare weight of 125mL-boiling flask

SW = Sample weight

#### 10.8. Data Evaluation

10.8.1. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs, to evaluate the data against the established criteria. Evaluations, that supersede those of DATA 051, are addressed in this document.

#### 10.9. Quality Control Sample

10.9.1. Refer to QUALITY 112, Quality Control, to determine the outcome of Quality Control (QC) samples. Guidelines, that supersede those of QUALITY 112, are addressed in this document.

## 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material                                 | Disposal Mechanism                                                                   |
|------------------------------------------|--------------------------------------------------------------------------------------|
| CS <sub>2</sub> :CCl <sub>4</sub>        | Dispose remainder in chlorinated waste container.                                    |
| Extracted Fertilizer, Cellulose thimbles | Allow any residual solvent to evaporate in hood or vented oven and dispose in trash. |
| Gloves & kimwipes                        | Allow any residual solvent to evaporate in hood and dispose in trash.                |
| Sulfur Waste                             | Dispose of in trash                                                                  |

## 12. REFERENCES

12.1. Sulfur in Fertilizers, Gravimetric Method. Method 980.02, AOAC 15<sup>th</sup> Edition.

12.2. FM-806, Secondary/Micronutrient, Sulfur – Free, 01/02/03 v.2

12.3. FSFL-SOP 521 Determination of Free (Elemental) Sulfur in Fertilizers by Gravimetric Method. 2006/01/17 v.1.0

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 521**

*Sample Preparation and Analysis of Free Sulfur [FM 806]*

Version: 1.0

Effective Date: 05/17/2022

Page 5 of 6

**13. VERSION HISTORY**

| <b>Version</b> | <b>Date</b> | <b>Description</b> | <b>Author</b>     | <b>Editor</b> |
|----------------|-------------|--------------------|-------------------|---------------|
| 1.0            | 04.11.2019  | Original Version   | William Meeks, Jr |               |

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 521**

*Sample Preparation and Analysis of Free Sulfur [FM 806]*

Version: 1.0

Effective Date: 05/17/2022

Page 6 of 6

**14. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer  
Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Wet Chemistry and Preparation Manager  
Date\_\_\_\_\_

\_\_\_\_\_  
William Meeks, Jr  
Author  
Date\_\_\_\_\_

# **ATTACHMENT**

## **BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

### METHOD 521 Attachment A Soxhlet Apparatus Assembly

Version: 1.0

Effective Date: 5/17/2022

Page 1 of 3

The following is a recommended setup of Soxhlet and Distillation apparatus. Variations of this setup are allowed in order to accommodate additional samples or variations in glassware or heating apparatus. All variations must be approved by a supervisor. Regardless of the configuration, care must be taken not to create a closed system.



#### Notes:

- The apparatus must be setup in a hood. Ensure the hood is in proper working condition.
- Do not create a closed system. Doing so may result in an explosion.
- Ensure the water is on and flowing from the bottom of the condenser to the top.
- Solvents such as Carbon Disulfide are extremely flammable. Ensure heating source is cool before removing/disassembling apparatus.
- Ensure that the receiver flask is large enough to accommodate total amount of solvent being distilled.
- Aluminum foil may be used to cover the water bath. This will help insulate the samples during distillation and maintain temperature for complete distillation.
- Store waste such as used cellulose thimbles, gloves, and Kimwipes in a box under the hood until residual solvent has evaporated.


METHOD 521 Attachment A  
**Soxhlet Apparatus Assembly**

Version: 1.0

Effective Date: 5/17/2022

Page 2 of 3

**Soxhlet Apparatus**



30 mm extraction condenser, Allihn, 24/40

Water flows into the bottom of the condenser and out of the top of the condenser.

30 mm Soxhlet extraction tube, 24/40 glass joint

Green Joint Clamp  
(Use on all 24/40 joints)

125 mL boiling flask, 24/40 glass joint

**ATTACHMENT**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 521 Attachment A  
**Soxhlet Apparatus Assembly**

Version: 1.0

Effective Date: 5/17/2022

Page 3 of 3

**Distillation Apparatus**

Distilling Adapter, 105°,  
24/40 glass joint, with  
Vacuum Take-Off

**(WARNING:** Vacuum  
Take-Off port must not  
be covered or closed.  
Doing so will make a  
closed system.)


Side Arm Adapter, 75°, 24/40  
glass joint

30 mm extraction condenser, Liebig,  
24/40 glass joint

Glass plug, 24/40 glass joint

Green Joint Clamp  
(Use on all 24/40 joints)

Water flows into the bottom of  
the condenser and out of the



Receiver flask, 24/40 glass joint  
(Large enough to receive total  
volume distilled.)

Glass Cow Distillation Receiver,  
24/40 glass joint

125 mL boiling flask, 24/40 glass joint

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 523

*Analysis of Soluble and Chelate Metals by AA- [FM830, FM831, FM832, FM833, FM840, FM841]*

Version: 1.0

Effective Date: 05/10/2023

Page 1 of 6

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for determination of soluble and chelated metals micronutrients: Iron, Manganese, Copper, and Zinc by Atomic Absorption Spectroscopy (AAS).

## 2. SCOPE

2.1. The procedure delineated below is applicable to Atomic Absorption Spectroscopy determination of soluble and chelated micronutrients: Iron, Manganese, Copper, and Zinc in fertilizer samples analyzed by the Bureau of Agricultural Environmental Laboratories. Quantitative determination of the amount of analyte present in the sample is achieved by measuring signal intensity at the wavelength of maximum absorption, which is specific to the element analyzed.

## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

3.2. Be careful with AAS Burner when handling it, it must be at room temperature.

3.3. Exhaust from combustion fumes and vapors from the AA flame housing can produce toxic vapors and may pose a health hazard if not adequately vented to the outside.

## 4. OUTLINE of PROCEDURE

- 4.1. 10.1 Sample Handling and Storage
- 4.2. 10.2 Standards
- 4.3. 10.3 Instrument Startup
- 4.4. 10.4 Sample Analysis
- 4.5. 10.5 Data Evaluation
- 4.6. 10.6 Quality Control Sample

## 5. INTERFERENCE

- 5.1. Polyammonium phosphate
- 5.2. Super phosphate

## 6. APPARATUS and MATERIAL

- 6.1. Perkin-Elmer PinAAcle 500 Atomic Absorption Spectrometer or equivalent
- 6.2. SC2-DX Autosampler
- 6.3. AA prep 2 Peristaltic Pump
- 6.4. Workstation Computer system capable of operating integration software which interfaces with the Atomic Absorption Spectrometer and acquires and stores all the required data.
- 6.5. Black-Black tubing (0.76mm ID, 0.86mm Wall, 72mm between Bridges Black\Black.
- 6.6. Volumetric flask Class A
- 6.7. Glass Pipettes Class A
- 6.8. Syringe with Luer lock Inlet at least for 20ml
- 6.9. Syringe Filter, 25mm 0.45µm aqueous base
- 6.10. 1000 ml Volumetric Flask, Class A

## 7 REAGENTS and MEDIA

- 7.1. Deionized Water (Type 1)
- 7.2. Isopropanol or Methanol
- 7.3. Nitric Acid, 69%

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 523

*Analysis of Soluble and Chelate Metals by AA- [FM830, FM831, FM832, FM833, FM840, FM841]*

Version: 1.0

Effective Date: 05/10/2023

Page 2 of 6

7.4 Triton X-100

7.5 Hydrochloric Acid

## 8 REAGENT PREPARATION GUIDELINES

### 8.1. 5% Nitric Acid

8.2.1 Directions: Add approximately 500 ml of DI Water (Type 2) in 1 Liter volumetric flask (class A). Add 7.4 ml of Nitric Acid (69%) concentrated. Gently swirl the volumetric. Bring up to volume with DI Water and mix well.

8.2. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015

8.3. The specific ingredient amounts can vary to accommodate the batch size.

8.4. Expiration Date: 12 months from date of preparation unless otherwise specified.

8.5. Store at room temperature

## 9 ASSOCIATED DOCUMENTS

9.1. ADMIN 023, Using the Laboratory Information Management System

9.2. ADMIN 023 Attachment M, Instrument Sequence Generation for Fertilizer

9.3. ADMIN 023 Attachment O, Automatic Result Entry for Fertilizer

9.4. DATA 050, Data Control

9.5. DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs

9.6. DATA 051 Attachment C, Fertilizer Continuing Calibration Verification (CCV)

9.7. DATA 051 Attachment A, Fertilizer Analysis Requirements

9.8. LABOP 140 Maintenance Attachment k, Operation, and Program of Atomic Absorption Spectroscopy by Flame

9.9. QUALITY 112, Quality Control

## 10 SPECIFIC PROCEDURES

### 10.1. Sample Handling and Storage

10.1.1. All extracts may be filtered. Filtered extracts have no specified hold time.

10.1.2. Sample extracts containing particulates should be filtered before analysis. Use a Syringe Filter, 25mm 0.45µm aqueous base.

10.1.3. Unfiltered extracts should be prepared and analyzed according to the table below:

| Micronutrient  | Unfiltered Extract Hold Times                                        |
|----------------|----------------------------------------------------------------------|
| Copper (Cu)    | Extracts can be held up to 24 hours                                  |
| Iron (Fe)      | not specified                                                        |
| Zinc (Zn)      | Extracts are to be analyzed within an hour from when buffer is added |
| Manganese (Mn) | Extracts must be analyzed on the same day as extraction              |

### 10.2. Standards

10.2.1. Primary standards (Iron, Manganese, Copper and Zinc) 1,000 ppm should be prepared from neat material or purchased from certified vendors that are NIST traceable.

10.2.2. Working standards are to be used for AAS calibration. They are made by diluting into DI Water (Type 1) to obtain the following concentrations.

#### AAS Standard Mix 1 (0.5ppm)

Add to a 1000 ml volumetric flask 500 ml of DI Water (Type 1) approximately. Pipet 0.5 ml of Iron, Manganese, Copper and Zinc of the primary standard and 20.0 ml of HCl. Bring up to volume and mix well.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 523

*Analysis of Soluble and Chelate Metals by AA- [FM830, FM831, FM832, FM833, FM840, FM841]*

Version: 1.0

Effective Date: 05/10/2023

Page 3 of 6

AAS Standard Mix 2 (5.0ppm Cu, Fe, Mn; 1.0ppm Zn)

Add to a 1000 ml volumetric flask 500 ml of DI Water (Type 1) approximately. Pipet 5.0 ml of Iron, Manganese and Copper. Pipet 1.0 ml of Zinc of the primary standard and 20.0 ml of HCl. Bring up to volume and mix well.

AAS Standard Mix 3 (8.0ppm Cu, Fe, Mn; 2.0ppm Zn)

Add to a 1000 ml volumetric flask 500 ml of DI Water (Type 1) approximately. Pipet 8.0 ml of Iron, Manganese and Copper. Pipet 2.0 ml of Zinc of the primary standard and 20.0 ml of HCl. Bring up to volume and mix well.

AAS Standard Mix 4 (10.0ppm Cu, Fe, Mn; 3.0ppm Zn)

Add to a 1000 ml volumetric flask 500 ml of DI Water (Type 1) approximately. Pipet 10.0 ml of Iron, Manganese and Copper. Pipet 3.0 ml of Zinc of the primary standard and 20.0 ml of HCl. Bring up to volume and mix well.

AAS Standard Mix 5 (20.0ppm Cu, Fe, Mn; 5.0ppm Zn)

Add to a 1000 ml volumetric flask 500 ml of DI Water (Type 1) approximately. Pipet 20.0 ml of Iron, Manganese and Copper. Pipet 5.0 ml of Zinc of the primary standard and 20.0 ml of HCl. Bring up to volume a mix well.

#### 10.2.3. Recommended Standards Level:

Soluble and Chelate for Iron, Manganese and Copper

| Standard Identification | Concentration |
|-------------------------|---------------|
| Standard 1              | 0.5 ppm       |
| Standard 2              | 5.0 ppm       |
| Standard 3              | 8.0 ppm       |
| Standard 4              | 10.0 ppm      |
| Standard 5              | 20.0 ppm      |

Soluble and Chelate Zinc

| Standard Identification | Concentration |
|-------------------------|---------------|
| Standard 1              | 0.5 ppm       |
| Standard 2              | 1.0 ppm       |
| Standard 3              | 2.0 ppm       |
| Standard 4              | 3.0 ppm       |
| Standard 5              | 5.0 ppm       |

#### 10.3. Instrument Startup

10.3.1. See LABOP 140 Maintenance, Operation and Program of AAS, Attachment K for details

10.3.2. Clean Burner with non-residual soap and water, do not scrape the slot flame burner exit, and dry. After inspecting, if there is still any residue, sonicate the burner for 20 minutes in a blend of methanol and acetone 50\50.

10.3.3. Set and align the burner into AAS analyzer

10.3.4. Turn on the AAS analyzer, pump and autosample

10.3.5. Turn on Fe, Mn, Cu and Zn multi lamp for 20 minutes

10.3.6. AAS spectrometer parameters used in soluble and chelate metal analysis:

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 523

*Analysis of Soluble and Chelate Metals by AA- [FM830, FM831, FM832, FM833, FM840, FM841]*

Version: 1.0

Effective Date: 05/10/2023

Page 4 of 6

| AAS (AA-01)<br>Perkin-Elmer<br>AAAnalyst 100 | Burner<br>Head<br>Alignment | Wavelength<br>(nm) | Slit<br>Width<br>(nm) | Acetylene<br>Flow<br>(L\Min) | Oxidant<br>Flow<br>(L\Min) |
|----------------------------------------------|-----------------------------|--------------------|-----------------------|------------------------------|----------------------------|
| Iron                                         | Straight                    | 248.3              | 0.2                   | 1.7                          | 4.0                        |
| Manganese                                    | Straight                    | 279.5              | 0.2                   | 1.7                          | 4.0                        |
| Copper                                       | Straight                    | 324.8              | 0.7                   | 1.7                          | 4.0                        |
| Zinc                                         | Straight                    | 213.9              | 0.7                   | 1.7                          | 4.0                        |

| AAS (AA-02)<br>Perkin-Elmer<br>Pinnacle 500 | Burner<br>Head<br>Alignment | Wavelength<br>(nm) | Slit<br>Width<br>(nm) | Acetylene<br>Flow<br>(L\Min) | Oxidant<br>Flow<br>(L\Min) |
|---------------------------------------------|-----------------------------|--------------------|-----------------------|------------------------------|----------------------------|
| Iron                                        | Straight                    | 248.3              | 0.2                   | 10.0                         | 2.50                       |
| Manganese                                   | Straight                    | 279.5              | 0.2                   | 10.0                         | 2.50                       |
| Copper                                      | Straight                    | 324.8              | 0.7                   | 10.0                         | 2.50                       |
| Zinc                                        | Straight                    | 213.9              | 0.7                   | 10.0                         | 2.50                       |

#### 10.4. Sample Analysis

10.4.1. Use the Sample Information File (Syngitix) to create the sequence.

10.4.2. An instrument blank must be injected at the beginning of a sample sequence to act as a correction.

10.4.3. After the initial calibration, there must be a continuing calibration check standard analyzed after every 10 samples and at the end of the sequence.

10.4.4. The Laboratory Control Samples (LCS) should be analyzed with their associated samples, not grouped together at the beginning of the sequence.

10.4.5. The sequence requirements are found in, DATA 051 Attachment A, Fertilizer Analysis Requirements.

10.4.6. Samples require the proper laboratory information management system (LIMS) codes to allow required uploading. Refer to ADMIN 023 Attachment M, Instrument Sequence Generation for Fertilizer, for additional information.

10.4.6.1. Example: AA99999-0AT

10.4.6.1.1. 0 corresponds to the sample weight (0, 1, 2.)

10.4.6.1.2. A corresponds to the different injections due to dilutions (A=1st injection, B=2nd injection, C=3rd injection.)

10.4.6.1.3. T or H corresponds to analysis type (T=Totals and H=Heavy Metals)

10.4.7. Dilution factors must be entered into the instrument sequence.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 523

*Analysis of Soluble and Chelate Metals by AA- [FM830, FM831, FM832, FM833, FM840, FM841]*

Version: 1.0

Effective Date: 05/10/2023

Page 5 of 6

#### 10.5. Data Evaluation

10.5.1. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs, to evaluate the data against the established criteria. Evaluations, that supersede those of DATA 051, are addressed in this document.

10.5.2. If the sample is excessive reprep using the excess method in PREP 462.

10.5.3. If liquid sample is deficient in Zinc, Manganese, or Copper and contains ammonium polyphosphate as a phosphate source, the sample should be re-prepped using the EDTA extraction in PREP 462.

#### 10.6. Quality Control Sample

10.6.1. Refer to QUALITY 112, Quality Control, to determine the outcome of Quality Control (QC) samples. Guidelines, that supersede those of QUALITY 112, are addressed in this document.

### 11 WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material             | Disposal Mechanism                         |
|----------------------|--------------------------------------------|
| Aqueous extracts     | Disposed of down sink while flushing water |
| AAS Instrument Waste | Disposed of down sink while flushing water |

### 12 REFERENCES

- 12.1. PinAAcle 500 users guide, external document 00100
- 12.2. AA Instrument Maintenance Checklist and log, external document 00101
- 12.3. AA Fast Flame Manual, external document 00102
- 12.4. Manganese- Soluble AA. FM-830, Rev. 2.0 (01/03/2003)
- 12.5. Copper- Soluble AA. FM-831, Rev. 2.0 (01/03/2003)
- 12.6. Iron- Soluble AA. FM-832, Rev. 2.0 (01/03/2003)
- 12.7. Zinc- Soluble AA. FM-833, Rev. 2.0 (01/03/2003)
- 12.8. Iron Chelate – Pure Material- AA. FM-840 Rev. 2.0 (12/23/02)
- 12.9. Chelated Secondary/Micronutrient Elements – AA. FM-841 Rev 1.0 (11/12/2001)

### 13 VERSION HISTORY

| Version | Date       | Description      | Author      | Editor |
|---------|------------|------------------|-------------|--------|
| 1.0     | 01.26.2022 | Original Version | Jack Andreu |        |

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 523**

*Analysis of Soluble and Chelate Metals by AA- [FM830, FM831, FM832, FM833, FM840, FM841]*

Version: 1.0

Effective Date: 05/10/2023

Page 6 of 6

**14 SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director  
Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer  
Date\_\_\_\_\_

\_\_\_\_\_  
Zachary Tower  
Hemp Analysis Manager  
Date\_\_\_\_\_

\_\_\_\_\_  
Jack Andreu  
Author  
Date\_\_\_\_\_

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 530

#### *Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FM415] [FM420] [FM435]*

Version: 1.0

Effective Date: 02/10/2020




Page 1 of 9

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the determination of ammoniacal nitrogen in fertilizer by direct ammonium distillation.

## 2. SCOPE

2.1. The procedure delineated below is applicable to the direct determination of ammoniacal nitrogen in fertilizer samples. Fertilizer samples are alkalinized converting ammonium ( $\text{NH}_4^+$ ) to the volatile ammonia ( $\text{NH}_3$ ). The alkalinization is followed by direct steam distillation into a pH-4.65 solution of 4% boric acid. The boric acid solution is then titrated back to pH 4.65 with sulfuric acid.



## 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

## 4. OUTLINE of PROCEDURE

- 4.1. System Startup
- 4.2. Ammoniacal Nitrogen Determination (FM420)
- 4.3. Ammoniacal Nitrogen Determination with Urea Present (FM415)
- 4.4. Ammoniacal Nitrogen Determination with Urea and Water-Insoluble Nitrogen Present (FM435)
- 4.5. Data Evaluation
- 4.6. Quality Control Sample

## 5. INTERFERENCE

5.1. Fertilizer samples containing urea can lead to false high results for ammoniacal nitrogen. The urea may breakdown to carbon dioxide and ammonia under alkaline conditions during steam distillation. This can be avoided by substituting the sodium hydroxide with an excess of magnesium oxide ( $\text{MgO}$ ) to make the sample mildly alkaline.

5.2. Fertilizer samples containing water-insoluble nitrogen can lead to false high results for ammoniacal nitrogen. This can be avoided by filtering out the water-insoluble nitrogen after performing a water extraction.

## 6. APPARATUS and MATERIAL

- 6.1. Buchi KjelMaster K-375 with pH electrode
- 6.2. 300mL sample tube
- 6.3. 500mL sample tube
- 6.4. Glass funnel
- 6.5. 50mL centrifuge tube
- 6.6. Watman No. 42 Filter Paper, or equivalent
- 6.7. Geno/Grinder
- 6.8. Analytical balance

## 7. REAGENTS and MEDIA

- 7.1. Boric Acid ( $\text{H}_3\text{BO}_3$ ), 4% w/v (Purchased)
- 7.2. Boric Acid ( $\text{H}_3\text{BO}_3$ ), solid

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 530**

*Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FM415]  
[FM420] [FM435]*

Version: 1.0

Effective Date: 02/10/2020

Page 2 of 9

- 7.3. Magnesium Oxide (MgO), solid, ACS grade
- 7.4. Sodium Hydroxide Solution (NaOH), 50% w/w for Kjeldahl Nitrogen Determination (Purchased)
- 7.5. Sulfuric Acid (H<sub>2</sub>SO<sub>4</sub>), concentrated, ACS grade
- 7.6. Sulfuric Acid, 0.5N certified
- 7.7. Water, Deionized (DI) Type I

## **8. REAGENT PREPARATION GUIDELINES**

- 8.1. Reagent 01: 4% Boric Acid, pH 4.65
  - 8.1.1. Directions: Dissolve 40g of boric acid in DI water to a final volume of 1000mL. (4% Boric Acid may also be purchased.)
  - 8.1.2. Adjust the pH of the 4% boric acid while stirring with ~10% H<sub>2</sub>SO<sub>4</sub> to a pH of 4.65 (use ~10% NaOH to adjust the pH). pH may be measured by selecting "System Preparation," select "Measuring pH or mV."
- 8.2. Reagent 02: ~10% H<sub>2</sub>SO<sub>4</sub>
  - 8.2.1. Directions: Add approximately 5mL conc. H<sub>2</sub>SO<sub>4</sub> to approximately 35mL of DI water. Bring to 50mL final volume.
- 8.3. Reagent 03: ~10% NaOH
  - 8.3.1. Directions: Add approximately 2.5mL 50% NaOH to approximately 35mL of DI water. Bring to 50mL final volume.
- 8.4. Reagent 04: 32% NaOH
  - 8.4.1. Directions: Add 1280mL of 50% NaOH to approximately 400mL of DI water. Mix and let cool. Bring to 2L final volume.
- 8.5. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015
- 8.6. The specific ingredient amounts can vary to accommodate the batch size.
- 8.7. Expiration Date: 12 months from date of preparation unless otherwise specified.
- 8.8. Store at room temperature

## **9. ASSOCIATED DOCUMENTS**

- 9.1. Form F00015, Reagent Preparation
- 9.2. Label L003, Reagent

## **10. SPECIFIC PROCEDURES**

- 10.1. System Startup
  - 10.1.1. Start chiller and allow to circulate for approximately 15 minutes prior to analysis.
  - 10.1.2. Check waste reservoir level. Empty the waste into the Kjeldahl waste container for disposal.
  - 10.1.3. Ensure the water reservoir, 4% Boric Acid reservoir, and H<sub>2</sub>SO<sub>4</sub> titrant bottle have adequate volume for analysis.
  - 10.1.4. Place instrument in Ready mode by selecting "Ready."
  - 10.1.5. Under "System Preparation," select "Preheating."
  - 10.1.6. pH Probe Calibration
    - 10.1.6.1. Select "System Preparation," select "Calibration pH electrode."
    - 10.1.6.2. Calibrate the pH probe using pH 4 buffer and pH 7 buffer.
    - 10.1.6.3. Acceptable slope is 95 to 105.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 530**

*Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FM415]  
[FM420] [FM435]*

Version: 1.0

Effective Date: 02/10/2020

Page 3 of 9

10.1.7. Priming

10.1.7.1. H<sub>2</sub>SO<sub>4</sub> Buret

10.1.7.1.1. Remove the receiving vessel and place a beaker under the titrant line.

10.1.7.1.2. Under "System Preparation," select "Buret functions" and "Start."

10.1.7.1.3. Repeat to ensure that all bubbles are removed from the buret and tubing.

10.1.7.1.4. Replace receiving vessel.

10.1.7.2. Water, 4% Boric Acid, NaOH Priming

10.1.7.2.1. Under "System Preparation," select "Priming"

**Note:** pH probe must be in the receiver while priming.

10.2. Ammoniacal Nitrogen Determination (FM420)

10.2.1. Blank

10.2.1.1. Add approximately 20mL of DI water to a 300mL sample tube.

10.2.1.2. Place the 300mL sample tube on the autosampler rack.

10.2.2. Lab Control Sample (LCS)

10.2.2.1. Weigh 0.2g of ammonium sulfate LCS into weighing vessel.

10.2.2.2. Transfer sample from the weighing vessel into a 300mL sample tube by rinsing the weighing vessel with DI water a minimum of 3 times.

10.2.2.3. Place the 300mL sample tube on the autosampler rack.

10.2.3. Samples

10.2.3.1. Weigh samples into weighing vessel using Table 1 as a guide.

| Table 1        |               |
|----------------|---------------|
| Ammoniacal N % | Sample Weight |
| <5%            | 1g            |
| >5%            | 0.2g          |

10.2.3.1. Transfer sample from the weighing vessel into a 300mL sample tube by rinsing the weighing vessel with DI water a minimum of 3 times.

10.2.3.2. Place the 300mL sample tube on the autosampler rack.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 530

#### *Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FM415] [FM420] [FM435]*

Version: 1.0

Effective Date: 02/10/2020

Page 4 of 9

#### 10.2.4. Method

| Table 2                 |                                           |                             |            |
|-------------------------|-------------------------------------------|-----------------------------|------------|
| 420 Ammoniacal          |                                           |                             |            |
| H <sub>2</sub> O volume | 50 mL                                     | Distillation mode           | Fixed time |
| NaOH volume             | 30 mL                                     | Distillation time           | 240 s      |
| Reaction time           | 5 s                                       | Determination mode          | Standard   |
| Titration type          | Boric acid titration                      | Receiving solution volume   | 60 mL      |
| Titration solution      | H <sub>2</sub> SO <sub>4</sub> 0.25 mol/L | Molarity                    | 0.25       |
| Valence factor          | 2                                         | Titer                       | 1          |
| Sensor type             | Potentiometric                            | Titration mode              | Standard   |
| Measuring mode          | Endpoint pH                               | Endpoint pH                 | 4.65       |
| Titration start volume  | volume 0.000 mL                           | Titration algorithm         | Normal     |
| Aspiration sample tube  | Yes                                       | Aspiration receiving vessel | Yes        |

#### 10.2.5. Sequence

10.2.5.1. Build sample sequence by selecting “Sequences” on the main home screen.

10.2.5.2. Name sequence.

10.2.5.3. “Preheating” and “Priming” steps may be added to the sequence as needed.

10.2.5.4. Add the appropriate sample rack to the sequence.

10.2.5.4.1. Select the rack position of blank. Choose “Blank” as the sample type. Enter the blank name and select “420 Ammoniacal” as the method.

**Note:** the blank is used for blank correction and must be run before the LCS and samples.

10.2.5.4.2. Select the rack position of the LCS. Choose “Reference substance” as the sample type and select the reference substance from the list. Enter the LCS name, weight, and method.

10.2.5.4.3. Select the rack position of the sample. Choose “Sample” as the sample type. Enter the sample name, weight, and method.

10.2.5.5. Add the “Cleaning” step to the sequence.

10.2.5.6. Add the “Dose H<sub>3</sub>BO<sub>3</sub>” step, unless the pH probe can be transferred to the probe storage solution immediately after the end of the sequence.

**Note:** The pH probe must not be stored dry.

#### 10.3. Ammoniacal Nitrogen Determination with Urea Present (FM415)

##### 10.3.1. Blank

10.3.1.1. Add approximately 20mL of DI water to a 300mL sample tube.

10.3.1.2. Place the 300mL sample tube on the sample rack.

##### 10.3.2. Lab Control Sample (LCS)

10.3.2.1. Weigh 0.2g of ammonium sulfate LCS into weighing vessel.

10.3.2.2. Transfer sample from the weighing vessel into a 300mL sample tube by rinsing the weighing vessel with DI water a minimum of 3 times.

10.3.2.3. Place the 300mL sample tube on the sample rack.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 530

#### *Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FM415] [FM420] [FM435]*

Version: 1.0

Effective Date: 02/10/2020

Page 5 of 9

#### 10.3.3. Samples

10.3.3.1. Weigh samples into weighing vessel using Table 1 as a guide.

10.3.3.2. Transfer sample from the weighing vessel into a 300mL sample tube by rinsing the weighing vessel with DI water a minimum of 3 times.

10.3.3.3. Place the 300mL sample tube on the sample rack.

#### 10.3.4. Method

| Table 3                      |                                           |                             |                       |
|------------------------------|-------------------------------------------|-----------------------------|-----------------------|
| 415-435 Ammoniacal with Urea |                                           |                             |                       |
| H <sub>2</sub> O volume      | 50 mL                                     | Distillation mode           | Automatic-IntelliDist |
| NaOH volume                  | 0 mL                                      | Distillation time           | 240 s                 |
| Reaction time                | 5 s                                       | Determination mode          | Standard              |
| Titration type               | Boric acid titration                      | Receiving solution volume   | 60 mL                 |
| Titration solution           | H <sub>2</sub> SO <sub>4</sub> 0.25 mol/L | Molarity                    | 0.250                 |
| Valence factor               | 2                                         | Titer                       | 1.0000                |
| Sensor type                  | Potentiometric                            | Titration mode              | Standard              |
| Measuring mode               | Endpoint pH                               | Endpoint pH                 | 4.65                  |
| Titration start volume       | volume 0.000 mL                           | Titration algorithm         | Normal                |
| Aspiration sample tube       | No                                        | Aspiration receiving vessel | Yes                   |

#### 10.3.5. Sample List

10.3.5.1. Build sample list by selecting "Sample List" on the main home screen.

10.3.5.2. Name Sample List.

10.3.5.2.1. Choose "Blank" as the sample type. Enter the blank name and select "415-435 Ammoniacal with Urea" as the method.

**Note:** the blank is used for blank correction and must be run before the LCS and samples.

10.3.5.2.2. Choose "Reference substance" as the sample type and select the reference substance from the list. Enter the LCS name, weight, and method.

10.3.5.2.3. Choose "Sample" as the sample type. Enter the sample name, weight, and method.

10.3.6. Add approximately 4g MgO (~1 scoop) directly to the 300mL sample tube. Immediately place the sample tube on distillation unit.

**Caution:** Sample tube will be very hot after distillation. Use care when removing the sample tube, allow distillation line to completely drain into sample tube by gently tapping the distillation line.

**Note:** Samples using MgO must not be aspirated.

#### 10.3.7. Shutdown

10.3.7.1. Transfer the pH probe to the probe storage solution immediately after the end of the sequence.

**Note:** The pH probe must not be stored dry.

10.3.7.2. Under "System Preparation," select "Cleaning."

#### 10.4. Ammoniacal Nitrogen Determination with Urea and Water-Insoluble Nitrogen Present (FM435)

##### 10.4.1. Blank

10.4.1.1. Add approximately 40mL of DI water to a 300mL sample tube.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 530

#### *Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FM415] [FM420] [FM435]*

Version: 1.0

Effective Date: 02/10/2020

Page 6 of 9

10.4.1.2. Place the 300mL sample tube on the sample rack.

#### 10.4.2. Lab Control Sample (LCS)

10.4.2.1. Weigh 0.2g of SRM 695 LCS into 50mL centrifuge tube.

10.4.2.2. Add approximately 40mL of DI water.

10.4.2.3. Agitate using Geno/Grinder for 1 minute at 1000rpm.

10.4.2.4. Filter sample into the 300mL sample tube using a folded filter placed in a glass funnel on the top of the 300mL sample tube. Wet with filter with DI water in insure a seal before adding the sample.

10.4.2.5. If sample is clumping, use 2 to 5mL of methanol to break up clumping.

10.4.2.6. Rinse the 50mL centrifuge tube a minimum of 3 times into the filter.

10.4.2.7. Remove filter paper and funnel from the 300mL sample tube.

#### 10.4.3. Samples

10.4.3.1. Weigh samples into 50mL centrifuge tube using Table 1 as a guide.

10.4.3.2. Add approximately 40mL of DI water.

10.4.3.3. Agitate using Geno/Grinder for 1 minute at 1000rpm.

10.4.3.4. Filter sample into the 300mL sample tube using a folded filter placed in a glass funnel on the top of the 300mL sample tube. Wet with filter with DI water in insure a seal before adding the sample.

10.4.3.5. If sample is clumping, use 2 to 5mL of methanol to break up clumping.

10.4.3.6. Rinse the 50mL centrifuge tube a minimum of 3 times into the filter.

10.4.3.7. Remove filter paper and funnel from the 300mL sample tube.

#### 10.4.4. Method

| <b>Table 3</b>                      |                                           |                             |                       |
|-------------------------------------|-------------------------------------------|-----------------------------|-----------------------|
| <b>415-435 Ammoniacal with Urea</b> |                                           |                             |                       |
| H <sub>2</sub> O volume             | 50 mL                                     | Distillation mode           | Automatic-IntelliDist |
| NaOH volume                         | 0 mL                                      | Distillation time           | 240 s                 |
| Reaction time                       | 5 s                                       | Determination mode          | Standard              |
| Titration type                      | Boric acid titration                      | Receiving solution volume   | 60 mL                 |
| Titration solution                  | H <sub>2</sub> SO <sub>4</sub> 0.25 mol/L | Molarity                    | 0.250                 |
| Valence factor                      | 2                                         | Titer                       | 1.0000                |
| Sensor type                         | Potentiometric                            | Titration mode              | Standard              |
| Measuring mode                      | Endpoint pH                               | Endpoint pH                 | 4.65                  |
| Titration start volume              | volume 0.000 mL                           | Titration algorithm         | Normal                |
| Aspiration sample tube              | No                                        | Aspiration receiving vessel | Yes                   |

#### 10.4.5. Sample List

10.4.5.1. Build sample list by selecting "Sample List" on the main home screen.

10.4.5.2. Name Sample List.

10.4.5.2.1. Choose "Blank" as the sample type. Enter the blank name and select "415-435 Ammoniacal with Urea" as the method.

**Note:** the blank is used for blank correction and must be run before the LCS and samples.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 530

#### *Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FM415] [FM420] [FM435]*

Version: 1.0

Effective Date: 02/10/2020

Page 7 of 9

10.4.5.2.2. Choose "Reference substance" as the sample type and select the reference substance from the list. Enter the LCS name, weight, and method.

10.4.5.2.3. Choose "Sample" as the sample type. Enter the sample name, weight, and method.

10.4.6. Add approximately 4g MgO (~1 scoop) directly to the 300mL sample tube. Immediately place the sample tube on distillation unit.

**Caution:** Sample tube will be very hot after distillation. Use care when removing the sample tube, allow distillation line to completely drain into sample tube by gently tapping the distillation line.

**Note:** Samples using MgO must not be aspirated.

10.4.7. Shutdown

10.4.7.1. Transfer the pH probe to the probe storage solution immediately after the end of the sequence.

**Note:** The pH probe must not be stored dry.

10.4.7.2. Under "System Preparation," select "Cleaning."

10.5. Data Evaluation

10.5.1. Calculation

$$\frac{[\text{mL of Titrant} - \text{mL of Titrant Blank}][0.25\text{mol of } H_2SO_4 \text{ Titrant}][2\text{mol } H_3BO_3][1\text{mol N}][14.0067\text{g N}][1\text{L } H_2SO_4]}{[\text{g sample weight}][1\text{ L of } H_2SO_4 \text{ Titrant}][1\text{mol } H_2SO_4][1\text{mol } H_3BO_3][1\text{mol N}][1000\text{ mL } H_2SO_4]} \times 100 = \%N$$

$$\frac{[\text{mL of Titrant} - \text{mL of Titrant for Blank}][0.5][1.40067]}{[\text{g sample weight}]} = \%N$$

10.5.2. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs, to evaluate the data against the established criteria. Evaluations, that supersede those of DATA 051, are addressed in this document.

10.6. Quality Control Sample

10.6.1. Refer to QUALITY 112, Quality Control, to determine the outcome of Quality Control (QC) samples. Guidelines, that supersede those of QUALITY 112, are addressed in this document.

## 11. WASTE MANAGEMENT

11.1. Material, marked for disposal, is removed via the following mechanism:

| Material        | Disposal Mechanism                         |
|-----------------|--------------------------------------------|
| Waste Reservoir | Neutralize to 6-8 pH and dispose down sink |

## 12. REFERENCES

- 12.1. Ammoniacal Nitrogen (Urea Present) – Direct. FM 415, January 29, 2003
- 12.2. Ammoniacal Nitrogen (Urea Not Present) – Direct. FM 420, February 5, 2003
- 12.3. Breakdown Nitrogen Preparation. FM 430, November 12, 2001
- 12.4. Ammoniacal Nitrogen (Urea Present) – Direct. FM 435, February 4, 2003
- 12.5. Selective nitrogen determination methods related to Kjeldahl, Buchi (Application Note No. 252/2016) November 2016.
- 12.6. AOAC, 920.03 *Nitrogen (Ammoniacal) in Fertilizer, Magnesium Oxide Method*, 2019

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 530**

*Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FM415]  
[FM420] [FM435]*

Version: 1.0

Effective Date: 02/10/2020

Page 8 of 9

**13. VERSION HISTORY**

| <b>Version</b> | <b>Date</b> | <b>Description</b> | <b>Author</b>     | <b>Editor</b> |
|----------------|-------------|--------------------|-------------------|---------------|
| 1.0            | 02.10.2020  | Original Version   | William Meeks, Jr |               |

***STANDARD OPERATING PROCEDURE***  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

METHOD 530

*Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FM415]  
[FM420] [FM435]*

Version: 1.0

Effective Date: 02/10/2020

Page 9 of 9

**14. SIGNATURE HISTORY**

\_\_\_\_\_  
Patricia Lucas  
Chief, Technical Director

Date\_\_\_\_\_

\_\_\_\_\_  
Teresa Rygiel  
Laboratory Director

Date\_\_\_\_\_

\_\_\_\_\_  
Quentin Cooper  
Quality Assurance Officer

Date\_\_\_\_\_

\_\_\_\_\_  
Taleya Williams  
Sample Prep and Wet Chemistry Manager

Date\_\_\_\_\_

\_\_\_\_\_  
William Meeks, Jr.  
MDTA Manager / Author

Date\_\_\_\_\_

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 531**

*Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer*  
[FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 1 of 15

## 1. PURPOSE

1.1. This standard operating procedure (SOP) provides instruction for the determination of urea and other water soluble organic nitrogen (WSON), water insoluble nitrogen (WIN), and total nitrogen (TN) in fertilizer by the Kjeldahl method. This SOP also includes instructions for the indirect determination of nitrate ( $\text{NO}_3$ ).

## 2. SCOPE

2.1. The procedure delineated below is applicable to the determination of urea and other water soluble organic nitrogen, water insoluble nitrogen, total nitrogen, and nitrate in fertilizer samples. Fertilizer samples are digested with a mixture of an inert salt, a catalyst, and sulfuric acid converting nitrogen to ammonium ( $\text{NH}_4^+$ ). The ammonium is then alkalized to convert to the volatile ammonia ( $\text{NH}_3$ ). The alkalization is followed by direct steam distillation into a pH-4.65 solution of 4% boric acid. The solution is then titrated back to pH 4.65 with sulfuric acid.

Digestion –  $\text{Organic N} + \text{H}_2\text{SO}_4 \rightarrow (\text{NH}_4)_2\text{SO}_4 + \text{H}_2\text{O} + \text{CO}_2$

Distillation –  $(\text{NH}_4)_2\text{SO}_4 + 2\text{NaOH} \rightleftharpoons 2\text{NH}_3 + \text{Na}_2\text{SO}_4 + 2\text{H}_2\text{O}$

Acid Trap –  $\text{B}(\text{OH})_3 + \text{NH}_3 + \text{H}_2\text{O} \rightleftharpoons \text{NH}_4^+ + \text{B}(\text{OH})_4^-$

Titration –  $2\text{B}(\text{OH})_4^- + \text{H}_2\text{SO}_4 \rightleftharpoons \text{SO}_4^{2-} + 2\text{B}(\text{OH})_3 + 2\text{H}_2\text{O}$

2.2. Potassium sulfate is used during the digestion to increase the boiling point of sulfuric acid from approximately 335 °C to approximately 354 °C, thereby increasing the rate of the reaction. Optimal acid to salt ratio is 2:1 (2 mL  $\text{H}_2\text{SO}_4$  to 1 g of potassium sulfate). A catalyst such as copper is utilized to increase the efficiency and speed of the digestion. Mercury is not used as a catalyst due to its toxicity.

2.3. For urea determination, AOAC Method 955.04 was used as a basis for this procedure and was modified for the use of an automated Kjeldahl system utilizing steam distillation and a boric acid titration. Copper was substituted for mercury as the catalyst, and the acid to salt ratio was optimized. AOAC Method 930.1 and 930.2 were incorporated to liberate nitrate interference by the addition of ferrous sulfate to the digestion.

2.4. Since Kjeldahl nitrogen digestions include Ammoniacal nitrogen, the Ammoniacal nitrogen results must be subtracted to determine urea results.

2.5. Samples containing WIN are extracted with water and filtered prior to the digestion step. AOAC Method 945.01 was used as a basis for this procedure. Mechanical agitation and centrifugation were adapted. A finer grade of filter paper is also utilized (Whatman No. 42 was substituted for Whatman No. 2). The aqueous extract is then digested as above for urea analysis.

2.6. For WIN determination, the filter with the remaining insoluble residue is digested. AOAC Method 945.01 was used as a basis for this procedure and was modified for the use of an automated Kjeldahl system utilizing steam distillation and a boric acid titration. Copper was substituted for mercury as the catalyst, and the acid to salt ratio was optimized. Hydrogen peroxide was incorporated during the early stage of digestion to mitigate excess foam.

2.7. For total nitrogen determination, salicylic acid and sulfuric acid are used to convert nitrate to nitrosalicylic acid. Sodium thiosulfate is then used to reduce the nitrosalicylic acid to its amino form. AOAC Method 955.04 D. was used as a basis for this procedure and was modified for the use of an automated Kjeldahl system utilizing steam distillation and a boric acid titration. Copper was substituted for mercury as the catalyst, and the acid to salt ratio was optimized.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES


### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 2 of 15



2.8. Nitrate can be calculated from subtracting the total Kjeldahl nitrogen (TKN) from the Kjeldahl total nitrogen (TN). TKN consists of the sum of nitrogen from organic nitrogen and ammoniacal forms.

2.9. Nitrogen sources not detectable with Kjeldahl include azide, azine, azo, hydrazine, nitrile, nitro, nitroso, oxime, and semi-carbazone.

### 3. HEALTH and SAFETY

3.1. A reference file of Safety Data Sheet (SDS) is available to all personnel involved in the chemical analysis. All essential personal protective equipment (PPE) including safety glasses, gloves and laboratory coat is required.

3.2. There is a risk of acid fume inhalation during digestion. The use of a scrubber during digestion and cooling will assist in removing acid fumes. A fume hood must be used during the operation of the scrubber and digester.

### 4. OUTLINE of PROCEDURE

#### 10.1 Sample Preparation

10.1.1 Preparation for Urea (FM425)

10.1.2 Preparation for Urea (FM440) containing WIN (FM445)

10.1.3 Preparation for Total Nitrogen (FM405)

#### 10.2 Digestion

10.2.1 Digestion of Urea (FM425/FM440)

10.2.2 Digestion of WIN (FM445)

10.2.3 Digestion of Total Nitrogen (FM405)

#### 10.3 Distillation, and Titration

10.3.1 System Startup

10.3.2 Sequence

10.3.3 Distillation and Titration of Urea (FM425/440)

10.3.4 Distillation and Titration of WIN (FM445)

10.3.5 Distillation and Titration of Total Nitrogen (FM405)

#### 10.4 Post-run Cleaning

#### 10.5 Data Export and Upload

#### 10.6 Data Evaluation

#### 10.7 Quality Control Sample

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 531**

*Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer*  
[FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025


Page 3 of 15

## 5. INTERFERENCE

5.1. A digestion temperature must not exceed 380 °C. Temperatures above 390 °C lead to the possibility of nitrogen gas (N<sub>2</sub>) formation resulting in low nitrogen results. The ratio of salt to sulfuric acid is critical in maintaining optimal temperature.

### 5.2. Interference for Urea Determination

5.2.1. Nitrate in fertilizer samples can cause low N recovery. Volatile N<sub>2</sub>O can be produced when nitrates oxidize a portion of the ammonia released from the digestion. This can be avoided by the addition of ferrous sulfate with sulfuric acid to reduce nitrates present in the sample to nitric oxide (NO<sup>-</sup>) gas.



5.2.2. WIN in fertilizer samples yields false high N results and must be filtered out for accurate determination of urea.

5.2.3. During digestion, black foam might form in the digestion tube. To achieve optimal reflux, efforts shall be made to keep the foam out of the condensation zone (5 cm below the constriction of the sample tube). This can be achieved by using antifoaming tablets and/or adjusting the temperature program to hold at a lower temperature, until the foaming has calmed, before continuing with the digestion ramp.

### 5.3. Interference for WIN Determination

5.3.1. It was found that, during sample digestion, the combination of water and a cellulose filter pad yields an excess of black foam that can rise into the suction module. To achieve optimal reflux, efforts shall be made to keep the foam out of the condensation zone (5 cm below the constriction of the sample tube). The addition of 30% H<sub>2</sub>O<sub>2</sub> early in the digestion was found to mitigate the foam.

### 5.4. Interference for Total Nitrogen Determination

5.4.1. The reaction of nitrate to nitrosalicylic acid is inhibited by the presence of water. Liquid fertilizers may be analyzed up to 1 g, for reproducible results.

## 6. APPARATUS and MATERIAL

- 6.1. Buchi KjelDigester K-449
- 6.2. Buchi Scrubber K-415 Triple Scrub
- 6.3. Buchi KjelMaster K-375 with pH electrode
- 6.4. Buchi KjelSampler K-376
- 6.5. 300-mL sample tube
- 6.6. 500-mL sample tube
- 6.7. Glass funnel
- 6.8. 50-mL centrifuge tube
- 6.9. Geno/Grinder
- 6.10. Glass rods
- 6.11. Analytical balance
- 6.12. Plastic weigh boats
- 6.13. Nitrogen-free weigh boats

## 7. REAGENTS and MEDIA

- 7.1. Antifoaming Tablets
- 7.2. Boric Acid (H<sub>3</sub>BO<sub>3</sub>), 4% w/v (Purchased)
- 7.3. Bromothymol blue
- 7.4. Ferrous Sulfate (FeSO<sub>4</sub>)

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 531**

*Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer*  
[FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 4 of 15

- 7.5. Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>), 30%, ACS grade
- 7.6. Kjeldahl Tablets, CT-50 5 g K<sub>2</sub>SO<sub>4</sub>, 0.15 g CuSO<sub>4</sub>
- 7.7. Salicylic Acid
- 7.8. Sodium carbonate (Na<sub>2</sub>CO<sub>3</sub>)
- 7.9. Sodium Hydroxide Solution (NaOH), 32% w/w for Kjeldahl Nitrogen Determination (Purchased)
- 7.10. Sodium Hydroxide Solution (NaOH), 50% w/w for Kjeldahl Nitrogen Determination (Purchased)
- 7.11. Sodium Thiosulfate
- 7.12. Sulfuric Acid (H<sub>2</sub>SO<sub>4</sub>), concentrated, ACS grade
- 7.13. Sulfuric Acid, 0.5N certified
- 7.14. Water, Deionized (DI), Type I

## 8. REAGENT PREPARATION GUIDELINES

8.1. **4% Boric Acid, pH 4.65**, Dissolve 40 g of boric acid in DI water to a final volume of 1000-mL. (4% Boric Acid may also be purchased.) Adjust the pH to 4.65 while stirring using ~10% H<sub>2</sub>SO<sub>4</sub> or ~10% NaOH. pH may be measured using the Buchi pH meter, by selecting "System Preparation," then "Measuring pH or mV."

**NOTE:** The pH of 4% Boric Acid that has been purchased must be adjusted to 4.65.

- 8.2. **~10% H<sub>2</sub>SO<sub>4</sub>**, Add approximately 5 mL conc. H<sub>2</sub>SO<sub>4</sub> to approximately 35 mL of DI water. Bring to 50-mL final volume.
- 8.3. **~10% NaOH**, Add approximately 10 mL 50% NaOH to approximately 35 mL of DI water. Bring to 50-mL final volume.
- 8.4. **32% NaOH**, Add 1280 mL of 50% NaOH to approximately 400 mL of DI water. Mix and let cool. Bring to 2-L final volume. (32% NaOH may be purchased)
- 8.5. **Scrubber Neutralization Solution**, Dissolve 600 g sodium carbonate in 3 L warm DI water. Add a spatula tip of Bromothymol blue as a color indicator.

Note: The fresh solution has a blue color while the neutralized or acidic solution has an orange yellowish color. The color transition occurs between pH 6.0 to 7.6.

- 8.6. Label the container appropriately, Label L003 indicates required information. Add relevant information to Reagent Preparation, Form F00015
- 8.7. The specific ingredient amounts can vary to accommodate the batch size.
- 8.8. Expiration Date: 12 months from date of preparation unless otherwise specified.
- 8.9. Store at room temperature

## 9. ASSOCIATED DOCUMENTS

- 9.1. LABOP 140 Attachment P, Maintenance, Operation, and Program for Büchi Kjeldahl Digestion Analyzer
- 9.2. Method 530, Automated Direct Distillation of Ammoniacal Nitrogen in Fertilizer [FN414] [FM420] [FM435]
- 9.3. External 00199 KjelDigester K-446/K-449 Operation Manual
- 9.4. External 00198 KjelMaster K-375 with KjelSampler K-376/K-377 Operation Manual

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 5 of 15

## 10. SPECIFIC PROCEDURES

### 10.1. Sample Preparation

#### 10.1.1. Preparation for Urea (FM425)

- 10.1.1.1. Weigh the sample into a plastic weigh boat, using Table 1 as a guide.
- 10.1.1.2. A reagent blank and a laboratory control sample (LCS) should also be digested and analyzed with each prep batch. (The blank is prepared with the same acid volume, water volume, and chemicals as the sample.)

| Table 1               |               |
|-----------------------|---------------|
| Urea & Ammoniacal N % | Sample Weight |
| <5%                   | 1 g           |
| >5%                   | 0.2 g         |

- 10.1.1.3. Transfer the sample into a 300-mL sample tube by rinsing the weigh boat with DI water.
- 10.1.1.4. Add approximately 2 g (~¼ teaspoon) of ferrous sulfate to the 300-mL sample tube.
- 10.1.1.5. Add one antifoaming tablet to the 300-mL sample tube.
- 10.1.1.6. Bring to a volume of approximately 100 mL with deionized water and mix.
- 10.1.1.7. Add two Kjeldahl Tablets and 20 mL sulfuric acid to the 300-mL sample tube and mix.
- 10.1.1.8. Use a glass boiling rod in each tube to prevent bumping or boil overs during digestion.
- 10.1.1.9. Cover any open positions on the digestion rack with glass caps.
- 10.1.1.10. Digest the sample as in step **10.2.1 Urea (FM425/FM440) Digestion**

#### 10.1.2. Preparation for Urea (FM440) containing WIN (FM445)

- 10.1.2.1. Weigh the sample into a 50-mL centrifuge tube using Table 2 as a guide.
- 10.1.2.2. A reagent blank and a laboratory control sample (LCS) should also be digested and analyzed with each prep batch. (The blank is prepared with the same acid volume, water volume, filter pad, and chemicals as the sample.)

| Table 2 (use the higher weight) |       |               |
|---------------------------------|-------|---------------|
| Sum of Urea & Ammoniacal N %    | WIN % | Sample Weight |
| <5%                             | <5%   | 1 g           |
| >5%                             | >5%   | 0.2 g         |

- 10.1.2.3. Add approximately 40 mL of DI water.
- 10.1.2.4. Agitate using Geno/Grinder for 1 minute at 1000rpm.
- 10.1.2.5. Centrifuge for 10 minutes at 4000 rcf.
- 10.1.2.6. Place a Whatman 42 filter in a glass funnel on the top of a 300-mL sample tube, for Urea **FM440**. Wet the filter with DI water to ensure a seal.
- 10.1.2.7. Decant the liquid portion of the sample through the filter into the 300-mL sample tube taking care not to disturb the solid sample pellet in the centrifuge tube.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 531**

*Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer*  
*[FM425] [FM440] [FM445] [FM405]*

Version: 2.0

Effective Date: 01/27/2025

Page 6 of 15

Note: A few milliliters may be left in the tube so not to transfer more solid sample to the filter than necessary.

10.1.2.8. Wash the sample pellet in the centrifuge tube 3 times following the steps below:

10.1.2.8.1. Add approximately 20 mL DI water.

10.1.2.8.2. Agitate using Geno/Grinder for 1 minute at 1000rpm.

10.1.2.8.3. Centrifuge for 10 minutes at 4000 RCF.

10.1.2.8.4. Decant the liquid portion of the sample through the filter into the 300-mL sample tube taking care not to disturb the solid sample pellet in the centrifuge tube.

10.1.2.9. Transfer the filter paper to a new 300-mL sample tube for WIN **FM445** and keep the solid sample pellet in the centrifuge tube for later transfer.

10.1.2.10. To the 300-mL sample tube containing the filtered sample **FM440**, add approximately 2 g (~¼ teaspoon) of ferrous sulfate.

10.1.2.11. Add one antifoaming tablet to the 300-mL sample tube.

10.1.2.12. Add two Kjeldahl Tablets and 20 mL sulfuric acid to the 300-mL sample tube and mix.

10.1.2.13. Use a glass boiling rod in each tube to prevent bumping or boil overs during digestion.

10.1.2.14. Cover any open positions on the digestion rack with glass caps.

10.1.2.15. Digest urea samples as in step **10.2.1 Digestion of Urea (FM425/FM440)**.

10.1.2.16. Transfer and wash the sample pellet into the 300-mL sample tube for WIN **FM445** repeating 5 times the following the steps below:

10.1.2.16.1. Add approximately 5 mL DI water to the centrifuge tube.

10.1.2.16.2. Shake or, if needed, Vortex to break up the sample pellet.

10.1.2.16.3. Transfer into the 300-mL sample tube for WIN.

10.1.2.16.4. Wash the 50-mL centrifuge tube into the 300-mL sample tube for WIN 4 times with approximately 5 mL DI water. Do not exceed ~30 mL total volume.

10.1.2.17. Add one antifoaming tablet to the 300-mL sample tube.

10.1.2.18. Add two Kjeldahl Tablets and 20 mL sulfuric acid to the 300-mL sample tube and mix by swirling.

10.1.2.19. Use a glass boiling rod in each tube to prevent bumping or boil overs during digestion.

10.1.2.20. Cover any open positions on the digestion rack with glass caps.

10.1.2.21. Digest the sample as in step **10.2.2 Digestion of WIN (FM445)**.

**10.1.3. Preparation for Total Nitrogen (FM405)**

10.1.3.1. Weigh the sample into a Nitrogen-free weigh boat, using Table 3 as a guide.

Liquid fertilizer samples are weighed directly into a 300-mL sample tube.

10.1.3.2. A reagent blank and a laboratory control sample (LCS) should also be digested and analyzed with each prep batch. (The blank is prepared with the same acid volume, water volume, weigh boat, and chemicals as the sample.)

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 7 of 15

| Table 3   |               |
|-----------|---------------|
| Total N % | Sample Weight |
| <5%       | 1 g           |
| >5%       | 0.2 g         |

- 10.1.3.1. Transfer the weigh boat containing the sample into a 300-mL sample tube.
- 10.1.3.2. Add one antifoaming tablet to the 300-mL sample tube.
- 10.1.3.3. Add one Kjeldahl Tablets to the 300-mL sample tube, and mix.
- 10.1.3.4. Weigh approximately 1 g of salicylic acid in a plastic cup. Add 20 mL sulfuric acid and mix by swirling occasionally, until completely dissolved.
- 10.1.3.5. Add the salicylic acid/sulfuric acid mixture to the 300-mL sample tube. Swirl to mix.
- 10.1.3.6. Let stand for at least 30 min at room temperature, swirling occasionally.
- 10.1.3.7. Add approximately 2.5 g sodium thiosulfate pentahydrate, and mix by swirling.
- 10.1.3.8. Let stand for 5 min., at room temperature. Swirl to mix.
- 10.1.3.9. Use a glass boiling rod in each tube to prevent bumping or boil overs during digestion.
- 10.1.3.10. Cover any open positions on the digestion rack with glass caps.
- 10.1.3.11. Digest the sample as in step **10.2.3 Digestion of Total Nitrogen (FM405)**.

## 10.2. Digestion

### 10.2.1. Digestion of Urea (FM425/FM440)

- 10.2.1.1. Ensure that the water is flowing through the condenser on the scrubber.
- 10.2.1.2. Ensure that the Scrubber Neutralization Solution is blue.
- 10.2.1.3. Digest the sample as in the table below:

**Note:** The digestion block can be preheated at the desired temperature using Method 0.

Table 4

| Urea Digestion (Method 1) |                  |            |
|---------------------------|------------------|------------|
| Step                      | Temperature (°C) | Time (min) |
| 1                         | 250              | 5          |
| 2                         | 280              | 20         |
| 3                         | 330              | 30         |
| 4                         | 420              | 65         |
| Cooling                   | -                | 35         |

10.2.1.4. **Due to the volume of water in each sample, the CONDENSATION TRAP MAY NEED TO BE EMPTIED during the digestion, depending on the sample set size.**

10.2.1.5. The digested sample is blue/green with a salt precipitate at the bottom. The salt must be re-dissolved into solution before analysis.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES


### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 8 of 15



10.2.1.6. After the sample cools, add approximately 100 mL of deionized water to the 300-mL sample tube.

10.2.1.7. Heat the sample until the salt re-dissolves, as in the table below, and allow to cool:

Table 5

| Re-dissolve (Method 2) |                  |            |
|------------------------|------------------|------------|
| Step                   | Temperature (°C) | Time (min) |
| 1                      | 280              | 0          |
| 2                      | 380              | 10         |
| Cooling                | -                | 20         |



Note: If the above re-dissolve step was not performed immediately after the digestion, precipitants may still be visually present. The sample should be immediately heated again with the above re-dissolve method.

10.2.1.8. Distill and titrate the sample as in step **10.3. Distillation and Titration**, using the appropriate instrument method.

#### 10.2.2. Digestion of WIN (FM445)

10.2.2.1. Ensure that the water is flowing through the condenser on the scrubber.

10.2.2.2. Ensure that the Scrubber Neutralization Solution is blue.

10.2.2.3. Heat samples for 5 min as in the table below:

**Note:** The digestion block can be preheated at the desired temperature using Method 0.

Table 6

| WIN Sample Heating (Method 3) |                  |            |
|-------------------------------|------------------|------------|
| Step                          | Temperature (°C) | Time (min) |
| 1                             | 300              | 5          |
| Cooling                       | -                | 20         |

10.2.2.4. **CAUTION:** Slowly add ~10 mL of 30% hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) in increments of 2.5 mL, 2.5 mL, and 5 mL, allowing ~5 seconds between dosage. (Some samples may foam as the peroxide is added.)

10.2.2.5. Digest the sample as in the table below:

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 9 of 15

Table 7

| WIN Digestion (Method 4) |                  |            |
|--------------------------|------------------|------------|
| Step                     | Temperature (°C) | Time (min) |
| 1                        | 310              | 10         |
| 2                        | 330              | 5          |
| 3                        | 420              | 70         |
| Cooling                  | -                | 35         |

10.2.2.6. The digested sample is blue/green.



10.2.2.7. If not distilling the sample immediately after digestion, add ~20 mL of water and swirl to mix, to prevent the sample from solidifying. If the sample has solidified, low heat may be applied to assistance in dissolving.

10.2.2.8. Distill and titrate the sample as in step **10.3. Distillation and Titration**, using the appropriate instrument method.

#### 10.2.3. Digestion of Total Nitrogen (FM405)

10.2.3.1. Ensure that the water is flowing through the condenser on the scrubber.

10.2.3.2. Ensure that the Scrubber Neutralization Solution is blue.

10.2.3.3. Digest the sample as in the table below:

**Note:** The digestion block can be preheated at the desired temperature using Method 0.

Table 8

| TN Digestion (Method 5) |                  |            |
|-------------------------|------------------|------------|
| Step                    | Temperature (°C) | Time (min) |
| 1                       | 300              | 0          |
| 2                       | 420              | 90         |
| Cooling                 | -                | 35         |

10.2.3.4. If the sample is not clear blue/green, digest for an additional 15 minutes at 420°C.

**STANDARD OPERATING PROCEDURE**  
**BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES**

**METHOD 531**

*Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer*  
[FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 10 of 15

10.2.3.5. If not distilling the sample immediately after digestion, add ~20 mL of water and swirl to mix, to prevent the sample from solidifying. If the sample has solidified, low heat may be applied to assistance in dissolving.

10.2.3.6. Distill and titrate the sample as in step **10.3. Distillation and Titration**, using the appropriate instrument method.

### **10.3. Distillation and Titration**

#### **10.3.1. System Startup**

10.3.1.1. Start the chiller and allow to circulate for approximately 15 minutes prior to analysis.

10.3.1.2. Check the waste reservoir level. Empty the waste into the Kjeldahl waste container, for disposal.

10.3.1.3. Ensure that the water reservoir, the 4% Boric Acid reservoir, and the H<sub>2</sub>SO<sub>4</sub> titrant bottle have adequate volume for analysis.

10.3.1.4. Place the instrument in Ready mode by selecting "Ready."

10.3.1.5. Under "System Preparation," select "Preheating."

#### **10.3.1.6. pH Probe Calibration**

10.3.1.6.1. Select "System Preparation," then "Calibration pH electrode."

10.3.1.6.2. Calibrate the pH probe using the pH-4 buffer and the pH-7 buffer.

10.3.1.6.3. Acceptable slope is 95% to 105%.

#### **10.3.1.7. Priming**

##### **10.3.1.7.1. H<sub>2</sub>SO<sub>4</sub> Buret**

10.3.1.7.1.1. Remove the receiving vessel and place a beaker under the titrant line.

10.3.1.7.1.2. Under "System Preparation," select "Buret functions," then "Start."

10.3.1.7.1.3. Repeat to ensure that all bubbles are removed from the buret and the tubing.

10.3.1.7.1.4. Replace the receiving vessel.

##### **10.3.1.7.2. Water, 4% Boric Acid, NaOH Priming**

10.3.1.7.2.1. Under "System Preparation," select "Priming"

**Note:** the pH probe must be in the receiver while priming.

#### **10.3.2. Sequence**

10.3.2.1. Build the sample sequence by selecting "Sequences" on the main home screen.

10.3.2.2. Name the sequence.

10.3.2.3. "Preheating" and "Priming" steps may be added to the sequence, as needed.

10.3.2.4. Add the appropriate sample rack to the sequence.

10.3.2.4.1. Select the rack position of the blank. Choose "Blank" as the sample type. Enter the blank name and select the appropriate method name (e.g., 425-440 Urea Kjeldahl, 445 WIN Kjeldahl, 405 TN Kjeldahl).

**Note:** the blank is used for blank correction and must be run before the LCS and samples.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 11 of 15

10.3.2.4.2. Select the rack position of the LCS. Choose "Reference substance" as the sample type, and select the reference substance from the list. Enter the LCS name, weight, and method.

10.3.2.4.3. Select the rack position of the sample. Choose "Sample" as the sample type. Enter the sample name, weight, and method.

10.3.2.5. Add the "Cleaning" step to the sequence.

10.3.2.6. Add the "Dose H<sub>3</sub>BO<sub>3</sub>" step, unless the pH probe can be transferred to the probe storage solution immediately after the end of the sequence.

**Note:** The pH probe must not be stored dry.

#### 10.3.3. Distillation and Titration of Urea (FM425, FM440)

Table 9

| 425-440 Urea Kjeldahl<br>Distillation Parameters |                                           |                             |        |
|--------------------------------------------------|-------------------------------------------|-----------------------------|--------|
| H <sub>2</sub> O volume                          | 0 mL                                      | Distillation time           | 240 s  |
| NaOH volume                                      | 90 mL                                     | Stirring Speed              | 5      |
| Reaction time                                    | 5 s                                       | Steam Output                | 100%   |
| Distillation mode                                | Fixed time                                |                             |        |
| Titration Parameters                             |                                           |                             |        |
| Titration type                                   | Boric acid titration                      | Molarity                    | 0.25   |
| Titration solution                               | H <sub>2</sub> SO <sub>4</sub> 0.25 mol/L | Titration mode              | Online |
| Sensor type                                      | Potentiometric                            | Endpoint pH                 | 4.65   |
| Measuring mode                                   | Endpoint pH                               | Stirring Speed              | 7      |
| Receiving solution volume                        | 60 mL                                     | Titration algorithm         | Normal |
| Titration start volume                           | volume 0.000 mL                           |                             |        |
| Determination Parameters                         |                                           |                             |        |
| Determination mode                               | Standard                                  |                             |        |
| Aspiration sample tube                           | Yes                                       | Aspiration receiving vessel | Yes    |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 12 of 15

#### 10.3.4. Distillation and Titration of WIN (FM445)

Table 10

| 445 WIN Kjeldahl<br>Distillation Parameters |                                           |                             |        |
|---------------------------------------------|-------------------------------------------|-----------------------------|--------|
| H <sub>2</sub> O volume                     | 60 mL                                     | Distillation time           | 240 s  |
| NaOH volume                                 | 90 mL                                     | Stirring Speed              | 5      |
| Reaction time                               | 5 s                                       | Steam Output                | 100%   |
| Distillation mode                           | Fixed time                                |                             |        |
| Titration Parameters                        |                                           |                             |        |
| Titration type                              | Boric acid titration                      | Molarity                    | 0.25   |
| Titration solution                          | H <sub>2</sub> SO <sub>4</sub> 0.25 mol/L | Titration mode              | Online |
| Sensor type                                 | Potentiometric                            | Endpoint pH                 | 4.65   |
| Measuring mode                              | Endpoint pH                               | Stirring Speed              | 7      |
| Receiving solution volume                   | 50 mL                                     | Titration algorithm         | Normal |
| Titration start volume                      | volume 0.000 mL                           |                             |        |
| Determination Parameters                    |                                           |                             |        |
| Determination mode                          | Standard                                  |                             |        |
| Aspiration sample tube                      | Yes                                       | Aspiration receiving vessel | Yes    |

#### 10.3.1. Distillation and Titration of Total Nitrogen (FM405)

Table 11

| 405 TN Kjeldahl<br>Distillation Parameters |                                           |                             |        |
|--------------------------------------------|-------------------------------------------|-----------------------------|--------|
| H <sub>2</sub> O volume                    | 60 mL                                     | Distillation time           | 240 s  |
| NaOH volume                                | 90 mL                                     | Stirring Speed              | 5      |
| Reaction time                              | 5 s                                       | Steam Output                | 100%   |
| Distillation mode                          | Fixed time                                |                             |        |
| Titration Parameters                       |                                           |                             |        |
| Titration type                             | Boric acid titration                      | Molarity                    | 0.25   |
| Titration solution                         | H <sub>2</sub> SO <sub>4</sub> 0.25 mol/L | Titration mode              | Online |
| Sensor type                                | Potentiometric                            | Endpoint pH                 | 4.65   |
| Measuring mode                             | Endpoint pH                               | Stirring Speed              | 7      |
| Receiving solution volume                  | 50 mL                                     | Titration algorithm         | Normal |
| Titration start volume                     | volume 0.000 mL                           |                             |        |
| Determination Parameters                   |                                           |                             |        |
| Determination mode                         | Standard                                  |                             |        |
| Aspiration sample tube                     | Yes                                       | Aspiration receiving vessel | Yes    |

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 13 of 15

#### 10.4. Post-run Cleaning

10.4.1. Wipe down any drips with a damp cloth.

10.4.2. Remove the suction module for cleaning. Place the suction module upside down in a sink, and flush with water. If major contaminants are visible, use a brush and detergent to clean the seals, the springs and the suction tubes. (The suction module may be cleaned in the dishwasher.)

10.4.3. Clean the drip tray with water.


10.4.4. Wash the heating block cover as needed.

10.4.5. Empty and rinse the condensation trap with water.

10.4.6. See the K-449 manual for additional cleaning instructions.

10.4.7. After the sample tubes are cleaned, inspect the glassware for scratches or chips. Discard any scratched or chipped glassware.


#### 10.5. Data Export and Upload



10.5.1. Select the samples to be exported.

10.5.2. Select Export

10.5.3. Select Print, then select report type "short"



# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]


Version: 2.0

Effective Date: 01/27/2025

Page 14 of 15

10.5.3.1. Copy the PDF report and the .kres files from the USB flash drive to P:\Pesticide Lab\BAEL\_INSTRUMENT\_DATA\AUTOMATED KJELDAHL

10.5.3.2. Open the file Buchi\_File\_Conversion.exe and follow the on-screen instructions to convert data to .csv file.



10.5.4. Upload the .csv file to the LIMS.

## 10.6. Data Evaluation

### 10.6.1. Titration Calculation

$$\frac{[\text{mL of Titrant} - \text{mL of Titrant Blank}][0.25 \text{ mol of } H_2SO_4 \text{ Titrant}][2 \text{ mol } H_3BO_3][1 \text{ mol N}][14.0067 \text{ g N}][1 \text{ L } H_2SO_4]}{[\text{g sample weight}][1 \text{ L of } H_2SO_4 \text{ Titrant}][1 \text{ mol } H_2SO_4][1 \text{ mol } H_3BO_3][1 \text{ mol N}][1000 \text{ mL } H_2SO_4]} \times 100 = \%N$$

$$\frac{[\text{mL of Titrant} - \text{mL of Titrant for Blank}][0.5][1.40067]}{[\text{g sample weight}]} = \%N$$

### 10.6.2. Urea Calculation

10.6.2.1. The Urea result can be calculated by subtracting the respective Ammoniacal (435 or 415 method) result from the respective Kjeldahl (440 or 425 method).

$$440 \%N (\text{Urea} + \text{Ammoniacal}) - 435 \%N (\text{Ammoniacal}) = \text{Urea}$$

or

$$425 \%N (\text{Urea} + \text{Ammoniacal}) - 415 \%N (\text{Ammoniacal}) = \text{Urea}$$

### 10.6.3. Nitrate Calculation

10.6.3.1. The Nitrate result can be calculated by subtracting the respective Kjeldahl (425 or 440 and 445 methods) result from the Total Kjeldahl (405 method).

$$405 \%N (\text{TN}) - 425 \%N = \text{Nitrate}$$

or

$$405 \%N (\text{TN}) - 440 \%N - 445 \%N = \text{Nitrate}$$

10.6.4. Refer to DATA 051, Data Evaluation for Feed, Fertilizer and Pesticide Residue Programs, to evaluate the data against the established criteria. Evaluations that supersede those of DATA 051 are addressed in this document.

## 10.7. Quality Control Sample

10.7.1. Refer to QUALITY 112, Quality Control, to determine the outcome of Quality Control (QC) samples. Guidelines that supersede those of QUALITY 112 are addressed in this document.

# STANDARD OPERATING PROCEDURE

## BUREAU OF AGRICULTURAL ENVIRONMENTAL LABORATORIES

### METHOD 531

#### *Automated Kjeldahl Digestion of WSON, WIN, and TN in Fertilizer* [FM425] [FM440] [FM445] [FM405]

Version: 2.0

Effective Date: 01/27/2025

Page 15 of 15

## 11. WASTE MANAGEMENT

11.1. Material marked for disposal is removed via the following mechanism:

| Material        | Disposal Mechanism                          |
|-----------------|---------------------------------------------|
| Waste Reservoir | Neutralize to 6-8 pH and dispose down sink. |
|                 |                                             |
|                 |                                             |

## 12. REFERENCES

- 12.1. AOAC, 930.01 *Nitrogen (Nitrate) in Fertilizer, Robertson Method*, 2019
- 12.2. AOAC, 930.02 *Nitrogen (Nitrate) in Fertilizer, Jones Modification of Robertson Method*, 2019
- 12.3. AOAC, 945.01 *Nitrogen (Water-Insoluble) in Fertilizer, Method I*, 2019
- 12.4. AOAC, 955.04 *Nitrogen (Total) in Fertilizer, Kjeldahl Method*, 2019
- 12.5. Buchi (2016). *Selective nitrogen determination methods related to Kjeldahl* (Application Note No. 252/2016). Retrieved from Buchi website:  
[https://static1.buchi.com/sites/default/files/2016\\_252\\_AN\\_Selective\\_nitrogen\\_determination\\_0.pdf](https://static1.buchi.com/sites/default/files/2016_252_AN_Selective_nitrogen_determination_0.pdf)
- 12.6. FDACS (2001). *Water Insoluble Nitrogen – Kjeldahl*. FM 445, November 1, 2001.
- 12.7. FDACS (2002). *Total Nitrogen – Kjeldahl*. FM 405, August 16, 2002.
- 12.8. FDACS (2003). *Other Water Soluble Nitrogen B Kjeldahl – Direct*. FM 440, February 4, 2003.
- 12.9. FDACS (2003). *Other Water Soluble Organic Nitrogen – Direct*. FM 425, February 5, 2003.
- 12.10. Labconco (1998). *A Guide to Kjeldahl Nitrogen Determination Methods and Apparatus* (2-41-9/98-JCON-5M-R3). Retrieved from Labconco website:  
<https://www.labconco.com/resources/brochures>

## 13. VERSION HISTORY

| Version | Date       | Description                                                                         | Author                                       | Editor           |
|---------|------------|-------------------------------------------------------------------------------------|----------------------------------------------|------------------|
| 1.0     | 01.16.2025 | Original Version                                                                    | William Meeks, Jr.,<br>Patrizia Lemma, Ph.D. |                  |
| 2.0     | 01.27.2025 | On item 8.3, the volume of 2.5 mL of sodium hydroxide solution was changed to 10 mL |                                              | William Meeks Jr |